

**Enzymatic Assay of CHONDROITINASE ABC  
(EC 4.2.2.4)**

**PRINCIPLE:**

Chondroitin Sulfate + H<sub>2</sub>O  $\xrightarrow{\text{Chondroitinase ABC}}$  Unsaturated Disaccharides

**CONDITIONS:** T = 37°C, pH = 8.0, A<sub>232nm</sub>, Light path = 1 cm

**METHOD:** Spectrophotometric Stop Rate Determination

**REAGENTS:**

- A. 250 mM Tris HCl and 300 mM Sodium Acetate Buffer with 0.05% (w/v) Bovine Serum Albumin, pH 8.0 at 37°C  
(Prepare 200 ml in deionized water using Trizma Base, Sigma Prod. No. T-1503, Sodium Acetate, Trihydrate, Sigma Prod. No. S-8625, and Albumin, Bovine, Sigma Prod. No. A-4503. Adjust to pH 8.0 at 37°C with 1 M HCl or 1 M NaOH.)
- B. 0.5% (w/v) Chondroitin Sulfate A Solution (Chon A)  
(Prepare 5 ml in Reagent A using Chondroitin Sulfate A, Sodium Salt, Sigma Prod. No. C-7571.)
- C. 0.5% (w/v) Chondroitin Sulfate B Solution (Chon B)  
(Prepare 5 ml in Reagent A using Chondroitin Sulfate B, Sodium Salt, Sigma Prod. No. C-2413.)
- D. 0.5% (w/v) Chondroitin Sulfate C Solution (Chon C)  
(Prepare 5 ml in Reagent A using Chondroitin Sulfate C, Sodium Salt, Sigma Prod. No. C-4384.)
- E. 0.01% (w/v) Bovine Serum Albumin Solution (BSA)  
(Prepare 20 ml in deionized water using Albumin Bovine Serum, Sigma Prod. No. A-4503.)
- F. 50 mM Potassium Chloride Solution, pH 1.8 at 25°C (KCl)  
(Prepare 50 ml in deionized water using Potassium Chloride, Sigma Prod. No. P-4504. Adjust to pH 1.8 at 25°C with 1 M HCl.)

**Enzymatic Assay of CHONDROITINASE ABC  
(EC 4.2.2.4)**

**REAGENTS:** (continued)

G. Chondroitinase ABC Enzyme Solution  
(Immediately before use, prepare a solution containing  
0.06 - 0.10 unit/ml of Chondroitinase ABC in cold  
Reagent E.)

**PROCEDURE:**

Pipette (in milliliters) the following reagents into  
suitable tubes:

	<u>Test 1</u>	<u>Test 2</u>	<u>Test 3</u>
Reagent G (Enzy Soln)	0.80	0.80	0.80

Equilibrate to 37°C for 2 minutes. Then add:

Reagent B (Chon A)	0.20	----	----
Reagent C (Chon B)	----	0.20	----
Reagent D (Chon C)	----	----	0.20

Mix by inversion. Incubate at 37°C for 21 minutes. At 0,  
3, 6, 9, 12, 15, 18, and 21 minutes ( $T_T$ ), transfer 0.10 ml  
from each Test to separate tubes containing 0.90 ml of  
Reagent F (KCl). The Blank for this assay is  $T_0$

Incubate each tube for an additional 10 minutes at 37°C.  
Centrifuge for 10 minutes and transfer the Test and Blank  
supernatants to suitable quartz cuvettes. Record the  
absorbance at 232 nm for each tube.

**CALCULATIONS:**

$T_T$  = Time intervals (minutes) for enzyme assay incubation

$$r A_{232nm} \text{ Test } (T_T) = A_{232nm} \text{ Test } (T_T) - A_{232nm} \text{ Blank } (T_0)$$

Plot the  $r A_{232nm} \text{ Test } (T_T)$  vs. Time. Determine the  $r A_{232nm}/\text{min}$   
Test from a linear portion of the graph.

**Enzymatic Assay of CHONDROITINASE ABC  
(EC 4.2.2.4)**

**CALCULATIONS:**

$$\text{Units/ml} = \frac{(\text{r } A_{232\text{nm}}/\text{min Test}) (1.0) (\text{df})}{(\text{E}_{\text{mM}}) (0.10) (0.8)}$$

df = Dilution factor

E<sub>mM</sub> = Millimolar extinction coefficient of unsaturated disaccharides: 5.1 for products from Chondroitin A and Chondroitin B, and 5.5 for products from Chondroitin C.

0.10 = Volume (in milliliter) of reaction mix used

T = Time (in minutes) of the assay as per the Unit Definition

1.0 = Total volume (in milliliter) of assay

0.8 = Volume (in milliliter) of enzyme used

$$\text{Unit/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

**UNIT DEFINITION:**

One unit will liberate 1.0 μmole of 2-acetamido-2-deoxy-3-O-(β-D-gluc-4-ene-pyranosyluronic acid)-4-O-sulfo-D-galactose from chondroitin sulfate A or 1.0 μmole of 2-acetamido-2-deoxy-3-O-(β-D-gluc-4-ene-pyranosyluronic acid)-6-O-sulfo-D-galactose from chondroitin sulfate C per minute at pH 8.0 at 37°C.

**FINAL ASSAY CONCENTRATION:**

In a 1.00 ml reaction mix, the final concentrations are 50 mM Tris, 60 mM sodium acetate, 0.02% (w/v) bovine serum albumin, 0.1% (w/v) chondroitin sulfate A, B or C, and 0.048 - 0.08 unit chondroitinase ABC.

**REFERENCE:**

Saito, H., Yamagata, T. and Suzuki, S. (1968) *J. Biol. Chem.* **243**, 1536-1542

Yamagata, T., Saito, H., Habuchi, O., and Suzuki, S. (1968) *J. Biol. Chem.* **243**, 1523-1535

**NOTES:**

1. This assay is based on the cited references.

**Enzymatic Assay of CHONDROITINASE ABC  
(EC 4.2.2.4)**

**NOTES:** (continued)

2. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

**This procedure is for informational purposes. For a current copy of Sigma's quality control procedure contact our Technical Service Department.**