

Acetic Anhydride

Product Specification

Acetic anhydride is an acylation reagent. Acylation reduces the polarity of amino, hydroxy, and thiol groups. Acylation may improve the stability of a compound by protecting unstable groups, and may increase volatility. Acetic anhydride can be used with a basic catalyst, such as pyridine. This combination promotes smooth reactions and has great solvent power. Pyridine acts as an acceptor for the acid by-product formed in the reaction. Pyridine may also react with acetic anhydride, however, forming N-acetyl-1,2-dihydro-2-pyridylacetic acid.

Features/Benefits

Acylation is an alternative to silylation, producing stable, volatile derivatives of alcohols, phenols, and amines for analysis by GC/FID. Acetylated compounds are more stable than corresponding silylated compounds.

Typical Procedures

These procedures are intended to be guidelines and may be adapted as necessary to meet the needs of a specific application. Always take proper safety precautions when using an acetylating reagent – consult MSDS for specific handling information.

Prepare a reagent blank (all components, solvents, etc., *except sample*), following the same procedure as used for the sample.

General Procedure

1. Dissolve 5mg sample in 5mL chloroform.
2. Add 0.5mL acetic anhydride and 1 mL acetic acid. Heat at 50°C for 2-16 hours.
3. Remove excess reagent by evaporating the mixture to dryness and redissolve the residue in chloroform for analysis by GC.


Alditol Formation

1. Dissolve 5mg sample in 5mL chloroform.
2. Add 1mL acetic anhydride:pyridine, 1:1. Heat at 100°C for 20 minutes.
3. Remove excess reagent by evaporating the mixture to dryness and redissolve the residue in ethyl acetate for analysis by GC.

Derivatization times vary widely, depending upon the specific compound(s) being derivatized. To determine when derivatization is complete, analyze aliquots of the sample at selected time intervals until no further increase in product peak(s) is observed. If derivatization is not complete under the procedure described here, the addition of a catalyst, use of another solvent, higher reaction temperature, longer reaction time, and/or higher reagent concentration should be evaluated.

Properties

Structure:



CAS Number:

108-24-7

Molecular Formula:

(CH₃CO)₂O

Formula Weight: 102.09


d: 1.080-1.085

n_D: 1.3901

Appearance:

Appearance: clear, colorless liquid

797-0243



Adapted from Blau and Halket (1).

Mechanism (1,2)

Acylation involves the introduction of an acyl group into a molecule that has a replaceable hydrogen atom (OH, NH, or SH group). Anhydride acylating reagents form acidic byproducts that must be removed prior to GC analysis, to prevent destructive effects on the phase in the column. Consequently, acylations with anhydride reagents normally are performed in pyridine, tetrahydrofuran, or another solvent capable of accepting the acid byproducts.

Toxicity – Hazards – Storage – Stability

Acetic anhydride is corrosive and a lachrymator. It can react vigorously with oxidizing materials, and will react violently on contact with water or steam. Store in a bottle or ampul at room temperature, in a dry, well ventilated area. Use only in a well ventilated area.

Recommended storage conditions for the unopened product are stated on the label. If you store an opened container or transfer the contents to another container for later reuse, validate that your storage conditions adequately protected the reagent.

Ordering Information

Description	Cat. No.
Acetic Anhydride	
10 x 2mL	33085
Microreaction Vessels with Hole Caps and Septa	
1mL, pk. of 12	33293
3mL, pk. of 12	33297
5mL, pk. of 12	33299
Books	
<i>Handbook of Derivatives for Chromatography</i> K. Blau and J. Halket	Z246220
<i>Handbook of Analytical Derivatization Reactions</i> D.R. Knapp	23561

References

1. K. Blau and J. Halket *Handbook of Derivatives for Chromatography* (2nd ed.) John Wiley & Sons, New York, 1993.
2. D.R. Knapp *Handbook of Analytical Derivatization Reactions* John Wiley & Sons, New York, 1979.

Additional Reading

C.J.W. Brooks and E.C. Horning, *Anal. Chem.* 36, 1540 (1964).
E.C. Horning, *et al.*, *Anal. Chem.* 36, 1546 (1964).

Contact our Technical Service Department (phone 800-359-3041 or 814-359-3041, FAX 800-359-3044 or 814-359-5468) for expert answers to your questions.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA • Sigma-Aldrich de Argentina, S.A. • Buenos Aires 1119 AUSTRALIA • Sigma-Aldrich Pty. Ltd. • Castle Hill NSW 2154 AUSTRIA • Sigma-Aldrich Handels GmbH • A-1110 Wien
BELGIUM • Sigma-Aldrich N.V./S.A. • B-2880 Bornem BRAZIL • Sigma-Aldrich Quimica Brasil Ltda. • 01239-010 São Paulo, SP CANADA • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC • Sigma-Aldrich s.r.o. • 186 00 Praha 8 DENMARK • Sigma-Aldrich Denmark A/S • DK-2665 Vallensbaek Strand FINLAND • Sigma-Aldrich Finland/YA-Kemia Oy • FIN-00700 Helsinki
FRANCE • Sigma-Aldrich Chimie • 38297 Saint-Quentin-Fallavier Cedex GERMANY • Sigma-Aldrich Chemie GmbH • D-82041 Deisenhofen GREECE • Sigma-Aldrich (o.m.) Ltd. • Ilioupoli 16346, Athens
HUNGARY • Sigma-Aldrich Kft. • H-1067 Budapest INDIA • Sigma-Aldrich Co. • Bangalore 560 048 IRELAND • Sigma-Aldrich Ireland Ltd. • Dublin 24 ISRAEL • Sigma Israel Chemicals Ltd. • Rehovot 76100
ITALY • Sigma-Aldrich s.r.l. • 20151 Milano JAPAN • Sigma-Aldrich Japan K.K. • Chuo-ku, Tokyo 103 KOREA • Sigma-Aldrich Korea • Seoul MALAYSIA • Sigma-Aldrich (M) Sdn. Bhd. • 58200 Kuala Lumpur
MEXICO • Sigma-Aldrich Química S.A. de C.V. • 50200 Toluca NETHERLANDS • Sigma-Aldrich Chemie BV • 3330 AA Zwijndrecht NORWAY • Sigma-Aldrich Norway • Torshov • N-0401 Oslo
POLAND • Sigma-Aldrich Sp. z o.o. • 61-663 Poznań PORTUGAL • Sigma-Aldrich Quimica, S.A. • Sintra 2710 RUSSIA • Sigma-Aldrich Russia • Moscow 103062 SINGAPORE • Sigma-Aldrich Pte. Ltd.
SOUTH AFRICA • Sigma-Aldrich (pty) Ltd. • Jet Park 1459 SPAIN • Sigma-Aldrich Quimica, S.A. • 28100 Alcobendas, Madrid SWEDEN • Sigma-Aldrich Sweden AB • 135 70 Stockholm
SWITZERLAND • Supelco • CH-9471 Buchs UNITED KINGDOM • Sigma-Aldrich Company Ltd. • Poole, Dorset BH12 4QH
UNITED STATES • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email: supelco@sial.com

H

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

BIW