

Technical Bulletin

Storing and Thawing Serum

Introduction

Animal serum is an excellent source of nutrients for cells in culture because it contains proteins, lipids, salts, vitamins, minerals, amino acids and other components necessary for growth. When stored and handled correctly, the performance characteristics of serum can be maintained for many years. Improperly storing and/or thawing serum products can decrease not only the immediate and long-term stability, but also their effectiveness.

Growth promotion data demonstrate that serum maintains its growth characteristics throughout its shelf life if stored correctly.

The table below shows stability data for Fetal Bovine Serum (FBS) for a five-year period. To assure stability, sera must be correctly stored and handled.

Storage

To effectively preserve the integrity of animal serum, it should be stored frozen and protected from light. The recommended storage temperature is -10 to -40 °C. At temperatures below -40 °C, the bottles may become brittle resulting in an increased risk of breakage.

Multiple thaw/freeze cycles should be avoided as they will hasten the degradation of serum nutrients and can induce the formation of insoluble precipitates. For this reason, serum should never be stored in "frost-free" freezers. These appliances occasionally warm themselves to avoid internal ice deposits and are detrimental to the clarity and stability of frozen serum products.

Long-term Storage of SERA: Test Results							
Product	Test	Time 0	Time 1 Year	Time 2 Year	Time 3 Year	Time 4 Year	Time 5 Year
Irradiated FBS	Multiple Passage	110	104	116	126	93	90
	Plating Efficiency	61	36	51	35	40	50
	Cloning Efficiency	35	43	31	21	25	23
FBS	Multiple Passage	113	103	99	95	103	102
	Plating Efficiency	59	58	57	74	32	42
	Cloning Efficiency	49	56	42	47	28	30

Multiple Passage is % of control growth, while Plating Efficiency (PE) and Cloning Efficiency (CE) are % PE and CE.



Recommended Thawing Procedure

- Remove the serum bottles from the freezer and allow them to acclimate to room temperature for approximately 10 minutes.
- Place each container in a 30 to 37 °C water bath or incubator. Excessive temperatures will degrade heat labile nutrients. If using a water bath, prevent the bottle caps from being completely submerged.
- 3. Gently swirl or shake the bottles every 10 15 minutes until the serum is completely thawed.

Cryoprecipitates

The method used to thaw serum is crucial to its optimum performance. The key to proper thawing is periodic agitation. If a bottle of serum is not periodically shaken or swirled as it thaws, salt and protein gradients will form throughout the liquid portion. Within these gradients are high concentrations of salts, proteins and lipids which can lead to the formation of crystalline or flocculent precipitates. These "cryoprecipitates" are not toxic to cell cultures, but they affect the appearance and consistency of each bottle of serum.

Small amounts of cryoprecipitates are not uncommon, even in serum that is thawed using the recommended procedure. This is normal and will not affect product performance.

If serum is thawed incorrectly, a much greater amount of cryoprecipitate will form and is often insoluble. Filtering serum to remove cryoprecipitates is not recommended and could result in the loss of nutrients, such as growth factors, mitogens and other proteins.

Additional Recommendation

After thawing, use the serum promptly. Liquid serum may be stored refrigerated (2 to 8 °C) up to four weeks. To avoid thaw/freeze cycles or long periods of refrigeration, it is recommended that any unused serum be immediately dispensed into small, useful aliquots and refrozen for future use.

For more information about this subject or other SAFC Biosciences' products and services, please call our Technical Services department.

United States

SAFC Biosciences, Inc. 13804 W. 107th Street Lenexa, Kansas 66215 USA

Phone: +1 913-469-5580

Toll free-USA: 1 800-255-6032

Fax: +1 913-469-5584 E-mail: info-na@sial.com

Europe

SAFC Biosciences Ltd.
Smeaton Road, West Portway
Andover, Hampshire SP10 3LF
UNITED KINGDOM

Phone: +44 (0)1264-333311
Fax +44 (0)1264-332412
E-mail info-eu@sial.com

Warranty, Limitation of Remedies

SAFC Biosciences warrants to the purchaser for a period of one year from date of delivery that this product conforms to its specifications. Other terms and conditions of this warranty are contained in SAFC Biosciences' written warranty, a copy of which is available upon request. ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE EXCLUDED. In no case will SAFC Biosciences be liable for any special, incidental, or consequential damages arising out of this product or the use of this product by the customer or any thrid party based upon breach of warranty, breach of contract, negligence, strict tort, or any other legal theory. SAFC Biosciences expressly disclaims any warranty against claims by any third party by way of infringement or the like. THIS PRODUCT IS INTENDED FOR PURPOSES DESCRIBED ONLY AND IS NOT INTENDED FOR ANY HUMAN OR THERAPEUTIC USE.

Additional Terms and Conditions are contained in the product Catalog, a copy of which is available upon request.

© 2010 SAFC Biosciences, Inc.

Issued February 2010 T010 0103 1005 1105 0306

> 73246-509942 1020