

Fast, Simple Separation of Dabsylated Amino Acids

Thirty-five DABS-AA can be separated in 25 minutes at room temperature on a 15cm x 4.6mm ID SUPEL COSIL LC-DABS column under HPLC conditions. DABS derivatization takes only 10 minutes at 70°C, permitting complete reaction of primary and secondary amino acids. These derivatives can be detected at visible wavelengths, at picomole concentrations. They also are stable at temperatures from -20°C to room temperature.

Key Words:

- DABS-AA • DABTH-AA • reversed phase HPLC
- derivatization

SUPEL COSIL™ LC-DABS and SUPEL COSIL LC-18-DB columns (15cm x 4.6mm ID, 3µm particles) will separate dimethylaminoazobenzene sulfonyl amino acids (DABS-AA) and dimethylaminoazobenzene thiohydantoin amino acids (DABTH-AA), respectively. Using reversed phase HPLC, these separations are fast, simple, and reproducible.

Precolumn derivatization of amino acids with DABS-chloride (Cl) takes only 10 minutes at 70°C and permits complete reaction of primary and secondary amino acids. DABS derivatives are also very stable at temperatures ranging from -20°C to room temperature.

DABS-AAs can be detected at visible light wavelengths. This permits HPLC analysis of amino acids at picomole concentrations and eliminates baseline noise that can occur when using UV wavelengths. The HPLC method discussed here takes only 25 minutes for a room temperature separation of approximately 35 DABS-AAs.

To develop this HPLC method for analyzing DABS-AAs, various reversed phase columns and several potassium dihydrogen phosphate and sodium acetate buffers were tested. Different pHs, flow rates, gradients, and organic solvents and modifiers also were evaluated. A 15cm x 4.6mm ID, 3µm particles SUPEL COSIL LC-DABS column, a 2cm x 4.6mm ID, 5µm particles Supelguard™ LC-18-T guard column, and a two-eluent mobile phase—consisting of 25mM potassium dihydrogen phosphate (pH 6.8) as solvent A, and acetonitrile:2-propanol (75:25) as solvent B—delivered the best analysis (Figure A1).

Along with the separation of 35 DABS-AAs, this method also resolves DABS derivatives of taurine, d-hydroxylysine, norleucine, cysteic acid, cystine, S-carboxymethylcysteine, and S-sulfocysteine. The separation of DABS-norleucine adds dimension to the method because this derivative can serve as an internal standard for checking different steps of the sample handling procedure. Resolution of cysteic acid, cystine, S-carboxy methylcysteine and

S-sulfocysteine indicates that the HPLC method will continue to work, even when different experimental conditions are used to treat a protein before hydrolysis.

Changing solvent B to contain acetonitrile and methanol in a 70:30 ratio delivers a chromatographic profile (Figure A2) similar to Figure A. This second analysis adds hydroxyproline, methionine sulfone, and methionine sulfoxide to the profile.

The resolution of hydroxyleucine, under conditions of the second HPLC procedure, makes this method useful for studying different aspects of metabolism in connective tissue. The method is also useful for evaluating distribution of hydroxyleucine in biological fluids. Resolution of methionine sulfone and methionine sulfoxide makes it possible to follow methionine oxidation of a protein *in vivo* under different biological conditions, or when treating it with performic acid.

In any amino acid analysis, tryptophan poses a special problem. Tryptophan is unstable under hydrochloric acid (HCl) hydrolysis conditions. Hydrolysis with methanesulfonic acid and derivatization with DABS-Cl prepare tryptophan for HPLC analysis, without adversely affecting the separation of other amino acids. But this procedure loses its effectiveness for samples with a sugar concentration of more than 20%. For analyzing tryptophan in protein samples with sugar concentrations greater than 20%, sodium hydroxide hydrolysis followed by DABS-Cl derivatization and HPLC analysis is a more effective procedure.

Resolution of DABTH-AAs on a 15cm x 4.6mm ID, 3µm particles SUPEL COSIL LC-18-DB column is fast and sensitive, permitting analyses of very small amounts of material. Figure B shows the separation of DABTH-AAs, including DABTH-norleucine, in only 18 minutes, by two HPLC methods. The gradient conditions remove excess reagent from the column during the analysis. The column is protected with a 2cm x 4.6mm ID, 5µm particles Supelguard LC-18-DB guard column. Although many analyses of DABTH derivatized amino acids can be performed on the same analytical column, you should change the guard column every 200-250 runs.

The 3µm particles in the SUPEL COSIL LC-18-DB column deliver chromatographic peaks nearly twice as high as those obtained from a column packed with 5µm particles. The columns also display excellent column-to-column reproducibility, indicating very reliable methodology. The consistently reproduced retention times for the DABTH-AAs confirm that these HPLC methods are excellent for routine analysis and identification of amino acids.


These methods offer exceptional sensitivity for detailed biological and structure-function studies of even trace concentrations of proteins and peptides.

Figure A. DABS Amino Acids by HPLC

Column: SUPELCOSIL LC-DABS, 15cm x 4.6mm ID, 3µm particles
Cat. No.: 59137
Mobile Phase: A = 25mM potassium dihydrogen phosphate (pH 6.8),
 B = acetonitrile:2-propanol, 75:25 (A1) or
 acetonitrile:methanol, 70:30 (A2)
Flow Rate: 2mL/min
Det.: 436nm UV (A1) or variable wavelength UV (A2)
Sample: 5µL DABS-derivatized amino acids

A1-System I

Gradient Program Time (min.)	%B
0-1	20
1-4	20-23
4-9	23
9-10	23-27
10-14	27
14-19	27-35
19-25	35-60
25-26	60-70
26-29	70
29-29.1	70-20
29.1-35.1	20



A2 — System II

Gradient Program Time (min.)	%B
0-1	25
1-5	25-28
5-12	28
12-13	28-30
13-19	30
19-30	30-40
30-39	40-65
39-40	65-75
40-44	75
44-44.1	75-25
44.1-53.1	25





Figure B. DABTH-Amino Acids by HPLC

Column: SUPELCOSIL LC-18-DB, 15cm x 4.6mm ID, 3µm particles
Cat. No.: 58993
Mobile Phase: gradient analysis, A = 35mM sodium acetate (pH 5.1),
 B = acetonitrile
Flow Rate: 2mL/min
Det.: 436nm UV
Sample: 5µL DABTH-derivatized amino acid standards,
 5-15 picomoles (A) or 0.5-1.5 picomoles (B)

Gradient Program Time (min.)	%B
0-1	35
1-2	35-44
2-5	44
5-10.5	44-54
10.5-24	54
24-25	54-85
25-31	85
31-31.1	85-35
31.1-37.1	35



Ordering Information:

Description	Cat. No.
SUPELCOSIL LC-DABS Column 15cm x 4.6mm ID x 3µm particles	59137
SUPELCOSIL LC-18-DB Column 15cm x 4.6mm ID x 3µm particles	58993
Supelguard Guard Column Kits includes guard column (2cm x 4.6mm ID x 5µm particles), holder, and connecting hardware	
Supelguard LC-18-T (for LC-DABS column)	59620
Supelguard LC-18-DB	59555
Replacement Supelguard Guard Columns, pk. 2	
Supelguard LC-18-T (for LC-DABS column)	59621
Supelguard LC-18-DB	59565

Acknowledgment

We gratefully acknowledge the contributions of Vilberto Stocchi, Istituto di Chimica Biologica, Università degli Studi di Urbino, Urbino, Italy.

Reference

1. Stocchi, V., *et al.*, Anal. Biochem., **179**, 107-117 (1989).
 Reference not available from Supelco.
 SUPELCOSIL and Supelguard are trademarks of Sigma-Aldrich Co.

Contact our Technical Service Department
 (phone 800-359-3041 or 814-359-3041, FAX 800-359-3044 or
 814-359-5468) for expert answers to your questions.

Note 124

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA • Sigma-Aldrich de Argentina, S.A. • Buenos Aires 1119 **AUSTRALIA** • Sigma-Aldrich Pty. Ltd. • Castle Hill NSW 2154 **AUSTRIA** • Sigma-Aldrich Handels GmbH • A-1110 Wien
BELGIUM • Sigma-Aldrich N.V./S.A. • B-2880 Bornem **BRAZIL** • Sigma-Aldrich Química Brasil Ltda. • 01239-010 São Paulo, SP **CANADA** • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC • Sigma-Aldrich s.r.o. • 186 00 Praha 8 **DENMARK** • Sigma-Aldrich Denmark A/S • DK-2665 Vallensbaek Strand **FINLAND** • Sigma-Aldrich Finland/YA-Kemia Oy • FIN-00700 Helsinki
FRANCE • Sigma-Aldrich Chimie • 38297 Saint-Quentin-Fallavier Cedex **GERMANY** • Sigma-Aldrich Chemie GmbH • D-82041 Deisenhofen **GREECE** • Sigma-Aldrich (o.m.) Ltd. • Ilioupoli 16346, Athens
HUNGARY • Sigma-Aldrich Kft. • H-1067 Budapest **INDIA** • Sigma-Aldrich Co. • Bangalore 560 048 **IRELAND** • Sigma-Aldrich Ireland Ltd. • Dublin 24 **ISRAEL** • Sigma Israel Chemicals Ltd. • Rehovot 76100
ITALY • Sigma-Aldrich S.r.l. • 20151 Milano **JAPAN** • Sigma-Aldrich Japan K.K. • Chuo-ku, Tokyo 103 **KOREA** • Sigma-Aldrich Korea • Seoul **MALAYSIA** • Sigma-Aldrich (M) Sdn. Bhd. • Selangor
MEXICO • Sigma-Aldrich Química S.A. de C.V. • 50200 Toluca **NETHERLANDS** • Sigma-Aldrich Chemie BV • 3330 AA Zwijndrecht **NORWAY** • Sigma-Aldrich Norway • Torshov • N-0401 Oslo
POLAND • Sigma-Aldrich Sp. z o.o. • 61-663 Poznań **PORTUGAL** • Sigma-Aldrich Química, S.A. • Sintra 2710 **RUSSIA** • Sigma-Aldrich Russia • Moscow 103062 **SINGAPORE** • Sigma-Aldrich Pte. Ltd.
SOUTH AFRICA • Sigma-Aldrich (pty) Ltd. • Jet Park 1459 **SPAIN** • Sigma-Aldrich Química, S.A. • 28100 Alcobendas, Madrid **SWEDEN** • Sigma-Aldrich Sweden AB • 135 70 Stockholm
SWITZERLAND • Supelco • CH-9471 Buchs **UNITED KINGDOM** • Sigma-Aldrich Company Ltd. • Poole, Dorset BH12 4QH
UNITED STATES • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email: supelco@sial.com

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

BHL