

SPPS-Based Synthesis of GLP-1 Analogs using Novabiochem® Solutions

Introduction


Glucagon-like peptide-1 (GLP-1) and its receptor agonists, such as semaglutide and liraglutide, play a key role in the management of type 2 diabetes and obesity by regulating glucose homeostasis and glycogen metabolism. With therapeutic indications and expanding research opportunities, GLP-1 analogs continue to gain significant attention, accompanied by a growing demand for GLP-1 based drugs. Recent innovations, emerging leads, and increasing patent activity in this space underscore the need for synthesis of next-generation GLP-1 receptor agonists (GLP-1RAs) with enhanced quality, safety, and efficacy.¹

Researchers are increasingly focused on addressing key challenges and optimizing early-stage peptide synthesis, particularly in the context of Structure-Activity Relationship (SAR) studies of GLP-1. Consequently, the need for high-quality materials and stringent standards becomes increasingly critical to ensure reliable and reproducible outcomes in GLP-1 research.

In this context, the Novabiochem® portfolio, offers high-quality, high-purity, and cost-effective resins, amide coupling agents, and protected amino acids for synthesis of GLP-1 analogs. These products support efficient, and reliable synthesis of GLP-1 peptides using Solid Phase Peptide Synthesis (SPPS), enabling consistently high success rates.

Challenges in GLP-1 Synthesis and Analog Development

Exenatide is recognized as the first GLP-1 analog synthesized using SPPS and the first to be approved for clinical use.² SPPS is a widely used method in GLP-1 research in which the peptide is assembled stepwise on an insoluble support. The process begins by anchoring the first amino acid to a solid resin, followed by alternating cycles of deprotection and coupling, where additional amino acids are sequentially added until the target



Molecular structure of semaglutide.

peptide is formed. The completed peptide is then cleaved from the support. This approach enables efficient, and rapid synthesis of high molecular weight and complex peptide sequences, including GLP-1 analogs.³

The use of orthogonal-protecting strategies, such as traditional Boc with benzyl-based side chain protection and advanced Fmoc with tBu protection, plays a pivotal role in successful SPPS. Fmoc with tBu provides greater flexibility for peptide chain modifications, thereby broadening the scope of SPPS. As a result, SPPS has evolved into a powerful, and versatile methodology that supports rapid, and efficient peptide synthesis, making it the method of choice for GLP-1 research.⁴

Despite these advancements, GLP-1 analog synthesis remains associated with several inherent challenges during drug discovery and development, largely due to the complexity of the synthetic process. Some of the typical challenges, underlying issues, and possible solution strategies are presented in **Table 1**.

Table 1. Typical challenges, underlying issues and solution strategies in GLP-1 peptide synthesis

Challenges	Underlying issues	Solution strategies
Synthesis complexity	- Increased risk of failure due to large amino acid chain length. (e.g., GLP-1; 30–31 amino acids) - Aggregation of peptides during synthesis - Challenging coupling steps requiring optimization	High-purity reagents and specialized coupling agents
Reproducibility concerns	Batch to batch variability in reagents affects the reproducible quality of final peptide products	High-quality reagents and with rigorous Quality Control (QC) ensure consistent and superior product quality
Racemization risks	Alterations in biological activity due to changes in racemization	Advanced coupling reagents with low racemization
Purification difficulties	- Difficulty in separating similar analogs - Need for high crude purity to simplify downstream processing	High coupling efficiency reduces impurity formation
Analog optimization	SAR studies require diverse modifications, optimizations and multiple variants	A broad product portfolio enables efficient exploration
Scaling challenges	Inconsistent performance across different scales (mg to g) during development	Reagents that perform reliably across different scales
Time pressure	Rapid turnaround required for analog library generation in competitive research landscape	Reliable reagents reduce troubleshooting time
Budget constraints	Inappropriate reagent selection increases cost and resource consumption	A broad portfolio of high-quality application-specific reagents reduces waste and failure in synthesis

Achieving high-purity peptide products is a key objective for peptide synthesis, as purity directly influences performance, reliability. The quality of starting materials plays a decisive role in determining the purity of final peptide product. High-grade resins, protected amino acids, and reagents form the foundation for robust, and reproducible peptide synthesis by minimizing synthesis related challenges.

Role of Novabiochem® Products in GLP-1 Analog Synthesis

The success of modern peptide therapeutics, particularly GLP-1RAs, has been strongly supported by advancements in SPPS. The use of high-purity peptide building blocks empowers SPPS workflows to achieve complex and previously inaccessible peptide targets, thereby ensuring reliability, yield, and overall synthesis success.⁵

The Novabiochem® portfolio encompasses a comprehensive range of innovative products and solutions for peptide synthesis, high-throughput organic chemistry, peptide labelling, and customized manufacturing. These products are designed to meet stringent quality and purity requirements, supporting consistent performance and high yields in SPPS across both standard and premium-grade peptide synthesis applications. This level of performance is ensured through rigorous specifications, defined by our expert R&D team and supported by advanced QC technologies, enabling high consistency, purity, and performance.

The key synthesis requirements of GLP-1 analogs and the corresponding significance of Novabiochem® products are summarized in **Table 2**.

Table 2. Synthesis requirements of GLP-1 analogs and significance of Novabiochem® products

Product type	Significance
Standard Fmoc- and Boc-Amino Acids	- Specialized Fmoc-amino acids help to overcome peptide synthesis challenges enabling reliable, consistent, and high-purity peptide synthesis. For example, Fmoc-Trp(Boc)-OH, Fmoc-Asn(Trt)-OH etc. - The 20 standard Fmoc-protected amino acid building blocks meet stringent specifications such as ◇ ≥99.0% HPLC purity ◇ ≥99.8% enantiomeric purity ◇ ≤0.2% free amine content (improves Fmoc stability during storage) ◇ ≤0.02% acetate content - Improve peptide yields and facilitate purification - Ensure consistent impurity profiles in crude peptides, enabling smoother scale-up and reducing Good Manufacturing Practices (GMP) related challenges - Boc-amino acids enable simultaneous peptide release and global deprotection, while preserving acid-sensitive residues and minimizing base-mediated side reactions
Unnatural Amino Acids	- Enable fine-tuning of peptide structure and function - Improve protease resistance, stability, and enable precise site-specific modifications - Broad portfolio ranging from modified lysines for late-stage conjugation to stereochemically unique residues like DL-isoserine

Product type	Significance
Resins	- Extensive resin portfolios - High-loading, low-swelling resins for large scale production of shorter peptides - Low-loading, high-swelling resins for synthesis of long or difficult sequences - Comprehensive portfolio ranging from traditional 4-Methylbenzhydrylamine (MBHA) and Rink amide resins to advanced preloaded and specialty supports - Unloaded resins for Fmoc SPPS - Preloaded resins for SPPS - Optimized for reliable SPPS, enabling high-purity peptide synthesis and efficient cleavage to support complex modifications - Specialized base resins such as NovaSyn® TG, NovaGel™, Polyethylene glycol-polyacrylamide (PEGA) resin for continuous flow of peptide synthesis
Amide Coupling Agents	- Tailored collection of coupling reagents, including C-free uronium and phosphonium types - Support high-yield and low-racemization peptide synthesis - Ensure compatibility with challenging sequences and reproducibility across batches
Peptide Labelling Agents	Enable site-specific incorporation of functional tags during SPPS supporting efficient monitoring, purification, characterization, quality assessment, and direct preparation of assay-ready peptides

Building on these capabilities, products are manufactured to meet the **MQ200-400** quality matrix, ensuring consistent performance and reliability of final peptide products. A broad range of supply options, from R&D to semi-bulk and full-scale production, supports synthesis across all development stages from a single source. Flexible solutions enable customization in both product specifications and volumes to address diverse application needs.

Recent Innovations in Synthesis of GLP-1 Analogs

Small-molecule GLP-1RAs, such as danuglipron and lotiglipron, have emerged as “next-in-class” candidates and have driven significant recent innovation, contributing to the development of next-generation GLP-1RAs. Notably, between 2021 to 2024, approximately 70 patents were published focusing on the design, synthesis, and functional activity of GLP-1RAs. This trend highlights the continued relevance and accelerating progress in the synthesis of novel GLP-1RAs, as illustrated in **Figure 1**.⁶

Building on this growing body of work, the methodologies used to develop GLP-1RAs show that subtle structural

modifications can improve half-maximal Effective Concentration (EC₅₀) and pharmacokinetic profiles, whereas more extensive redesigns are often required to deliver broader improvements in efficacy and ADME (Absorption, Distribution, Metabolism, Excretion) properties. These insights highlight the scientific complexity and the inherently non-trivial nature of designing clinically viable GLP-1RAs.⁶

These innovations have reinforced the importance of Structure–Activity Relationship (SAR) for the GLP-1 receptor. Accordingly, pharmaceutical companies and research institutions are increasingly prioritizing optimized peptide synthesis approaches, with a strong emphasis on structural modifications that influence SAR. Even minor chemical modifications and transformations to lead scaffolds can result in improved binding capacity, therapeutic activity and overall yield.

Conclusion

The rapid evolution of GLP-1 research continues to reshape the therapeutic innovation in metabolic diseases, driven largely by advances in SPPS. However, the development of GLP-1 analogs remains challenging due to structural complexity, sequence dependent constraints, which demand high-purity materials, high-quality reagents, optimized workflows, and reliable synthesis strategies. More than 160 drug candidates are currently in development, including multi-agonist and once-weekly therapies targeting multiple pathways to improve efficacy.⁷ In summary, the GLP-1 drug landscape is poised for significant growth, driven by their strong clinical efficacy, expanding indications, and a robust pipeline of innovative therapies.

Against this backdrop, the Novabiochem® portfolio serves as a trusted partner for peptide synthesis, by delivering premium-grade resins, amino acids, and reagents designed for superior purity, consistency, and efficiency. Backed by proven innovation in peptide chemistry, these solutions enable the acceleration of GLP-1 research and confidently support the development of next-generation therapeutics.

[Click here](#) to explore the comprehensive range of Novabiochem® products for GLP-1 analog synthesis.




Figure 1. Graphical distribution of published patents dedicated to GLP-1RAs during 2018–2024.⁶

Image Credit: [DOI](#)

References

1. Nauck, M. A.; Quast, D. R.; Wefers, J.; Meier, J. J. GLP-1 receptor agonists in the treatment of type 2 diabetes-state of the art. *Mol. Metab.* **2021**, *46*, 101-102. DOI: [10.1016/j.molmet.2020.101102](https://doi.org/10.1016/j.molmet.2020.101102)
2. Barnett, A. Exenatide. *Expert Opin. Pharmacother.* **2007**, *8* (15), 2593-2608. DOI: [10.1517/14656566.8.15.2593](https://doi.org/10.1517/14656566.8.15.2593)
3. Merrifield, R. B. Solid phase peptide synthesis. I. The synthesis of a tetrapeptide. *J. Am. Chem. Soc.* **1963**, *85* (14), 2149-2154. DOI: [10.1021/ja00897a025](https://doi.org/10.1021/ja00897a025)
4. Amblard, M.; Fehrentz, J.-A.; Martinez, J.; Subra, G. Methods and protocols of modern solid phase peptide synthesis. *Mol. Biotechnol.* **2006**, *33* (3), 239-254. DOI: [10.1385/MB:33:3:239](https://doi.org/10.1385/MB:33:3:239)
5. Behrendt, R.; White, P.; Offer, J. Advances in Fmoc solid phase peptide synthesis. *J. Pept. Sci.* **2016**, *22* (1), 4-27. DOI: [10.1002/psc.2836](https://doi.org/10.1002/psc.2836)
6. Tolkacheva, E. V.; Saliev, A. Y.; Salakhov, T. L.; Balakin, K. V.; Ivanov, R. A.; Chernyshov, V. V. Next-in-class GLP-1R danuglipron- and lotiglipron-like agonists: a patent review (2020-2024). *Curr. Med. Chem.* **2025**, *32*. DOI: [10.2174/0109298673366258250710101146](https://doi.org/10.2174/0109298673366258250710101146)
7. Zheng, Z.; Zong, Y.; Ma, Y.; Tian, Y.; Pang, Y.; Zhang, C.; Gao, J. lucagon-like peptide-1 receptor: mechanisms and advances in therapy. *Signal Transduct. Target. Ther.* **2024**, *9* (1), 1-29. DOI: [10.1038/s41392-024-01931-z](https://doi.org/10.1038/s41392-024-01931-z)



To place an order or receive technical assistance:
SigmaAldrich.com/support



For local contact information:
SigmaAldrich.com/offices

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany
SigmaAldrich.com

We have built a unique collection of life science brands with unrivalled experience in supporting your scientific advancements.

Millipore® Sigma-Aldrich® Supelco® Milli-Q® SAFC® BioReliance®

© 2026 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, BioReliance, Millipore, Milli-Q, SAFC, Sigma-Aldrich, Supelco, Novabiochem, NovaSyn and NovaGel are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MK_AN14341EN Ver. 1.0 70708 06/2026

