


Basics of Chiral HPLC and High Throughput Chiral Methods Development Screening

Presented by Tom Beesley
Sigma-Aldrich/Supelco

T408174

The Field of Stereochemistry

1. All isomers have same chemical formula but differ in the arrangement of certain chemical groups in space.
2. Many types of isomers have different physical and chemical properties and can be separated by conventional phases like C18.
3. Enantiomers have same physical and chemical properties and can be recognized or separated only in chiral environment.





The Inclusion Complex



The basis for many chiral separations, especially in the reversed phase mode is a phenomena called inclusion complexing. First described for the polyglucose structures, cyclodextrins, it has been identified as a mechanism for the macrocyclic glycopeptides as well as the cellulose and amylose CSPs. An understanding of this phenomena is then imperative to an understanding of how these phases separate.

Structure and dimensions of the most common cyclodextrin molecules.



α -cyclodextrin




β -cyclodextrin




γ -cyclodextrin



Schematic of one glucose unit in a cyclodextrin molecule



Inclusion Complex Schematic




Inclusion complexing effectiveness is dependent upon the position and strength of substituents on benzene or fused ring analytes. The elution order is typically meta<ortho<para. The bulky character of the meta and ortho limits full inclusion. The linear nature of the para structure results in full inclusion. Selectivity for this isomer is generally highest when strong hydrogen bonding groups are present.

Acetonitrile accelerates release from the cavity. Methanol reduces hydrogen bonding effects to the cyclodextrin hydroxyl groups. Combinations of acetonitrile and methanol in water are useful in optimizing a separation.


Examples of Inclusion Phenomena

Positional Isomers

Xylenes




Peak 1 - 5.74 min. (ortho)
Peak 2 - 6.38 min. (meta)
Peak 3 - 6.83 min. (para)




CYCLOBOND I 2000
Mobile Phase: 30/70: CH₃CN/H₂O

Cresols




Peak 1 - 6.62 min. (ortho)
Peak 2 - 7.31 min. (meta)
Peak 3 - 9.82 min. (para)




CYCLOBOND I 2000
Mobile Phase: 40/60: CH₃OH/H₂O


Structural Isomers



Corticosterone



Cortisone




Hydrocortisone

Peak 1 - 6.29 min.

Peak 2 - 8.02 min.

Peak 3 - 9.18 min.




CYCLOBOND I 2000

Mobile Phase: 40/60: CH₃CN/H₂O


Geometric Isomers

cis/trans Stilbene



Peak 1 - 5.60 min.

Peak 2 - 7.14 min.




CYCLOBOND I 2000

Mobile Phase: 70/30: MeOH/H₂O


Enantiomers

Remacemide



Peak 1 - 5.74 min.

Peak 2 - 6.79 min.



CYCLOBOND I 2000


Mobile Phase:

3.5/96.5: CH₃CN/PO₄(0.1M) pH 3.5


Models for Identification of Potential Racemates

A. Enantiomers based on **CARBON**

Type 1: **Stereogenic Center**




Ibuprofen



Peak 1 - 5.77 min.

Peak 2 - 6.47 min.




CHIROBIOTIC V


10/90: THF/20mM Na Citrate, pH 6.3

Models for Identification of Potential Racemates


Type 2: **Axis of Symmetry**



1,1'-Binaphthyl-2,2'-diylhydrogenphosphate



Peak 1 - 9.12 min.
Peak 2 - 10.03 min.




CHIROBIOTIC V
30/70: ACN/0.1% TEAA, pH 4.1

Models for Identification of Potential Racemates

B. Enantiomers based on **SULFUR**


Possibilities include sulfoxide, sulfoximide, sulfinate and sulfonium ion.

Methyl phenyl sulfoxide



Peak 1 - 10.86 min.

Peak 2 - 14.49 min.




CHIROBIOTIC TAG
40/60: EtOH/Hex

Models for Identification of Potential Racemates


C. Enantiomers based on **PHOSPHORUS**

Possibilities include phosphine, phosphine oxide, phosphinate and phosphonium ion.


Phosphine Oxide



Ifosfamide



Peak 1 - 7.53 min.
Peak 2 - 8.48 min.




CHIROBIOTIC T
10/90: THF/H₂O


Models for Identification of Potential Racemates

D. Enantiomers based on **NITROGEN**

Possibilities include amine oxide and ammonium ion.



Peak 1 - 8.31 min.
Peak 2 - 11.40 min.





**(+) and (-)-Methylallylphenylbenzyl
ammonium iodide**

CHIROBIOTIC V
100/0.02/0.01: MeOH/HOAc/TEA

MESO COMPOUND:

A compound whose functional groups are superimposable on their mirror images even though they contain chiral centers. It is optically inactive.

2,6-bis-(1-phenyl-ethyl)-4-methylaniline




Peak 1 – 6.19 (meso)
Peak 2 – 7.96
Peak 3 – 8.47

CYCLOBOND I 2000 RSP
250x4.6mm


30/70: CH₃CN/0.1% TEAA, pH 6.5
1.0 mL/min.

DIASTEREOMERIC/ENANTIOMERIC STRUCTURES:




Note: Conventional reversed phase would separate this analyte into 2 peaks (diastereomers). Only a chiral stationary phase would resolve all 4 enantiomers.

- Peak 1 – 11.54 min.
- Peak 2 – 12.29 min.
- Peak 3 – 16.57 min.
- Peak 4 – 19.35 min.



Labetalol



CHIROBIOTIC V
100%: MeOH/0.1% ATFA, v/w

Principles Governing Chiral Separation

Concept: formation of a **diastereomeric complex** in a chromatographic equilibrium such that the nonchiral interactions are at minimum strength and the differential chiral interaction is at maximum strength. Identifying those points of interaction between the stationary phase and the racemate guides you in the choice of CSPs and the best conditions under which to operate.

Non-chiral interactions generally anchor a molecule and, therefore, assist in the formation of the diastereometric complex. Both **dipole interactions** driven in normal phase and polar organic solvents; **inclusion complexation** driven in reversed phase modes and **ion exchange** driven in polar ionic mode are the first significant areas to address the potential of an appropriate chiral stationary phase.

Types of Interaction

All forces in the chiral interaction process **do not** have to be attractive, they can be attractive as well as repulsive. Commonly described forces for chiral recognition are listed as follows:

- $\pi-\pi$
- Hydrogen bonding
 - a. Hydrogen donor site
 - b. Hydrogen acceptor site
- Inclusion complexation
- Steric hindrance
- Dipole-dipole
- Ionic interaction

Modern Chiral Stationary Phases

Polymeric

Synthetic

- Methacrylate
- Polycyclic amine-3

Focus – normal phase,
low cost,
reproducible, high
capacity

Astec-Supelco
Phases

Natural

- Cellulose
- Amylose
- Proteins

Small molecule ligands

- Copper complex-2
- π -complex
- Crown ether
- Cyclodextrin-12
- Macrocyclic glycopeptides-6

These 3 cover applications that
are possible on 140 CSPs that
have been created worldwide in
this area.

Astec-Supelco Chiral Phases

Polycyclic Amine Phases

P-CAP, P-CAP-DP

Copper Complex Phases

CLC-D, CLC-L

Macrocyclic Glycopeptide Phases

CHIROBIOTIC V2, CHIROBIOTIC V,
CHIROBIOTIC T, CHIROBIOTIC T2,
CHIROBIOTIC TAG, CHIROBIOTIC R

Astec-Supelco Chiral Phases con't


Cyclodextrin Phases

CYCLOBOND I 2000

CYCLOBOND I 2000 AC, CYCLOBOND I 2000 SP,
CYCLOBOND I 2000 RSP, CYCLOBOND I 2000 HP-RSP,
CYCLOBOND I 2000 SN, CYCLOBOND I 2000 RN,
CYCLOBOND I 2000 DM, CYCLOBOND I 2000 DMP,
CYCLOBOND I 2000 DNP,
CYCLOBOND II, CYCLOBOND II AC

Structure of Vancomycin CSP

CHIROBIOTIC V




Structure of Teicoplanin CSP

**Teicoplanin,
CHIROBIOTIC T**

sugar and alkyl chain


→ **Key sites**



ionic site

**Teicoplanin Aglycone,
CHIROBIOTIC TAG**

ionic site




ionic site

CHIROBIOTIC V2 and CHIROBIOTIC T2


- Extensions of the CHIROBIOTIC V2 and T2 created by changing the silica porosity; position and type of several linkages and the chain length used to anchor the ligand.
- These changes have **enhanced** the **selectivity** and **capacity** of these phases for a large number of analyte classes
- The CHIROBIOTIC V2 has statistically a greater number of chiral separations than the V while the CHIROBIOTIC T performs better than the T2.

Proposed Structure of Ristocetin A CSP

CHIROBIOTIC R



Bonded Derivatized Cyclodextrins



* Signifies stereogenic center

R =	Suffix	CD Type
- OCH ₃	DM (methylated)	β-CD
- COCH ₃	AC (acetylated)	β-CD γ-CD
$\begin{array}{c} \text{OH} \\ \\ -\text{CH}_2\text{CHCH}_3 \\ * \end{array}$	SP or RSP/ HP-RSP (hydroxypropyl ether)	β-CD
$\begin{array}{c} \text{CH}_3 \\ \\ -\text{CONHCH} \\ * \end{array} \text{---} \text{Naphthyl}$	RN or SN (naphthylethyl carbamate)	β-CD
$\begin{array}{c} \text{CH}_3 \\ \\ -\text{CONH} \text{---} \text{C}_6\text{H}_3 \text{---} \text{CH}_3 \end{array}$	DMP (3,5-dimethylphenyl carbamate)	β-CD
$\begin{array}{c} \text{O}_2\text{N} \\ \\ \text{---} \text{C}_6\text{H}_3 \text{---} \text{CF}_3 \\ \\ \text{O}_2\text{N} \end{array}$	DNP (2,6-dinitro-4-trifluoromethyl phenyl ether)	β-CD

Most Productive CYCLOBOND Phases

- CYCLOBOND I 2000
- CYCLOBOND I 2000 HP-RSP (highest hit rate)


Primarily basic chiral compounds have been resolved on the RSP and to a lesser extent both neutrals and acidics. All successful separations have been in the reversed phase mode with the organic component in the range of 5-40%.

- CYCLOBOND I 2000 DMP (second most productive CSP in this line) This is a π -basic phase.
- CYCLOBOND I 2000 DNP (new addition)

This is a π -acidic phase that has demonstrated separations not previously possible on any cyclodextrin phase.

CYCLOBOND

two different enantioselective retention mechanisms



A: Polar organic mode (surface)

B: Reversed phase mode (inclusion)

Mobile phases types and mechanisms

Normal Phase

Reversed Phase

Polar Organic

Polar Ionic

Normal phase and types of interaction

Composition:

Hydrocarbon solvent: heptane + polar alcohol: IPA or EtOH. Addition of MtBE a useful additive to enhance steric effects


Dominant interactions:

π - π interaction, hydrogen bonding

Type of CSPs:



CYCLOBOND I 2000 DMP, SN, RN, CHIROBIOTIC V2, T, R and TAG

Normal Phase Solvents Mephenytoin



Ref: Anal. Chem., Vol. 66 (9),
p 1473-184, May 1994.
Armstrong, Tang, S. Chen, Zhou,
Bagwill and J.R. Chen.

Gradient Separation in the Normal Phase



Derivative: 3,4-Dinitrobenzoyl
CYCLOBOND I 2000 SN
Hexane/EtOH Gradient

Reversed Phase

Composition:

Organic solvent: ACN, MeOH or THF + aqueous buffer: TEAA, NH₄OAc

Dominant interactions:


Inclusion complexation, hydrogen bonding

Type of CSPs:

CYCLOBOND I 2000, RSP, AC, DMP, DNP
CHIROBIOTIC V2, T, R and TAG

Reversed Phase Separation

5-Methyl-5-phenyl hydantoin



Ref: Anal. Chem.,
Vol.66(9), p. 1473-
1484, May 1994.

Reversed Phase Solvents

Factors Influencing a Separation

- pH: CHIROBIOTICS 3.0-7.0
- Organic modifier: MeOH, ACN
- Buffer type and concentration
- Flow rate: 0.2 to 1.0 mL/min
- Temperature: 0 to 50°C

Polar Organic Phase

Composition:

ACN + MeOH + HOAc + TEA

Dominant interactions:

Hydrogen bonding, dipole-dipole

Type of CSPs:

Cyclodextrins only –
CYCLOBOND I 2000, HP-RSP, AC, DMP, DNP

New Polar Ionic Mode: Ionizable analytes; CHIROBIOTIC phases only!

Composition:

MeOH + HOAc + TEA

Dominant interactions:

Ionic interaction, hydrogen bonding

Type of CSPs:

Macrocyclic glycopeptides only -
CHIROBIOTIC V2, T, R and TAG

Polar Ionic Mode: *CHIROBIOTIC Phases Only*

BENEFITS:

- Faster, more efficient separations, low pressure, long column life
- Faster method development, simple optimization, broad selectivity
- Best use ammonia and formic acid or acetic acid modifiers (0.01-1.0 parts) in 100 parts methanol for LC/MS/MS compatibility.
- Analyte salts easily disassociated
- Complimentary to normal phase separations on polysaccharide CSPs
- Very useful in preparative purifications to replace hexane/ethanol.
 - Lower boiling point than heptane or hexane, higher evaporation rate
 - Less toxic
 - Higher evaporation rate

POLAR IONIC MODE

Mobile Phase Components

- **Methanol/Acid/Base Components:**

- Methanol, anhydrous
- Acid: anhydrous TFA, acetic acid, formic acid
- Base: TEA, DEA, NH₃
- Alternatives: ammonium formate or acetate

- **Ratio of Acid to Base**




- Controls only selectivity
- Rate: 4:1 to 1:4, most typical 2:1

- **Concentration of Acid + Base**




- Controls retention and selectivity
- Concentration range: 0.001 to 1.0 part per 100, most typical 0.10.

Salt Effects in the Polar Ionic Mode




Atrolactic acid, CHIROBIOTIC T2



Mianserin CHIROBIOTIC V2




Importance of Acid/Base Ratios In PIM

Example	CHIROBIOTIC V
Mobile Phase	MeOH/Acid/Base
100/0.1/0.1 Peak 1 – 6.2 min. Peak 2 – 7.4 min. Ratio: 1:1	
100/0.15/0.05 Peak 1 – 10.4 min. Peak 2 – 14.5 min. Ratio: 3:1	
100/0.05/0.15 Peak 1 – 3.4 min. Peak 2 – 3.6 min. Ratio: 1:3	

Note:
short
retention
times


Mianserin



Nitrogen on Mianserin group is
NH(+), COOH on V is (-)

Nitrogen on Mianserin
group as free amine, but
COOH on V is fully charged:
weak ionic interaction

Comparison of Four Basic Mobile Phase Types for CHIROBIOTIC CSPs




Temperature Effects

- Higher efficiency
- Shorter analysis time

Arginine

CHIROBIOTIC T,
25°C

30/70 MeOH/10 mM
NH₄OAc (pH4.0)
1.0 mL/min, ELSD



Arginine

CHIROBIOTIC T,
45°C


30/70 MeOH/10 mM
NH₄OAc (pH4.0)
1.0 mL/min, ELSD



Diacetyl-cysteine (N¹⁵)

CHIROBIOTIC T,
25°C


100/0.1% (v/w)
MeOH/NH₄OAc
1.0 mL/min, ELSD



Diacetyl-cysteine (N¹⁵)

CHIROBIOTIC T,
55°C

100/0.1% (v/w)
MeOH/NH₄OAc
1.0 mL/min, ELSD



Application of Temperature Effect in Loading

Study of Diacetyl-Cysteine (N¹⁵)


CHIROBIOTIC T, 250x4.6mm

Mobile Phase: MeOH/NH₄OAc (0.1 % v/w)

Detection: UV@230nm

Less Time!

Less mobile phase additive!



Practical Guides for CSP Screening



Aim Method Development Process

- To quickly identify a suitable column for selectivity, in the minimum number of experiments



Method Development Approaches

- Single column with multi-mobile phases (mainly CHIROBIOTIC phases)
- Multi-column switching with minimal number of simple mobile phases

Screening Strategy



Choose a small set of CSPs that:

- are broad-based to increase chances of success for a wide range of molecular types
- offer selectivity in a wide range of mobile phases for increased selectivity possibilities and sample solubility
- are complementary to each other, minimum overlap in selectivity

Complementary Separations

CHIROBIOTIC T versus CHIROBIOTIC R

Naproxen



T

R

Mobile Phase: 30/70: MeOH/0.1% TEAA, pH 4.1

Screening Strategy


Current industry trend is for generic screening methods with:

- a simple set of columns combining CHIROBIOTIC and Chiralcel/pak* CSP's or CHIROBIOTICS and CYCLOBOND
- a minimum number of solvents covering all five mobile phase types

*Chiracel and Chiralpak are the trademarks of Daicel.

Normal Phase vs Polar Ionic Mode[©]

Metoprolol




CHIRACEL OD[®]

Peak 1 – 11.9 min.

Peak 2 – 18.2 min.

Normal Phase

20/80/0.1: IPA/Hex/DEA




CHIROBIOTIC T[®]

Peak 1 – 15.36 min

Peak 2 – 17.11 min


Polar Ionic Mode

MeOH/0.1% ATFA



Normal Phase vs Polar Ionic Mode[©]

Alprenolol




CHIRACEL OD

Peak 1- 12.4 min.

Peak 2 - 16.4 min.

Normal Phase

20/80/0.1: IPA/Hex/TFA




CHIROBIOTIC V

Peak 1 – 7.69 min.

Peak 2 – 8.33 min.


Polar Ionic Mode

100/0.01/0.01: MeOH/HOAc/TEA




Cellulosic vs CHIROBIOTIC CSP's : RP vs PIM

Tolperisone




CHIRALPAK® AD-RH Peak 1 – 5.23 min
Peak 2 – 6.06 min

Reversed Phase 60/40: ACN/20mM Borate



CHIROBIOTIC V2 Peak 1 – 8.80 min
Peak 2 – 10.69 min

Polar Ionic Mode 100/0.1%: MeOH/ATFA



CHIRALPAK® is the registered trademark of Daicel Chemical Industries, Ltd

Purpose of Assay

Generic screening develops options for a variety of applications

- Quick analytical method (possibly suitable for later optimization and validation)
- Suitable for possible small scale prep
- Trace analysis of unwanted isomer
- Impurity profiling
- LC-MS method

Complementary Method Development

Cellulose & Amylose (CHIRALCEL/CHIRALPAK)	Macrocyclic Glycopeptides (CHIROBIOTIC)
Broad applicability over wide range of compound types	Broad applicability over wide range of compound types
Traditionally NP, but now also RP, POM (different columns)	Designed to provide selectivity in PIM, RP, NP (same column)
Chiral interaction sites via different chemistries & helical structure	Large number of chiral interaction sites
Most useful are AD, OD, AS, OJ	Most useful are V2, T, R, TAG
Coated phases	Chemically bonded

Complementary Method Development

CHIRALCEL/PAK:

- Compound must be in **neutral form** - interaction always non-ionic
- Separate samples into acids, bases and neutrals (neutrals can be screened with either acids or bases)

CHIROBIOTIC:

- Compound can be **ionised**, or can be a salt – ionic interactions are a key mechanism
- Same mobile phase screens are used for all samples, but can choose selective screening for acid, bases or neutrals

Note: Functional group on or near stereogenic center dictates whether analyte is acid or base

1. COLUMN INSTALLATION

CHIROBIOTIC™ columns are shipped in methanol. Before starting to use a new column, wash with 20 mL HPLC grade methanol at 1 mL/min. The column test standard, 5-methyl-5-phenylhydantoin, can be injected at this stage.

CYCLOBOND™ columns are shipped in IPA and should be washed with 30 mL HPLC grade water at 0.8 mL/min before starting the method development screen.

2. MOBILE PHASE CHOICE

No.	Mobile Phase	Composition (% v)
REVERSED PHASE MODE:		
1	MeOH/20mM NH ₄ OAc, pH 4.0	20/80
2	MeOH/20mM NH ₄ OAc, pH 6.0	20/80
3	ACN/20mM NH ₄ OAc, pH 4.0	30/70
4	ACN/20mM NH ₄ OAc, pH 6.0	30/70
POLAR IONIC MODE®:		
5	MeOH/HOAc/TEA*	100/0.1/0.1
POLAR ORGANIC MODE:		
6	ACN/MeOH/HOAc/TEA	95/5/0.3/0.2
CHIROBIOTIC PHASES ONLY:		
	If not progressing to normal phase, wash with MeOH at this stage to test and store the column.	100% MeOH
NORMAL PHASE MODE:		
7	EtOH/Hexane (or heptane, isohexane)	30/70
8	Washing cycle	100% EtOH
9	Column storage: CHIROBIOTIC	100% MeOH
	CYCLOBOND	100% IPA

* Use salts (NH₄O₂CCF₃ for bases, NH₄OAc for acids) when developing methods for prep.

3. COLUMN CHOICE AND RUN TABLE

Select your choice of columns from the list below. For a 6-column switching system, we recommend CHIROBIOTIC V2, T, R and CYCLOBOND I 2000, DNP and HP-RSP.

No.	Column Type (250x4.6mm)	1	2	3	4	5	6	7	8
I	CHIROBIOTIC V2	y		y		y		y	y
II	CHIROBIOTIC T	y	y		y	y		y	y
III	CHIROBIOTIC R		y	y	y	y		y	y
IV	CHIROBIOTIC TAG	y	y		y	y		y	y
V	CYCLOBOND I 2000	y		y			y		y
VI	CYCLOBOND I 2000 DNP	y		y			y	y	y
VII	CYCLOBOND I 2000 DMP	y		y			y	y	y
VIII	CYCLOBOND I 2000 HP-RSP	y		y	y		y		y

4. RUN CONDITIONS

Flow Rate:	1.0 mL/min.
Equilibration Time:	25 minutes
Run Time:	25 minutes
Temperature:	Ambient
Detector:	UV - 230nm
Sample:	1 mg/mL in MeOH

Notes

- The recommended protocol assumes the use of 250 x 4.6mm columns. For 100 x 4.6mm columns, use the same conditions at 0.5 mL/min.
- It is permissible to run straight from the reversed phase to the polar ionic mode®, and from the polar ionic mode to normal phase without an intermediate solvent wash.
- If any screening run results in a retention time less than 5 minutes, reduce the strength of the mobile phase and re-run. Aim for retention times from 5 to 20 minutes. In reversed phase mode reduce organic component, in polar ionic mode® or polar organic mode reduce acid/base concentration. Retention times can be later reduced in the optimization process.
- If a separation occurs in the polar ionic mode®, for a neutral molecule, change to 100% organic solvent (i.e. MeOH, EtOH or ACN).
- If the compound does not elute in reversed phase, increase the organic content to 40%. In the polar ionic mode®, increase the acid/base concentration up to 1.0/1.0. In the polar organic mode for the CYCLOBOND columns, increase the MeOH concentration up to 10%.

5. OPTIMIZATION PROCEDURES

Polar ionic mode® (CHIROBIOTIC phases only)	- Test alternative acid/base ratios (generally higher acid for basic molecules, higher base for acidic molecules) - To change acid/base to a volatile salt, use ammonium trifluoroacetate for basic compounds and ammonium acetate for acidic compounds at a concentration of 0.1wt%, adjust accordingly. Ammonium formate may be used as a compromise for both acidic and basic compounds.
Polar organic mode	- Eliminate MeOH - Test alternative acid/base ratios
Reversed phase mode	- Test smaller pH changes - Change organic to THF, ACN, MeOH - Change buffer type and buffer concentration - Change temperature
Normal phase mode	Change EtOH concentration

6. OPTIMIZING FOR MS DETECTION

CHIROBIOTIC: Use salts, as in Step 5, when using the polar ionic or polar organic modes.

CYCLOBOND: Use NH₄OH to replace TEA in polar organic mode, lower concentration by 50 to 75%.
Both Phases: Use ammonium acetate or ammonium formate when using in reversed phase.

7. RETESTING YOUR METHOD DEVELOPMENT COLUMNS

To ensure the selectivity performance of CHIROBIOTIC columns, periodically test with 5-methyl-5-phenylhydantoin in 100% MeOH. For testing CYCLOBOND columns, please refer to your CYCLOBOND Handbook.

ADVANCED SEPARATION TECHNOLOGIES

World Headquarters: 37 Leslie Court, Post Office Box 297, Whippany, NJ 07981 USA Tel: (973) 428-9080 Fax: (973) 428-0152
E-mail: info@astecusa.com www.astecusa.com

UK and Ireland Sales Office: 1 Blake Street, Congleton, Cheshire CW12 4DS UK Tel: +44 (0) 1260 276276 Fax: +44 (0) 1260 290067
E-mail: info@asteceuro.com www.asteceuro.com

©2006 Advanced Separation Technologies Inc. All rights reserved. Printed in USA.

Chiral Screening Protocol

1. Mobile Phase Composition:

No.	Mobile Phase	Composition (v/v)
1	MeOH/20mM NH ₄ OAc, pH 4.0	20/80
2	MeOH/20mM NH ₄ OAc, pH 6.0	20/80
3	ACN/20mM NH ₄ OAc, pH 4.0	30/70
4	ACN/20mM NH ₄ OAc, pH 6.0	30/70
5	MeOH/HOAc/TEA (PIM)*	100/0.1/0.1
6	ACN/MeOH/HOAc/TEA (POM)**	95/5/0.3/0.2
7	EtOH/Heptane	30/70
8	EtOH	100

*PIM: Polar Ionic Mode, CHIROBIOTIC

** POM: Polar Organic Mode, CYCLOBOND


Chiral Screening Protocol

2. Run Table:


No	Column Type	1	2	3	4	5	6	7	8
I	CHIROBIOTIC V2	y		y		y		y	y
II	CHIROBIOTIC T	y	y		y	y		y	y
III	CHIROBIOTIC R		y		y	y		y	y
VI	CHIROBIOTIC TAG	y	y		y	y		y	y
V	CYCLOBOND I 2000	y		y			y		y
VI	CYCLOBOND DNP	y		y			y	y	y
VII	CYCLOBOND DMP	y		y			y	y	y
VIII	CYCLOBOND HP-RSP	y		y	y		y		y

Note: Each CSP favors certain molecular types.

Typical Screening Results




Optimized CHIROBIOTIC V




MeOH/NH₄TFA; 100/0.02 w%

Chiralyzer® Optical Rotation




Dorzolamide *HCl Method Development








- Processed through method development screen:
- 3 CHIROBIOTIC phases: V2,T,TAG
- 3 CYCLOBOND phases: B-CD,DNP,HP-RSP.



**Dorzolamide
Hydrochloride
4S,6S**



**Dorzolamide
Hydrochloride
Related Compound A
4R,6R**


Spectrum	Column	mode	elution	File
	CHIROBIOTIC TAG	RP	No Retention	C:\asci\AC D_Chrom1 373.cdf
	CHIROBIOTIC TAG	PIM	Separation	C:\asci\AC D_Chrom1 384.cdf
	CHIROBIOTIC V2	RP	Separation	C:\asci\AC D_Chrom1 397.cdf
	CHIROBIOTIC V2	PIM	Separation	C:\asci\AC D_Chrom1 404.cdf
	CHIROBIOTIC T	RP	No Separation	C:\asci\AC D_Chrom1 416.cdf
	CHIROBIOTIC T	PIM	No Separation	C:\asci\AC D_Chrom1 423.cdf
	Cyclobond I 2000	RP	No Retention	C:\asci\AC D_Chrom1 435.cdf
	Cyclobond I 2000	POM	Unknown	C:\asci\AC D_Chrom1 442.cdf
	Cyclobond 2000 HP-RSP	RP	No Retention	C:\asci\AC D_Chrom1 461.cdf
	Cyclobond 2000 HP-RSP	POM	Unknown	C:\asci\AC D_Chrom1 485.cdf
	Cyclobond 2000 DNP	RP	No Separation	C:\asci\AC D_Chrom1 497.cdf
	Cyclobond 2000 DNP	POM	Unknown	C:\asci\AC D_Chrom1 504.cdf

Chromatographic Results

- PIM mode best!

Peak 1 (R_t1): 19.81 min.

Peak 2 (R_t2): 22.71 min.



CHIROBIOTIC V2, 250x4.6mm I.D., 5 μ particles (15024AST)

100:5 MeOH:H₂O, 3.81mM NH₄TFA (or 0.05% NH₄TFA)

22°C, 0.3mL/min., UV-254nm

Inj. Vol. 10 μ L @1.0 mg/mL in MeOH

Racemic Switches: Optimized Conditions

CHIROBIOTIC Phases

Final Results

Compounds	Mobile Phase	Column	k _{1'}	α	R _s
Albuterol	100/0.1w%, MeOH/NH ₄ TFA	T2	2.5	1.3	2.8
Amlodipine	100/0.1w%, MeOH/NH ₄ TFA	V2	3.2	1.1	1.5
Bupivacaine	100/0.1w%, MeOH/NH ₄ TFA	V2	0.8	1.4	2.5
Citalopram	100/0.1w%, MeOH/NH ₄ TFA	V	4.8	1.1	1.5
Clenbuterol	100/0.1w%, MeOH/NH ₄ TFA	T2	2.1	1.2	2.5
Fluoxetine	100/0.1w%, MeOH/NH ₄ TFA	V2	1.5	1.2	2.5
Formoterol	100/0.6/0.4, MeOH/HOAc/TEA	T2	2.6	1.2	1.5

Racemic Switches: Optimized Conditions

CHIROBIOTIC Phases con't

Final Results

Compounds	Mobile Phase	Column	k _{1'}	α	Rs
Ibuprofen	10/90, THF/20mM NaCitrate, pH6.3	V	0.8	1.24	1.5
Isoproterenol	100/0.1/0.1, MeOH/HOAc/TEA	T2	5.5	1.33	3.3
Ketoprofen	20/80, MeOH/0.1% TEAA, pH 6.0	R	2.6	1.60	3.0
Lercanidipine	100/0.02w%, MeOH/ NH ₄ TFA	V	1.0	1.18	1.5
Lorazepam	100 % MeOH	T	0.4	4.01	11
Methylphenidate	100/0.1w%, MeOH/NH ₄ TFA	V2	2.4	1.53	4.0
Metoprolol	100/0.1w%, MeOH/NH ₄ TFA	T	4.5	1.14	2.0
Mianserin	100/0.1w%, MeOH/NH ₄ Formate	V2	0.5	1.86	2.8

Racemic Switches: Optimized Conditions

CHIROBIOTIC Phases con't

Final Results

Compounds	Mobile Phase	Column	k _{1'}	α	Rs
Mosapride	30/70, ACN/0.1% TEAA, pH 4.1	V	3.0	1.32	3.3
Naproxen	10/90, THF/NaCitrate, 20mM, pH 6.3	V	1.4	1.22	1.5
Nicardipine	100/0.02w%, MeOH/NH ₄ TFA	V	0.7	1.74	4.5
Oxamniquine	100/0.1w%, MeOH/NH ₄ TFA	V2	2.7	1.23	2.2
Oxazepam	100% MeOH	T	0.5	4.31	12
cis-Permethrin	99.5/0.5, Hex/EtOH	TAG	4.6	1.22	2.5
Propranolol	100/0.1w%, MeOH/NH ₄ TFA	T	4.2	1.14	1.9
Pseudo-ephedrine	100/0.1w%, MeOH/NH ₄ TFA	T2	2.7	1.13	1.7

Racemic Switches: Optimized Conditions

CHIROBIOTIC Phases con't

Final Results

Compounds	Mobile Phase	Column	k ₁ '	α	Rs
Sotalol	100/0.2/0.1, MeOH/HOAc/TEA	T	3.5	1.12	1.5
Terbutaline	100/0.1w%, MeOH/NH ₄ TFA	T2	2.2	1.92	7.0
Thalidomide	100% MeOH	V	0.7	3.44	7.0
Tolperisone	100/0.1w%, MeOH/NH ₄ TFA	V2	1.8	2.33	2.7
Trimipramine	100/0.1w%, MeOH/NH ₄ TFA	V2	1.7	1.28	2.3
Warfarin	40/60, EtOH/0.1% TEAA, pH 4.1	V	1.3	1.45	3.0

CYCLOBOND Screening Hits Leading to Baseline Separations

Compounds	Columns*	Compounds	Columns*
Chlorpheniramine	β -CD	Sertraline	HP-RSP/ β -CD
Chlorthalidone	HP-RSP/ β -CD	Thioridazine	β -CD
Epinephrine	HP-RSP	Thalidomide	DNP
Methadone	HP-RSP	Tramadol	DMP
Miconazole	HP-RSP	Warfarin	β -CD
Nefopam	DNP/HP-RSP	Zopiclone	β -CD
Omeprazole	DMP		
Oxazepam	DNP		
Propranolol	β -CD		
Pseudo-ephedrine	β -CD		

***CYCLOBOND I 2000: β -CD/HP-RSP/DNP/DMP – 40 % hit rate!**

Note: Compounds from racemic switch list.


Final Results

Racemic Switches: Optimized Conditions CYCLOBOND Phases


Compounds	Mobile Phase	Column	k _{1'}	α	Rs
Chlorpheniramine	10/90, ACN/0.1% TEAA, pH4.1	β -CD	4.7	1.14	1.5
Chlorthalidone	15/85, ACN/10mM NH ₄ OAc, pH4.0	HP-RSP	1.4	1.22	1.6
Epinephrine	25/75, ACN/20mM NH ₄ OAc, pH4.1	HP-RSP	7.1	1.16	2.0
Methadone	20/80, ACN/ 10mM NH ₄ OAc, pH 3.6	HP-RSP	1.5	1.24	1.7
Miconazole	25/75, ACN/10mM NH ₄ OAc, pH 4.0	HP-RSP	4.1	1.18	1.6
Nefopam	25/75, ACN/20mMNH ₄ OAc,pH 4.1	DNP	0.5	1.52	2.2
Omeprazole	100/0.4/0.1, ACN/HOAc/NH ₄ OH	DMP	3.3	1.41	2.4
Sertraline	30/70, ACN/20mM NH ₄ OAc, pH4.1	HP-RSP	5.4	1.19	1.7
Thioridazine	20/80, ACN/0.1% TEAA, pH4.1	β -CD	6.9	1.21	2.5
Tramadol	20/80, ACN/20mM NH ₄ OAc, pH5.5	DMP	1.6	2.20	2.8
Zopiclone	95/5/0.3/0.2,ACN/MeOH/HOAc/TEA	β -CD	1.0	1.13	1.5

Options from Results of Dual Screen


Omeprozole (Nexium)




CHIROBIOTIC R (NP)



CHIROBIOTIC R (RP)



CYCLOBOND I 2000 DMP (POM)



Statistics

- Of forty racemic switch drugs tested in this study:
 - The success rate for CHIROBIOTIC phases is 80%.
 - The success rate for CYCLOBOND phases is 40%.
 - Eight compounds can be separated by both types (20%).
 - Resolution range:
 - CYCLOBOND Phases: 1.5 to 2.8
 - CHIROBIOTIC Phases: 1.5 to 11.0


Best Positive Screening Results

Column: **CHIROBIOTIC V2, 250x4.6mm**

Mobile Phase: 40/60, ACN/10 mM NH₄OAc, pH 3.8

Flow Rate: 0.8mL/min

UV: 380 nm




Column: **CHIROBIOTIC V2, 150x4.6mm**

Mobile Phase: 50/50, MeOH/5 mM NH₄OAc, pH 3.5

Flow Rate: 0.8mL/min

UV: 380 nm

Optimized Method for LC/MS application



Method Development Screen

- Monitor column performance *regularly*
- Store columns in correct solvents (free of additives):
 - CHIROBIOTIC: 100% MeOH
 - CYCLOBOND: 100% 2-PrOH
 - CHIRALCEL/PAK: Hexane/2-PrOH, 90/10

Isocratic or Gradient?

- Gradients suitable for NP, unsuitable for RP
- Recent papers suggest that gradient method development not much faster and has no greater success rate
- In PIM, ammonium formate gradient possible from 0.01% to 0.5% in methanol to replace HOAc which has a very broad window

Results of Complementary Generic Screening (1)

Previous studies reported:

- Screening* four polysaccharides over 5 (NP, POM) mobile phases provided selectivity for 87% of a set of 53 compounds: the same study tested three CHIROBIOTIC phases over just 2 (RP, PIM) mobile phases with selectivities of 65%.
- It was noted that, together, they provided a 96% success rate, with increased selectivity for certain families of compounds on the CHIROBIOTIC CSPs confirming the complementary effect.


* M E Andersson et al, J Chrom A, 1005 (2003) 83

Sales Tools Available

- Product Guide
- Handbooks (4)
- New Product Bulletins
- New Applications now available
- Seminar Program (CD)
- Method Development program
- Website: www.sigma-aldrich.com
- E newsletter

Chiral Analytical Services From Supelco

- Chiral column screening; HPLC & GC
- Method optimization
- Small-scale purification
- Achiral method development services as well



Mission

- Provide fast, efficient and effective analytical & purification services
- Provide consistent, informative and friendly communication – before, during and following each study
- Provided timely assistance with legal and technical issues
- Present a streamlined, easy to use system
- Offer fair pricing structure

Supelco Chiral Services Laboratory



Supelco Chiral Services Laboratory



Key Features

- Fully confidential
- Fast service – goal is within 48 hours of sample receipt for screening projects, time variable for method development and purification
- Full report given upon completion of project
- Elution order determination by measurement of optical rotation
- Purified material to the customer's specification

GC Chiral Column Screening

- Manual exploration of 3-4 GC column chemistries
- Samples that require derivatization are verified by GC-MS




HPLC Method Optimization

- Certain customers may be interested in obtaining enantiomerically-pure materials
- Following successful chiral screening or successful reproduction of methods provided by customer, a method optimization and loading study is conducted to establish the most efficient means of chromatographically purifying the materials
- The output of a loading study is approximately 100 mg of purified material and the methodology utilized to obtain such material

Output



595 North Harrison Rd.
Bellefonte, Pennsylvania 16823
(814) 359-3041



Conditions:

Column: CHIROBIOTIC V2, 25 cm x 4.6 mm I.D., 5 μ m particles (15024AST)
Mobile Phase: 5:100, H₂O:0.05w% NH₄TFA in MeOH
Temperature: 22 °C
Flow Rate: 0.3 mL/min
Detection: UV at 254 nm
Injection Volume: 10 μ L
Sample: 1.0 mg/mL in MeOH

Peak 1 retention time (R_{t1}): 19.81 min.

Peak 2 retention time (R_{t2}): 22.71 min.

Recommendation for Optimization/Conclusions:

1. The best resolution of the chiral enantiomers Dorzolamide Hydrochloride and Dorzolamide Hydrochloride Related Compound A was observed on the CHIROBIOTIC V2. The aforementioned method was used for confirmation and optimization purposes.
2. The buffer salt ammonium trifluoroacetate was employed to sharpen the peaks, and thus, increase resolution.

We are committed to the success of our Customers, Employees and Shareholders through leadership in Life Science, High Technology and Service.
The SIGMA-ALDRICH Family

- Client specified information
- Detailed report
- Full conditions
- Recommendations

} Project-dependent recommendations by our experts for advice in optimization, scale-up, etc.

CONCLUSIONS

- Chirobiotic phases are complimentary to cellulose and amylose CSPs
- Polar ionic mode is best mobile phase for ionic racemates
- Reversed phase can be run on all columns for any type of analyte
- Mobile phase can be run sequentially on all Supelco phases.
- Cyclobond phases have solved chiral separation problems not possible on any other CSP especially the Cyclobond HP-RSP.
- Generic screening has proven to be the most efficient methodology for chiral selectivity can be streamlined to six columns and 3-4 mobile phases.

The Supelco Chiralyser

Optical rotation is a very useful tool when dealing with chiral entities. Up to now the major problem has been the sensitivity and ease of use of the commercially available devices. The Astec Chiralyser has come a long way in solving those problems. Two available scientific tools have made this a reality.

1. The availability of single wavelength LED's which allowed the choice of the best, most sensitive wavelength (430 nm) as a light source. In addition, from the LED technology you get long life > 100,000 hours.
2. Use of the Faraday effect that allows nulling of the parallel magnetic field electronically so that the entire sample in the cell can be read instantaneously.

What are the benefits of a Chiralyser in the area of chiral separations?

1. Validates peaks as (-) and (+). Eliminates peaks that are achiral contaminants.
2. Identifies a reversal of elution order. A common occurrence between CSP's and utilizing different mobile phase conditions.
3. Identifies enantiomeric pairs. A useful determination as the introduction of an achiral guard column can easily complete resolution if diastereoisomers are overlapping.
4. Detects analytes with no UV unlike circular dichroism.
5. System is easy to set up and economical to use, no complicated software.

Below is a list of compounds that have been run in our lab. It demonstrates the reversal of elution order problem and the variety of chiral molecules that have been assayed.

Compound	Column	Mobile Phase	First peak
5,5 Hydantoin	All Chirobiotics	MeOH	(-)
Oxazepam	V/T/R	MeOH	(+)
*N-Amine	V/V2	100/0.1w%,MeOH/ATFA	R(+)
*N-Amine	V	30/70,MeOH/TEAA,4.1	S(-)
Methodone	V2	100/0.1w% MeOH/ATFA	R(-)
Methodone	HP-RSP	20/80, ACN/NH ₄ OAc,3.6	R(-)
*Propranolol	T/T2	100/0.1w% MeOH/ATFA	S(-)
*Propranolol	TAG	100/0.1w% MeOH/ATFA	R(+)
Terbutaline	T/T2/TAG/V	100/0.1w% MeOH/ATFA	(-)
Metoprolol	T	100/0.1w%,MeOH/ATFA	R(+)
*Bendroflumethiazide	V	10/90,THF/NH ₄ NO ₃	(+)
*Bendroflumethiazide	T	30/70, MeOH/H ₂ O	(-)
Naproxen	V	10/90, THF/NaCitrate	S(+)
Naproxen	R	20/80, MeOH/TEAA,5.5	S(+)
Ibuprofen	V	10/90, THF/NaCitrate	R(-)
Albuterol	V/T	100/0.1w%,MeOH/ATFA	(-)

* Reversed elution order

*Mianserin	V	100/0.1w%,MeOH/ATFA	(-)
*Mianserin	T	100/0.1w%,MeOH/ATFA	(+)
Mandelic acid	T/R	30/70, MeOH/pH 4.5	S(+)
Mandelic acid	T/R	100/0.1w%,MeOH/ATFA	S(+)
Bupvacaine	V/V2	100/0.1w%,MeOH/ATFA	S(-)
Nicardipine	V/V2/T	100/0.1w%,MeOH/ATFA	(-)
*Ritalin	V/V2	100/0.1w%,MeOH/ATFA	(-)
*Ritalin	T2	100/0.1w%,MeOH/ATFA	(+)
Omeprazole	R	40/60, EtOH/Heptane	R(+)
Omeprazole	R	30/70, MeOH/HOAc,4.1	R(+)
Amphetamine	V2	100/0.05w%,MeOH/ATFA	S(+)
Methamphetamine	V2	100/0.05w%,MeOH/ATFA	S(+)
Dextro/Levorphanol	V/V2	100/0.1w%,MeOH/ATFA	D(+)
Dextro/Levo methorphan	V/V2	100/0.1w%,MeOH/ATFA	D(+)
Butorphanol	T/T2	100/0.1w%,MeOH/ATFA	D(+)
Ketoprofen	R	100/0.02w%,MeOH/HOAc	S(+)
*a-Me a-Ph succinimide	R	20/80,EtOH/Hex	(-)
*a-Me a-Ph succinimide	TAG/T/V	20/80,EtOH/Hex	(+)
*Warfarin	B-CD	100/0.3/0.2, ACN/HOAc/TEA	S(-)
*Warfarin	V	30/70, ACN/NH ₄ OAc,4.1	R(+)

* Reversed elution order