

2007 VOLUME 2 NUMBER 4



INNOVATION @ WORK

ABC TRANSPORTERS AND CANCER DRUG RESISTANCE

NF-_KB AND INFLAMMATION

DNA DAMAGE AND REPAIR


NEW PRODUCTS

Cultured HeLa cells in metaphase stage of mitosis.

Cellular Mechanisms and Cancer



sigma.com



NEW from Sigma Your comprehensive gene search tool

Your Favorite Gene[™] is a new search tool that matches your gene of interest to thousands of research products available from Sigma.

Flexible search options:

- Use a wide variety of terms to locate genes
- Refine or expand your search by pathways, function & more
- Retrieve detailed gene information

Unparalleled product coverage matched to your gene of interest:

- 150,000 shRNAs
- 4,000 antibodies, proteins & kits

Integrated links to public databases:

- Compendia Bioscience
- Human Protein Reference Database

- 1,000 bioactive small molecules
- 725,000 siRNAs coming soon!
- NCBI Entrez Gene
- Cytoscape[™]





Introduction

The medical definition of cancer appears simple and straightforward. According to the U.S. National Cancer Institute, cancer is "(a) term for diseases in which abnormal cells divide without control."¹ Behind this basic definition is a complex and unpredictable spectrum of over 100 types of cancer.

The human aspect of cancer cannot be completely separated from the scientific research. The World Health Organization recognizes cancer as a leading cause of death worldwide, and emphasizes prevention and early detection as crucial to reducing the global burden of the disease.² The American Association for Cancer Research established its Scientist-Survivor program in 1999 to encourage communication between patients, patient advocates, and leading scientists in the field.³ Newspapers, magazines, and other media sources announce breakthrough discoveries to the public, increasing awareness that cancer is not an individual disease but a collection that has no single initiating event or defined evolution.

Cancer results from a cascade of abnormal cell reactions. When a cellular mechanism goes wrong, the resulting damage, if not repaired, may contribute to a cell's evolution into malignancy. Because cancers begin with a single cell, cancer investigators use genomics, proteomics, and signaling techniques to determine and evaluate cellular changes and contributing cause and effect. Discoveries such as the correlation between the human papilloma virus and cervical cancer encourage the scientific community to seek similar breakthroughs for other cancer types. The understanding of cellular mutations and signaling pathways involved in mutagenesis and abnormal cell function has been used to screen potential new drugs with more efficacy and/or less toxicity.

It's impossible to comprehensively review current cancer research; the amount of information is enormous and the rate of discovery is increasing. In the preface to his book "The Biology of Cancer", Robert Weinberg writes "...we are deluged with a vast amount of genetic, biochemical, and cell biological information about cancer development, far more almost than any human mind can assimilate and comprehend."⁴ For this issue of BioFiles we have selected three aspects of cancer biology to review.

- The exploitation of ABC transporter proteins by cancer cells to export chemotherapeutic drugs
- The activation of NF– κ B in response to inflammation and its role in cancer progression
- Genetic damage, mutagenesis, and cellular repair processes

Innovation @ Work - Innovative products from Sigma for genomics and proteomics studies, including gene silencing, protein expression profiling, phosphopeptide enrichment, and whole genome amplification, are also highlighted.

References:

1. World Health Organization, Cancer Fact Sheet No. 297, Feb. 2006, www.who.int/mediacentre/factsheets/fs297/en/index.html

- 2. National Cancer Institute, U.S. National Institutes of Health, www.cancer.gov
- 3. American Association for Cancer Research, www.aacr.org

4. The Biology of Cancer, Robert A. Weinberg, Garland Science, Taylor & Francis Group, LLC, New York, NY (2007)

Cover image: Cultured HeLa cells labeled with anti-tubulin antibody and propidium iodide (to label the DNA) in metaphase stage of mitosis. Photographed by Dr. K. G. Murti, St. Jude Children's Research Hospital, Memphis, TN.

Technical content: Libby Yunger, Ph. D., Chloe McClanahan, B. Sc., Vicki Caligur, B. Sc



2007 Volume 2 Number 4

Table of Contents

Cellular Mechanisms and Cancer

Role of ABC Transporters and Multi-Drug Resistance Reversal Using

RNA Interference	7
MDR-substrate Anticancer Agents	9
ABC Transporter Membrane Proteins	10
Antibodies to ABC Transporters	.11
Compounds for MDR Detection	13

NF-_KB and Inflammation

nflammation and Cancer:	
Гhe NF-кВ Connection	14
Natural Product Inhibitors of	
NF-ĸB Activation	16
Antibodies and Kits	18

DNA Damage and Repair

ONA Damage and Repair Mechanisms	.20
DNA Alkylating Agents	. 24
DNA Crosslinking Agents	. 26
DNA Synthesis Inhibitors	. 27
DNA Repair Enzymes and Antibodies	

New Products

Г

Antitumor Agents	30
Chemopreventative Agents	
Enzymes	
Multi-Drug Resistance	32
Tumor Growth Regulation	



MISSION® TRC ShRNA

2

MISSION[®] shRNA Human Tumor Suppressor Gene Panel


Silencing Tumor Suppressor Genes Using MISSION® TRC shRNA

RNAi methodologies have been utilized in small- to large-scale screening projects, allowing researchers to perform gene-silencing experiments in timeframes and target cells not previously possible. While siRNA screens have become fairly common, large-scale screening with shRNA is still evolving. Advantages of shRNA-based experiments include long-term knockdown and viral delivery to nontransfectable cell types.

The MISSION TRC shRNA Human Tumor Suppressor Gene Family Set is a panel of lentiviral-based shRNAs consisting of approximately 75 gene sets, each represented by 3–5 individual shRNA clones. The shRNA constructs included in the tumor suppressor set allow you to identify genes involved in cell survival by assisting in making a cell more sensitive, or more resistant, to an added chemotherapeutic agent through the knockdown of gene expression and resulting function.

Scientists at Sigma have utilized the MISSION TRC shRNA Human Tumor Suppressor Gene Family Set in order to address biologically relevant questions while simultaneously developing screening strategies that can be used by researchers in the field. The tumor suppressor gene family set delivers a potential pharmacogenomics solution, allowing you to profile the gene expression of a cancer cell line and compare it against the response of the cell line to chemotherapeutic treatment.¹

In our study of paclitaxel resistance by the human lung cancer cell line, A549, various transcripts were down-regulated in the lung cancer cells using the tumor suppressor set. Several genes that are known to correlate with paclitaxel resistance, including *p53* and *VHL*, were identified by the panel. In addition, some genes, when down-regulated by these shRNAs, led to increased sensitivity of the cells to paclitaxel treatment.



Features of the pLKO.1-Puro vector allow for transient or stable transfection of the shRNA as well as production of lentiviral particles. Unlike adenovirus or murine-based MML or MSCV retroviral systems, lentiviral-based particles permit efficient infection and integration of the specific shRNA construct into differentiated and nondividing cells, overcoming low transfection and integration difficulties when using these cell lines, pLKO.1 allows for long-term knockdown and phenotypic observation and transduction of difficult or sensitive cell lines.

Features and Benefits

- Rapidly validate and compare results in multiple model systems – spend more time analyzing your data and less time gathering it
- Simplified, ready to transduce lentiviral format – straightforward shRNA protocols let you progress to new cell line studies faster
- shRNA enables you to down-regulate genes long-term – easily transfer your experiments from in vitro to in vivo

1. Ji, D., Deeds, S.L., and Weinstein, E.J., shRNAs Targeting Tumor Suppressor Genes, *Gen. Eng. News*, **27**, 26-27 (2007).

Ordering Information


Cat. No.	Product Name	Pack Size
SH0531	MISSION TRC shRNA Human Tumor Suppressors	1 set
SHP001	MISSION Lentiviral Packaging Mix	0.25 mL 1.7 mL
H9268	Hexadimethrine bromide	5 g 10 g 50 g 100 g
P9620	Puromycin Ready Made Solution, 10 mg/mL in water	10 mL
N4888	Nutrient Mixture F-12 Ham with sodium bicarbonate, without L-glutamine, sterile-filtered liquid, cell culture tested	500 mL 6 × 500 mL 24 × 500 mL
QR0200	Quantitative RT-PCR ReadyMix™	1 kit
T7402	Paclitaxel from <i>Taxus</i> brevifolia	1 mg 5 mg 25 mg

For more information, please visit our Web site at sigma.com/rnai.




Panorama[®] Antibody Arrays

How the Panorama Antibody Arrays Work



Induction with Retinoic Acid in Mouse F9 Cells



Using the Panorama Antibody Microarray – Cell Signaling Kit, proteins from mouse F9 cells are observed to be upregulated or down-regulated in expected fashion following stimulation with retinoic acid. F9 cells were treated for 96 hours with all-trans-retinoic acid (10^{-7} M). Extracts were prepared from untreated and treated cells using Extraction/ Labeling Buffer and labeled with Cy3 and Cy5, respectively. A mixture containing equal amounts of each labeled extract (5 µg/mL) was incubated on the array as described in the kit protocol. The same slide is shown at the two fluorescence emission wavelengths for Cy3 and Cy5. (Note that the blue background for Cy5 represents the normal and unavoidable interference of the nitrocellulose membrane.) Changes in expression level were confirmed by immunoblot, as indicated.

Protein Expression Profiling

Features

- Rapid, multiplex analysis of protein expression
- Profile 224 protein targets in half a day
- Antibodies react with human, mouse and rat species
- Comprehensive kits include buffers, extraction reagents, analysis software, and labware

High Throughput Protein Expression Analysis

Panorama Antibody Microarrays are designed to identify the changes in expression of multiple targets in a biological sample with a single assay. During cellular processes, mRNA undergoes a number of alternative processing steps prior to and following translation so there is often poor correlation between mRNA expression profiles and protein expression.¹ Antibody arrays provide an effective solution for analyzing biologically-relevant events at the proteome level. Highly specific antibodies are arrayed in duplicate on nitrocellulose slides compatible with most DNA array scanners. Cy™3/Cy5-based detection chemistries

may be used with cell or tissue extracts. Array antibodies are available separately for additional downstream analyses.

Applications

- Sample response to external treatments (e.g., starvation, addition of growth factors, serum, drugs/biomolecules)
- Comparison of different cell or tissue samples (clinical specimens, tissue, heart, brain, lung, etc.)
- Determine off-target effects in gene silencing (RNAi) studies

Components

Panorama Antibody Slides	QuadriPERM [®] Cell Culture Vessel
Extraction/Labeling Buffer	Protease Inhibitor Cocktail
Phosphatase Inhibitor Cocktail II	Benzonase® Ultrapure
Array Incubation Buffer	Collection Tubes, Polypropylene
Phosphate Buffered Saline, pH 7.4, with TWEEN® 20 (Washing Buffer)	SigmaSpin [™] Post-Reaction Clean-Up Columns
Panorama Antibody List on diskette including analysis software	

1. Gygi, S.P., et al., *Mol. Cell Bio.*, **19**, 1720-1730 (1999)

Ordering Information

Cat. No.	Product Name	Pack Size
CSAA1	Panorama Antibody Microarray - Cell Signaling	1 kit
GRAA2	Panorama Antibody Microarray - Gene Regulations	1 kit
MPAA3	Panorama Antibody Microarray - MAPK & PKC Pathways	1 kit
PPAA4	Panorama Antibody Microarray - p53 Pathways	1 kit
XP725	Panorama Antibody Microarray - XPRESS Profiler725	1 kit

For more information, please visit our Web site at sigma.com/arrays.

3





4

GenomePlex[®] Single Cell Whole Genome Amplification Kit


Whole Genome Amplification from a Single Cell

The GenomePlex Single Cell Whole Genome Amplification Kit (Cat. No. WGA4) is designed to amplify the genome of a single cell. This rapid and straightforward method provides million-fold amplification yielding microgram quantities of genomic DNA from a single cell. Traditional single cell whole genome amplification methods yield insufficient quantities with significantly biased representation. In contrast, the GenomePlex technology provides enhanced amplification efficiency by using WGA DNA polymerase


The Single Cell WGA kit includes all of the reagents necessary for cell lysis and subsequent whole genome amplification. Single cells can be isolated by fluorescence-activated cell sorting (FACS). laser capture microdissection (LCM), dilution, or any other applicable method. Single Cell WGA has been successfully applied to Comparative Genome Hybridization (CGH), STR analysis of amniocentesis samples, and genomic analysis of in vitro fertilized embrvos.

Human Control Chromosome 3

Human Kidney Tumor Chromosome 3



Human leukemia U937 cells were isolated using flow cytometric analysis and sorting (FACS), lysed, and amplified using the GenomePlex Single Cell Whole Genome Amplification Kit. The DNA was then purified with the GenElute™ PCR Cleanup Kit. An estimated million-fold amplification from the WGA process resulted in a final yield ranging from 5.4-6.2 μ g. The Single Cell WGA Kit produces consistent yield and size range as visualized by a 1% agarose gel



DNA from normal and tumorigenic human kidney cells were amplified using the GenomePlex Single Cell WGA Kit. Amplified material was hybridized to metaphase BAC arrays to determine chromosomal karyotype. Data that falls outside of the red line indicates chromsomal loss, while data that continues past the green line suggests chromosomal amplification. As expected the control sample demonstrated a balanced chromosomal copy number. Chromosome 3 for the amplified kidney tumor single cell displayed under representation as depicted by the red bar. These results match previous microarray work using an amplification method that took three days. GenomePlex technology accurately amplifies genomic material down to single cell resolution. Data is courtesy of Dr. Michael Speicher from the Institute of Human Genetics, TU Munich.

Cat. No.	Product Name	Pack Size
WGA2	GenomePlex Complete Whole Genome Amplification (WGA) Kit	10 reactions 50 reactions
WGA4	GenomePlex Single Cell Whole Genome Amplification Kit	10 reactions 50 reactions
WGA3	GenomePlex WGA Reamplification Kit	50 reactions
WGA5	GenomePlex Tissue Whole Genome Amplification Kit	10 reactions 50 reactions
NA1020	GenElute PCR Clean-Up Kit	1 kit
G1N10	GenElute Mammalian Genomic DNA Miniprep Kit	1 kit
S4438	SYBR® Green JumpStart™ Taq ReadyMix™	100 reactions 500 reactions

To see Sigma's complete line of Whole Genome Amplification products, visit our Web site at sigma.com/wga.


Ordering Information

GenomePlex[®] Complete Whole Genome Amplification Kit

Amplification of Genome-Representative DNA from Limited Starting Material


The GenomePlex Complete Whole Genome Amplification Kit (**Cat. No. WGA2**) contains everything required for whole genome amplification including an optimized enzyme, WGA DNA Polymerase. The WGA DNA Polymerase provides increased accuracy in amplification, as evidenced by producing no amplicon in the negative control reactions. WGA has been used in a variety of applications, and is suitable for use with purified genomic DNA from a variety of sources including blood cards, whole blood, buccal swabs, tissue, soil, plant, and serum. GenomePlex WGA uses nanogram quantities of starting genomic DNA, which after PCR yields on average 10 µg of amplified DNA. After purification, the WGA product can be analyzed in a manner similar to any genomic or chromosomal DNA sample. A number of downstream applications may be performed including TaqMan® assays, CGH analysis, SNP analysis, and sequencing.

Achieve Robust Amplification Representative of the Original Input Genome



WGA was performed on increasing concentrations of human genomic DNA. Amplification product can be detected on an agarose gel with as little as 10 pg of input DNA. Optimal performance is found with 1 to 10 ng of starting material. Increasing the amount of input DNA to 100 ng is not recommended.

Lane M:	DNA Marker	Lanes 7,8:	100 pg DNA
Lanes 1,2:	no template Lanes	Lanes 9,10:	1 ng DNA
Lanes 3,4:	1 pg DNA Lane	Lanes 11,12:	10 ng DNA
Lanes 5,6:	10 pg DNA Lanes	Lanes 13,14:	100 ng DNA



GenomePlex WGA was performed on genomic DNA isolated from HT29 colon carcinoma cells and from a healthy human male. 2.5 µg of WGA product was labeled with Cy[™]3 or Cy5 dye using the Genomic DNA Labeling Kit PLUS (Agilent). The entire labeled sample was loaded onto an Agilent Human Genome CGH Microarray 105A. Specific activities were between 28 and 43 for all samples, and always within 50% of samples being compared. The dye swaps (**A & B**) demonstrate that there was no bias in the DNA labeling and the aberrations detected are consistent with the HT-29 karyotype.


Looking for a high-throughput system for the rapid and highly representative amplification of genomic DNA from trace amounts of starting material? Visit **sigma.com/wgaautomation** for automated WGA protocols and methods.

NEW! Technical note on Agilent Array CGH with WGA, visit sigma.com/wgacgh.

5



PhosphoProfile[™] Phosphopeptide Enrichment Kit



Workflow highlighting the alternate use of trypsin spin columns (Cat. No. TT0010) for digestion followed by selective enrichment of phosphopeptides from sample.



The Complete Solution Kit 24 samples up to 25 nmoles phosphopeptide each.

Features and Benefits

- Mass spec compatible save time, reduce sample handling and potential loss
- Proteomics Grade Trypsin for clean and complete digests
- Phos-Select[™] Gallium Spin Columns – high capacity, novel Ga⁺³ silica media for fast, unbiased capture and recovery
- Controls validate your process for confidence and reporting
- Buffers enzyme reaction, binding, washing and elution formulations optimized for robust performance
- Consumable equipment included matched equipment means no risk of sample loss, additional purchases, or waste

Phosphoproteomics — Phosphopeptide Enrichment

Matrix Assisted Laser Desorption/Ionization, Time of Flight (MALDI-TOF) or Electrospray Ionization (ESI) Mass Spectrometry of phosphopeptides from tryptic protein digests are powerful tools for characterization and identification of phosphorylation sites. A combination of Iow intrinsic abundance, inefficient ionization, and/or signal suppression of most phosphopeptides may limit or even prevent detection, unless the phosphopeptide content is significantly enriched prior to analysis. This kit conveniently includes all materials needed to enrich phosphopeptides from digested samples in a robust and unbiased manner (see binding comparison). The phosphopeptide capture matrix is a novel Ga+³ chelate silica based on a proprietary nitriloacetic acid (NTA) analog. The silica beads are approximately 20 microns in diameter with a pore size of 1,000 Angstroms. The matrix is packed into spin columns for easy, microscale affinity capture of phosphopeptides.

Binding Comparison: Demonstrating enrichment of phosphopeptides for IMAC technologies.

Standard phosphopeptides representing the three most common sites of phosphorylation (phosphoserine, phosphothreonine, and phosphotyrosine) were used. The lyophilized solids were first dissolved in water, and an approximately equimolar mixture of the peptides was formulated. Each phosphopeptide was added to a BSA digest at a weight ratio of ~1.7% to produce a total phosphopeptide content of ~5% by weight. Quantitation results are given in the table below. Note that competitor A, B, and C technologies were biased in selecting Peptide 2, while Sigma's technology bound and eluted the peptides in approximately the same ratio as applied to the column.

	Recovery of phosphopeptide standards				
Kit	1	2	3	Total	Specificity*
Sigma	59%	52%	74%	59%	50%
Competitor A	6%	19%	11%	13%	28%
Competitor B	39%	56%	17%	42%	28%
Competitor C	37%	65%	37%	46%	25%

Performance summary of the IMAC technologies tested within this study.

*Specificity was measured as a percentage of the total HPLC peak area corresponding to phosphorylated peptides that appeared in the elution fractions.

Ordering Information

Cat. No.	Product Name	Pack Size
PP0410	PhosphoProfile I Phosphopeptide Enrichment Kit	1 kit
TT0010	Trypsin Spin Columns Proteomics Grade	10 each

For more information, please visit our Web site at sigma.com/pep.

6

Role of ABC Transporters and Multi-Drug Resistance Reversal Using RNA Interference (RNAi)*

Chemotherapy is the treatment of choice against many types of cancer. However, over time chemotherapeutic drugs can become less effective due to the development of resistance that involves a group of membrane proteins. These multi-drug transporters expel cytotoxic molecules from the cell, thus keeping intracellular drug concentrations below the cell-killing threshold. These transporters belong to the superfamily of ATP Binding Cassette (ABC) proteins that are present in all organisms from bacteria to humans. The transport activity of ABC proteins has an important effect on the efficacy of pharmaceuticals by modulating the absorption, distribution, and excretion of these xenobiotics.

ABC transporter proteins are located in the plasma membrane of cells and in the membranes of cellular organelles where they mediate the transport of various substrate molecules. These substrates exhibit a wide variety of chemical structures. Most ABC proteins are active transporters, which utilize the energy generated by ATP hydrolysis; however, some ABC transporters form transmembrane channels.

Cancer Multi-Drug Resistance - The Players

Numerous clinical data revealed that the multi-drug resistance phenotype in tumors is associated with the overexpression of certain ABC transporters, termed multi-drug resistance (MDR) proteins. P-glycoprotein (P-gp, MDR1, ABCB1) was the first discovered ABC transporter¹⁻³ and is likely to be responsible for the most widely observed mechanism in clinical multi-drug resistance.⁴⁻⁷ Soon after the cloning and characterization of MDR1, it became evident that other efflux pumps also play significant roles in transport-associated drug resistance. Two other ABC transporters have definitively demonstrated participation in the multi-drug resistance of tumors: the <u>multi-drug resistance</u> protein 1 (MRP1, ABCC1), and the <u>mitoxantrone resistance</u> protein (MXR/BCRP, ABCG2).⁷⁻¹¹

Basic Mechanism of Cancer Multi-Drug Resistance and Substrate Specificity of MDR-ABC Transporters

The generally accepted mechanism of multi-drug resistance is that the MDR proteins actively expel the cytotoxic drugs from cells, maintaining the drug concentration within the cells below the toxic level. The drug efflux mediated by these primary active transporters is driven by the energy of ATP hydrolysis. Tumors with MDR protein overexpression (e.g., hepatomas, lung or colon carcinomas) often show primary (or intrinsic) resistance to chemotherapy treatment. In addition, chemotherapy itself may induce the overexpression of these proteins, resulting in the multi-drug resistant clones becoming less sensitive to treatment (secondary drug resistance).¹²⁻¹⁴

The most intriguing characteristic distinguishing the MDR proteins from other mammalian transporters is their broad substrate specificity. Unlike other selective (classical) transport proteins, multidrug transporters recognize and handle a large number of structurally diverse, mainly hydrophobic compounds, which explains cross-resistance to several chemically unrelated compounds, a characteristic feature of the multi-drug resistance phenotype.⁴⁻⁷ In addition to their overlapping substrate specificity, each transporter can handle unique compounds. The following table of MDR-substrate anticancer agents provides a selection of anticancer agents available from Sigma and identifies the key ABC transporter(s) responsible for each agent's cellular efflux.

Circumvention of Cancer Multi-Drug Resistance

Prevention of clinical multi-drug resistance should significantly improve therapeutic response in a large number of cancer patients. The initial search for pharmacological modulators of MDR transporters yielded two generations of compounds having poor clinical response profiles. Therefore, there has been a shift to structure-based drug design to synthesize modulator compounds characterized by a high affinity to MDR transporters.^{15,16} Additionally, research that utilizes siRNA and shRNA-mediated RNAi-based gene silencing methodology has recently delivered promising results.

MDR Modifying Agents

MDR modifying agents, which competitively or non-competitively inhibit the MDR proteins, may increase the cytotoxic action of multi-drug resistant related drugs by preventing the active efflux of these drugs from the target cells. The co-application of an "MDR-modulating" compound in combination with chemotherapy would be expected to significantly improve the cancer cure rate. The first-generation modulators consisted of compounds that were already in clinical use. Calcium channel blockers, quinine derivatives, calmodulin inhibitors, and the immunosuppressive agent cyclosporin A, were all shown to interact with the MDR transporters *in vitro* and *in vivo*. These modulators were not specifically developed for MDR protein inhibition, and their inherent low affinity for MDR transporters resulted in a high toxicity profile and were never shown to inhibit P-gp in patients.¹⁵⁻¹⁸

Most of the second-generation modulators were derivatives of the first-generation compounds that retained MDR modulatory effects, but with reduced activity toward other physiological targets. Prominent examples of this group are R-verapamil, biricodar (VX-710), and valspodar (PSC-833). These modulators were shown to inhibit P-gp in patients, but further study revealed significant pharmokinetic interaction with several anticancer drugs, which delayed excretion of the anticancer agent, resulting in toxicity requiring reduction of anticancer drug doses.¹⁵⁻¹⁸

The third-generation MDR modulators are designed to interact with specific MDR transporters7, 11-13, 19, 20 with high affinity and with efficacy at nanomolar concentrations. Development of this class of MDR modulators employed combinatorial chemistry to produce potent and selective inhibitors. Examples are the small hydrophobic peptide derivatives named reversins, which were shown to have a strong inhibitory effect on P-gp/MDR1-mediated drug efflux without any toxic effect in the control cells.²¹ Ongoing clinical trials using third-generation MDR modulators for specific cancer types include: elacridar (GF120918), tariguidar (XR9576), zosuquidar (LY335979), laniquidar (R101933), and ONT-093. Still the shortcomings of earlier generation modulators continue to exist.^{17, 22} Other approaches to prevent the expression or function of multi-drug transporters are being considered, including the use of MDR protein targeted antibodies, the use of carriers that deliver these drugs selectively to tumor tissues, and the use of RNA interference

琳SIGMA

ABC Transporters and Cancer Drug Resistance

Multi-Drug Resistance Reversal Using RNAi

The stable reversal of MDR protein-mediated drug efflux by RNAi technology has been demonstrated in vitro for MDR1, MXR, MRP2, and MDR3. One of the early multi-drug resistance studies using RNAi technology reported a complete suppression of MDR1 expression on the mRNA and protein level in human gastric carcinoma cells.²³ A subsequent study further demonstrated inhibition of both MDR1 and MDR3 expression in conjunction with the reversal of paclitaxel resistance in human ovarian cancer cells. Treatment of ovarian cancer cell lines with either chemically synthesized siRNAs or transfection with specific vectors that express targeted siRNAs resulted in decreased mRNA and protein levels. In this study, MTT assays of siRNA-treated cells demonstrated 7 to 12.4-fold reduction of paclitaxel resistance in the lines treated with the synthesized siRNA of MDR1 and 4.7 to 7.3-fold reduction of paclitaxel resistance in the cell lines transfected with siRNA of MDR1 expressing vectors.²⁴ A more recent study surprisingly showed that the MDR1 phenotype in human hepatoma cells was completely reversed by using two transfected clones.²⁵ Aside from the more frequently studied MDR1 phenotype, reversal of the drug-resistant MXR and MRP2 phenotype using both siRNA and shRNA-mediated approaches was also demonstrated in human carcinoma cells.26,27

In a pre-clinical study the ablation of MDR1 in cells stably transduced with shRNA was functionally confirmed by increased sensitivity of MDR1-transfected cells toward the cytotoxic drugs vincristine, paclitaxel, and doxorubicin as well as by transport of 99mTc-sestamibi. In the same study, shRNA-mediated down-regulation of MDR activity in tumor implants in living animals was followed by direct noninvasive bioluminescence imaging using the fluorophore coelenterazine, a known MDR1 transport substrate. Additionally, a MDR1-firefly luciferase (MDR1-FLuc) fusion construct was used to document the effect of shRNA delivered in vivo on MDR1-FLuc protein levels with D-luciferin bioluminescence imaging.²⁸ A similar study validated selective MRP2 gene function inhibition: after the intravenous delivery of siRNA effectors into mice, researchers observed a significantly reduced calcein excretion rate and resultant siRNA accumulation in the kidney.29

RNAi is proving to be a powerful laboratory tool for better understanding the multi-drug resistance genotype and phenotype. Its future therapeutic utility in suppressing gene expression in cancer patients will likely be dependent on the availability of effective RNAi delivery systems. Lessons can be learned from the history of gene therapy and antisense technologies. These technologies ultimately failed to produce successful clinical outcomes due to potentially harmful and inefficient delivery systems. The use of RNAi in complex genetic diseases, such as cancer, will not see a quick and straightforward transition from research to clinical success, but with time the promise of viable RNAi therapies may be realized.³⁰ Additionally, innovative technologies combined with new directions in the study of ABC transporters will lead to an understanding of whether or not ABC transporters are important molecular targets for anticancer drug development.

ma-aldrich.com/biofiles

sign

- *Technical content provided by: Balázs Sarkadi¹, Gergely Szakács^{1,2} and András Váradi³ 1. National Medical Center, Institute of Haematology and Immunology, Budapest Hungary
- 2. National Cancer Institute, NIH, Bethesda, MD, USA
- 3. Institute of Enzymology, Hungarian Academy of Sciences, Budapest, Hungary

References:

- 1. Juliano, R.L. and Ling, V., Biochim. Biophys. Acta, 455, 152-162 (1976).
- 2. Chen, D., et al., Cell, 47, 381-389 (1986)
- 3. Ueda, K., et al., Proc. Natl. Acad. Sci. USA, 84, 3004-3008 (1987).
- 4. Endicott, J.A. and Ling, V., Annu. Rev. Biochem., 58, 137-171 (1989).
- 5. Higgins, C.F., Ann. Rev. Cell Biol., 8, 67-113 (1992).
- 6. Gottesman, M.M. and Pastan, I., Annu. Rev. Biochem., 62, 385-427
- (1993)
- 7. Gottesman, M.M., et al., Nat. Rev. Cancer, 2, 48-58 (2002).
- 8. Cole, S.P.C., et al., Science, 258, 1650-1654 (1992)
- 9. Borst, P., et al., J. Natl. Cancer Inst., 92, 1295-1302 (2000).
- Deeley, R.G. and Cole, S.P.C., Sem. Cancer Biol. 8, 193-204 (1997).
 Litman, T., et al., Cell. Mol. Life Sci., 58, 931-959 (2001).
- 12. List, A.F., et al., Blood, 98, 3212-3220 (2001)
- 13. Ambudkar, S.V., et al., Annu. Rev. Pharmacol. Toxicol., 39, 361-398 (1999).
- . 14. Sarkadi, B., et al., *Semin. Cancer Biol.*, **8**, 171-182 (1997)
- 15. Perez-Tomas, R., Curr. Med. Chem., 13(16), 1859-76 (2006).
- 16. McDevitt, C.A. and Callaghan, R., Pharmacol. Ther., 113(2), 429-441 (2007)

17. B. A. Teicher (ed.), Humana Press, Totowa, NJ, Cancer Drug Resistance, (2006), 2.15

- 18. Leonard, G.D., et al., Curr. Opin. Investig. Drugs, 3, 1652-1659 (2002). 19. Krishna, R., et al., Eur. J. Pharm. Sci., 11, 265-283 (2000).
- 20. Leslie, E.M., et al., Toxicology, 5, 3-23 (2001)
- 21. Sharom, F.J., et al., Biochem. Pharmacol., 58, 571-586 (1999).
- 22. Nobili, S., et al., Curr. Drug Targets, 7(7), 861-879 (2006).
- 23. Stege, A., et al., Cancer Gene Ther., 11(11), 699-706 (2004)
- 24. Zhenfeng, D., et al., Mol. Cancer Ther., 3(7), 833-838 (2004)
- 25. Chen, X.P., et al., World J. Gastroenterol., 12(21), 3332-3337 (2006).

26. Priebsch, A., et al., Oligonucleotides, 16(3), 263-274 (2006). 27. Materna, V., et al., Biochem. Biophys. Res. Commun., 348(1), 153-157 (2006)

28. Pichler, A., et al., Clin. Cancer Res., 11(12), 4487-4494 (2005).

29. van de Water, F.M., et al., Drug Metab. Dispos., 34(8), 1393-1397 (2006)

30. Downward, J., BMJ, 328(7450), 1245-1248 (2004).

MISSION". Accelerate, Innovate, Create

For more information on siRNA available through Sigma-Genosys and the comprehensive MISSION TRC shRNA library, targeting 15,000 human and 15,000 mouse genes, visit our Web site at sigma.com/rnai.

To obtain more information on specific ABC transporter shRNA clones, email RNAi@sial.com.

Drug Resistance

ABC Transporters and Cancer

SIGMA

Compound Actinomycin D	MDR1	¹ Cat. No. A1410	from <i>Streptomyces</i> sp., ~98% (HPLC)	2 mg
				2 Hiy
				5 mg
				10 mg 25 mg
				100 mg
		A4262	from <i>Streptomyces</i> sp., ~95% (HPLC)	2 mg
				5 mg
				10 mg 25 mg
Carmustine	MRP1	C0400	≥98%	25 mg
Jannustine		00400	27070	100 mg
Chlorambucil	MRP1	C0253		1 g
				10 g
Colchicine	MDR1	C9754	~95% (HPLC)	100 mg
				500 mg 1 g
				5g
Cyclophosphamide	MRP1	C0768	Bulk package, 97.0-103.0% (HPLC)	100 mg
monohydrate				500 mg
				1 g 5 g
		C7397	ISOPAC ^{®2}	3
Daunorubicin	MDR1	D8809	Meets USP testing specifications	1 mg
hydrochloride	MRP1	20007		10 mg
	MXR			50 mg
cis-Diammineplatinum(II)	MRP1	P4394		25 mg
dichloride (cisplatin)				100 mg 250 mg
,oropratinity				1 g
Doxorubicin HCI	MDR1, MRP1, MXR	D1515	~98% (TLC)	10 mg
Etoposide	MDR1	E1383	Synthetic, ≥98%	25 mg
	MRP1			100 mg
Folinic acid calcium salt	MRP1	F7878	>00%	250 mg
(Leucovorin)	IVIRP I	F/0/0	≥90%	100 mg 500 mg
, , , ,				1 g
Hydroxyurea	MRP1	H8627	≥98% (TLC)	1 g
				5 g 10 g
				10 g 25 g
				100 g
Methotrexate hydrate	MRP1	A6770	>98% (HPLC)	10 mg
				25 mg
				100 mg 500 mg
				1 g
Mitomycin C	MDR1	M0503	from Streptomyces caespitosus	2 mg
			Vial contains 2 mg mitomycin C and 48 mg NaCl.	5 X 2 mg
				10 X 2 mg 50 X 2 mg
Mitoxantrone	MDR1	M6545	≥97% (HPLC)	10 mg
dihydrochloride	MXR		· /	50 mg
Paclitaxel	MDR1	T7402	from <i>Taxus brevifolia</i> , ≥95% (HPLC)	1 mg
				5 mg 25 mg
Tamoxifen	MDR1	T5648	≥99%	23 mg
	MRP1	10010	_,,,,,	5 g
Vinblastine sulfate salt	MDR1	V1377	≥97% (TLC)	1 mg
	MRP1			5 mg
				10 mg 25 mg
				50 mg
Vincristine sulfate salt	MDR1	V8879	95.0-105.0% (HPLC)	1 mg
	MRP1			5 mg 25 mg

MDR-substrate Anticancer Agents available from Sigma

The ABC transporters indicated are the proteins which are responsible for the cellular efflux of the specific anticancer agent.
 ISOPAC[®] products are packaged in 100 mL serum bottles with a butyl rubber stopper and aluminum tear seal. Injecting any compatible solvent permits preparation of any desired strength solution without exposure.

SIGMA

sigma-aldrich.com/biofiles

ABC Transporters and Cancer Drug Resistance

ABC Transporter Membrane Proteins

MDR1 human

ABCB1; Pgp

The MDR1 protein is involved in cancer drug resistance and in the transport of hydrophobic drugs and xenobiotics in the bowel, kidney, liver, and the blood-brain barrier. Drugs interacting with this protein may be useful for the reversal of cancer drug resistance or increasing the absorption or brain entry of various pharmacological agents.

Membrane preparation, for ATPase, recombinant, expressed in Sf9 cells

Supplied as isolated Sf9 cell membranes containing human MDR1 (Pgp) suspended in TMEP solution.

Distributed for SOLVO Biotechnology, Inc.

-70°C DRY ICE

500 ul

500 µL

500 uL

Mdr1b from rat

The MDR1 protein is involved in cancer drug resistance and in the transport of hydrophobic drugs and xenobiotics in the bowel, kidney, liver, and the blood-brain barrier. In rodents, there are two MDR1 genes, mdr1a and mdr1b, while in human, there is a single MDR1 gene. Based on function and tissue distribution in rodents, the equivalent of the human MDR1 gene product (PgP) is the product of the rodent mdr1b gene. There have been no reported significant differences in function, substrate specificity, or substrate affinity between these two proteins.

$\label{eq:membrane} \mbox{Membrane preparation, for ATPase, recombinant, expressed in Sf9 cells.}$

Supplied as isolated Sf9 cell membranes containing rat Mdr1b suspended in TMEP solution.

Distributed for SOLVO Biotechnology, Inc.

-70°C DRY ICE

M9319

MRP2 human

MRP2 (ABCC2) is an organic anion transporter found in the liver, kidney, and gut epithelium apical membranes. The transport of glucuronate conjugates plays a role in the detoxification of endogenous and xenobiotic substances, and may cause multidrug resistance (MDR) in tumor cells.

Membrane preparation, for Vesicular Transport, recombinant, expressed in Sf9 cells

Supplied as isolated Sf9 cell membranes containing human MRP2 suspended in TMEP solution.

Distributed for SOLVO Biotechnology, Inc.

-70°C DRY ICE

M9069

Mrp2 from rat

MRP2 (ABCC2) is an organic anion transporter found in liver, kidney, and gut epithelium apical membranes. The transport of glucuronate conjugates plays a role in the detoxification of endogenous and xenobiotic substances, and may cause multidrug resistance (MDR) in tumor cells. The rat Mrp2 transporter shows 72.3% sequence identity and 85.6% sequence similarity with human MRP2. Both transporters are expressed on the canalicular membrane of the liver and are known to be responsible for the transport of some organic molecules and their conjugates to the bile.

Membrane preparation, for ATPase, recombinant, expressed in Sf9 cells

Supplied as isolated Sf9 cell membranes containing rat Mpr2 suspended in TMEP solution. Distributed for SOLVO Biotechnology, Inc.

500 µL

MXR human

MXR membrane vesicles are purified from recombinant baculovirus transduced Sf9 cells or selected, MXR over-expressing mammalian cells.

Distributed for SOLVO Biotechnology, Inc.

Membrane preparation, wild type variant, for Vesicular Transport

Supplied as isolated mammalian cell membranes containing human MXR (wild type variant) suspended in TMEP solution. The MXR transporter can be produced in sufficient quantity by selected, MXR over-expressing mammalian cell lines.

-70°℃ DRY ICE

500 µL

Membrane preparation, wild type variant, for Vesicular Transport, recombinant, expressed in Sf9 cells

Supplied as isolated Sf9 cell membranes containing wild type human MXR suspended in TMEP solution.

The MXR transporter can be expressed in Sf9 insect cells using the baculoviral expression system, yielding high protein levels (up to 5% of total membrane protein) in the cell membrane of infected cells.

-70°C DRY ICE

M9444 500 μL

MXR Control

Membrane preparation, from Sf9 cells

Control for ATPase and vesicular transport assays.

Supplied as isolated Sf9 cell membranes containing a defective MXR gene suspended in TMEP solution.

Distributed for SOLVO Biotechnology, Inc.

-70°C DRY ICE

M9944

500 µL

MXR Control

Membrane preparation, mammalian

Control for ATPase and vesicular transport assays.

Supplied as isolated mammalian cell membranes (not selected for transport expression) suspended in TMEP solution. Distributed for SOLVO Biotechnology, Inc.

isti buteu toi 30140 biotechnology,

-70°C DRY ICE

C3992



Drug Resistance

ABC Transporters and Cancer

Monoclonal Anti-Breast Cancer Resistance Protein antibody produced in mouse

Anti-ABCG2: Anti-BCRP

250 µg/mL, clone BXP-21, tissue culture supernatant

Antibodies to ABC Transporters

Immunogen: fusion protein containing human BCRP (amino acids 271-396) and maltose-binding protein. Reacts with an internal epitope of BCRP. Does not cross-react with the human MDR1, MRP1, MRP2 gene products.

Supplied in serum-free medium containing 0.7% bovine serum albumin and 0.1% sodium azide.

Species reactivity: human

Antigen mol wt ~70 kDa

Application(s)

Immunocytochemistry.	1:20-1:50 using acetone-fixed
	cytospin preparations
Immunohistochemistry	(frozen sections) 1:20 using acetone-fixed
	frozen sections
Immunohistochemistry	(formalin-fixed, paraffin-embedded sections)
	suitable using pretreated human tissue
Immunoblotting	suitable
Isotype	IgG2a
-20°C DRY ICF	

L BITTIGE	
B7059	1 mL

Anti-Breast Cancer Resistance Protein antibody produced in rabbit

Anti-BCRP

Affinity isolated antibody

Immunogen: synthetic peptide corresponding to amino acid residues 150-167 of human breast cancer resistance protein with C-terminal added cysteine, conjugated to KLH.

Solution in 0.01 M phosphate buffered saline containing 15 mM sodium azide.

Species reactivity: human, mouse

Antigen mol wt ~70 kDa

Application(s)

Immunoblotting......2.5-5 µg/mL using whole extract of human term placenta or mouse kidney and a chemiluminescent detection reagent

Immunohistochemistry (formalin-fixed, paraffin-embedded

sections)	20-40 µg/mL using heat-retrieved tissue
	sections from human term placenta by indirect
	immunoperoxidase staining of syncytiotrophoblasts.

-20°C DRY ICE	
B7185	200 µL

Monoclonal Anti-P-Glycoprotein (MDR) antibody produced in mouse

Clone F4, ascites fluid

Immunogen: mixture of human and hamster drug-resistant whole cells and crude plasma membranes.

The antibody recognizes an epitope located in the amino terminal half of P-glycoprotein (Pgp), at the third extracellular loop of the molecule. The epitope is resistant to formalin fixation and periodate oxidation.¹ The antibody detects specifically human MDR1 P-glycoprotein, but does not appear to recognize the human MDR3 product¹, nor the mouse mdr1a, mdr1b or the mdr3 P-glycoprotein.²

Contains 15 mM sodium azide

Species reactivity: human, hamster

Antigen mol wt 170-180 kDa

Application(s)

Immunoblottingsuitable
Radioimmunoassay suitable using cell-surface RIA
Immunoprecipitation suitable
Immunohistochemistry (frozen sections)suitable
Immunocytochemistry suitable
Flow cytometrysuitable
Indirect ELISAsuitable
Immunohistochemistry (formalin-fixed, paraffin-embedded sections)
1:500 using human kidney sections
Isotype IgG1
Lit. cited: 1. Chu, T.M., et al., Hybridoma 12, 417 (1993)
2. Chu, T.M., et al., Biochem. Biophys. Res. Commun. 203, 506 (1994)
-20°C DRY ICF

P7965 0.2 mL

Monoclonal Anti-MDR3 P-Glycoprotein antibody produced in mouse

250 μ g/mL, clone P₃II-26, tissue culture supernatant

Immunogen: MDR3 P-gp (amino acids 629-692) GST fusion protein.

Reacts with an internal epitope of MDR3. Does not cross-react with human MDR1 P-gp.

Supplied in serum-free medium containing 0.7% bovine serum albumin and 0.1% sodium azide.

Species reactivity: human

Application(s)

M7317	1 mL
-20°C DRY ICE	
Isotype	lgG2b
Immunoblotting	
	Not suitable
Immunohistochemistry (formalin-fixed, paraffi	n-embedded sections)
	acetone-fixed sections
Immunohistochemistry (frozen sections)	1:20 using
	cytospin preparations
Immunocytochemistry1:20	

SIGMA

ABC Transporters and Cancer Drug Resistance

Monoclonal Anti-MRP1 antibody produced in mouse

Anti-Multidrug Resistance Associated Protein 1

~2 mg/mL, clone QCRL-1, purified immunoglobulin

Immunogen: non-denatured membrane preparations of H69AR human small cell lung cancer cell line, which highly expresses MRP1. The epitope resides within amino acids 918-924 of human MRP1. Does not cross-react with the human MDR1 and MDR3. mouse MRP1, and human MRP2, MRP3, MRP4, MRP5, and MRP6 gene products. The antibody may cross-react with canine MRP1.

Solution in 0.01 M phosphate buffered saline, pH 7.4, containing 1% BSA and 15 mM sodium azide.

Species reactivity: human

Antigen mol wt ~190 kDa

Application(s)

Flow cytometry	0.5-2 µg/mL using human H69AR cells
Immunoblotting	suitable
Immunoprecipitation	suitable
Immunohistochemistry	suitable
Immunocytochemistry	suitable
lsotype	lgG1
-20°C DRY ICE	
M9067	200 µL

Anti-MRP2 antibody produced in rabbit

Anti-ABCC2; Anti-cMOAT; Anti-cMRP; Anti-Multidrug Resistance Associated Protein 2

Affinity isolated antibody

Immunogen: synthetic C-terminal peptide corresponding to amino acids 1528-1545 of human MRP2 conjugated to KLH.

Additional lower bands including an approx. 175 kDa band representing an immature unglycosylated form may be detected in various extract preparations.

Solution in 0.01 M phosphate buffered saline, pH 7.4, containing 1% bovine serum albumin and 15 mM sodium azide.

Species reactivity: human, rat

Antigen mol wt ~190 kDa

Application(s)

Immunoblotting	1:1,000 using whole extract of
-	HepG2 human hepatoblastoma cells
Immunoblotting	.1:1,000 using 293T cells expressing
-	recombinant human MRP2
Indirect immunofluorescence	1:100 using
	paraformaldehyde-fixed HepG2 cells
Indirect immunofluorescence 1:100 using rat liver frozen sections	
-20°C DRY ICE	
M8316	200 µL

SIGMA

Monoclonal Anti-MRP2 antibody produced in mouse

Anti-ABCC2; Anti-cMOAT; Anti-cMRP; Anti-Multidrug Resistance Associated Protein 2

~1.5 mg/mL, clone CPR96, purified immunoglobulin

Immunogen: synthetic C-terminal peptide corresponding to amino acid residues 1528-1545 of human MRP2 with N-terminal added cysteine conjugated to KLH.

Solution in 0.01 M phosphate buffered saline, pH 7.4, containing 15 mM sodium azide.

Species reactivity: human

Antigen mol wt ~180 kDa

Application(s)

Immunoblotting	1-2 µg/mL using cell extracts of
-	293T cells transfected with human MRP2
Indirect ELISA	suitable
Array	suitable
lsotype	lgG1
-20°C DRY ICE	
M3692	200 µL

Anti-MRP3 antibody produced in rabbit

Anti-Multidrug Resistance Associated Protein 3

Affinity isolated antibody

Immunogen: synthetic C-terminal peptide of human MRP3 (amino acids 1507-1527) with N-terminal cysteine conjugated to KLH. The sequence in mouse and rat differs by three residues.

Solution in 0.01 M phosphate buffered saline, pH 7.4, containing 15 mM sodium azide.

Species reactivity: human, rat

Antigen mol wt ~180 kDa

Application(s)

Immunoblotting	2.5-5 μg/mL
US	ing whole extract of cultured human
	colon carcinoma HCT-116 cells
Indirect immunofluorescence	. 15-30 µg/mL using frozen sections
	of rat skin (outer root sheaths
	of hair follicles)
-20°C DRY ICE	
M0318	200 µl

Monoclonal Anti-MRP5 antibody produced in rat

Anti-Multidrug Resistance Associated Protein 5

250 µg/mL, clone M₅I-1, tissue culture supernatant

Immunogen: human MRP5 (amino acids 82-168) bacterial fusion protein. Reacts with an internal epitope of MRP5. Does not crossreact with the human MDR1, MRP1, MRP2, MRP3 gene products. Supplied in serum-free culture medium containing 0.7% bovine serum albumin and 0.1% sodium azide.

Species reactivity: human

Antigen mol wt 190-200 kDa

Application(s)

M6067	1 mL
-20°C DRY ICE	
Isotype	lgG2a
Immunoblotting	1:20-1:50
	acetone-fixed frozen sections
Immunohistochemistry (frozen sections)	
	preparations
Immunocytochemistry 1:20-1:50) using acetone-fixed cytospin
Application(3)	

Compounds for MDR Detection

bisBenzimide H 33342 trihydrochloride

bisBenzimide; 2'-(4-Ethoxyphenyl)-5-(4-methyl-1-piperazinyl)-2,5'-bi-1Hbenzimidazole trihydrochloride; HOE 33342; Hoechst 33342 C₂₇H₂₈N₆O · 3HCl · xH₂O FW 561.93 (Anh) [23491-52-3]

Membrane-permeable, fluorescent DNA stain with low cytotoxicity that intercalates in A-T regions of DNA.

Useful for staining DNA, chromosomes and nuclei. May be used for fluorescence microscopy or flow cytometry.

Excitation max. = 346 nm

Emission max. = 460 nm

≥95% (TLC)

-20°C EC No. 2456911 BRN 1234011

B2261

25 mg, 100 mg, 500 mg, 1 g

5 g, 10 g, 25 g

100 µL

Calcein

Bis[N,N-bis(carboxymethyl)aminomethyl]fluorescein; Fluoresceinbis(methyliminodiacetic acid); Fluorexon

C30H26N2O13 FW 622.53 [1461-15-0]

Used for the fluorometric determination of calcium and EDTA titration of calcium in the presence of magnesium.

-20°C EC No. 2159571 BRN 8181068

C0875

Calcein AM solution

C46H46N2O23 FW 994.86 [148504-34-1]

 \geq 4 mM in DMSO, \geq 90% (TLC)

Fluorescent cell permeable derivative of calcein.

-20°C

C1359

Coelenterazine, native

CLZN; 3,2-Dihydro-2-(p-hydroxybenzyl)-6-(p-hydroxyphenyl)-8-benzylimidazo-lo[1,2-a]pyrazin-3-one

C₂₆H₂₁N₃O₃ FW 423.46 [55779-48-1]

Luminophore of the aequorin complex which is oxidized by oxygen to illuminate at 465 nm when Ca2+ binds to the complex; used to measure Ca²⁺ concentration in cells with high sensitivity and large dynamic range

-20°C

C2230 50 µg

Luciferase from Photinus pyralis (firefly)

[61970-00-1] E.C. 1.13.12.7

Lyophilized powder, 15-30 × 10⁶ light units/mg protein Arsenate free.

Lyophilized powder approximately 20% protein; balance is primarily NaCl, HEPES buffer salts, and carbohydrate.

Mol wt 120 kDa (two non-identical subunits, each containing four free thiol groups, one of which is necessary for activity)

Composition Protein ~20% (E^{1%}₂₈₀)

Chromatographically prepared and crystallized.

Sold on the basis of protein content

One light unit produces a biometer peak height equivalent to 0.02 µCi of ¹⁴C in PPO/POPOP cocktail. Light units measured in 50 µl assay mixture containing 5 pmol ATP and 7.5 nmol luciferin in Tris-glycine buffer, pH 7.6, at 25 °C.

Sensitivity≤1 femtomole	ATP
(using 0.2 µg of luciferase and suitably sense	itive
liquid scintillation counters or luminome	ters)
ATPase	tein
Nucleoside diphosphokinase	tein
-20°C EC NO. EINECS	
L9506 1 mg, 2 mg, 10	mg

D-Luciferin

4,5-Dihydro-2-(6-hydroxy-2-benzothiazolyl)-4-thiazolecarboxylic acid; Firefly luciferin; (S)-2-(6-Hydroxy-2-benzothiazolyl)-2-thiazoline-4-carboxylic acid $\rm C_{11}H_8N_{2}O_3S_2$ FW 280.32 [2591-17-5] EC No. 2199813 BRN 30484

Synthetic, SigmaUltra, ≥99% (HPLC)

Highly purified grade with reduced levels of inhibitor relative to L9504

L6152	1 mg, 10 mg
Synthetic	
L9504	1 mg, 5 mg, 10 mg, 50 mg, 100 mg

Rhodamine 123

2-(6-Amino-3-imino-3H-xanthen-9-yl)benzoic acid methyl ester; Rh 123 C₂₁H₁₇CIN₂O₃ FW 380.82 [62669-70-9]

Fluorescent dye most commonly used in flow cytometry as functional reporter for P-glycoprotein (P-gp). Functional assays for MDR proteins are better prognostic indicators in cancer therapy than levels of MDR protein expression. Rh 123 can be used in multiparameter analyses without fluorescence interference in combination with common protein labeling dyes such as PE-Cy™5 and AMCA (7-amino-4-methylcoumarin-3-acetic acid). Recent reports indicate Rh 123 may also be a substrate of MRP1. Used as a laser dye and for selective cell growth effects.

EC No. 2636878 BRN 6030951

R8004

5 mg, 25 mg, 100 mg

G ma-aldrich.com/biofiles SIGMA

NF-κB and Inflammation

Inflammation and Cancer: the NF-kB Connection

Chronic inflammation is an underlying factor in the development and progression of many of the chronic diseases of aging, such as arthritis, atherosclerosis, diabetes, and cancer. Oxidative cellular stress induced by environmental factors, such as cigarette smoke, UV or ionizing radiation, bacterial or viral infection, or any number of oxidizing xenobiotic compounds, triggers a wide range of cellular responses, some of which are proinflammatory and proapoptotic, while others protect the cell against apoptosis and enhance cellular adhesion, cell proliferation, and angiogenesis. The inappropriate induction or constitutive activation of these protective responses in mutated or damaged cells appears to be a major factor in the transformation and proliferation of cancer cells.¹⁻⁵ Two nuclear transcription factors that are involved in mediating the cellular responses to oxidative cell stress and proinflammatory stimuli are activator protein-1 (AP-1) and nuclear factor- κ B (NF- κ B). The role of these transcription factors on cancer initiation and progression has been studied in cell culture and *in vivo* models.^{5,6} The activation of two AP-1 components, c-Jun and c-Fos, by JNK and by ERK1/2 or p38 MAPK, respectively, is involved in the malignant transformation of cells stimulated by the tumor promotor phorbol 12-myristate 13-acetate (PMA). The proinflammatory and antiapoptotic response to tumor promotion is primarily mediated through activation of NF- κ B by the IKK family of serine/threonine kinases. The following discussion will focus on the NF- κ B pathway as a target for cancer chemotherapy and chemoprevention.




Figure 1. Many diverse stimuli utilize intracellular signaling pathways to activate NF-kB, a nuclear transcription factor that regulates proinflammatory and cell survival pathways.

SIGMA

A

NF-κB Transcription Factors

NF-κB refers to a family of transcription factors that has been highly conserved through evolution and is present in the cytoplasm of all cells. NF-κB has been called a "stress sensor" because its activity is induced by a wide variety of stimuli,⁷ including tumor necrosis factor (TNF-α), PMA and other tumor promoters, cigarette smoke extract (CSE), lipopolysaccharide (LPS), oxidants, and pathogenic bacteria. The NF-κB family comprises five members: p50 (NF-κB1), p52 (NF-κB2), RelA (p65), RelB, and c-Rel. p50 and p52 are cleaved from inactive precursor proteins, p105 and p100, respectively, prior to translocation to the nucleus. NF-κB family members are characterized by having:

- · a Rel homology domain that binds to DNA
- a dimerization domain
- the ability to bind to the intracellular inhibitor complex, $I\kappa B$

The most widely studied NF- κ B heterodimers are p50/p65 and p50/c-Rel (both associated with the classical or canonical pathway) and p52/RelB (alternative pathway). The classical pathway is activated by inflammatory cytokines, bacterial and viral infections, and oxidative stimuli and induces gene expression responsible for the antiapoptotic actions of NF- κ B. The alternative pathway is primarily involved in B cell survival.^{2,7-9} The classical pathway is illustrated in **Figure 1**.

Cytoplasmic NF-KB is sequestered as an inactive complex with its regulatory subunit, IkB. The most abundant member of the IkB family of proteins is $I\kappa B\alpha$. Phosphorylation of two conserved serine residues in the N-terminal domain of the NF- κ B/I κ B complex induces the rapid dissociation and polyubiguitination of IkB followed by its degradation by the 26S proteasome. Activated NF- κ B translocates to the nucleus where specific subunit lysines are acetylated by SRC-1 and p300 histone acetyltransferases. Acetylation promotes DNA binding and NF-kB-induced gene transcription.¹⁰ A list of more than 200 proteins that are regulated by NF-KB is given by Ahn and Aggarwal.⁷ Many of the genes regulated by NF-kB code for inflammatory cytokines and proteins that mediate cell survival, cellular adhesion, cell cycle activation, cell proliferation, angiogenesis, and oncogenesis. However, not all the actions of NF-kB promote cell survival. Activation of NF-kB also appears to be essential for p53-induced apoptosis in response to oxidative stress or to the anticancer agents, doxorubicin and etoposide.^{2,11}

Activation of NF-κB via Phosphorylation

Phosphorylation of the NF-kB/lkB complex is catalyzed by IKK, a protein complex that contains two homologous kinase subunits (IKK α and IKK β) and a regulatory subunit IKK γ /NEMO.¹⁷ Activation of the IKK complex can be initiated by any one of several intracellular phosphorylation pathways, including NF-kB-inducing kinase (TRAF/NIK),^{8,12,13} MEK1,^{11,14} ERK5,¹⁵ and PI3K/Akt.¹⁶ T cell, B cell and lysophosphatidic acid receptors activate a kinase cascade that results in the activation of the IKK complex by Bcl10, Malt1 and CARMA-1.9,17 Acetylation of serines and threonines in the activation loop of the IKK α and IKK β subunits can prevent phosphorylation and activation of the IKK complex.¹⁸ In addition, regulation of the NF-kB/kB complex can also occur independently of IKK activation or inhibition. Both the PI3K/Akt and JAK/STAT/Pim kinase pathways activate NF-kB by phosphorylating Cot, a serinethreonine kinase that can induce the proteasomal degradation of IKB.19 Upregulation or over-expression of mitogen activated protein kinase phosphatase-5 (MKP5) can decrease cytokine-induced phosphorylation of NF-kB/lkB and of p38 MAP kinase.²⁰

There are two mechanisms by which the NF- κ B-induced gene transcription is terminated. Genes coding for I κ B complex proteins are upregulated by NF- κ B. Newly formed I κ B α subunits can enter the nucleus where they bind to and inactivate NF- κ B, and the NF- κ B/I κ B α complex is exported back to the cellular cytoplasmic compartment.⁹ Alternatively, the translational action of DNA-bound NF- κ B can be terminated by deacetylation of the p65 subunit at Lys³¹⁰ by the histone deacetylases SIRT1 or HDAC1.^{10,21} Inhibition of p300 histone acetyltransferase or the overexpression or activation of SIRT1 has been shown to inhibit NF- κ B-mediated gene expression.^{10,22} Conversely, decreased histone deacetylase level or activity has been shown to increase the expression of inflammatory cytokines, presumably through enhancement of NF- κ B-mediated gene transcription.²³

NF-κB and its Relationship to Disease

Altered regulation of NF- κ B activity is observed in many genetically-linked diseases and chronic diseases of aging, including cancer.7,9 NF-KB activation has been linked to inflammation-driven tumor promotion and progression. In addition, many solid and hematopoietic cancers express constitutively active NF-KB that contributes to the pathogenesis of the disease by inducing factors that promote proliferation, invasiveness, angiogenesis, and resistance to chemotherapeutic agents and radiation. Researchers have hypothesized that inhibition of NF-KB activation or transcriptional activity may delay cancer onset or may be used as an adjunct to more traditional chemotherapeutics. Many compounds, including many phytochemicals and micronutrients having putative chemopreventive properties, inhibit NF-kB activation or constitutive NF-kB activity in cellular or *in vivo* models of cancer.^{2,5,7,24} The anticancer activities of some of these phytochemicals are summarized in the table on page 16.

Many of the inducers of NF-kB also regulate other intracellular pathways that mediate cell cycle arrest or apoptosis. One hypothesis holds that inhibition of prosurvival pathways mediated via NF-kB allows expression of proapoptotic mechanisms that may also be mediated by the various cell signaling pathways. For example, TNF- α is an inflammatory cytokine that mediates a broad spectrum of biological actions via activation of the TNFR1 receptor and, depending on the cellular environment, can promote either cell survival or programmed cell death. On the one hand, TNF-α/TNFR1-induced activation of TRAF2 initiates kinase cascades that lead to phosphorylation and activation of AP-1 and NF-ĸB, thereby promoting gene expression, cell survival and proliferation. Conversely, through its interactions with TRADD or FADD, TNF- α / TNFR1 initiates signaling pathways that activate caspase 8 and the proteolytic cascade that ends in apoptotic cell death. Furthermore, some of the genes induced by NF- κ B activation, such as Gadd45 β and XIAP, are inhibitors of prolonged JNK activation. JNK activation facilitates mitogen-induced and oxidative stress-induced apoptosis. Thus, in the presence of NF-kB inhibitors, TNFR1 activation would favor the proapoptotic, anticancer actions of TNF- α .^{2,13,25}

It should be noted that many phosphorylation pathways, e.g., PI3K-Akt, JNK, and ERK, are involved in both prosurvival and proapoptotic cellular processes. Therefore, the physiological response to procarcinogenic and anticarcinogenic xenobiotics represents the sum of all the specific intracellular signaling pathways that are upregulated and down-regulated in response to these stimuli. While chemopreventive natural products may tip the balance toward cell death in some damaged or transformed cells, their actions may promote carcinogenesis or drug resistance in other cells or organisms.

SIGMA

$\text{NF-}\kappa\text{B}$ and Inflammation

Natural Product Inhibitors of NF-κB Activation

Compound	Inducer	Mechanism of NF-κB Inhibition	Chemopreventive/ Chemotherapeutic Actions	Ref.	Cat. No.	Pack Size
Betulinic acid	TNF- α , IL-1 β , PMA, H ₂ O ₂ , okadaic acid, cigarette smoke	Inhibits activation of IKK α .	Downregulates COX-2 and MMP-9.	26, 27	855057	100 mg 500 mg
Caffeic acid phenethyl ester (CAPE)		Suppresses p65 nuclear translocalization; inhibits MAPK.	Decreases HIF-1α expression, induces HO-1 expression, reduces iNOS and MMP-9 expression.	5, 27, 28	C8221	1 ς
Curcumin	TNF- α , PMA, H ₂ O ₂	Inhibits phosphorylation of Akt and activation of ΙΚΚα; upregulates MKP5 expression, inhibits p38 activation.	Also inhibits STAT3 and AP-1 activation and activates PPAR- γ ; suppresses expression of c-Fos and c-Jun, inhibits COX-2 expression, antioxidant; blocks tumor	5, 20, 29, 30, 31	C7727 ≥94% (curcuminoid content), ≥80% (Curcumin)	500 mg
			initiation, promotion, invasion, angiogenesis, and metastasis.		28260 Mixture of curcumin, desmethoxycurcumin and bisdesmethoxycurcumin, ≥95.0% (TLC)	10 g 50 g
Embelin	TNF- α , LPS, IL-1 β , PMA, H ₂ O ₂ , okadaic acid, cigarette smoke	Inhibits IKKα-induced phosphorylation of ΙκΒα.	Enhances apoptosis induced by cytokines and chemotherapeutic agents.	32	E1406	10 mg 50 mg
Epigallocatechin gallate	PMA	Inhibits IKK activation	Also inhibits activation of AP-1, inhibits PI3K-Akt and ERK1/2 pathways, suppresses COX-2	5, 33	50299 ≥97.0% (HPLC)	1 mç
			induction, blocks <i>H. pylori</i> - induced glycosylation of TLR-4.		E4143 ≥95%, from green tea	50 mg
Evodiamine	TNF- α , LPS, IL-1 β , PMA, cigarette smoke	Inhibits phosphorylation of Akt, association of Akt with ΙΚΚα, and activation of ΙΚΚα.	Inhibits COX-2 induction, enhances apoptosis induced by cytokines and chemotherapeutic agents, inhibits TNF-induced invasive activity. Inhibits both constitutive and induced NF- κ B; decreases HIF-1 α expression.	28, 34,35	E3531	250 mg
Genistein	⁶⁰ Co Radiation, docetaxel, cisplatin		Potentiates radiation-induced apoptosis; upregulates p21 ^{WAF1/Cip1} .	36	G6649 Synthetic, ≥98% (HPLC)	5 mg 25 mg 100 mg
					G6776 from <i>Glycine max</i> (soybean), ~98% (HPLC)	5 mg 10 mg
					91955 BioChemika, ≥98.0% (HPLC)	25 mg 100 mg
6-Gingerol	TNF- α , IL-1 β	Upregulates MKP5 expression, inhibits p38 activation.	Anti-inflammatory	20	G1046	n/a
Honokiol	TNF-α, P. acnes	Inhibits activation of ΙΚΚα and phos- phorylation of Akt; operates downstream of ΜΕΚΚ-1.	Downregulates genes involved in antiapoptosis, proliferation, invasion, and angiogenesis; downregulates COX-2.	37, 38, 39	H4914	10 mg 25 mg
Indole-3-carbinol	TNF- α , LPS, IL-1 β , PMA, H ₂ O ₂ , okadaic acid, cigarette smoke	Inhibits activation of ΙΚΚα.	Inhibits constitutive and induced NF-κB activation; binds to ERα, AhR, and AR and suppresses invasion and migration of human breast cancer cells.	40	17256	1 g 5 g 25 g
Magnolol	TNF-α, <i>P. acnes</i>	Inhibits activation of IKKα and phosphorylation of Akt, operates downstream of MEKK-1.	Downregulates COX-2	37,38	M3445	10 mg
Partenolide		Inhibits IKK		27	P0667	5 mg 25 mg
Piceatannol	LPS	Inhibits activation of IKK α and IKK β , inhibits phosphorylation of Akt and Raf-1; inhibits Syk and p56 ^{lck} activity.	Inhibits iNOS and COX-2 induction and transcription of proinflammatory cytokine genes	41,42	P0453	5 mg 25 mg

sigma-aldrich.com/biofiles



16

Natural Product Inhibitors of NF-kB Activation (continued)

Compound	Inducer	Mechanism of NF-κB Inhibition	Chemopreventive/ Chemotherapeutic Actions	Ref.	Cat. No.	Pack Size
Silibinin	TNF-α	Inhibits IKK $lpha$ activity	Inhibits constitutive and induced NF-κB, inhibits PI3K-Akt and ERK1/2 pathways, inhibits c-Jun and c-Fos activation; inhibits tumor invasion, angiogenesis.	33, 47, 48, 49, 50	S0417	1 g 10 g
Tetrandrine	CD28, ConA, PHA, PMA + ionomycin	Inhibits IKK α and β activity, inhibits activation of JNK, p38 MAPK, ERK.	Antiinflammatory; also blocks AP-1 activation	51	T2695	1 g
Xanthohumol	TNF-α		Antiangiogenic; suppresses iNOS expression.	52	X0379	5 mg

Abbreviations: PMA - phorbol myristate acetate; HO-1 - Heme oxygenase 1; MKP5 - mitogen activated protein kinase phosphatase-5;

LPS – lipopolysaccharide; TNR-4 – Toll-like receptor-4; P. acnes – Propionibacterium acnes; ER- α – estogen receptor α ; AhR – aryl hydrocarbon receptor; AR - androgen receptor; Con A - concanavalin A; PHA - Phasolus vulgaris lectin

References:

- 1. Kundu, J.K., and Surh, Y.J., Mutat. Res., 591, 123-146 (2005).
- 2. Luo, J.L., et al., J. Clin. Invest., 115, 2625-2632 (2005).
- 3. Nickoloff, B.J., et al., J. Invest. Dermatol., 124, x-xiv (2005)
- 4. Schottenfeld, D., and Beebe-Dimmer, J., CA Cancer J. Clin., 56, 69-83 (2006).
- 5. Surh, Y.J., et al., J. Nutr., 135, 2993S-3001S (2005).
- 6. Hou, D.X., et al., J. Biomed. Biotechnol., 5, 321-325 (2004)
- 7. Ahn, K.S., and Aggarwal, B.B., Ann. NY Acad. Sci., 1056, 218-233 (2005).
- 8. Ramakrishnan, P., et al., Immunity, 21, 477-489 (2004)
- 9. Yates, L. L., and Górecki, D.C., Acta Biochim. Pol., 53, 651-662 (2006).
- 10. Yeung, F., et al., EMBO J., 23, 2369-2380 (2004).
- 11. Armstrong, M.B., et al., Neoplasia, 8, 967-977 (2006).
- 12. Hu, W.H., et al., J. Biol. Chem., 280, 29233-29241 (2005).
- 13. van Horssen, R., et al., Oncologist, 11, 397-408 (2006).
- 14. Zhi, L., et al., J. Leukoc. Biol., e-pub Feb. 8, 2007.
- 15. Garaude, J., et al., J. Immunol., 177, 7607-7617 (2006)
- 16. Agarwal, A., et al., Oncogene, 24, 1021-1031 (2005).
- 17. Klemm, S., et al., Proc. Natl. Acad. Sci. USA, 104, 134-138 (2007)
- 18. Mittal, R., et al., Proc. Natl. Acad. Sci. USA, 103, 18574-18579 (2006).
- 19. Amaravadi, R., and Thompson, C.B., J. Clin. Invest., 115, 2618-2624 (2005)
- 20. Nonn, L., et al., Carcinogensis, e-pub Dec. 6, 2006.
- 21. Kim, Y.K., et al., Biochem. Biophys. Res. Commun., 347, 1088-1093 (2006)
- 22. Kaeberlein, M., et al., J. Biol. Chem., 280, 17038-17045 (2005).

23. Yang, S.R., et al., Am. J. Physiol. Lung Cell. Mol. Physiol., 291, L46-L57 (2006)

24. Amiri, K.I., and Richmond, A., Cancer Metastasis Rev., 24, 301-313 (2005).

- 25. Papa, S., et al., J. Cell Sci., 117, 5197-5208 (2004).
- 26. Takada, Y., and Aggarwal, B.B., J. Immunol., 171, 3278-3286 (2003).
- 27. Zhou, Y., et al., Endocr. Relat. Cancer, 12, S37-S46 (2005).
- 28. Sun, H.L., et al., Oncogene, e-pub Jan. 8, 2007.
- 29. Aggarwal, S., et al., Mol. Pharmacol., 69, 195-206 (2006).
- 30. Shishodia, S., et al., Ann. NY Acad. Sci., 1056, 206-217 (2006). 31. Thangapazham, R.L., et al., AAPS J., 8, E443-E449 (2006).
- 32. Ahn, K.S., et al., Mol. Pharmacol., 71, 209-219 (2007).
- 33. Singh, R.P., and Agarwal, R., Endocr. Relat. Cancer, 13, 751-778 (2006).

- 34. Choi, Y.H., et al., Arch. Pharm. Res, 29, 293-297 (2006).
- 35. Takada, Y., et al., J. Biol. Chem., 280, 17203-17212 (2005).
- 36. Raffoul, J.J., et al., BMC Cancer, 6, 107 (2006).
- 37. Ahn, K.S., et al., Mol. Cancer Res., 4, 621-633 (2006).
- 38. Lee, J., et al., Planta Med., 71, 338-343 (2005).
- 39. Tse, A.K., et al., Biochem. Pharmacol., 70, 1443-1457 (2005).
- 40 Takada, Y., et al., Blood, 106, 641-649 (2005).
- 41. Islam, S., et al., Microbiol. Immunol., 48, 729-736 (2004).
- Jin, C.Y., et al., Pharmacol. Res., 54, 461-467 (2006). 42.
- 43. Birrell, M.A., et al., FASEB J., 19, 840-841 (2005).
- 44. Borra, M.T., et al., J. Biol. Chem., 280, 17187-17195 (2005).
- 45. Lagouge, M., et al., Cell, 127, 1109-1122 (2006).
- 46. Meng, Y., et al., World J. Gastroenterol., 11, 525-528 (2005).
- 47. Chen, P.N., et al., Chem. Biol. Interact., 156, 141-150 (2005).
- 48. Dhanalakshmi, S., et al., Oncogene, 21, 1759-1767 (2002).
- 49. Singh, R.P., et al., Oncogene, 24, 1188-1202 (2005)
- 50. Yoo, H.G., et al., Int. J. Mol. Med., 13, 81-86 (2004)
- 51. Ho, L.J., et al., Brit. J. Pharmacol., 143, 919-927 (2004).
- 52. Albini, A., et al., FASEB J., 20, 527-529 (2006).

The 2008-2009 Sigma Life Science Catalog

Your work is unique...innovative...groundbreaking.

The new 2008-2009 Sigma Life Science Catalog is the perfect match. Reserve you copy today (available January 2008).

To reserve a copy visit us online at sigma.com/sigmacat1

 Over 30.000 Products · 2.000 New Products · New! Life Science Web Tool Guide

SIGMA

g

NF-κB and Inflammation

Antibodies and Kits

Anti-IKKα (699-715) antibody produced in rabbit

Anti-IkB Kinase o

Affinity isolated antibody

Immunogen: synthetic peptide corresponding to amino acids 699-715 of human IkB Kinase α (IKK α) (C2)^{1,2}. This sequence contains one amino acid substitution with the mouse sequence. Solution in phosphate buffered saline containing 0.02% sodium azide

Species reactivity: human

Antigen mol wt 85 kDa

Application(s)

2. Regnier, C.H., et al., Cell, 90, 373-383 (1997) -20°C DRY ICF 0.1 mg

17778

Anti-IKKα (716-734) antibody produced in rabbit

Anti-IkB Kinase o

Affinity isolated antibody

Immunogen: synthetic peptide corresponding to amino acids 716-734 of human IkB Kinase α (IKK α) (C1).

Solution in phosphate buffered saline containing 0.02% sodium azide

Species reactivity: human

Antigen mol wt 85 kDa

Application(s)

Immunoblotting...... 0.5 µg/mL using human HeLa cell extract -20°C DRY ICE 17903

0.1 ma

0.2 mL

Anti-IKKβ antibody produced in rabbit

Anti-IKK2; Anti-I κ B Kinase β

Affinity isolated antibody

Immunogen: synthetic peptide located near the C-terminal region of human IKKβ (amino acids 720-735).

Solution in 0.01 M phosphate buffered saline, pH 7.4, containing 1% bovine serum albumin and 15 mM sodium azide.

Antigen mol wt 85 kDa

Application(s)

Immunoblotting	1:500 using a whole extract
	of human epitheloid carcinoma HeLa cells
-20°C DRY ICE	

Anti-IKKy/NEMO (N-Terminal) antibody produced in rabbit

Anti-IKKAP1

19767

~0.5 mg/mL, affinity isolated antibody

Immunogen: synthetic peptide corresponding to amino acids 400-416 of human IKKγ. Solution in phosphate buffered saline containing 0.02% sodium azide Species reactivity: human, mouse, rat

Antigen mol wt ~52 kDa

Application(s)

Immunoblotting	0.5-1 µg/mL
-20°C DRY ICE	
15032	100 µg

Monoclonal Anti-NF-KB antibody produced in mouse

clone NF-12, ascites fluid

Immunogen: recombinant mouse NF-kB p65 fragment (C-terminal, 151 amino acids).

Contains 15 mM sodium azide

Species reactivity: mouse, human

Antigen mol wt 65 kDa

Application(s)

N8523	0.2 mL
-20°C DRY ICE	
Isotype	lgG1
Array	suitable
	from a mouse fibroblast cell line
Immunoblotting	1:1,000 using a whole cell extract
Immunocytochemistry	suitable

NF-κB (p50) human

Full length, 453 amino acid protein from spliced cDNA containing the p105 precursor polypeptide. For use in gel shift and footprinting assays, as well as in vitro transcriptional activation.

recombinant, expressed in Escherichia coli, buffered aqueous glycerol solution

Solution in 50 mM NaCl, 5 mM DTT, 0.5 mM PMSF, 20 mM HEPES, 10 µM zinc acetate, 0.1% NP-40, and 10% glycerol. Unit Definition : One gel shift unit is defined as the amount of p50 required to shift 0.38 pmol of NF-κB oligonucleotide. -70°C DRY ICE

N9909

Histone Deacetylase Assay Kit, Fluorometric

50 un

Histone deacetylases (HDACs) are enzymes that remove acetyl groups from histone proteins. HDACs are important regulators of gene expression and genome function. They are implicated in a number of human disease states, namely several cancers, neurological disorders, and aging. Therefore, HDAC are relevant key targets for therapeutic intervention. Moreover, histone deacetylase inhibitors have been shown to serve as antitumor agents.

The Histone Deacetylase Assay Kit provides a simple method for the detection of HDAC activity based on a two-step enzymatic reaction. The substrate for this enzymatic assay is a substituted peptide with an acetylated lysine residue and a bound fluorescent group. The first step of the reaction is deacetylation of the acetylated lysine side chain by the HDAC containing sample (HeLa cell extract, purified enzyme, etc.). The second step is the cleavage of the deacetylated substrate by the developer solution and the release of the free highly fluorescent group. The measured fluorescence is directly proportional to the deacetylation activity of the sample.

The Histone Deacetylase Assay Kit has been tested on HeLa, HEK 293T, NIH 3T3, and U 973 cell extracts. The kit includes all the reagents required for the fast and easy measurement of HDAC activity in cell or nuclear extracts, or with purified enzyme preparations. In addition, the kit provides HeLa cell lysate as a source of HDAC activity for inhibitor screening or as a positive control, a HDAC inhibitor (Trichostatin A), and a standard to enable activity quantitation.

Sufficient for 100 assays (96 well plates)

-20°C DRY ICF

CS1010	1 kit

SIGMA

ma-aldrich.com/biofiles

sign

Sirt1 human

SIR2a; SIR2L1; Sirtuin1

Sirtuins are a family of NAD⁺ dependent deacetylases that remove an acetyl group from the ϵ -amino group of lysine residues. The proteins within this family are named after the first protein discovered, from yeast, called Sir2 (Silent Information Regulator 2). The proteins are conserved from bacteria to higher eukaryotes. In humans, there are seven Sir2 family members (SIRT1 to SIRT7). SIRT1 plays a pivotal role in the regulation of cellular differentiation, metabolism, cell cycle, apoptosis and regulation of p53. Several targets for SIRT1 were identified, among them Lys³⁸² of p53.¹ Using RNA interference, additional targets were identified. It was demonstrated that reduced levels of human SIRT1 led to increased acetylation of Histone H4-Lys¹⁶, H4-Lys²⁰, and Histone H3-Lys⁹ as well as histone H1-Lys²⁶.²

Recombinant, expressed in *Escherichia coli*, N-terminal histidine tagged, >90% (SDS-PAGE)

Solution containing 50 mM Tris, pH 7.4, 100 mM NaCl, 1 mM DTT, protease inhibitors (Cat. No. P8340) 1:200, and 10% glycerol (w/v).

Lit cited: 1. Vaziri, H., et al., *Cell* **107**, 149-159 (2001) 2. Vaquero, A., et al., *Mol. Cell.* **16**, 93-105 (2004)

CO/	14

Monoclonal Anti-Sirt1 antibody produced in mouse

Anti-mouse homologue of yeast Sir2; Anti-Sir2a; Anti-Sirtuin1

~2 mg/mL, clone SIR11, purified immunoglobulin

Immunogen: synthetic peptide corresponding to amino acids 722-737 at the C-terminus of mouse Sir2 with C-terminal added lysine.

Solution in 0.01 M phosphate buffered saline, pH 7.4, containing 15 mM sodium azide.

Species reactivity: mouse

Antigen mol wt ~105 kDa

Application(s)

Immunoblotting	1-2 µg/mL using total cell extract of C-2 cells
Indirect ELISA	suitable
Array	suitable
Isotype	IgG1
-20°C DRY ICE	
S5196	200 µL

NEW SIRT1 Assay Kit

The assay procedure is based on a two-step enzymatic reaction. The first step is deacetylation by SIRT1 of a substrate that contains an acetylated lysine side chain. The second step is the cleavage of the deacetylated substrate by the Developing Solution and the release of a highly fluorescent group. The measured fluorescence is directly proportional to the deacetylation activity of the enzyme in the sample. The kit offers all the reagents required for the fast and easy measurement of purified SIRT1 activity and for screening of inhibitors/activators. Moreover, the kit contains an inhibitor (nicotinamide) and an activator (resveratrol) as negative and positive controls, respectively.

Sufficient for 100 assays

Components

Assay Buffer SIRT1 Substrate (Fluorometric) Standard (non-acetylated) 20 mM Developer Solution Nicotinamide Solution (inhibitor) Sirt1 human NAD* Solution Resveratrol Solution (activator) Terror DRY ICE	100 μL 100 μL 1.5 mL 100 μL 150 μg 1 mL
CS1040	1 kit

Anti-Sir2 (AS-16) antibody produced in rabbit

Affinity isolated antibody

Immunogen: synthetic peptide corresponding to mouse Sir2 sequence (amino acids 722-737) with N-terminal cysteine added conjugated to KLH. This sequence is 62% homologous to the corresponding human sequence.

Solution in 0.01 M phosphate buffered saline, pH 7.4 containing 1% bovine serum albumin and 15 mM sodium azide

Antigen mol wt ~110 kDa

Species reactivity: mouse

Application(s)

150 µg

Array	suitable
Immunoblotting	
0	mouse NIH 3T3 cells
Immunoprecipitation	
Indirect immunofluorescence	suitable using
	mouse NIH 3T3 cells
-20°C DRY ICE	
S5313	0.2 mL

Splitomicin

≥98% (HPLC)

Sir2p (silent information regulator) and HDAC inhibitor.

2-8°C

5 mg

19

gma-aldrich.com/biofiles

DNA Damage and Repair Mechanisms

Damage to cellular DNA is involved in mutagenesis and the development of cancer. The DNA in a human cell undergoes several thousand to a million damaging events per day, generated by both external (exogenous) and internal metabolic (endogenous) processes. Changes to the cellular genome can generate errors in the transcription of DNA and ensuing translation into proteins necessary for signaling and cellular function. Genomic mutations can also be carried over into daughter generations of cells if the mutation is not repaired prior to mitosis.

Once cells lose their ability to effectively repair damaged DNA, there are three possible responses (see **Figure 1**).

- The cell may become senescent, i.e., irreversibly dormant. In 2005, multiple laboratories reported that senescence could occur in cancer cells *in vivo* as well as *in vitro*, stopping mitosis and preventing the cell from evolving further.^{1.4}
- The cell may become apoptotic. Sufficient DNA damage may trigger an apoptotic signaling cascade, forcing the cell into programmed cell death.
- 3. The cell may become malignant, i.e., develop immortal characteristics and begin uncontrolled division.




Figure 1. The pathway of cellular DNA damage and repair that leads to senescence, apoptosis, or cancer.

To compensate for the degree and types of DNA damage that occur, cells have developed multiple repair processes including mismatch, base excision, and nucleotide excision repair mechanisms, with little process redundancy. Cells may have evolved to proceed into apoptosis or senescence if overwhelming damage occurs rather than expend energy to effectively repair the damage. The rate at which a cell is able to make repairs is contingent on factors including cell type and cell age.

Sources of DNA Damage

For many years, exogenous sources of damage have been thought to be the primary cause of DNA mutations leading to cancer. However, Jackson and Loeb proposed that endogenous sources of DNA damage also contribute significantly to mutations that lead to malignancy.⁵ Both environmental and cellular sources can result in similar types of DNA damage.

DNA can be attacked by physical and chemical mutagens. Physical mutagens are primarily radiation sources, including UV (200-300 nm wavelength) radiation from the sun. UV radiation produces covalent bonds that crosslink adjacent pyrimidine (cytosine and thymine) bases in the DNA strand. Ionizing radiation (X-rays) initiates DNA mutations by generating free radicals within the cell that create reactive oxygen species (ROS) and result in single-strand and double-strand breaks in the double helix. Chemical mutagens can attach alkyl groups covalently to DNA bases; nitrogen mustard compounds that can methylate or ethylate the DNA base are examples of DNA alkylating agents. Procarcinogens are chemically inert precursors that are metabolically converted into highly reactive carcinogens. These carcinogens can react with DNA by forming DNA adducts, i.e., chemical entities attached to DNA. Benzo/a/pyrene, a polyaromatic heterocycle, is not itself carcinogenic. It undergoes two sequential oxidation reactions mediated by cytochrome P450 enzymes, which results in benzo/a/pyrenediol epoxide (BPDE), the carcinogenic metabolite that is able to form a covalent DNA adduct (see Figure 2)




Figure 2. Benzo[a]pyrene is oxidized by P450 enzymes to create the highly carcinogenic benzo[a]pyrenediolepoxide.

SIGMA

A

DNA damage can also result from endogenous metabolic and biochemical reactions, some of which are not well understood.⁶ Hydrolysis reactions can partially or completely cleave the nucleotide base from the DNA strand. The chemical bond connecting a purine base (adenine or guanine) to the deoxyribosyl phosphate chain can spontaneously break in the process known as depurination. An estimated 10,000 depurination events occur per day in a mammalian cell.⁷ Depyrimidination (loss of pyrimidine base from thymine or cytosine) also occurs, but at a rate 20 to 100-fold lower than depurination.

Deamination occurs within the cell with the loss of amine groups from adenine, guanine, and cytosine rings, resulting in hypoxanthine, xanthine, and uracil, respectively. DNA repair enzymes are able to recognize and correct these unnatural bases. However, an uncorrected uracil base may be misread as a thymine during subsequent DNA replication and generate a $C \rightarrow T$ point mutation.

DNA methylation, a specific form of alkylation, occurs within the cell due to a reaction with S-adenosylmethionine (SAM). SAM is an intracellular metabolic intermediate that contains a highly reactive methyl group. In mammalian cells, methylation occurs at the 5-position of the cytosine ring of a cytidine base (C) that is 5' to a guanosine base (G), i.e., sequence CpG. A significant source of mutation error is the spontaneous deamination of the 5-methyl-cytosine product of methyl-ation. Loss of the amine group results in a thymine base, which is not detected by DNA repair enzymes as an unnatural base. The resulting substitution is retained in DNA replication, creating a C \rightarrow T point mutation (see **Figure 3**).




Figure 3. The 2-phase mutation of cytosine results in thymine, creating a C→T point mutation.

Normal metabolic processes generate reactive oxygen species (ROS), which modify bases by oxidation. Both purine and pyrimidine bases are subject to oxidation. The most common mutation is guanine oxidized to 8-oxo-7,8-dihydroguanine, resulting in the nucleotide 8-oxo-deoxyguanosine (8-oxo-dG). The 8-oxo-dG is capable of base pairing with deoxyadenosine, instead of pairing with deoxycytotidine as expected. If this error is not detected and corrected by mismatch repair enzymes, the DNA subsequently replicated will contain a C \rightarrow A point mutation. ROS may also cause depurination, depyrimidination, and single-strand or double strand breaks in the DNA.

Other genomic mutations may be introduced during DNA replication in the S phase of the cell cycle. Polymerases that duplicate template DNA have a small but significant error rate, and may incorporate an incorrect nucleotide based on Watson-Crick pairing versus the template DNA. Chemically altered nucleotide precursors may be incorporated into the generated DNA by the polymerase, instead of normal bases. In addition, polymerases are prone to "stuttering" when copying sections of DNA that contain a large number of repeating nucleotides or repeating sequences (microsatellite regions). This enzymatic "stuttering" is due to a strand slippage, when the template and replicated strands of DNA slip out of proper alignment. As a result, the polymerase fails to insert the correct number of nucleotides indicated by the template DNA, resulting in too few or too many nucleotides in the daughter strand.

Single strand and double strand cleavage of the DNA may occur. Single strand breaks may result from damage to the deoxyribose moiety of the DNA deoxyribosylphosphate chain. Breaks also result as an intermediate step of the base excision repair pathway after the removal of deoxyribose phosphate by AP-endonuclease 1.⁸ When a single strand break occurs, both the nucleotide base and the deoxyribose backbone are lost from the DNA structure. Double strand cleavage most often occurs when the cell is passing through S-phase, as the DNA may be more susceptible to breakage while it is unraveling for use as a template for replication.

Mechanisms of DNA Repair

While the cell is able to evolve into either an apoptotic or senescent state, these actions are performed as a last resort. For each type of DNA damage, the cell has evolved a specific method of repairing the damage or eliminating the damaging compound.

O⁶-Methylguanine DNA methyltransferase (MGMT; DNA alkyltransferase) cleaves both methyl and ethyl adducts from guanine bases on the DNA structure. The reaction is not a catalytic (enzymatic) reaction but is stoichiometric (chemical), consuming one molecule of MGMT for each adduct removed. Cells that have been engineered to overexpress MGMT are more resistant to cancer, likely because they are able to negate a larger amount of alkylating damage. A recent study by Niture, et al., reports an increase in MGMT expression by use of cysteine/glutathione enhancing drugs and natural antioxidants.⁹

DNA polymerases such as polymerase- δ contain proofreading activities and are primarily involved in replication error repair. When an error is detected, these polymerases halt the process of DNA replication, work backward to remove nucleotides from the daughter DNA chain until it is apparent that the improper nucleotide is gone, and then reinitiate the forward replication process. Mice with a point mutation in both copies of the *Pold*1 gene demonstrated a loss of proofreading activity by DNA polymerase- δ and developed epithelial cancers at a significantly higher rate than did mice with wild-type genes or with a single copy mutation.¹⁰

A group of proteins known as mismatch excision repair (MMR) enzymes is capable of correcting errors of replication not detected by the proofreading activities of DNA polymerase. MMR enzymes excise an incorrect nucleotide from the daughter DNA and repair the strand using W-C pairing and the parent DNA strand as the correct template.¹¹ This is especially crucial for errors generated during the replication of microsatellite regions, as the proofread-ing activity of DNA polymerase does not detect these errors. To a lesser degree, MMR enzymes also correct a variety of base pair anomalies resulting from DNA oxidation or alkylation. These mutations include modified base pairs containing O⁶-methylguanine and 8-oxoguanine, and carcinogen and cisplatin adducts.^{12,13} Mutations in the human mismatch excision repair genes *MSH2* and *MLH1* are associated with hereditary non-polyposis colorectal cancer (HNPCC) syndrome.¹⁴

SIGMA

Base Excision Repair and Nucleotide Excision Repair

Base excision repair (BER) involves multiple enzymes to excise and replace a single damaged nucleotide base. The base modifications primarily repaired by BER enzymes are those damaged by endogenous oxidation and hydrolysis. A DNA glycosylase cleaves the bond between the nucleotide base and ribose, leaving the ribose phosphate chain of the DNA intact but resulting in an apurinic or apyrimidinic (AP) site. 8-Oxoguanine DNA glycosylase I (Ogg1) removes 7,8-dihydro-8-oxoguanine (8-oxoG), one of the base mutations generated by reactive oxygen species. Polymorphism in the human *OGG1* gene is associated with the risk of various cancers such as lung and prostate cancer. Uracil DNA glycosylase, another BER enzyme, excises the uracil that is the product of cytosine deamination, thereby preventing the subsequent $C \rightarrow T$ point mutation.¹⁵ N-Methylpurine DNA glycosylase (MPG) is able to remove a variety of modified purine bases.¹⁶

The AP sites in the DNA that result from the action of BER enzymes, as well as those that result from depyrimidination and depurination actions, are repaired by the action of AP-endonucle-ase 1 (APE1). APE1 cleaves the phosphodiester chain 5' to the AP site. The DNA strand then contains a 3'-hydroxyl group and a 5'-abasic deoxyribose phosphate. DNA polymerase β (Pol β) inserts the correct nucleotide based on the corresponding W-C pairing

and removes the deoxyribose phosphate through its associated AP-lyase activity. The presence of X-ray repair cross-complementing group 1 (XRCC1) is necessary to form a heterodimer with DNA ligase III (LIG3). XRCC1 acts as a scaffold protein to present a non-reactive binding site for Polβ, and bring the Polβ and LIG3 enzymes together at the site of repair.¹⁷ Poly(ADP-ribose) polymerase (PARP-1) interacts with XRCC1 and Polβ and is a necessary component of the BER pathway.^{18,19} The final step in the repair is performed by LIG3, which connects the deoxyribose of the replacement nucleotide to the deoxyribosylphosphate backbone. This pathway has been named "short-patch BER".²⁰

An alternative pathway called "long-patch BER" replaces a strand of nucleotides with a minimum length of 2 nucleotides. Repair lengths of 10 to 12 nucleotides have been reported.^{21,22} Long-patch BER requires the presence of proliferation cell nuclear antigen (PCNA), which acts as a scaffold protein for the restructuring enzymes.²³ Other DNA polymerases, possibly Polð and Pole, ²⁴ are used to generate an oligonucleotide flap. The existing nucleotide sequence is removed by flap endonuclease-1 (FEN1). The oligonucleotide is then ligated to the DNA by DNA ligase I (LIG1), sealing the break and completing the repair.¹⁷ The process used to determine the selection of short-patch versus long patch BER pathways is still under investigation (see **Figure 4**).²⁵




Figure 4. Schematic of both short-patch and long-patch BER pathways.

SIGMA

While BER may replace multiple nucleotides via the long-patch pathway, the initiating event for both short-patch and long-patch BER is damage to a single nucleotide, resulting in minimal impact on the structure of the DNA double helix. Nucleotide excision repair (NER) repairs damage to a nucleotide strand containing at least 2 bases and creating a structural distortion of the DNA. NER acts to repair single strand breaks in addition to serial damage from exogenous sources such as bulky DNA adducts and UV radiation.²⁶ The same pathway may be used to repair damage from oxidative stress.²⁷ Over 20 proteins are involved in the NER pathway in mammalian cells, including XPA, XPC-hHR23B, replication protein A (RPA), transcription factor TFIIH, XPB and XPD DNA helicases, ERCC1-XPF and XPG, Polo, Pole, PCNA, and replication factor C.²⁸ Overexpression of the excision repair cross-complementing (ERCC1) gene has been associated with cisplatin resistance by non-small-cell lung cancer cells²⁹ and corresponds to enhanced DNA repair capacity.³⁰ Global genomic NER (GGR) repairs damage throughout the genome, while a specific NER pathway called Transcription Coupled Repair (TCR) repairs genes during active RNA polymerase transcription.³¹


Repair of Double-Strand Breaks

Double-strand breaks in DNA can result in loss and rearrangement of genomic sequences. These breaks are repaired by either nonhomologous end-joining (NHEJ) or by homologous recombination (HR), also called recombinational repair or template-assisted repair. The HR pathway is activated when the cell is in late S/G₂ phase and the template has recently been duplicated. This mechanism requires the presence of an identical or nearly identical sequence linked to the damaged DNA region via the centromere for use as a repair template. Double-stranded breaks repaired by this mechanism are usually caused by the replication machinery attempting to synthesize across a single-strand break or unrepaired lesion, resulting in the collapse of the replication fork.

Non-homologous end-joining (NHEJ) is used at other points of the cell cycle when sister chromatids are not available for use as HR templates. When these breaks occur, the cell has not yet replicated the region of DNA that contains the break, so unlike the HR pathway, there is no corresponding template strand available. In NHEJ, the Ku heterodimeric protein positions the two ends of the broken DNA strands for repair without an available template, losing sequence information in the process. Multiple enzymes are involved in the rejoining process, including DNA ligase IV, XRCC4, and DNA-dependent protein kinase (DNA-PK).^{32,33} NHEJ is inherently mutagenic as it relies on chance pairings, called microhomologies, between the single-stranded tails of the two DNA fragments to be joined (see Figure 5). In higher eukaryotes, DNA-PK is required for NHEJ repair, both via the primary mechanism and via an alternative back-up mechanism (D-NHEJ).3

Future Applications

While DNA damage is a key factor in the development and evolution of cancer cells, continued damage is used as part of clinical treatments for cancer, forcing malignant cells into apoptosis or senescence. Many chemotheraputic drugs such as bleomycin, mitomycin, and cisplatin, are effective because they cause further DNA damage in cancer cells that replicate at a faster rate than surrounding tissue. Cellular DNA repair mechanisms are a doubleedged sword; by reducing mutations that may lead to cancer, these processes strive for genomic integrity, but the same mechanisms in malignant cells allow those cells to survive additional DNA damage and continue uncontrolled growth. In order to block this survival mechanism within cancer cells, clinical trials are now being performed using inhibitors to specific DNA repair enzymes, including MGMT, PARP, and DNA-PK.35-38



Loss of base pairs present in original wild-type sequence

Figure 5. General mechanism of NHEJ repair of double-strand breaks in DNA.

References:

- Collado, M., et al., Nature, 436, 642 (2005).
- 2 Chen, Z., et al., Nature, 436, 725 (2005).
- 3. Michaloglou, C., et al., Nature, 436, 720 (2005).
- 4. Braig, M., et al., Nature, 436, 660 (2005).
- Jackson, A.L., and Loeb, L.A., Mutat. Res., 477, 7-21, (2001). 5
- De Bont, R., and van Larebeke, N., Mutagenesis, 19, 169-85 (2004). 6.
- Lindahl, T. and Nyberg, B., Biochemistry, 11, 3610-3618 (1972).
- 8. Brem., R. and Hall, J., Nucleic Acids Res., 33, 2512-2520 (2005).
- 9. Niture S.K., et al, Carcinogenesis, Epub ahead of print (Aug. 31, 2006)
- 10. Goldsby, et al., Proc. Natl. Acad. Sci. USA, 99, 15560 (2002).
- 11. Yang, W., Mutat. Res., 460, 245-56 (2000).
- lyer, R.R., et al., (Review), Chem. Rev., 106, 302-23 (2006). 12
- 13. Modrich. P. J. Biol. Chem., 281, 30305-30309 (2006).
- 14. Muller, A. and Fishel, R., Cancer Invest., 20, 102-9 (2002).
- 15. Lindahl, T., Proc. Natl. Acad. Sci. USA, 71, 3649-3653 (1974).
- 16. Singer, B., and Hang, B., Chem. Res. Toxicol., 10, 713-32 (1997).
- 17. Lindahl, T. and Wood, R.D., Science, 286, 1897-1905 (1999)
- 18. Caldecott, K.W., et al., Nucleic Acids Res., 24, 4387-4394 (1996).
- Dantzer, F., et al., Biochemistry, 39, 7559-7569 (2000). 19
- Srivastava, D.K., et al., J. Biol. Chem., 273, 21203-9 (1998) 20
- 21. Ranalli, T.A., et al., J. Biol. Chem., 277, 41714-41724 (2002).
- 22 Sattler, U., et al., EMBO Reports, 4, 363-367 (2003).
- 23. Fortini, P., et al., *Biochemistry*, **37**, 3575-3580 (1998)
- 24. Klungland, A. and Lindahl, T., EMBO Journal, 16, 3341-8 (1997) 25. Sung, J.S. and Demple, B., (Review), FEBS Journal, 273, 1620-1629
- (2006)
- 26. Balajee, A.S. and Bohr, V.A., Gene, 250, 15-30 (2000)
- 27 Gros, L., et al., Oncogene, 21, 8905-8925 (2002).
- 28 You, J.S., J. Biol. Chem., 278, 7476-7485 (2003).
- Rosell, R., et al., Cancer Control, 10, 297-305 (2003)
- Vogel, U., et al., Mutat. Res., 461, 197-210 (2000). 30
- 31. Hanawalt, P.C., Oncogene, 21, 8949-8956 (2002).
- 32. Critchlow, S.E. and Jackson, S.P., Trends Biochem. Sci., 23, 394-398 (1998)
- 33. Wang, H., et al., Nucleic Acids Res., 31, 5377-88 (2003).
- 34. Perrault, R., et al., J. Cell. Biochem., 92, 781-794 (2004).
- 35. Sanchez-Perez, I., Clin. Transl. Oncol., 8, 642-6 (2006).
- 36. Madhusudan, S. and Hickson, I.D., Trends Mol. Med., 11, 503-11 (2005).
- 37. Plummer, E.R., Curr. Opin. Pharmacol., 6, 364-8 (2006).
- 38. Sabharwal, A. and Middleton, M.R., Curr. Opin. Pharmacol., 6, 355-63 (2006).

DNA Damage and Repair

ma-aldrich.com/biofiles

<u>0</u>

IGMA

g

DNA Alkylating Agents

Compound	Synonyms	Biological Action	Cat. No.	Description	Pack Size
Angelicin		Upon long wavelength UV irradiation, forms monoadduct with double-stranded DNA; inhibits DNA and RNA synthesis and cell replication in Ehrlich ascites tumor cells.	A0956		10 mg 25 mg
Azoxymethane	AOM	Carcinogen that induces O ⁶ -methylguanine adducts in DNA leading to G to A transitions. Induces colon tumors in laboratory animals; used to study the mechanism of cancer progression and chemoprevention.	A2853	13.4 M, ≥90% (GC), contains ≤15% ethanol and acetic acid	25 mg 100 mg
Benzo <i>[a]</i> pyrene	Benzo <i>[def]</i> chrysene; 3,4-Benzopyrene; 3,4- Benzpyrene	A polycyclic aromatic hydrocarbon (PAH), environmental pollutant and carcinogen. Forms depurinating adducts with DNA that dissociate to form abasic lesions.	B1760	≥96% (HPLC)	1 g
Bleomycin sulfate from <i>Streptomyces</i> verticillus	Blenoxane; Bleo; Blexane	Anticancer antibiotic that binds to DNA, inhibits DNA synthesis and causes DNA scissions at specific base sequences. Needs to bind oxygen and a metal ion such as copper or iron to cleave DNA. Metallobleomycin conjugates also cleave tRNA and DNA-RNA hybrids.	B5507	Activity: 1.2-1.7 units/mg solid A mixture of bleomycin sulfate salts Cu <0.1%	15 units
			B2434	Biotechnology Performance Certified Cell culture tested Potency: 1.5-2.0 units/mg	20 mg
Chlorambucil	4-[Bis(2-chloroethyl)amino] benzenebutyric acid	Anticancer agent that alkylates DNA and induces apoptosis by a p53-dependent mechanism.		See C0 2	253 on page 9
Daunorubicin hydrochloride	Daunomycin hydrochloride	Anthracycline anticancer antibiotic. Induces damage of DNA by intercalating into DNA and inhibiting topoisomerase II.		See D8	809 on page 9
3,3'-Dichlorobenzidine dihydrochloride	DCB	Industrial carcinogen that induces liver and bladder tumors in laboratory animals.	D9886		25 g
7,12-Dimethyl- benz[a]anthracene	9,10-Dimethyl-1,2- benzanthracene; DMBA	Carcinogen that is activated to form covalent adducts with DNA. Used to induce colon, skin and mammary tumors in mice.	D3254	≥95%	100 mg 1 g 5 g
N,N'-Dimethylhydrazine dihydrochloride	1,2-Dimethylhydrazine dihydrochloride; DMH; N,N-Dimethylhydrazine (sym.) dihydrochloride	Methylating agent that induces O ⁶ -methyl- guanine adducts in DNA; induces colon cancer in rats and mice.	D161802	≥99%	1 g 10 g
Doxorubicin hydrochloride	Adriamycin hydrochloride; DOX; Hydroxydaunorubicin hydrochloride	Anthracycline anticancer antibiotic. Induces damage of DNA by intercalating into DNA and inhibiting topoisomerase II.	See D1515 on page		5 15 on page 9
L-Ethionine	L-2-Amino-4-(ethylthio)- butyric acid	An analog of methionine that interferes with the normal methylation of DNA and other methylation pathways, and induces pancreatic toxicity and liver cancer. However, in animal studies, it has also been shown to act synergistically with methionine-depletion to block the growth and metastasis of methionine- dependent tumors.	E1260	~99% (TLC)	250 mg 1 g
N-(2-Fluorenyl)acetamide	2-AAF; 2-Acetamido- fluorene; N-Acetyl-2- aminofluorene	A genotoxic carcinogen that forms adducts with DNA following its N-hydroxylation by CYP1A2; used to model liver carcinogenesis in rat.	A7015	≥90% (TLC)	5 g 25 g
Hexestrol	meso-3,4-Bis(4-hydroxy- phenyl)hexane; 4,4'-(1,2- Diethylethylene)diphenol; Dihydrodiethylstilbestrol	Nonsteroidal synthetic estrogen; reduced analog of diethylstilbestrol. Hexestrol is metabolized to its quinone, which forms depurinating adducts with DNA leading to oncogenic mutations and cancer initiation.	H7753	≥98%	1 g 5 g
Idarubicin hydrochloride	DMDR; Idamycin; IMI-30	Anthracycline anticancer antibiotic. Groove- binding intercalator of DNA with higher capacity than daunorubicin; topoisomerase II inhibitor.	l1656		10 mg
Methoxyamine hydrochloride	Methoxylamine hydrochloride; O-Methylhydroxylamine hydrochloride	Binds to abasic sites on DNA and acts as an inhibitor of base excision repair.	226904	98%	1 g 5 g 25 g
Methyl methanesulfonate	5	Randomly methylates DNA in both the major (7-MeG) and minor (3-MeA and 3-MeG) grooves.	M4016	~99% (GC)	1 g 5 g 25 g

DNA Damage and Repair

sigma-aldrich.com/biofiles

DNA Alkylating Agents (continued)

Compound	Synonyms	Biological Action	Cat. No.	Description	Pack Size
Mitoxantrone dihydrochloride	1,4-Dihydroxy-5,8- bis[[2-[(2-hydroxyethyl)- amino]ethyl]-amino]- 9,10-anthracenedione dihydrochloride	DNA intercalating drug. Inhibits DNA synthesis.			5 45 on page 9
N-Nitrosodiethylamine	Diethylnitrosamine; NDEA	Ethylating agent that is carcinogenic in all animal species tested. Main target organs are the nasal	N0258	ISOPAC ^{®1}	1 g
		cavity, trachea, lung, esophagus and liver.	N0756	Bulk Package	10 mL 25 mL
N-Nitroso-N-ethylurea	ENU; N-Ethyl-N-nitroso- urea	DNA alkylating agent that is carcinogenic in many animal species. Induces benign and	N3385	ISOPAC®1	1 g
		malignant tumors in numerous organs, including nervous system, stomach, esophagus, pancreas, respiratory tract, intestine, lymphoreticular tissues, skin and kidney.	N8509	Bulk Package Contains aqueous acetic acid as stabilizer Weight on dry basis	5 g 25 g
N-Nitroso-N-methylurea	MNU; N-Methyl-N-nitroso- urea	DNA alkylating agent that is carcinogenic in many animal species. Induces benign and	N1517	ISOPAC ^{®1}	1 g
		malignant tumors in numerous organs, including nervous system, stomach, esophagus, pancreas, respiratory tract, intestine, lymphoreticular tissues, skin and kidney.	N4766	Bulk Package Contains aqueous acetic acid as stabilizer Weight on dry basis	25 g 100 g
1-Nitrosopiperidine	NPIP	Carcinogen that induces benign and malignant tumors of the respiratory tract, stomach and esophagus in rodents. Potent nasal carcinogen in rats.	N6007		5 mL
1-Nitrosopyrrolidine	N-Nitrosopyrrolidine; NPYR	Carcinogen that forms DNA adducts that primarily result in A:T to G:C mutations. Induces hepatocellular carcinomas and lung adenomas in mice.	158240	99%	10 g
Ochratoxin A	N-[(3R)-(5-Chloro-8- hydroxy-3-methyl-1-oxo- 7-isochromanyl)carbonyl]- L-phenylalanine	Mycotoxin found in food that is nephrotoxic and carcinogenic in the kidney. Its metabolites are genotoxic and form adducts with DNA in kidney cells from several species.	O1877	from Aspergillus ochraceus (Aspergillus oryzae), ≥98% (TLC) Benzene-Free	1 mg 5 mg 25 mg
2,3,4,5,6-Pentafluoro- benzyl bromide	α-Bromo-2,3,4,5,6-penta- fluorotoluene	Used to derivatize <i>N</i> -7-substituted guanine adducts of DNA for determination by GC- electron capture mass spectrometry		99%	1 g 5 g 25 g
Rebeccamycin		Rebeccamycin intercalates into the DNA and inhibits topoisomerase I by stabilizing the covalent topoisomerase I-DNA intermediate (the cleavable complex) resulting in increased DNA strand breaks.	R4900	≥98% (HPLC), from Saccharothrix aerocolonigenes	250 µg
Safrole	5-Allyl-1,3-benzodioxole; 4-Allyl-1,2-methylene-di- oxybenzene	Genotoxic hepatocarcinogen; metabolites of safrole form adducts with DNA and induce chromasomal aberrations and sister chromatid exchanges.	S9652	≥97%	50 mL 500 mL

 Image: Line Control of the control strength solution without exposure.

SIGMA

DNA Crosslinking Agents

Compound	Ind Synonyms DNA Crosslinking Action (Synonyms DNA Crosslinking Action Cat. No. De				Pack Size
4'-Aminomethyltrioxsalen hydrochloride	4'-Aminomethyl- 4,5',8-trimethylpsoralen hydrochloride	Used to inactivate DNA and RNA viruses, including HIV-1, by nucleic acid crosslinking followed by UV irradiation.	A4330	≥98%	5 mg		
Carboplatin	<i>cis</i> -Diammine(1,1-cyclo- butanedicarboxylato) platinum	A platinum-based antineoplastic drug that forms intrastrand crosslinks with neighboring guanine residues in DNA. Tumors acquire resistance to these drugs through the loss of DNA-mismatch repair (MMR) activity.	C2538	≥98% (HPLC)	100 mg 250 mg		
Carmustine	BCNU; 1,3-Bis(2-chloro- ethyl)-1-nitrosourea	DNA alkylating agent causing DNA interstrand crosslinks.		See C0	400 on page 9		
CB 1954	5-(1-Aziridinyl)-2,4-di- nitrobenzamide	Anticancer prodrug that is activated by nitroreductase to form a potent DNA crosslinking agent.			50 mg		
Cyclophosphamide	Cytoxan	Cytotoxic nitrogen mustard derivative used in cancer chemotherapy. Cyclophosphamide crosslinks DNA, causes strand breakage, and induces mutations.		See C0768, C7397 on page			
cis-Diammineplatinum(II) dichloride	Cisplatin	Potent platinum-based antineoplastic agent. Forms cytotoxic adducts with the DNA dinucleotide d(pGpG), inducing intrastrand crosslinks.		See P4394 on page 9			
Melphalan	4-[Bis(2-chloro- ethyl)amino]-L-phenyl- alanine; L-PAM; L-Phenyl- alanine mustard	Antineoplastic agent that forms DNA intrastrand crosslinks by bifunctional alkylation of 5'-GGC sequences.	M2011	M2011 ≥95% (TLC) 100 250			
Mitomycin C	MMC	Following enzymatic reduction, activated mitomycin C forms covalent adducts and crosslinks with CpG sequences in DNA.	See M0503 on page 9				
Nimustine hydrochloride	N'-[(4-Amino-2-methyl- 5-pyrimidinyl)methyl]- N-(2-chloroethyl)-N- nitrosourea hydrochloride; ACNU	Antineoplastic compound that is a bifunctional DNA alkylating agent.	N8659 1		1 g		
Oxaliplatin	[SP-4-2-(1R- <i>trans</i>)]- (1,2-Cyclohexanediamine- N,N')[ethanedioata (2–)-O,O'] platinum	Platinum-based antitumor agent that forms interstrand crosslinks with neighboring guanine residues in DNA.	09512		5 mg		
Psoralen	7H-Furo[3,2- <i>g</i>] benzopyran-7-one; Furo[3,2- <i>g</i>]coumarin	When activated by UV light, psoralen induces interstrand crosslinks in DNA. Photochemical reagent for the investigation of nucleic acid structure and function.	erstrand crosslinks in DNA. Photochemical agent for the investigation of nucleic acid		10 mg 25 mg 100 mg		
Trioxsalen	TMP; 4,5',8-Trimethyl- psoralen; Trisoralen	Photochemical crosslinker of DNA has been used as a probe for nucleic acid structure and function.	T6137	≥98% (HPLC)	100 mg 500 mg 1 g		



DNA Synthesis Inhibitors

Compound	Synonyms	Inhibition Specificity	Cat. No.	Description	Pack Size
Aminopterin	4-Aminofolic acid; 4-Amino-PGA; 4-Amino- pteroyl-L-glutamic acid	Blocks DNA synthesis by blocking the production of tetrahydrofolate cofactors required for the synthesis of thymidine. In the cell, it is converted by folylpolyglutamate synthase to a high molecular weight polyglutamate metabolite that binds to dihydrofolate reductase and inhibits its activity. Aminopterin is more potent but more toxic than methotrexate.		~98% (TLC)	50 mg 100 mg 1 g
Methotrexate hydrate	(+)-Amethopterin hydrate; 4-Amino-10-methylfolic acid; Methylaminopterin; MTX			See Ad	5770 on page 9
arabinofuranoside binofuranosyl)cytosine; replication by fo Arabinosylcytosine; topoisomerase I		Ara-C incorporates into DNA and inhibits DNA replication by forming cleavage complexes with topoisomerase I resulting in DNA fragmentation; does not inhibit RNA synthesis.	C1768	Free base ≥90% (HPLC)	100 mg 500 mg 1 g 5 g
	Cytosine arabinoside		C6645	Hydrochloride ≥99% (HPLC)	25 mg 100 mg 500 mg 1 g 5 g
5-Fluoro-5'-deoxyuridine	5'dFUrd; Doxifluridine	5-Fluoro-5'-deoxyuridine is a prodrug that is converted to 5-fluorouracil by uridine phosphorylase.	F8791	≥98% (TLC)	25 mg 100 mg 250 mg
5-Fluorouracil	2,4-Dihydroxy-5-fluoro- pyrimidine; 5-Fluoro- 2,4(1 <i>H</i> ,3 <i>H</i>)-pyrimidine- dione; 5-FU	5-Fuorouracil (5-FU) is a potent antitumor agent that inhibits thymidylate synthetase thus depleting intracellular dTTP pools. 5-FU is metabolized to ribonucleotides and deoxyribonucleotides, which can be incorporated	F6627	≥99% (TLC)	1 g 5 g 10 g
		into RNA and DNA. Treatment of cells with 5-FU leads to an accumulation of cells in S-phase and induces p53 dependent apoptosis.	47576	BioChemika, ≥99.0% (HPLC)	1 g 5 g 25 g
Ganciclovir		Prodrug nucleoside analog that is phosphorylated following expression of a viral suicide gene encoding thymidine kinase. The phosphorylated active analog is incorporated into the DNA of replicating malignant and normal eukaryotic cells, causing cell death. The cell cycle is irreversibly arrested at the G ₂ M checkpoint.	G2536	≥99% (HPLC)	100 mg



DNA Repair Enzymes and Antibodies

Monoclonal Anti-AP Endonuclease antibody produced in mouse

Anti-APE1; Anti-APEN; Anti-APE/Ref1; Anti-APEX1; Anti-APEX nuclease; Anti-AP lysase; Anti-Apurinic/Apyrimidinic Endonuclease; Anti-HAP1

~2 mg/mL, clone APEREF, purified immunoglobulin

Immunogen: recombinant human AP endonuclease.

Solution in 0.01 M phosphate buffered saline, pH 7.4, containing 15 mM sodium azide.

Species reactivity: human, dog, rat, mouse

Antigen mol wt ~37 kDa

Application(s):

A2105	200 µL
-20°C DRY ICE	
lsotype	IgG1
Array	suitable
Immunocytochemistry	suitable
Indirect ELISA	suitable
Immunoblotting 0.5-1 µg/mL using total cell externation	ract of Raji cells

Monoclonal Anti-DNA Ligase I antibody produced in mouse

~1 mg/mL, clone 1A9, purified immunoglobulin

Immunogen: bovine DNA ligase I.

Solution in phosphate buffered saline, pH 7.4, containing 0.08% sodium azide.

Species reactivity: human, bovine

Antigen mol wt ~125 kDa

Application(s):

	IgG1 100 μα
0	1-5 µg/mL using HeLa or DiFi cells

Monoclonal Anti-Excision Repair Cross Complementing Protein-1 antibody produced in mouse

Anti-ERCC-1

~1 mg/mL, clone 8F1, purified immunoglobulin

Immunogen: recombinant full length human ERCC1 protein. Solution in phosphate buffered saline pH 7.4, containing 0.08% sodium azide.

Species reactivity: human, rat Antigen mol wt 33-36 kDa

Application(s):

Immunoblotting	
Immunoprecipitation	2 µg using 1 mg protein lysate
Immunohistochemistry	
Isotype	lgG1
-20°C WET ICE	_
F8903	100 µg

Methylpurine DNA Glycosylase human

AAG; MPG

MPG is a base excision repair (BER) protein that removes the mutated N-methylpurine nucleotide from alkylated DNA, creating an apurinic/apyrimidinic site. MPG recognizes and excises 3-methyladenine, 7-methylguanine, 3-methylguanine, N-6-ethanoladenine, and hypoxanthine.

≥80% (SDS-PAGE), recombinant, expressed in *Escherichia coli* (as a MBP fusion protein), activity: ≥5000 units/mg protein

Solution in 50% glycerol (w/v) containing 20 mM Tris-HCI, pH 7.4, 1 mM DTT, 1 mM EDTA, 250 mM NaCl, 0.25% CHAPS, 1% protease inhibitor cocktail (v/v) (Cat. No. P8340).

Apparent mol wt ~70 kDa

Unit definition: One unit will cleave 50% of 0.5 pmole of double stranded DNA oligomer substrate containing a hypoxanthine site lesion in 10 min. at 37 $^\circ$ C.

-70°C DRY ICE

25 µg

Monoclonal Anti-Methylpurine DNA N-Glycosylase antibody produced in mouse

Anti-AAG; Anti-3-Alkyladenine DNA glycosylase; Anti-Alkylpurine DNA N-glycosylase; Anti-APNG; Anti-Methyladenine DNA Glycosylase; Anti-MPG

~2 mg/mL, clone MPG80, purified immunoglobulin Immunogen: recombinant human MPG protein.

Solution in 0.01 M phosphate buffered saline, pH 7.4, containing 15 mM sodium azide.

Species reactivity: human

Antigen mol wt ~32 kDa

Application(s):

M6195	25 µl, 200 µl
-20°C DRY ICE	
Isotype	lgG1
Immunocytochemistry	suitable
Indirect ELISA	suitable
	of 293T cells expressing MPG
Immunoblotting	0.5-1 µg/mL using total cell extract

Monoclonal Anti-Mismatch Repair Protein 2 antibody produced in mouse

Anti-MSH2

~1 mg/mL, clone 2MSH01, tissue culture supernatant

Immunogen: recombinant human mismatch repair protein 2 (MSH2).

Contains 15 mM sodium azide.

Species reactivity: human

Antigen mol wt 102 kDa

Application(s):

Immunohistochemistry (formalin-fixed, paraffin-embedded sections)	1:25
Isotype	
M6315	500 µL

28

Ogg1 Assay Kit

The assay is based on Ogg1 glycosylase activity that recognizes and removes the mutated base (8-oxo-G). The substrate is a 23 oligonucleotide containing 8-oxo-dG at its 11th base, labeled with ³²P at its 5' end, and annealed to its complementary strand (containing dC at the opposite base position to the 8-oxo-dG). Upon cleavage of the substrate by the Ogg1 enzyme, the oligonucleotide strands are run on a denaturing gel and a 10 base fragment (labeled cleavage product) is revealed in addition to the original 23 base oligonucleotide band. The detection is performed by autoradiography.

The Ogg1 Assay Kit:

- Provides a tool for measurement of Ogg1 enzyme activity
- Contains control Ogg1 enzyme enabling screening of Ogg1 modulators
- · Can be used with cell lysates and purified enzyme preparation

Sufficient for 180 assays

Components

Ogg1 Substrate 8-oxo-G strand	1 vial
Ogg1 Substrate complementary strand	1 vial
Reaction Buffer	10 X 1.5 mL
Stop Solution	1mL
Ogg1 (control enzyme)	10 µg
Putrescine	1 mL
SigmaSpin [™] Post-Reaction Purification Columns	2 each
CS0710	1 kit

CS0710

Ogg1 from mouse

8-Oxoguanine-DNA glycosylase

≥90% (SDS-PAGE), recombinant, expressed in Escherichia coli, activity: >20,000 units/mg protein

Solution in 50% (w/v) glycerol containing 20 mM Tris, pH 7.5, 200 mM sodium chloride, 1 mM EDTA, and 1 mM DTT.

Apparent mol wt 38.8 kDa (345 amino acid protein)

Unit definition: One unit is the amount of protein that cleaves 50% of 0.5 pmol double-stranded, 8-oxoguanine-mutated DNA oligomer substrate in 10 min at 37 °C.

-20°C WET ICE	
O2135	10 µg

Uracil DNA Glycosylase from Escherichia coli

DNA Uracil Glycosylase; UDG; Uracil N-glycosylase [59088-21-0]

UDG catalyzes the removal of uracil residues from both single- and double-stranded DNA, but not RNA. This reaction leaves the DNA sugar-phosphodiester backbone intact. The resulting DNA is not suitable for use as a hybridization target or as a template for DNA polymerases.

Solution in 30 mM Tris-HCl, pH 7.5, 1 mM DTT, 0.05% (w/v) Tween® 20, 1 mM EDTA, 150 mM NaCl, 50% (v/v) glycerol.

Unit definition: One unit catalyzes the release of 1 nmol of free uracil from ³H-poly(dU) in 1 hr at 37°C.

Concentration.....1 unit/µL Purchase of this product conveys the licensed right under U.S. Patent No. 5,035,996 and foreign equivalents to use this product only in internal research conducted by the purchaser. No rights under the aforementioned patents are conveyed which permit the resale, transfer, or use for purposes other than research. Rights under the aforementioned patents, for purposes other than internal research, may be obtained by contacting Life Technologies, Inc. -20°C DRY ICE

U1257

100 units



Sigma-Aldrich cytochemistry assays are used in the diagnosis of acute myeloid leukemia and in the staining of inflammatory tissue exudates for characterizing neutrophils, lymphocytes and macrophages. The prepackaged kits are a cost effective way to bring testing in-house with no capital equipment requirements

Classic methods include:

- Naphthol AS-D Chloroacetate Esterase
- Tartrate Resistant Acid Phosphatase
- α-Naphthyl Acetate Esterase
- Mveloperoxidase
- Nitroblue Tetrazolium

Cytochemistry

- Periodic Acid Schiff
- Alkaline Phosphatase
- α-Naphthyl Butyrate Esterase

From traditional hematoxylin and eosin stains to special stain kits designed to visualize specific tissue features Sigma-Aldrich offers a comprehensive solution of quality histology reagents you need and the support you demand.

Please visit sigma-aldrich.com/handh for all your histology needs



SIGMA

g

ma-aldrich.com/biofiles

New Products

Sigma is committed to continually adding new products for use in cancer research and cell biology, as well as new technology plat-forms for cellular studies. We welcome your recommendations for new products. Please tell us your suggestions by visiting **Sigma.com/product_suggestions**.

Antitumor Agents

Cephalomannine

C₄₅H₅₃NO₁₄ FW 831.90 [71610-00-9]

Antitumor; antiproliferative. Promotes the formation of highly stable microtubules that resist depolymerization but does not induce cytokine secretion.

≥97% (HPLC) C4991

1 mg, 5 mg

100 µg

10 mg

Chrysomycin A

C28H28O9 FW 508.52 [82196-88-1]

Antibiotic from *Streptomyces* sp. Inhibits the catalytic activity of human topoisomerase II. Exhibits antitumor activity against human cell lines K562, HT29, MCF7, PC6, and MKN28.

≥98% (HPLC)

2-8°C

C6616			

Chrysomycin B

 $C_{27}H_{28}O_9 \; FW \; 496.51 \;\; [83852\text{-}56\text{-}6]$

Antibiotic from *Streptomyces* sp. Inhibits the catalytic activity of human topoisomerase II. Exhibits antitumor activity against human cell lines K562, HT29, MCF7, PC6, and MKN28.

≥98% (HPLC)

2-8°C	
C6491	100 µg

CIL-102

1-[4-(Furo[2,3-b]quinolin-4-ylamino)phenyl]ethanone

C₁₉H₁₄N₂O₂ FW 302.33 [479077-76-4]

Tubulin polymerization inhibitor: apoptosis inducer.

US Patent No. 6,750,223. Sold for non-human research purposes under a non-exclusive license agreement with Kaohsiung Medical University.

≥95% (HPLC)

-20°C WET ICE

C5

17-Dimethylaminoethylamino-17-demethoxygeldanamycin

17-DMAG

 $C_{32}H_{48}N_4O_8\ FW\ 616.75$

17-DMAG is a more potent water soluble analog of geldanamycin. 17-DMAG has shown more antitumor activity than 17-AAG. Inhibits cancer growth and promotes apoptosis in multiple cell lines.

-20°C	
D5193	1 mg

Finasteride

 $\begin{array}{l} \text{N-}\textit{tert}\text{-}\text{Butyl-3-oxo-4-aza-}5\alpha\text{-}\textit{androst-1-en-}17\beta\text{-}\textit{carboxamide};\\ \text{N-}(2\text{-}\text{Methyl-2-propyl)-3-oxo-4-aza-}5\alpha\text{-}\textit{androst-1-ene-}17\beta\text{-}\textit{carboxamide};\\ \text{MK-906 } C_{23}\text{H}_{36}\text{N}_2\text{O}_2 \ \ \text{FW 372.54 [98319-26-7]}\\ \text{Selective } 5\alpha\text{-}\textit{reductase inhibitor}; \ \text{antiandrogen.} \end{array}$

≥97% (HPLC)

-20°C

F1293

ICRF-193

meso-4,4'-(3,2-butanediyl)-bis(2,6-piperazinedione) $C_{12}H_{18}N_4O_4\;\; FW\;\; 282.30\;\; [21416-68-2]$

ICRF-193 induces a G₂ checkpoint that is associated with an ATR dependent inhibition of polo-like kinase 1 (Plk1) activity and a decrease in cyclin B1 phosphorylation. Induces apoptosis in several cell lines including K562 and Molt-4 cells. ICRF-193 is a topoisomerase II inhibitor that targets topoisomerase II- β to a greater extent than it targets topoisomerase II- α and does not cause DNA damage.

≥**9**5%

-20°C

14659

Kazusamycin A

C33H48O7 FW 556.73 [92090-94-3]

Kazusamycin A, an unsaturated branched chain fatty acid with a terminal lactone ring, is a hydroxy analog of Leptomycin B. It has significant *in vitro* cytotoxic activity against various human and mouse tumor lines encompassing a wide range of tissue types. Kazusamycin A exhibits *in vivo* antitumor activity against experimental murine tumors. It inhibits nuclear export and Rev translocation, a regulatory gene product in the HIV genome, at nanomolar concentrations.

from Streptomyces sp., ≥95% (HPLC)

70% methanol solution.

-20°C DRY ICE

K1764

Linomide

LS-2616; N-Phenylmethyl-1,2-dihydro-4-hydroxyl-1-methyl-2-oxoquinoline-3-carboxamide; Quinoline-3-carboxamide

 $C_{18}H_{16}N_2O_3 \ FW \ 308.33 \ [84088-42-6]$

Immunomodulator; antiangiogenic. Upregulates IL-10 and decreases CXC chemokine production. Upregulates mRNA of IL-23, IL-12p35, and IFN γ . Also, antagonizes hepatic effects of TNF α by blocking recruitment of inflammatory cell infiltration. ≥98% (HPLC)

27070 (i

2-8°C

10 mg

1 µg

100 mg

1 mg, 5 mg

Myoseverin B

2,6-Bis(4-methoxybenzylamino)-9-cyclohexylpurine C $_{\rm 27}H_{\rm 32}N_6O_2$ FW 472.58 [361431-27-8] Tubulin polymerization inhibitor.

≥97% (HPLC)

-20°C DRY ICE	
M3316	

1 mg, 5 mg

New Products

ma-aldrich.com/biofile

sign

Neocarzinostatin from *Streptomyces carzinostaticus*

Holoneocarzinostatin; NCS; NSC-69856; Zinostatin [9014-02-2]

Protein-small molecule complex composed of an enediyne chromo-phore tightly bound to a 113-amino acid single chain protein. Chemotherapeutic agent which inhibits DNA synthesis, possesses antitumor activity, and induces apoptosis.

Supplied as 100 μg (~9 nmoles) in a solution containing 20 mM MES buffer, pH 5.5.

≥90% (SDS-PAGE), ~0.5 mg/mL

2-8°C WET ICE		
N9162		

Nutlin-3

(±)-4-[4,5-Bis(4-chlorophenyl)-2-(2-isopropoxy-4-methoxyphenyl)-4,5-dihydroimidazole-1-carbonyl]-piperazin-2-one

 $C_{30}H_{30}\check{CI}_2N_4O_4 \ FW \ 581.49 \ [548472-68-0]$

Mdm2 (mouse double minute 2) antagonist, p53 pathway activator, apoptosis inducer.

Sold under license from Hoffman-La Roche, Inc. US patent 6,734,302

≥98% (HPLC)

-20°C WET ICE	
N6287	1 mg, 5 mg

p53 human

 $\mathsf{p53}$ gene is highly conserved and expressed in normal tissues. It is the most commonly mutated gene in human cancer and more than

500 gene mutations have been described in various types of malignancies, hematologic as well as solid tumors. Intact p53 function is essential for the maintenance of the non-tumorogenic phenotype of cells. Thus, p53 plays a vital role in suppressing the development of cancer. Useful for the study of postranslational modification of p53, gel-shift assays, or protein-protein interaction studies.

Supplied at a concentration of 0.2 mg/ml in PBS, pH 7.0.

recombinant, expressed in Escherichia coli

>90% (SDS-PAGE)

-20°C WET ICE	
P6249	50 µg
>95% (SDS-PAGE)	
-20°C WET ICE	
P6374	20 µg

PTACH

Cpd 51; S-[6-(4-Phenyl-2-thiazolylcarbamoyl)hexyl] thioisobutyrate $C_{20}H_{26}N_2O_2S_2\;\;FW\;\;390.56\;\;[848354-66-5]$

HDAC inhibitor; more potent than the majority of HDAC inhibitors except for SAHA.

≥98% (HPLC)

2-8°C

P5874

RG108

N-Phthalyl-1-tryptophan C₁₉H₁₄N₂O₄ FW 334.33 [48208-26-0]

DNA methyltransferase (DMNT) inhibitor. Reactivates tumor suppressor gene expression (p16, SFRP1, secreted frizzled related protein-1, and TIMP-3) in tumor cells by DNA demethylation. Inhibits human tumor cell line (HCT116, NALM-6) proliferation and increases doubling time in culture. ≥98% (HPLC)

≥ 70 70

R8279	10 mg

SKI II

100 µg

4-[[4-(4-Chlorophenyl)-2-thiazolyl]amino]phenol C₁₅H₁₁ClN₂OS FW 302.78 [312636-16-1]

Sphingosine kinase (SK) plays a pivotal role in regulating tumor growth and SK can act as an oncogene. Expression of SK RNA is significantly elevated in a variety of solid tumors, compared with normal tissue from the same patient. A number of novel inhibitors of human SK were identified, and several representative compounds were characterized in detail. These compounds demonstrated activity at sub- to micromolar concentrations, making them more potent than any other reported SK inhibitor, and were selective toward SK compared with a panel of human lipid and protein kinases. Kinetic studies revealed that the compounds were not competitive inhibitors of the ATP-binding site of SK. SKI-II inhibitor is orally bioavailable, detected in the blood for at least 8 h, and showed a significant inhibition of tumor growth in mice with IC₅₀ = 0.5 µM; SKI II does not act at ATP-binding site. Displays no inhibition of ERK2, PI 3-kinase, or PKCa at concentrations up to 60 µM. SKI II induces apoptosis and inhibits proliferation in several other tumor cell lines in vitro (IC₅₀ = $0.9-4.6 \mu$ M).

≥98% (HPLC)

S5696

10 mg, 50 mg

Sulochrin

$C_{17}H_{16}O_7 \ FW \ 332.30 \ [519-57-3]$

Antibiotic from *Aspergillus* and *Penicillium* sp. Fungal metabolite; VEGF inhibitor and antiangiogenic that inhibits the VEGF-induced tube formation of human umbilical vein endothelial cells.

≥98% (HPLC)

2-8°C \$4570

2 mg

1 mg

sigma-aldrich.com/biofiles

New Products

SIGMA

New Products

Chemopreventitive Agents

2-Amino-N-quinolin-8-yl-benzenesulfonamide

OBS C15H13N3O2S FW 299.35 [16082-64-7] Inhibitor of cell cycle at G₂ phase; apoptosis inducer.

≥98% (HPLC)

A3105

Borrelidin

Borrelidine; 2-(7-Cyano-8,16-dihydroxy-9,11,13,15-tetramethyl-18-oxoacyclooctadeca-4,6-dien-2-yl)-cyclopentanecarboxylic acid; Cyclopentanecarboxylic acid; NSC 216128; Treponemycin C₂₈H₄₃NO₆ FW 489.64 [7184-60-3]

Potent angiogenesis inhibitor that induces apoptosis of the capillary tube-forming cells. Also, displays antimalarial activity against drug-resistant Plasmodia. Antimicrobial and selective threonyl t-RNA synthetase inhibitor.

from Streptomyces parvulus, ≥98% (HPLC)

-20°C WET ICE

B3061

Carnosol

C₂₀H₂₆O₄ FW 330.42 [5957-80-2]

A phenolic diterpene with antioxidant and anticarcinogenic activities

from Rosemarinus officinalis L

-20°C DRY ICE	
C9617	5 mg

Etidronate disodium

Dihydrogen (1-hydroxyethylidene)bisphosphonate disodium C₂H₆Na₂O₇P₂ FW 249.99 [7414-83-7]

Bisphosphonate antiresorptive agent. Less potent inhibitor of farnesyl diphosphate synthase (IC₅₀ = 80 μ M) as compared to the nitrogen containing bisphosphonates

2-8°C WET ICE	
P5248	10 mg, 100 mg

Ifosfamide

N,3-Bis(2-chloroethyl)tetrahydro-2H-1,3,2-oxazaphosphorin-2-amine-2-oxide; Ifex

C7H15CI2N2O2P FW 261.09 [3778-73-2]

Ifosfamide is a nitrogen mustard compound that is a structural isomer of cyclophosphamide. Ifosfamide is a prodrug that must be transformed by cytochrome P450 to the biologically active component. It is used as an antineoplastic agent in cancer chemotherapy, but ifosfamide is more likely to cause renal toxicity than cyclophosphamide

>98%

2-8°C	WET	ICE	
14909			

Lomustine

1-(2-Chloroethyl)-3-cyclohexyl-1-nitrosourea

C₉H₁₆CIN₃O₂ FW 233.70 [13010-47-4]

Antineoplastic agent with cellular DNA effects. Lomustine induces p53 expression in A2870 cells.

≥**98%**

-20°C

L5918

Enzymes

5 mg

1 mg

Breast Tumor Kinase Active human

Brk; Protein Tyrosine Kinase 6; PTK6

recombinant, expressed in Escherichia coli, ≥85% (SDS-PAGE)

Breast Tumor Kinase (Brk, PTK6) is a member of the non-receptor tyrosine kinases (PTKs). Brk is expressed in breast tumors. Brk promotes cell motility and invasion and functions as a mediator of EGF-induced migration and invasion.

Protein apparent mol wt ~80 kDa

-70°C DRY ICE B9810 5 µg

Multi-Drug Resistance

Fumitremorgin C

FTC C22H25N3O3 FW 379.45 [118974-02-0]

from Neosartorya fischeri, >98% (HPLC and TLC)

Fumitremorgin C (FTC) is a fungal toxin of the diketopiperazines family of compounds. In mammalian cells, FTC is tremorgenic and causes cell cycle arrest. Fumitremorgin C has been shown to reverse resistance to doxorubicin, mitoxantrane, and topotecan in non-Pgp (P-glycoprotein), non-MRP (multidrug resistance protein) multidrug-resistance (MDR) cells. This reversal of resistance is associated with an increase in drug accumulation. Fumitremorgin C is a specific, selective, and potent inhibitor at micromolar concentrations of the breast cancer resistant protein (BCRP/ABCG2), an ABC transporter associated with chemotherapy resistance. FTC, in combination with mitoxantrone, can be used for the detection of ABCG2 functional activity in several cell lines.

2-8°C WFT ICF F9054

250 µg

mg

Tumor Growth Regulation

(E/Z)-4-HTA monohydrate

(E,Z)-2-[4-[1-(p-Hydroxyphenyl)-2-phenyl]-1-butenyl]phenoxyacetic acid; 4-Hydroxytamoxifen acid C24H22O4 · H2O FW 392.44 Non-steroidal estrogen.

≥90% (HPLC)

_/0/	• (,
2-8°C		

2-80	
H5039	2 mg, 10

PD 173074

N-[2-[4-(Diethylamino)butyl]amino-6- (3,5-dimethoxyphenyl)pyrido[2,3-d]pyrimidin-7-yl]-N'-(1,1-dimethylethyl)urea C₂₈H₄₁N₇O₃ FW 523.67 [219580-11-7]

Fibroblast growth factor receptor 3 (FGFR3) inhibitor: $IC_{50} = 5 \text{ nM}$ in inhibition of FGFR3 autophosphorylation; arrested the G0/G1 phase of the cell cycle of FGFR3-expressing cells. 100-fold more selective for FGFR3 than for VEGF receptors, IGF-1 receptors, and MAPKs.

≥98% (HPLC)

2-8°C
P2499

1 g, 5 g

100 ma

5 mg, 25 mg

Reveromycin A

C36H52O11 FW 660.79 [134615-37-5]

Polyketide antibiotic from Streptomyces sp. Epidermal growth factor (EGF) inhibitor; apoptosis inducer; G1 phase cell cycle inhibitor having antiproliferative behavior against human cell lines KB and K562 as well as antifungal activity.

≥98% (HPLC)

-20°C WET ICE R0654

32

com/biofile

ma-aldrich.

ص ص

s i

PathFinder



INNOVATION @ WORK

Discover Your Path to Innovation

On your path to innovation, Sigma is with you every step of the way. PathFinder is an online collection of interactive, interconnected maps showing biological signaling and metabolic pathways. For you to explore the relationships between different pathway elements, individual components are linked with related high-quality products.

You know your destination. PathFinder will get you there.

With Sigma's broad range of products, you will discover that we offer everything from small molecules to antibodies and enzymes. You will also be able to use PathFinder to locate qPCR components and siRNAs for gene knockdown. A valuable resource, PathFinder provides fast and accurate information on numerous levels, all in one place – and all linked to the important products that are key to the success of your research. In addition to products and services, you'll have immediate access to these helpful tools:

- Specific workflow analysis
- Detailed product descriptions
- In-depth technical information
- Relevant technical articles

Learn how PathFinder can help you discover your path to innovation by visiting us at: sigma.com/pathfinder

Accelerating Customers' success through Leadership in Life Science, High Technology and Service SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI 63178 • USA



Argentina

SIGMA-ALDRICH DE ARGENTINA S.A. Free Tel: 0810 888 7446 Tel: (+54) 11 4556 1472 Fax: (+54) 11 4552 1698

Australia

SIGMA-ALDRICH PTY ITD Free Tel: 1800 800 097 Free Fax: 1800 800 096 Tel: (+61) 2 9841 0555 Fax: (+61) 2 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH Tel: (+43) 1 605 81 10 Fax: (+43) 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA. Free Tel: 0800 14747 Free Fax: 0800 14745 Tel: (+32) 3 899 13 01 Fax: (+32) 3 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA Free Tel: 0800 701 7425 Tel: (+55) 11 3732 3100 Fax: (+55) 11 5522 9895

Canada

SIGMA-ALDRICH CANADA LTD. Free Tel: 1800 565 1400 Free Fax: 1800 265 3858 Tel: (+1) 905 829 9500 Fax: (+1) 905 829 9292

China

SIGMA-ALDRICH (SHANGHAI) TRADING CO. LTD. Free Tel: 800 819 3336 Tel: (+86) 21 6141 5566 Fax: (+86) 21 6141 5567

Czech Republic

SIGMA-ALDRICH S.R.O. Tel: (+420) 246 003 200 Fax: (+420) 246 003 291

World Headquarters

3050 Spruce St., St. Louis, MO 63103 (314) 771-5765 sigma-aldrich.com

Denmark

SIGMA-ALDRICH DENMARK A/S Tel: (+45) 43 56 59 10 Fax: (+45) 43 56 59 05

Finland SIGMA-ALDRICH FINLAND OY Tel: (+358) 9 350 9250

Fax: (+358) 9 350 92555 France SIGMA-ALDRICH CHIMIE S.à.r.I. Free Tel: 0800 211 408 Free Fax: 0800 031 052

Tel: (+33) 474 82 28 00 Fax: (+33) 474 95 68 08

Germany SIGMA-ALDRICH CHEMIE GmbH Free Tel: 0800 51 55 000 Free Fax: 0800 64 90 000 Tel: (+49) 89 6513 0 Fax: (+49) 89 6513 1160

Greece SIGMA-ALDRICH (O.M.) LTD. Tel: (+30) 210 994 8010 Fax: (+30) 210 994 3831

Hungary SIGMA-ALDRICH Kft Ingyenes zöld telefon: 06 80 355 355 Ingyenes zöld fax: 06 80 344 344 Tel: (+36) 1 235 9055 Fax: (+36) 1 235 9050

India

SIGMA-ALDRICH CHEMICALS PRIVATE LIMITED Telephone Bangalore: (+91) 80 6621 9600 New Delhi: (+91) 11 4165 4255 Mumbai: (+91) 22 2570 2364 Hyderabad: (+91) 40 6684 5488 Fax

Bangalore: (+91) 80 6621 9650 New Delhi: (+91) 11 4165 4266 Mumbai: (+91) 22 2579 7589 Hyderabad: (+91) 40 6684 5466

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052 Technical Service (800) 325-5832 • sigma-aldrich.com/techservice Development/Bulk Manufacturing Inquiries SAFC[®] (800) 244-1173

The SIGMA-ALDRICH Group

SIGMA ALDRICH

© 2007 Sigma-Aldrich Co. All rights reserved. SIGMA, ♣, SAFC, SAFC; SIGMA-ALDRICH, ♣, ISOTEC, ALDRICH, ֎, FLUKA, ♠, and SUPELCO are trademarks belonging to Sigma-Aldrich Co. and its affiliate Sigma-Aldrich Biotechnology, LP. Redel-de Haen^{*} trademark under license from Riedel-de Haen GmbH. MISSION, Panorama, Isopac: are registered trademarks, Your Favorite Gene, Phosphorofile, Sigma-Aldrich Biotechnology LF. Redel-de Haen^{*} trademark of Sigma-Aldrich Biotechnology LF. GenomePiex is a registered trademark of Rubicon Genomics. QuadriPERM is a registered trademark of Heraeus Instruments GmbH. Benzonaes is a registered trademark of Merck KGaA, Darmstadt, Germany. TWEEN is a registered trademark of Inderma. SYBE is a registered trademark of Oflocular Probes, Inc. TaqMan is a registered trademark of Amersham Biosciences Limited. Cytoscape is a trademark of Cytoscape Consortium Corporation. Cy[™] is a registered trademark of Amersham Biosciences. Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing sign.



3050 Spruce Street • St. Louis, MO 63103 USA

Ireland

SIGMA-ALDRICH IRFLAND LTD. Free Tel: 1800 200 888 Free Fax: 1800 600 222 Tel: (+353) 1 404 1900 Fax: (+353) 1 404 1910

Israel

SIGMA-ALDRICH ISRAEL LTD. Free Tel: 1 800 70 2222 Tel: (+972) 8 948 4100 Fax: (+972) 8 948 4200

Italy SIGMA-ALDRICH S.r.I. Numero Verde: 800 827018 Tel: (+39) 02 3341 7310 Fax: (+39) 02 3801 0737

Japan SIGMA-ALDRICH JAPAN K.K. Tokyo Tel: (+81) 3 5796 7300 Tokyo Fax: (+81) 3 5796 7315

Korea SIGMA-ALDRICH KOREA Free Tel: (+82) 80 023 7111 Free Fax: (+82) 80 023 8111 Tel: (+82) 31 329 9000 Fax: (+82) 31 329 9090

Malaysia SIGMA-ALDRICH (M) SDN. BHD Tel: (+60) 3 5635 3321 Fax: (+60) 3 5635 4116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V. Free Tel: 01 800 007 5300 Free Fax: 01 800 712 9920 Tel: 52 722 276 1600 Fax: 52 722 276 1601

The Netherlands

SIGMA-ALDRICH CHEMIE BV Free Tel: 0800 022 9088 Free Fax: 0800 022 9089 Tel: (+31) 78 620 5411 Fax: (+31) 78 620 5421

New Zealand

Free Tel: 0800 936 666

Free Fax: 0800 937 777 Tel: (+61) 2 9841 0555 Fax: (+61) 2 9841 0500 Norway

SIGMA-ALDRICH NEW ZEALAND LTD.

SIGMA-ALDRICH NORWAY AS Tel: (+47) 23 17 60 60 Fax: (+47) 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o. Tel: (+48) 61 829 01 00 Fax: (+48) 61 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A. Free Tel: 800 202 180 Free Fax: 800 202 178 Tel: (+351) 21 924 2555 Fax: (+351) 21 924 2610

Russia

SIGMA-ALDRICH RUS, LLC Tel: +7 (495) 621 6037 Fax: +7 (495) 621 5923

Singapore

SIGMA-ALDRICH PTE. LTD. Tel: (+65) 6779 1200 Fax: (+65) 6779 1822

South Africa

SIGMA-ALDRICH SOUTH AFRICA (PTY) LTD. Free Tel: 0800 1100 75 Free Fax: 0800 1100 79 Tel: (+27) 11 979 1188 Fax: (+27) 11 979 1119

ISOTEC

Spain

SIGMA-ALDRICH QUÍMICA, S.A. Free Tel: 900 101 376 Free Fax: 900 102 028 Tel: (+34) 91 661 99 77 Fax: (+34) 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB Tel: (+46) 8 742 4200 Fax: (+46) 8 742 4243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH Free Tel: 0800 80 00 80 Free Fax: 0800 80 00 81 Tel: (+41) 81 755 2828 Fax: (+41) 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD. Free Tel: 0800 717 181 Free Fax: 0800 378 785 Tel: (+44) 1747 833 000 Fax: (+44) 1747 833 313 SAFC (UK) Free Tel: 0800 71 71 17

United States

SIGMA-ALDRICH P.O. Box 14508 St. Louis, Missouri 63178 Toll-Free: 800 325 3010 Toll-Free Fax: 800 325 5052 Call Collect: (+1) 314 771 5750 Tel: (+1) 314 771 5765 Fax: (+1) 314 771 5757

Internet sigma-aldrich.com

SIGMA-ALDRICH

Accelerating Customers' Success through Leadership in Life Science, High Technology and Service

02047-504243 0057

SAFC **#Fluka** SUPELCO Riedel-de Haën*