

**Impurity Assay of α -D-GLUCOSE 1-PHOSPHATE
present in α -D-GALACTOSE 1-PHOSPHATE**

PRINCIPLE:

α -D-Glucose 1-Phosphate $\xrightarrow{\text{Phosphoglucomutase}}$ α -D-Glucose 6-Phosphate

α -D-Glucose 6-Phosphate + β -NADP $\xrightarrow{\text{G-6-PDH}}$ 6-Phospho-D-gluconate + β -NADPH

Abbreviations used:

β -NADP = β -Nicotinamide Adenine Dinucleotide Phosphate, Oxidized Form

β -NADPH = β -Nicotinamide Adenine Dinucleotide Phosphate, Reduced Form

G-6-PDH = Glucose-6-Phosphate Dehydrogenase

CONDITIONS: T = 25°C, pH = 7.4, $A_{340\text{nm}}$, Light path = 1 cm

METHOD: Continuous Spectrophotometric Rate Determination

REAGENTS:

- A. 250 mM Glycylglycine Buffer, pH 7.4 at 30°C
(Prepare 50 ml in deionized water using Gly-Gly, Free Base, Sigma Prod. No. G-1002. Adjust to pH 7.4 at 25°C with 1 M HCl.)
- B. 10 mg/ml Galactose 1-Phosphate Solution, Potassium (Gal-1-P)¹
(Prepare 2 ml in deionized water using α -D-Galactose 1-Phosphate, Dipotassium Salt, Sigma Prod. No. G-0380.)
- C. 20 mM β -Nicotinamide Adenine Dinucleotide Phosphate Solution (β -NADP)
(Prepare 1 ml in deionized water using β -Nicotinamide Adenine Dinucleotide Phosphate, Sodium Salt, Sigma Prod. No. N-0505 or dissolve the contents of one 30 mg vial of β -Nicotinamide Adenine Dinucleotide Phosphate, Sodium Salt, Sigma Stock No. 240-330, in the appropriate volume of deionized water. **PREPARE FRESH.**)

**Impurity Assay of α -D-GLUCOSE 1-PHOSPHATE
present in α -D-GALACTOSE 1-PHOSPHATE**

REAGENTS: (continued)

- D. 1.2 mM Glucose 1,6-Diphosphate Solution¹ (G 1,6-P)
(Prepare 1 ml in deionized water using α -D-Glucose 1,6-Diphosphate, Cyclohexylammonium Salt, Hydrate Sigma Prod. No. G-7137.)
- E. 300 mM Magnesium Chloride Solution (MgCl_2)
(Prepare 1 ml in deionized water using Magnesium Chloride, Hexahydrate, Sigma Prod. No. M-0250.)
- F. Glucose-6-Phosphate Dehydrogenase Solution (G-6-PDH)
(Immediately before use, prepare a solution containing 25 units/ml in cold Reagent A using Glucose-6-Phosphate Dehydrogenase, Sigma Prod. No. G-7877.)
- G. Phosphoglucomutase Enzyme Solution (PGLUM)
(Immediately before use, prepare a solution containing 25 units/ml of Phosphoglucomutase, Sigma Prod. No. P-3397, in cold Reagent A.)

PROCEDURE:

Pipette (in milliliters) the following reagents into suitable cuvettes:

	<u>Test</u>	<u>Blank</u>
Reagent A (Buffer)	1.40	1.40
Deionized Water	-----	1.00
Reagent B (Gla 1-P)	1.00	-----
Reagent D (G 1,6-P)	0.05	0.05
Reagent C (β -NADP)	0.10	0.10
Reagent E (MgCl_2)	0.25	0.25
Reagent F (G-6-PDH)	0.10	0.10

Mix by inversion and equilibrate to 25°C, using a suitably thermostatted spectrophotometer. Record the initial $A_{340\text{nm}}$ for both the Test and Blank versus air. Then add:

Reagent G (PGLUM)	0.10	0.10
-------------------	------	------

Immediately mix by inversion and allow the reaction to proceed until $\Delta\text{Abs}_{340}/\text{minute}$ is ≤ 0.0020 , (and this rate is maintained for a minimum of five minutes). Record the final $A_{340\text{nm}}$ for both the Test and Blank versus air.

**Impurity Assay of α -D-GLUCOSE 1-PHOSPHATE
present in α -D-GALACTOSE 1-PHOSPHATE**

CALCULATIONS:

$$\Delta A = A_f - A_i$$

A_f = Final Absorbance

A_i = Initial Absorbance

$$\text{Micromoles of Glu-1-P/1 ml of Gal-1-P} = \frac{(\Delta A_{\text{Test}} - \Delta A_{\text{Blank}})(3)}{6.22}$$

Glu-1-P = α -D-Glucose-1-Phosphate

Gal-1-P = α -D-Galactose-1-Phosphate

3 = Total volume (in milliliters) of assay

6.22 = Millimolar extinction coefficient of β -NADPH at 340 nm

$$\text{Molar \% Glu-1-P} = \frac{(\mu\text{moles of Glu-1-P/1 ml of Gal-1-P})(100)}{\mu\text{moles of Gal-1-P/1 ml of Gal-1-P}}$$

$$\frac{\mu\text{moles Glu-1-P}}{1.0 \text{ ml of Gal-1-P}} = \frac{\text{mgs X } 1000}{\text{Enzymatic MW}}$$

FINAL ASSAY CONCENTRATION:

In a 3.00 ml reaction mix, the final concentrations are 117 mM glycylglycine, 0.67 mM β -nicotinamide adenine dinucleotide phosphate, 0.02 mM glucose 1,6-diphosphate, 25 mM magnesium chloride, 2.5 units glucose 6-phosphate dehydrogenase and 2.5 units phosphoglucomutase and varying concentrations of α -D-galactose 1-phosphate.

REFERENCE:

Bergmeyer, H.A., and Michal, G. (1974) in *Methods of Enzymatic Analysis* (Bergmeyer, H.U. ed.) Volume 3, 2nd ed., 1233-1237, Academic Press, Inc., New York, NY

NOTES:

1. If the test material is α -Galactose-1-Phosphate, Barium Salt, Sigma Product No. G-1300, prepare in 30 mM EDTA, using Ethylenediaminetetraacetic Acid, Tetrasodium, Hydrate, Sigma Stock No. ED4SS.
2. Glucose 1,6-Diphosphate is required in order to obtain maximum activity.
3. Glucose-6-Phosphate Dehydrogenase unit definition: One unit will oxidize 1.0 μ mole of D-glucose 6-phosphate to 6-phospho-D-gluconate per minute in the presence of β -NADP at pH 7.4 at 25°C.
4. Phosphoglucomutase Unit Definition: One unit will convert 1.0 μ mole of α -D-glucose

1-phosphate to α -D-glucose 6-phosphate per minute at pH 7.4 at 30°C.

NOTES: (continued)

5. This assay is based on the cited references.
6. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

This procedure is for informational purposes. For a current copy of Sigma's quality control procedure contact our Technical Service Department.