


Biofiles

SIGMA-ALDRICH®

SIGMA Life Science Where *bio* begins™



Enzymes and Reagents for Alternative Energy

Metabolic strategies in alternative energy

Enzyme-based fuel cells

Lignocellulosic depolymerization

Enzymes for biofuel research

Reagents for ethanol and biodiesel analysis

Biofilesonline

Your gateway to Biochemicals and Reagents for Life Science Research



BioFiles Online allows you to:

- Easily navigate the content of the current BioFiles issue
- Access any issue of BioFiles
- Subscribe for email notifications of future eBioFiles issues

Register today for upcoming issues and eBioFiles announcements at sigma.com/biofiles



Greener Products

Sigma-Aldrich is committed to doing our part to minimize our footprint on the environment. As a life science and high technology company, we recognize your need for high quality chemicals comes first. In cases where it becomes possible to consider alternatives, we have made it easier for you to find greener options through our new "Greener Alternatives and Technologies" product lists.

Products Including:

- Quality Environmentally-Friendly Enzymes
 - Ionic Liquids
 - Greener Alternatives: Solvents
 - Greener Alternatives: Reagents

sigma-aldrich.com/green



Life Science Labware Center

The Life Science Labware Center allows easy browsing of labware products, including disposables, books, and equipment

Features Include:

- Corning® cell culture selection guide
 - Products by top manufacturers
 - New products by research areas
 - New lab start-up programs

Start using the Life Science Labware center at sigma-aldrich.com/lslabware

Biofilescontents

Introduction	3
Advancements in Enzyme-Based Fuel Cell Batteries, Sensors and Emissions Reduction Strategies	8
Selected Enzymes for Biofuel Cell Research	10
Enzymes and Reagents for Ethanol Research	12
Enzymes for Lignocellulosic Ethanol Research	12
Enzymes for Starch Hydrolysis	20
Ethanol Analysis	23
Enzymes for BioDiesel Research	24
Lipase	24
Phospholipase C	24
BioDiesel Analysis	25
Metabolic Standards for Mevalonate and Isoprenoid Pathway Analysis	30

Technical content: Robert Gates, M. S. M.

Introduction

Robert Gates

Market Segment Manager

robert.gates@sial.com



Leveraging nature's metabolic strategies for tomorrow's energy resources.

Although most of the world's population understands that carbon dioxide is the predominant end product of energy production from combustible fuels, most of us forget that the "raw material" for most of our combustible fuel sources is also atmospheric carbon dioxide. Photosynthetic carbon dioxide fixation is most often the first step in the synthesis of both our renewable fuels and non-renewable fossil fuels. Even much of our hydrogen production is an indirect product of carbon dioxide fixation as it is often derived from fossil fuels.

The basic differences between non-renewable fossil fuels and our "new generation" renewable fuels is time, energy input and mass. Our non-renewable fuels, such as petroleum and coal required thousands of years of collection of solar energy to power carbon fixation for the biosynthesis of enormous amounts of biomass. That biomass then underwent millions of years of anaerobic decomposition to yield today's fossil fuels. In our attempt to circumvent or speed up this process, we are finding natural "renewable analogs" to some of our fossil fuels, i.e., ethanol and butanol to replace or supplement gasoline and fatty acid esters to replace diesel.

However, we have not advanced far enough in our manufacturing technology to divest ourselves of the need for biological systems to do the complicated work for us. Current processes for ethanol production rely on plants to fix carbon dioxide, synthesize a "common currency" molecule such as

glucose, and then needlessly polymerize this glucose into starch or cellulose. Then we often rely on microbial and enzymatic hydrolysis of these biopolymers to give back glucose and then ferment it to yield ethanol. Sugarcane ethanol employs a slightly simpler path via sucrose.

Our challenge for the future may be to find the "straightest metabolic line" between carbon dioxide and ethanol, butanol, fatty acids or even hydrogen. The ideal route to ethanol production might be to devise an enzymatic or cellular system to fix carbon dioxide and take the direct route from ribose-1,5-bis-phosphate to pyruvate to ethanol and bypass glucose biosynthesis, polymerization and depolymerization. We have included the Nicholson/IUBMB Metabolic Pathway chart on the following pages illustrating the specific pathways related to carbon dioxide fixation and glucose, cellulose, pectin, starch, lignin and ethanol production highlighted in yellow.

Enzymes and their metabolic pathways enter the energy equation from other directions. Microbial strains are being discovered that can ferment cellulosic biomass to yield hydrogen. Some photosynthetic microbes, such as algae, are capable of splitting water to produce hydrogen utilizing solar energy. Hydrogenases are being targeted not only for hydrogen production but also for the oxidation of H₂ in fuel cells. Immobilized enzyme-based fuel cell batteries can utilize multi-enzyme systems that mimic *in vivo* glycolytic and Krebs Cycle pathways. Immobilized enzymes are also being studied for the sequestration of carbon dioxide emissions.


Engineering organisms and enzyme systems to maximize biofuel production could require insertion of exogenous enzymatic pathways from various organisms into a new host. In addition, alternative metabolic pathways such as the mevalonate and isoprenoid pathways are being studied as sources of biofuel.

This issue of BioFiles is devoted to how Sigma-Aldrich enzymes and reagents are involved in many of these areas of research. As summarized above and illustrated

below, many of these areas require an understanding of metabolic pathways and enzymatic processes.

References

- (1) S. Atsumi, J. C. Liao, Metabolic Engineering for Advanced Biofuels Production from *Escherichia coli*, *Curr. Opin. Biotechnol.*, **19**, 414–419 (2008).
- (2) H. Alper, G. Stephanopoulos, Engineering for biofuels: exploiting innate microbial capacity or importing biosynthetic potential?, *Nature Reviews Microbiology*, **7**, 715–723 (2009).
- (3) P. Peralta-Yahya, J. D. Keasling, Advanced biofuel production in microbes, *Biotechnol. J.*, **5**, 147–162 (2010).
- (4) J. M. Clomburg, R. Gonzalez, Biofuel production in *Escherichia coli*: the role of metabolic engineering and synthetic biology, *Appl. Microbiol. Biotechnol.*, **86**, 419–434 (2010).
- (5) D. Klein-Marcuschamer, V. G. Yadav, A. Ghaderi, G. Stephanopoulos, De Novo Metabolic Engineering and the Promise of Synthetic DNA, *Adv. Biochem. Eng. Biotechnol.* in press (2010).
- (6) J. L. Fortman, S. Chhabra, A. Mukhopadhyay, H. Chou, T. S. Lee, E. Steen, J. D. Keasling, Biofuel alternatives to ethanol: pumping the microbial well, *Trends in Biotechnol.*, **26**, 375–381 (2010).
- (7) D. B. Levin, R. Islam, N. Cicic, R. Sparling, Hydrogen production by *Clostridium thermocellum* 27405 from cellulosic biomass substrates, *Int. J. Hydrogen Energy*, **31**, 1496–1503 (2006).
- (8) M. L. Ghirardi, S. Kosourov, A. Tsygankov, A. Rubin, M. Seibert, Cyclic Photobiological Algal H₂ Production, *Proceedings of the 2002 U. S. DOE Hydrogen Program Review*
- (9) A. A. Karyakin, S. V. Morozov, E. E. Karyakina, N. A. Zorin V. V. Pereygin, S. Cosnie, Hydrogenase Electrodes for Fuel Cells, *International Hydrogenases Conference 2004*.



The New Enzyme Explorer

Redesigned Format and Expanded Resources.

The **Enzyme Explorer Indices** provide paths to find more than 3,000:

- Enzyme/proteins
- Inhibitors
- Substrates
- Cofactors

The **Enzyme Explorer Product Features**

include enzyme and protein related products for Cell Biology, Diagnostics and industrial scale manufacturing.

The **Enzyme Explorer Learning Center**

includes:

- IUBMB Metabolic Pathways
- Proteolytic Enzyme Resources
- Enzyme Assay Library
- Carbohydrate Analysis Enzyme Resource
- Product Guides and Protocols
- The New Protease Finder

Allow the **New Protease Finder** to help you locate both the endo and exoproteases needed for precision protein cleavage.





- Request Enzyme Literature
- Proteolytic Enzymes
- Protease Inhibition
- Diagnostic and Metabolic Research
- Detection Substrates
- Metabolomics and Dietary Research
- Cell Dissociation and Lysis
- Plasma Proteins
- Carbohydrate Analysis
- Proteomics and Protein Expression



sigma-aldrich.com/enzymeexplorer

Metabolic Origins of Fuel





© 2003 International Union of Biochemistry and
Molecular Biology www.iubmb.org

22nd Edition Designed by Donald E. Nicholson, D.Sc.,
The University of Leeds, England – and Sigma-Aldrich Cat. No. M3782





Advancements in Enzyme-Based Fuel Cell Batteries, Sensors and Emissions Reduction Strategies

Dr. Shelley Minteer, Ph. D.

Sigma® Life Science and Aldrich Materials Science merge to help researchers completely rethink tomorrow's immobilized enzyme technologies.

Today's fuel cells and batteries lack the ability to produce electrical power for the long periods of time we desire and commonly incorporate toxic or precious metals to produce and store their energy.

Enzyme-based fuel cells have potentially higher energy density than traditional batteries using much more environmentally compatible materials. While the idea of enzyme-based biofuel cells is not a new concept, recent developments have captured the interests of the research community. Dr. Shelley Minteer at Saint Louis University Department of Chemistry is one of the pioneers in this area. Dr. Minteer and her team have developed a new enzyme immobilization technology that has begun to revolutionize the way we think about using enzymes for long-term catalytic applications.

Dr. Minteer's enzyme immobilization strategy encapsulates and stabilizes Sigma enzymes in modified Nafion® polymer matrices from Aldrich. The efficacy of Dr. Minteer's immobilized enzymes can be measured in years compared to days for other biofuel cell technologies. Dr. Minteer has successfully demonstrated the use of several metabolic fuels such as glycerol, fatty acids, ethanol and other alcohols, pyruvate and glucose to provide the fuel source, while a wide variety of enzymes and coupled multienzyme systems provide the catalytic power to convert the fuel to electrical current.


In addition to synthetic polymer supports, Dr. Minteer and her collaborator, Dr. Michael Cooney of the University of Hawaii, have immobilized dehydrogenases on macroporous chitosan scaffolds. The scaffolds were fabricated from solutions of native and hydrophobically modified chitosan polymers through the process of thermally induced phase separation. The hydrophobically modified chitosan is proposed to possess amphiphilic micelles into which the enzyme can be encapsulated and retained.

One of Dr. Minteer's basic biofuel cell batteries employs an ethanol/alcohol dehydrogenase system to generate a proton gradient across the conductive membrane. They have also refined these dehydrogenase systems by replacing traditional NAD(H) cofactor-dependent dehydrogenases with pyrroloquinoline quinone(PQQ)-dependent dehydrogenases to improve efficiency. Other multienzyme coupled systems use classical glycolytic and Krebs cycle enzyme systems and corresponding metabolic substrates to generate power.


Dr. Minteer's lab has also utilized mitochondria to generate current by metabolism of classic cellular energy sources. The resulting electrons are shuffled through the electron transport chain where they reduce oxygen at cytochrome c oxidase and combine with the protons from oxidation to form water. During this metabolic process it is possible for cytochrome c to undergo direct electron transfer to a carbon electrode and complete oxygen reduction at a platinum cathode.

Taking advantage of the fact that nitroaromatics can decouple the inhibition of pyruvate metabolism, mitochondria-modified electrodes have also been used to electrochemically sense the presence of nitroaromatic explosive compounds. These oligomycin inhibited mitochondria-modified electrodes were employed as the bioanode in a pyruvate/air biofuel cell. This concept could be used as a traditional electrochemical sensor or as a self-powered sensor.

Dr. Minteer's Nafion-based immobilization technology is also addressing the reduction of fluegas emissions from coal-fired power plants using carbonic anhydrase to sequester CO₂.



Enzymes immobilized in Nafion® membrane are encapsulated in the polymer resulting in increased stability while retaining active site availability and optimized orientation.



An example of one of Dr Minteer's fuel cells; utilizing pyruvate as the initial fuel source, this fuel cell incorporates a multi-enzyme coupled system that generates a proton gradient source of energy analogous to that produced *in vivo* in the Krebs Cycle.

References

- (1) S.D. Minteer, B.Y. Liaw, and M.J. Cooney, "Enzyme-based biofuel cells," *Current Opinions in Biotechnology*, 2007, **18**, 228–234.
- (2) C.M. Moore, N.L. Akers, and S.D. Minteer, "Improving the Environment for Immobilized Dehydrogenase Enzymes by Modifying Nafion with Tetraalkylammonium Bromides," *Biomacromolecules*, **5**, 2004, 1241–1247.
- (3) M.J. Cooney, M. Windmeisser, B.Y. Liaw, C. Lau, T. Klitzbach, and S.D. Minteer, "Design of Chitosan Gel Pore Structure: Towards Enzyme Catalyzed Flow-Through Electrodes," *Journal of Materials Chemistry*, 2008, **18**, 667–674.
- (4) P. Atanassov, C. Aplett, S. Banta, S. Brozik, S. Calabrese Barton, M. Cooney, B.Y. Liaw, S. Mukerjee, and S.D. Minteer, "Enzymatic Biofuel Cells," *Electrochemical Society Interface*, **16**, 2007, 28–31.
- (5) T.L. Klitzbach, M. Watt, Y. Ansari, and S.D. Minteer, "Improving the microenvironment for enzyme immobilization at electrodes by hydrophobically modifying chitosan and Nafion polymers," *Journal of Membrane Science*, 2008, **311**, 81–88.
- (6) V. Svoboda, M. Cooney, B.Y. Liaw, S. Minteer, E. Piles, D. Lehnert, S.C. Barton, R. Rincon, and P. Atanassov, "Standardized characterization of biocatalytic electrodes," *Electroanalysis*, 2008, **20**, 1099–1109.
- (7) B.L. Treu, R. Arechederra, and S.D. Minteer, "Bioelectrocatalysis of Ethanol via PQQ-Dependent Dehydrogenases Utilizing Carbon Nanomaterial Supports," *Journal of Nanoscience and Nanotechnology*, 2009, **9**, 2374–2380.
- (8) B.L. Treu and S.D. Minteer, "Isolation and Purification of PQQ-dependent Lactate Dehydrogenase from Gluconobacter and Use for Direct Electron Transfer at Carbon and Gold Electrodes," *Bioelectrochemistry*, 2008, **74**, 71–77.
- (9) D. Sokic-Lazic and S.D. Minteer, "Citric acid cycle biomimetic on a carbon electrode," *Biosensors and Bioelectronics*, 2008, **24**, 945–950.
- (10) A.E. Blackwell, M.J. Moehlenbrock, J.R. Worsham, and S.D. Minteer, "Comparison of Electropolymerized Thiazine Dyes as an Electrocatalyst in Enzymatic Biofuel Cells and Self Powered Sensors," *Journal of Nanoscience and Nanotechnology*, 2009, **9**, 1719–1726.
- (11) M. Cooney, V. Svoboda, C. Lau, G. Martin, and S.D. Minteer, "Enzymatic Biofuel Cells," *Energy and Environmental Science*, 2008, **1**, 320–337.
- (12) M.J. Cooney, J. Petermann, C. Lau, S.D. Minteer, "Characterization and Evaluation of Hydrophobically Modified Chitosan Scaffolds: towards design of enzyme immobilized flow-through electrodes," *Carbohydrate Polymers*, 2009, **75**, 428–435.
- (13) G. Martin, S. Minteer, and M.J. Cooney, "Spatial distribution of malate dehydrogenase in amphiphilic chitosan scaffolds," *Applied Materials and Interfaces*, 2009, **1**(2) 367–372.
- (14) C. Hettige, S.D. Minteer, and S. Calabrese Barton, "Simulation of Multi-Step Enzyme Electrodes," *ECS Transactions*, 2008, **13**/21, 99–109.
- (15) K. Sjoholm, M.J. Cooney, and S.D. Minteer, "Effects of Degree of Deacetylation on Enzyme Immobilization in Hydrophobically Modified Chitosan," *Carbohydrate Polymers*, 2009, **77**, 420–424.
- (16) G. Martin, S.D. Minteer, and M.J. Cooney, "Fluorescence characterization of chemical microenvironments in hydrophobically modified chitosan," *Carbohydrate Polymers*, 2009, **77**, 695–702.
- (17) D. Sokic-Lazic and S.D. Minteer, "Pyruvate/air enzymatic biofuel cell capable of complete oxidation," *Electrochemical and Solid-State Letters*, 2009, **12**(9), F26–F28.
- (18) M.C. Beilke, T.L. Klitzbach, B.L. Treu, D. Sokic-Lazic, J. Wildrick, E.R. Amend, L.M. Gebhart, R.L. Arechederra, M.N. Germain, M.J. Moehlenbrock, Sudhanshu, and S.D. Minteer, "Enzymatic Biofuel Cells," in *Micro-Fuel Cells: Principles and Applications*, Elsevier, 2009, 179–241.
- (19) M.J. Moehlenbrock, R.L. Arechederra, K.H. Sjoholm, and S.D. Minteer, "Analytical Techniques for Characterizing Enzymatic Biofuel Cells," *Analytical Chemistry*, 2009, **81**(23), 9538–9545.
- (20) P. Addo, R. Arechederra, and S.D. Minteer, "Evaluating Enzyme Cascades for Methanol/Air Biofuel Cells Based on NAD-dependent Enzymes," *Electroanalysis*, 2010, **22**(7–8), 807–812.
- (21) R. Rincon, M. Germain, K. Artyushkova, M. Mojica, M. Germaine, P. Atanassov, and S.D. Minteer, "Structure and Electrochemical Properties of Electrocatalysts for NADH Oxidation," *Electroanalysis*, 2010, **22**(7–8), 799–806.
- (22) W. Gellert, M. Kesmez, J. Schumacher, N. Akers, and S.D. Minteer, "Biofuel Cells for Portable Power," *Electroanalysis*, 2010, **22**(7–8), 727–731.
- (23) D. Sokic-Lazic, R.L. Arechederra, B.L. Treu, and S.D. Minteer, "Oxidation of Biofuels: Fuel Diversity and Effectiveness of Fuel Oxidation through Multiple Enzyme Cascades," *Electroanalysis*, 2010, **22**(7–8), 757–764.
- (24) C. Lau, G. Martin, S.D. Minteer, and M.J. Cooney, "Development of a Chitosan Scaffold Electrode for Fuel Cell Applications," *Electroanalysis*, 2010, **22**(7–8), 793–798.
- (25) R. Arechederra, B.L. Treu, and S.D. Minteer, "Development of Glycerol/Oxygen Biofuel Cells," *Journal of Power Sources*, **173**, 2007, 156–161.
- (26) R. Arechederra and S.D. Minteer, "Organelle-Based Biofuel Cells: Immobilized Mitochondria at Carbon Paper Electrodes," *Electrochimica Acta*, 2008, **53**, 6698–6703.
- (27) M. Germain, R. Arechederra, and S.D. Minteer, "Nitroaromatic Actuation of Mitochondrial Bioelectrocatalysis for Self Powered Explosive Sensors," *Journal of the American Chemical Society*, 2008, **130**(46), 15272–15273.
- (28) R. Arechederra, K. Boehm, and S.D. Minteer, "Mitochondrial Bioelectrocatalysis for Biofuel Cell Applications," *Electrochimica Acta*, 2009, **54**, 7268–7273.



Selected Enzymes for Biofuel Cell Research

Enzyme-based biofuel cells are incorporating a wide variety of enzymes. Sigma-Aldrich manufactures hundreds of enzymes, many of which may have potential applications in biofuel cells. For a complete list of enzymes, data sheets, assay protocols, and other enzyme-related technical resources visit the web-based Sigma Enzyme Explorer at sigma-aldrich.com/enzymeexplorer. The Enzyme Explorer also features the IUBMB on-line metabolic pathway chart with over 500 links to enzymes, substrates and cofactors.

The following products are only a small sample of the enzymes available from Sigma-Aldrich of potential use in biofuel cell research.

▼ Alcohol Dehydrogenase

ADH

Alcohol Dehydrogenase from *Parvibaculum lavamentivorans*

► recombinant, expressed in *Escherichia coli*, activity: ≥40.0 units/mg

One unit corresponds to the amount of enzyme which reduces 1 μmole ethyl acetate per minute at pH 7.0 and 30 °C.

75449-1ML	1 mL
75449-5ML	5 mL

Alcohol Dehydrogenase from *Saccharomyces cerevisiae*

Alcohol:NAD⁺ oxidoreductase; Alcohol Dehydrogenase from yeast [9031-72-5]

A 141 kDa tetramer containing 4 equal subunits. The active site of each subunit contains a zinc atom. Each active site also contains 2 reactive sulphhydryl groups and a histidine residue.

Isoelectric point: 5.4–5.8

Optimal pH: 8.6–9.0

Substrates: Yeast ADH is most active with ethanol and its activity decreases as the size of the alcohol increases or decreases. Branched chain alcohols and secondary alcohols also have very low activity.

K_M (ethanol) = 2.1 × 10⁻³ M

K_M (methanol) = 1.3 × 10⁻¹ M

K_M (isopropanol) = 1.4 × 10⁻¹ M

Inhibitors: Compounds that react with free sulphydryls, including N-alkylmaleimides and iodoacetamide.

Zinc chelator inhibitors, including 1,10-phenanthroline, 8-hydroxyquinoline, 2,2'-dipyridyl, and thiourea. Substrate analogue inhibitors, including β-NAD analogs, purine and pyrimidine derivatives, chloroethanol, and fluoroethanol.

Extinction Coefficient: E^{1%} = 14.6 (water, 280 nm)

One unit will convert 1.0 μmole of ethanol to acetaldehyde per min at pH 8.8 at 25 °C.

► lyophilized powder (contains buffer salts), activity: ≥300 units/mg protein

Solids containing ≤ 2% citrate buffer salts suitable for conventional determination of β-NAD, β-NADH, ethanol and acetaldehyde

Contains bound β-NAD and β-NADH and is not suitable for the recycling microassay of β-NAD and β-NADH. If you require ADH for this purpose, see Cat. No. A3263.

A7011-7.5KU	7500 units
A7011-15KU	15000 units
A7011-30KU	30000 units
A7011-75KU	75000 units
A7011-150KU	150000 units
A7011-300KU	300000 units

► activity: ≥300 units/mg protein

Solids containing <2% citrate buffer salts suitable for recycling micro-assay of β-NAD and β-NADH

A3263-7.5KU	7500 units
A3263-15KU	15000 units
A3263-30KU	30000 units
A3263-75KU	75000 units
A3263-150KU	150000 units

Alcohol Dehydrogenase, *Deinococcus radiodurans*, recombinant from *E. coli*

ADH
[9031-72-5]

► activity: ≥10000 units/mL

One unit corresponds to the amount of enzyme which reduces 1 μmole ethyl-2-oxo-4-phenylbutyrate per minute at pH 7.0 and 37 °C

16892-1ML	1 mL
16892-5ML	5 mL

Alcohol Dehydrogenase, NADP⁺ dependent from *Thermoanaerobium brockii*

Alcohol:NADP⁺ oxidoreductase; TBADH
[9028-12-0]

An extremely thermostable enzyme with broad substrate specificity for alcohols, ketones and acetaldehyde.

One unit will oxidize 1.0 μmole of 2-propanol to acetone per min at pH 7.8 at 40 °C in the presence of NADP⁺.

► lyophilized powder, activity: 5–15 units/mg protein

Contains phosphate buffer salts and dithioerythritol

A8435-100UN	100 units
A8435-250UN	250 units

► lyophilized powder, activity: 30–90 units/mg protein

Purified

Contains sodium citrate and dithioerythritol.

A9287-100UN	100 units
--------------------	-----------

Alcohol Dehydrogenase ▲

Bilirubin Oxidase from *Myrothecium verrucaria*

Bilirubin:oxygen oxidoreductase
[80619-01-8]

► lyophilized powder, activity: 15–65 units/mg protein

One unit will oxidize 1.0 μmole of bilirubin per min at pH 8.4 at 37 °C.

B0390-25UN	25 units
B0390-100UN	100 units

Cytochrome c Oxidase from bovine heart

Ferrocytochrome-c:oxygen oxidoreductase, Complex IV; EC 1.9.3.1
[9001-16-5]

Cytochrome c oxidase is the principal terminal oxidase of high oxygen affinity in the aerobic metabolism of all animals, plants, yeasts and some bacteria.

Supplied as a solution in 25 mM Tris-HCl buffer, pH 7.8, 5 mM EDTA, and 39 mM sodium dodecyl maltoside.

activity: ≥20 units/mg protein

One unit of cytochrome c oxidase will oxidize 1.0 μmole of ferrocytochrome c per min. at 25 °C at pH 6.0.

C5499-250UG 250 μg

D-Fructose Dehydrogenase from *Gluconobacter industrius*

D-Fructose:(acceptor) 5-oxidoreductase
[37250-85-4]

► lyophilized powder, activity: 400–1,200 units/mg protein

Lyophilized powder containing citrate-phosphate buffer salts, TRITON® X-100, and stabilizer

One unit will convert 1.0 μmole D-fructose to 5-ketofructose per min at pH 4.5 at 37 °C.

F5152-250UN 250 units

Glucose Dehydrogenase from calf liver

β-D-Glucose: NAD[P]⁺ 1-oxidoreductase
[9028-53-9]

► activity: ~15 units/g solid

Crude powder

One unit will oxidize 1.0 μmole of β-D-glucose to D-glucono-δ-lactone per min at pH 7.6 at 25 °C.

G5625 Inquire

Glucose Dehydrogenase from *Pseudomonas sp.*

[9028-53-9]

► powder, activity: ≥200 units/mg

One unit corresponds to the amount of enzyme which will oxidizes 1 μmole β-D-glucose to D-glucono-δ-lactone per minute at pH 8.0 and 37 °C

19359-10MG-F 10 mg

Glucose Dehydrogenase from *Thermoplasma acidophilum*

β-D-Glucose:NAD[P]⁺ 1-oxidoreductase
[9028-53-9]

► recombinant, expressed in *Escherichia coli*, ammonium sulfate suspension, activity: ≥100 units/mg protein (glucose substrate), activity: ≥50 units/mg protein (galactose substrate)

A thermostable enzyme with resistance to organic solvents. Enzyme also shows substantial activity with D-galactose as substrate.

Suspension in 2.8 M (NH₄)₂SO₄, 50 mM phosphate buffer, pH 7.0

One unit will oxidize 1.0 μmole of D-glucose to D-glucono-δ-lactone per min at pH 7.0 at 37 °C in the presence of NADP⁺.

G5909-100UN	100 units
G5909-500UN	500 units

► lyophilized powder, activity: 100,000–250,000 units/g solid (without added oxygen)

May contain traces of amylase, maltase, glycogenase, invertase, and galactose oxidase.

Protein determined by biuret

G7141-10KU	10000 units
G7141-50KU	50000 units
G7141-250KU	250000 units
G7141-1MU	1000000 units
G7141-2.5MU	2500000 units

Laccase

[80498-15-3]

Laccase from *Rhus vernicifera*

Benzenediol:oxygen oxidoreductase

► crude acetone powder, activity: ≥50 units/mg solid

One unit will produce a ΔA₅₃₀ of 0.001 per min at pH 6.5 at 30 °C in a 3 mL reaction volume using syringaldazine as substrate.

L2157-10KU 10000 units

Laccase from *Trametes versicolor*

former nomenclature: *Coriolus versicolor*

► powder, activity: >20 units/mg

One unit corresponds to the amount of enzyme which converts 1 μmole of catechol per minute at pH 4.5 and 25 °C

53739-100MG-F	100 mg
53739-1G-F	1 g

► powder, activity: ≥0.5 units/mg

One unit corresponds to the amount of enzyme which converts 1 μmole of catechol per minute at pH 6.0 and 25 °C

38429-1G	1 g
38429-10G	10 g

Laccase ▲



Enzymes and Reagents for Ethanol Research

Enzymes for Lignocellulosic Ethanol Research


Enzymes for Lignocellulosic Ethanol Research

Sigma-Aldrich is committed to helping unlock today's challenge of releasing the energy of glucose that nature has cleverly locked into lignocellulosic biomass. To the untrained eye, starch and cellulose structures appear almost identical. Yet the barriers between cellulosic biomass hydrolysis and ethanol production appear have proven to be orders of magnitude more challenging than with starch derived biomass.

Although many of the biochemical aspects of these challenges are well understood, they remain an enigma of nature. The two main, and closely related hindrances to the enzymatic hydrolysis of cellulosic polymers are:

How do we physically dismantle biomass on a molecular level to yield cellulosic polymers available to hydrolysis by cellulase and other glycolytic enzymes? How do we

penetrate the barrier posed by the closely integrated network of lignins surrounding cellulose microfibrils?



Visit the Enzyme Explorer Assay Library


Features over 600 detailed procedures for measuring enzyme activity and related metabolites.

The Library is the result of over ten years of in-house work on analytical procedures developed by Sigma-Aldrich® scientists.

Visit sigma.com/enzymeexplorer



Cellulases and Related Enzyme Activities




Cellulase catalyzes the endohydrolysis of β -(1→4)-D-glucosidic linkages in cellulose, lichenin and cereal β -D-glucans


Most cellulase preparations contain various classes of cellulases and related activities.

- The endo- β -(1→4)-D-glucanases (EGs) cleave cellulose in its internal non-terminal regions, yielding oligosaccharides.
- The two forms of cellobiohydrolases (CBH) are considered exo- β -(1→4)-D-glucanases, releasing cellobiose disaccharides from the cellulose chain reducing end (CBHI), and non-reducing end (CBHII).
- β -Glucosidase is an exohydrolase that yields glucose monosaccharides from soluble oligosaccharides such as cellobiose.

Other activities that can be found in cellulase preparations include xylanase, hemicellulase, and laminaranase.



Cereal β -Glucan is a polymer of β -(1→4)-D-glycopyranosyl units occupying as predominantly cellobiose and cellobiotetraose separated by single β -(1→3)-D-glycopyranosyl units. Crosslinking can occur within the consecutive cellobiose regions.



Hemicelluloses are a group of plant-derived heteropolysaccharides associated with cellulose and lignin. The most common hemicelluloses are xylan, glucuronoxylan, arabinoxylan, glucomannan and xyloglucan. In angiosperms, the principal hemicellulose component, xylan, is a polymer of β -(1→4)-D-xylopyranose. In arabinoxylan, branching occurs at the C2 & C3 positions with a-L-arabinofuranose. Glucuronoxylan, also found in angiosperms, has the xylan backbone with 4-O methylglucuronic acid branching. In addition, arabinose branching as well as acetylation may be present. Gymnosperms contain glucomannans comprised primarily of D-mannosyl and D-glucosyl residues..

- Hemicellulase preparations are typically a mixture of glycolytic enzymes containing xylanase, mananase and other activities.
- Xylanase catalyzes the endohydrolysis of β -(1→4)-D-xylosidic linkages in xylans yielding various β -(1→4)-D-xylooligosaccharides.

Cellobiase from *Aspergillus niger*

Cellobiase enzyme preparation obtained by submerged fermentation of an *Aspergillus niger* microorganism. The cellobiase hydrolyzes cellobiose to glucose.

- Novozyme 188 liquid, activity: ≥250 units/g

A product of Novozyme Corp.

C6105-50ML	50 mL
C6105-250ML	250 mL

▼ Cellulase

[9012-54-8]

Cellulase from *Aspergillus* sp.

Carezyme 1000L®

activity: ≥1000 units/g

produced by submerged fermentation of a genetically modified *Aspergillus* microorganism

A product of Novozyme Corp.

C2605-50ML	50 mL
C2605-250ML	250 mL

Cellulase from *Trichoderma longibrachiatum*

- powder, activity: ≥1.0 unit/mg solid

An acid cellulase with xylanase, pectinase, mannanase, xyloglucanase, laminarase, β-glucosidase, β-xylosidase, α-L-arabinofuranosidase, amylase, and protease activities

One unit will liberate 1.0 μmole of glucose from cellulose in one hour at pH 5.0 at 37 °C (2 hr incubation time).

C9748-100G	100 g
-------------------	-------

Cellulase from *Trichoderma reesei*

ATCC 26921

1,4-(1,3;1,4)-β-D-Glucan 4-glucano-hydrolase

- lyophilized powder, activity: ≥1 unit/mg solid

One unit will liberate 1.0 μmole of glucose from cellulose in one hr at pH 5.0 at 37 °C (2 hr incubation time).

C8546-2.5KU	2500 units
C8546-5KU	5000 units
C8546-10KU	10000 units

- Celluclast® 1.5L aqueous solution, activity: ≥700 units/g

Produced by submerged fermentation of a selected strain of the fungus *Trichoderma reesei* and catalyzes the breakdown of cellulose into glucose, cellobiose, and higher glucose polymers.

A product of Novozyme Corp.

C2730-50ML	50 mL
-------------------	-------

Cellulase from *Trichoderma viride*

1,4-(1,3;1,4)-β-D-Glucan 4-glucano-hydrolase

One unit will liberate 1.0 μmole of glucose from cellulose in one hour at pH 5.0 at 37 °C (2 hr incubation time).

- crude powder, activity: 3–10 units/mg solid

C9422-5KU	5000 units
C9422-10KU	10000 units

- Onozuka RS powder, activity: ≥5,000 units/g solid

Manufactured by Yakult

C0615-1G	1 g
-----------------	-----

Cellulase ▲**Driselase® from *Basidiomycetes* sp.**

[85186-71-6]

Crude powder containing laminarinase, xylanase and cellulase.

- powder

D9515-1G	1 g
D9515-5G	5 g
D9515-25G	25 g

▼ β-Glucanase**β-Glucanase from *Aspergillus niger***

[9074-98-0]

- powder, activity: ~1 units/mg

One unit corresponds to the amount of enzyme which will release 1 μmole of reducing sugar equivalents (expressed as glucose) per minute at pH 5.0 and 55 °C, using β-D-glucan (**Cat. No. 49102**) as substrate

49101-100MG	100 mg
49101-500MG	500 mg

β-Glucanase from *Bacillus subtilis*

[9074-98-0]

- powder, activity: ~1 units/mg

One unit corresponds to the amount of enzyme which will release 1 μmol of reducing sugar equivalents (expressed as glucose) per minute of pH 6.0 and 55 °C, using β-D-glucan (**Cat. No. 49102**) as substrate

49106-100MG	100 mg
--------------------	--------

β-Glucanase from *Trichoderma longibrachiatum*

A mixture composed mainly of β-1→3 / 1→4-glucanase, xylanase, and cellulase activities. β-glucosidase, β-xylosidase, α-L-arabinofuranosidase,

amylase, and protease activities are also present.

activity: ≥1.0 units/mg solid

One unit will liberate 1.0 μmole of glucose from cellulose in one hr at pH 5.0 at 37 °C (2 hr incubation time).

G4423-100G

100 g

β-Glucanase ▲**β-(1→3)-D-Glucanase from *Helix pomatia***

[9044-93-3]

- powder, activity: 0.5–1.5 units/mg

One unit corresponds to the amount of enzyme which liberates 1 μmole of glucose from laminarin (**Cat. No. 61340**) per minute at pH 5.0 and 37 °C

Improved filterability of wines by enzymic decomposition of carbohydrate-containing colloids¹; Induction of hydrolases as a defense reaction against pathogens, review²

Lit cited: 1. Wucherpfennig, K., and Dietrich, H., *Weinwirtschaft* **118**, 598 (1982); 2. Boller, T., *UCLA Symp. Mol. Cell Biol.*, New Ser. **22**, 247 (1985)

49103-10MG

10 mg

49103-50MG

50 mg

β-Glucosidase from almonds

β-D-Glucoside glucohydrolase

[9001-22-3]

- lyophilized powder, activity: ≥2 units/mg solid

Crude

One unit will liberate 1.0 μmole of glucose from salicin per min at pH 5.0 at 37 °C.

G0395-2.5KU

2500 units

G0395-5KU

5000 units

G0395-50KU

50000 units

- lyophilized powder, activity: 20–40 units/mg solid

Chromatographically purified

One unit will liberate 1.0 μmole of glucose from salicin per min at pH 5.0 at 37 °C.

G4511-100UN

100 units

G4511-250UN

250 units

G4511-1KU

1000 units

Hemicellulase from *Aspergillus niger*

[9025-56-3]

- powder, activity: 0.3–3.0 unit/mg solid (using a β-galactose dehydrogenase system and locust bean gum as substrate)

An undefined mixture of glycolytic enzymes usually containing xylanase, mannanase and other activities.

One unit will produce a relative fluidity change of 1 per 5 minutes using locust bean gum as substrate at pH 4.5 at 40 °C

H2125-150KU

150000 units

▼ Laminarinase

Endo-1,3(4)- β -D-Glucan 3(4)-glucanohydrolase
[62213-14-3]

Laminarinase from *Penicillium* sp.

1,3-(1,3;1,4)- β -D-Glucan 3(4)-glucanohydrolase

- lyophilized powder, activity: 5–10 units/mg protein

Lyophilized powder containing acetate buffer salts

One unit will liberate 1.0 mg of reducing sugar (measured as glucose) from laminarin per min at pH 5.0 at 37 °C.

L9259-25UN 25 units

Laminarinase from *Trichoderma* sp.

1,3-[1,3;1,4]- β -D-Glucan 3(4)-glucanohydrolase

View more information on enzymes for complex carbohydrate analysis at sigma-aldrich.com/enzymeexplorer

- powder, activity: 100–400 units/g solid

Contains chitinase activity.

One unit will liberate 1.0 mg of reducing sugar (measured as glucose) from laminarin per min at pH 5.0 at 37 °C.

L5272-5UN 5 units
L5272-25UN 25 units

Laminarinase ▲

Online Metabolomics Resource

Use the IUBMB–Sigma-Aldrich Interactive Metabolic Pathway Chart to find the Metabolite Standards you need.

The Metabolic Pathways Map contains over 500 hyperlinks to Sigma product listings. Just click on the metabolite name or the enzyme's E.C. number to access product information.

You can access the chart at
sigma-aldrich.com/metpath

Xylanase from *Trichoderma viride*

endo-1,4- β -Xylanase; 1,4- β -D-Xylanxylohydrolase [9025-57-4]

- lyophilized powder, activity: 100–300 units/mg protein

Contains sorbitol and sodium acetate buffer salts

One unit will liberate 1 μ mole of reducing sugar measured as xylose equivalents from xylan (Cat. No. X0627) per min at pH 4.5 at 30 °C.

X3876-250UN	250 units
X3876-1KU	1000 units

Xylanase ▲**endo-1,4- β -Xylanase from *Trichoderma longibrachiatum***


1,4- β -D-Xylanxylohydrolase

Primary activity is an acid-neutral endo-1,4- β -D-xylanase, additional activities include β -glucanase, cellulase, pectinase, mannanase, xyloglucanase, laminarase, β -glucosidase, β -xylosidase, α -L-arabinofuranosidase, amylase, and protease.

activity: ≥1.0 units/mg solid


One unit will liberate 1 μ mole of reducing sugar measured as xylose equivalents from xylan (Cat. No. X0627) per min at pH 4.5 at 30 °C.

X2629-100G	100 g
-------------------	-------




Enzymes for Pectin Hydrolysis


Pectins are complex-branched heteropolysaccharides primarily containing an α -(1 \rightarrow 4) polygalacturonic acid backbone which can be randomly acetylated and methylated. Three different pectins have been isolated from plant cell walls.



Homogalacturonan is composed of a poly- α -(1 \rightarrow 4)-galacturonic acid backbone with random partial methylation and acetylation.



Rhamnogalacturonan I is composed of alternating α -(1 \rightarrow 2)-L-rhamnosyl- α -(1 \rightarrow 4)-D-galacturonosyl backbone with two types of branching composed of a ribofuranose or galactose oligomers.



Rhamnogalacturonan II is composed of poly- α -(1 \rightarrow 4)-D-galacturonic acid backbone with random partial methylation, acetylation and four types of branching.

▼ Pectinase**Pectinase from *Aspergillus aculeatus***

- Pectinex® Ultra SPL aqueous solution, activity: ≥9,500 units/mL

Highly active pectolytic enzyme preparation produced by a selected strain of *Aspergillus aculeatus*. A product of Novozyme Corp.

P2611-50ML	50 mL
P2611-250ML	250 mL

Pectinase from *Aspergillus niger*

Polygalacturonase solution from *Aspergillus niger*; Poly-(1,4- α -D-galacturonide) glycanohydrolase [9032-75-1]

Used in plant protoplast preparation to digest cell wall prior to organelle isolation.

- aqueous glycerol solution, activity: ≥5 units/mg protein (Lowry)

Solution in 40% glycerol

One unit will liberate 1.0 μ mole of galacturonic acid from polygalacturonic acid per min at pH 4.0 at 25 °C.

P4716-5KU	5000 units
P4716-10KU	10000 units
P4716-25KU	25000 units
P4716-100KU	100000 units

Pectinase from *Aspergillus niger*

- Pectinex 3XL® aqueous solution

Pectolytic enzyme preparation produced from a selected strain of *Aspergillus niger*: contains mainly pectintranseliminase, polygalacturonase, and pectinesterase and small amounts of hemicellulases and cellulases.

A product of Novozyme Corp.

P2736-50ML	50 mL
P2736-250ML	250 mL

Pectinase from *Rhizopus* sp.

Poly-(1,4- α -D-galacturonide) glycanohydrolase; Poly-galacturonase [9032-75-1]

Used in plant protoplast preparation to digest cell wall prior to organelle isolation.

- Macerozyme R-10 powder, activity: 400–800 units/g solid

A source of pectinase activity, also containing cellulase and hemicellulase activities.

One unit will liberate 1.0 μ mole of galacturonic acid from polygalacturonic acid per min at pH 4.0 at 25 °C.

P2401-500UN	500 units
P2401-1KU	1000 units
P2401-5KU	5000 units

Pectinase ▲

Pectolyase from *Aspergillus japonicus*

Reported to contain two types of pectinase, endopolypalacturonase (EC3.2.1.15), endopectin lyase (EC4.2.2.10) and a maceration stimulating factor.

Used in plant protoplast preparation to digest cell wall prior to organelle isolation.

- lyophilized powder, activity: ≥0.3 units/mg solid

One unit will liberate 1.0 μ mole of galacturonic acid from polygalacturonic acid per min at pH 5.5 at 25 °C.

P3026-100MG	100 mg
P3026-250MG	250 mg
P3026-1G	1 g

Feeling the Heat? Then Capture It!



ALDRICH
Chemistry

Aldrich Materials Science — We focus on materials so you can focus on results.

Materials for Renewable and Alternative Energy Solutions

- Solar Energy Conversion (Organic, Crystalline and Thin Film Photovoltaics)
- Hydrogen Economy
- Supercapacitors/Advanced Batteries
- Alternative Fuels



Check out the Sigma-Aldrich Renewable & Alternative Energy Web Portal: sigma-aldrich.com/renewable

sigma-aldrich.com

SIGMA-ALDRICH®

Enzymes for Lignin Depolymerization

To some, lignin appears to be a foreign substance resembling a randomly modified polystyrene polymer, not resembling any of the more familiar biopolymers composed of carbohydrates, amino acids or nucleotides. However, next to cellulose, it is the second most abundant biopolymer on earth comprising as much as 30% of the dry mass of wood. Unlike other biopolymers, enzymatic degradation is oxidative versus hydrolytic and is limited to a select group of fungi and bacteria. Lignins appear to be heterogeneous and random with regards to degree of polymerization, branching, and monomer composition and sequence.



▼ Laccase

[80498-15-3]

Laccase from *Rhus vernicifera*

Benzenediol:oxygen oxidoreductase

- crude acetone powder, activity: ≥ 50 units/mg solid

One unit will produce a ΔA_{530} of 0.001 per min at pH 6.5 at 30 °C in a 3 mL reaction volume using syringaldazine as substrate.

L2157-10KU	10000 units
-------------------	-------------

Laccase from *Trametes versicolor*former nomenclature: *Coriolus versicolor*

- powder, activity: >20 units/mg

One unit corresponds to the amount of enzyme which converts 1 μ mole of catechol per minute at pH 4.5 and 25 °C

53739-100MG-F	100 mg
53739-1G-F	1 g

- powder, activity: ≥ 0.5 units/mg

One unit corresponds to the amount of enzyme which converts 1 μ mole of catechol per minute at pH 6.0 and 25 °C

38429-1G	1 g
38429-10G	10 g

Laccase, *Coriolus versicolor*, CLEA

- Laccase, *Coriolus versicolor*, cross-linked enzyme aggregate activity: ≥ 0.3 units/mg

One unit corresponds to the amount of enzyme which oxidizes 1 μ mol ABTS per minute at pH 4.5 and 25°C

38837-100MG	100 mg
--------------------	--------

Laccase ▲**Lignin Peroxidase**

LiP; Peroxidase, lignin; Ligninase
[42613-30-9]

- powder, activity: >0.1 units/mg

One unit corresponds to the amount of enzyme, which oxidizes 1 μ mole 3,4-dimethoxybenzyl alcohol per minute at pH 3.0 and 30 °C

42603-10MG-F	10 mg
---------------------	-------

Manganese Peroxidase from *Nematoloma frowardii*

[114995-15-2]

41563

Inquire

Manganese peroxidase from white-rot fungus (*Phanerochaete chrysosporium*)

MnP; Peroxidase, manganese; Manganese-dependent lignin peroxidase
[114995-15-2]


- powder, activity: ≥ 20 U/g

only partially soluble in water or buffer

One unit corresponds to the amount of enzyme, which oxidizes 1 μ mole Mn²⁺ per minute to Mn³⁺ at pH 4.5 and 25 °C

This product is supplied as a mixture of isozymes, mol. mass (40–65 kDa).

93014-10MG-F	10 mg
---------------------	-------




Lignin peroxidase catalyzes the oxidative cleavage of carbon-carbon and ether bonds in lignin-related compounds.



Manganese Peroxidase

Manganese peroxidase is a hemoprotein involved in the oxidative degradation of lignin in white rot basidiomycetes.

Enzymes for Starch Hydrolysis




Enzymatic hydrolysis of starch-based biomass to yield monomeric glucose is often a multi-enzyme process and can involve α -amylase, β -amylase, amyloglucosidase, isoamylase and α -glucosidase.

- α -Amylase catalyzes the endohydrolysis of α -(1 \rightarrow 4)-D-glucosidic linkages in polysaccharides containing three or more α -(1 \rightarrow 4)-linked D-glucose units.
- β -Amylase catalyzes the exohydrolysis of α -(1 \rightarrow 4)-D-glucosidic linkages in polysaccharides resulting in the successive liberation of maltose units from the non-reducing ends of the chains.
- Amyloglucosidase (glucoamylase) catalyzes the hydrolysis of terminal α -(1 \rightarrow 4)-D-glucose residues successively from the non-reducing ends of maltooligo- and polysaccharides with release of β -D-glucose. Most forms of the enzyme can rapidly hydrolyze α -(1 \rightarrow 6)-D-glucosidic bonds when the next bond in the sequence is

1 \rightarrow 4- and some preparations of this enzyme hydrolyze 1 \rightarrow 6- and 1 \rightarrow 3-D-glucosidic bonds in other polysaccharides.

- Isoamylase catalyzes the hydrolysis of α -(1 \rightarrow 6)-D-glucosidic branch linkages in amylopectin and β -limit dextrans.

- α -Glucosidase catalyzes the hydrolysis of terminal 1 \rightarrow 4-linked α -D-glucose residues successively from the non-reducing ends of maltooligo- and to a lesser extent polysaccharides with release of β -D-glucose. Most forms of the enzyme can slowly hydrolyze α -(1 \rightarrow 6)-D-glucosidic bonds.



▼ α -Amylase

α -Amylase from *Aspergillus oryzae*

1,4- α -D-Glucan-glucanohydrolase
[9001-19-8]

- lyophilized powder, activity: 150–250 units/mg protein (biuret)

Crude

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 6.9 at 20 °C.

A6211 Inquire

- Fungamyl® 800 L aqueous solution, activity: ≥800 FAU/g

A product of Novozyme Corp.

A8220-50ML	50 mL
A8220-250ML	250 mL

► Taka-Diastase from *Aspergillus oryzae*; Taka-Amylase A powder, activity: ~1.5 units/mg (~0.2 U acc. to Willstätter)

One unit corresponds to the amount of enzyme which liberates 1 μ mol maltose per minute at pH 6.0 and 25 °C (starch acc. to Zulkowsky, **Fluka No. 85642**, as substrate).

Molecular weight¹; Product distribution in the hydrolysis of amylose²; Application in (selective) hydrolysis/condensation of glycosidic bonds³; Suitable for use in thiamine determination⁴.

Lit cited: 1. T. Takagi, *J. Biochem.* **89**, 363 (1981); 2. H. Kondo et al., *J. Biochem.* **87**, 1053 (1980); 3. K. Fujita et al., *Bull. Chem. Soc. Jpn.* **62**, 3150 (1989); 4. P.W. Defibaugh, *J. Assoc. Off. Anal. Chem.* **70**, 514 (1987)

86250-100G	100 g
86250-500G	500 g

- spray-dried (powder), activity: ≥150 units/mg protein (biuret)

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 6.9 at 20 °C.

A9857-250KU	250000 units
A9857-1MU	1000000 units
A9857-5MU	5000000 units

α -Amylase from *Bacillus* sp.

- ammonium sulfate suspension, activity: 40–90 units/mL packed gel

Suspension in 2.0 M $(\text{NH}_4)_2\text{SO}_4$, pH 7.0

One unit will liberate 1.0 μ mole of maltose from starch per min at pH 6.9 at 30 °C. (Each unit is equivalent to 0.7 of mg maltose liberated per 3 min. at 20 °C)

A0909 Inquire

α -Amylase from *Bacillus amyloliquefaciens*

1,4- α -D-Glucan Glucano-hydrolase
[9000-85-5]

An endoamylase that randomly hydrolyzes α (1→4)-glycosidic linkages in amylose and amylopectin. The breakdown products are oligosaccharides and dextrans of varying chain length. This enzyme is active at high temperatures (70–90 °C).

One unit is the amount of enzyme which dextrinizes 5.26 g of dry starch per hour under standard conditions.

- BAN™ 240L liquid, activity: ≥250 units/g
A product of Novozyme Corp.

A7595-50ML	50 mL
A7595-250ML	250 mL

α -Amylase from *Bacillus licheniformis*

1,4- α -D-Glucan-glucanohydrolase
[9000-85-5]

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 6.9 at 20 °C.

- Termamyl® 120 saline solution, activity: 500–1,000 units/mg protein (biuret)

Aqueous solution containing approx. 15% sodium chloride and 25% sucrose.

Reported to be heat stable at temperatures as high as ~90 °C.

A product of Novozyme Corp.

A3403-500KU	500000 units
A3403-1MU	1000000 units
A3403-5MU	5000000 units

- lyophilized powder, activity: 500–1,500 units/mg protein

Heat-stable formulation

Lyophilized powder containing potassium phosphate

A4551-100MG	100 mg
A4551-1G	1 g

α -Amylase from *Bacillus subtilis*

[9000-90-2]

- powder, activity: ~380 units/mg

One unit is the amount of enzyme which liberates 1 μ mole of maltose per minute at pH 6.9 and 25 °C (using **Cat. No. 85642** as substrate)

Heat stability of bacterial α -amylases¹; Action pattern on sweet potato starch, amylose and amylopectin²; Action on native wheat starch.³

Lit cited: 1. J.E. Anderson et al., *J. Food Sci.* **48**, 1622 (1983); 2. P.L. Chang Rupp, S.J. Schwartz, *J. Food Biochem.* **12**, 191 (1988); 3. P. Colonna et al., *Biotechnol. Bioeng.* **31**, 895 (1988).

10069-250MG	250 mg
10069-1G	1 g

- powder, activity: ~50 units/mg

One unit corresponds to the amount of enzyme which liberates 1 μ mole maltose per minute at pH 6.9 at 25 °C (starch acc. to Zulkowsky, **Cat. No. 85642**, as substrate).

Applications in (selective) hydrolysis/condensations of glycosidic bonds.^{1,2}

Lit cited: 1. Y.Takasaki, *Agric. Biol. Chem.* **49**, 1091 (1985); 2. C.S. Su, C.P. Yang, *J. Chem. Technol. Biotechnol.* **48**, 313 (1990)

10070-10G	10 g
10070-50G	50 g

α -Amylase from barley malt

1,4- α -D-Glucan-glucanohydrolase
[9000-90-2]

- powder

β -amylase activity: ≥1 unit/mg solid

α -amylase activity: ≥1 unit/mg solid

Package size based on α -amylase activity

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 6.9 at 20 °C.

A2771 Inquire

α -Amylase from human pancreas

[9000-90-2]

- lyophilized powder, activity: ≥100 units/mg protein

Lyophilized from Tris buffer containing NaCl and CaCl_2 .

Prepared by modified method of Levitzki et al.

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 6.9 at 20 °C.

A9972-100UG	100 μ g
--------------------	-------------

α -Amylase from human saliva

1,4- α -D-Glucan-glucanohydrolase
[9000-90-2]

- lyophilized powder, activity: 300–1,500 units/mg protein

Lyophilized powder containing $(\text{NH}_4)_2\text{SO}_4$ and sodium citrate.

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 6.9 at 20 °C.

A1031-1KU	1000 units
A1031-5KU	5000 units

- lyophilized powder, activity: 1,000–3,000 units/mg protein

Lyophilized powder containing $(\text{NH}_4)_2\text{SO}_4$ and sodium citrate

Chromatographically purified

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 6.9 at 20 °C.

A0521-100UN	100 units
A0521-500UN	500 units
A0521-2.5KU	2500 units

α -Amylase from porcine pancreas

Molecular mass: 51–54 kDa.

α -Amylase isolated from porcine pancreas is a glycoprotein. It is a single polypeptide chain of ~475 residues containing two SH groups and four disulfide bridges and a tightly bound Ca^{2+} necessary for stability. Chloride ions are necessary for activity and stability. The pH range for activity is 5.5 to 8.0, with the pH optimum at 7.

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 6.9 at 20 °C.

► activity: ≥10 units/mg solid

Contains lactose

Package size based on α -amylase activity

A3176-500KU	500000 units
A3176-1MU	1000000 units
A3176-2.5MU	2500000 units
A3176-5MU	5000000 units
A3176-10MU	10000000 units

► DFP Treated, saline suspension, activity: 700–1400 units/mg protein ($E^{1\%}_{280}$)

Suspension in 2.9 M NaCl solution containing 3 mM CaCl_2

DFP treated. 2x crystallized

View more information on enzymes for complex carbohydrate analysis at sigma-aldrich.com/enzymeexplorer

A6255-10MG	10 mg
A6255-25MG	25 mg
A6255-100MG	100 mg

► DFP treated, ammonium sulfate suspension, activity: ≥500 units/mg protein ($E^{1\%}_{280}$)

Suspension in 3.2 M $(\text{NH}_4)_2\text{SO}_4$, pH 6.1

A2643	Inquire
--------------	---------

► PMSF treated, saline suspension, activity: 700–1400 units/mg protein ($E^{1\%}_{280}$)

Suspension in 2.9 M NaCl solution containing 3 mM CaCl_2

2x crystallized

A4268-25MG	25 mg
A4268-100MG	100 mg

 α -Amylase ▲**▼ β -Amylase**

1,4- α -D-Glucan maltohydrolase
[9000-91-3]

 β -Amylase from barley

► activity: 20–80 units/mg protein (biuret)

Crude

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 4.8 at 20 °C.

A7130-10KU	10000 units
A7130-50KU	50000 units
A7130-250KU	250000 units

 β -Amylase from sweet potato► ammonium sulfate suspension, activity: ≥750 units/mg protein ($E^{1\%}_{280}$)

Crystalline suspension in 2.3 M $(\text{NH}_4)_2\text{SO}_4$

One unit will liberate 1.0 mg of maltose from starch in 3 min at pH 4.8 at 20 °C.

View more information on enzymes for complex carbohydrate analysis at sigma-aldrich.com/enzymeexplorer

A7005-10KU	10000 units
A7005-25KU	25000 units
A7005-50KU	50000 units
A7005-100KU	100000 units

 β -Amylase ▲ **α -Amylase, heat-stable**

α -Amylase from *Bacillus licheniformis*; 1,4- α -D-Glucan-glucanohydrolase
[9000-85-5]

► solution, for use in Total Dietary Fiber Assay, TDF-100A

A3306-10ML	10 mL
-------------------	-------

Amylase, Maltogenic from *Bacillus* sp.

Maltogenic Amylase; Glucan 1,4- α -maltohydrolase
Novamyl 1000BG

A2986-10G	10 g
------------------	------

▼ Amyloglucosidase

Exo-1,4- α -glucosidase; 1,4- α -D-Glucan glucohydrolase;
Glucoamylase
[9032-08-0]

Amyloglucosidase from *Aspergillus niger*

► lyophilized, powder, activity: ~70 units/mg

For glycogen determination in whole yeast cells; Synthesis of hetero-oligosaccharides by glucoamylase in reverse.

One unit corresponds to the amount of enzyme which liberates 1 μ mole of glucose per minute at pH 4.8 and 60 °C (starch acc. to Zulkowsky, Cat. No. 85642, as substrate).

10115-1G-F	1 g
10115-5G-F	5 g

► lyophilized powder, activity: 30–60 units/mg protein (biuret)

Lyophilized powder containing less than 0.02% glucose

One unit will liberate 1.0 mg of glucose from starch in 3 min at pH 4.5 at 55 °C.

A7420-5MG	5 mg
A7420-25MG	25 mg
A7420-100MG	100 mg

► powder, activity: ~120 units/mg

One unit corresponds to the amount of enzyme which liberates 1 μ mol glucose per minute at pH 4.8 and 60 °C (starch, Fluka No. 85642, as substrate).

Synthesis of hetero-oligosaccharides by glucoamylase in reverse¹.

Lit cited: 1. R.A. Rastall et al., *Biotechnol. Lett.* **13**, 501 (1991)

10113-1G	1 g
10113-5G	5 g

Amyloglucosidase from *Rhizopus* sp.

► activity: ≥5,000 units/g solid

One unit will liberate 1.0 mg of glucose from starch in 3 min at pH 4.5 at 55 °C.

Using soluble starch, almost theoretical yields of glucose are obtained.

A7255	Inquire
--------------	---------

Amyloglucosidase ▲**▼ α -Glucosidase**

α -D-Glucosidase glucohydrolase
[9001-42-7]

Protein determined by biuret.

 α -Glucosidase from *Bacillus stearothermophilus*

Maltase

► Maltase lyophilized powder, activity: ≥50 units/mg protein

Lyophilized powder containing potassium phosphate buffer salt

One unit will liberate 1.0 μ mole of D-glucose from p-nitrophenyl- α -D-glucoside per min at pH 6.8 at 37 °C.

Activity using maltose as substrate at pH 6.0 at 25 deg C is ~2X greater than that obtained using p-nitrophenyl- α -D-glucoside as substrate at pH 6.8 at 37 °C.

G3651-250UN	250 units
--------------------	-----------

 α -Glucosidase from rice

Maltase

► ammonium sulfate suspension, activity:

40–80 units/mg protein

Suspension in 2.8 M $(\text{NH}_4)_2\text{SO}_4$ solution

One unit will convert 1.0 μ mole of maltose to 2.0 μ moles of D-glucose per min at pH 4.0 at 37 °C.

G9259-100UN	100 units
--------------------	-----------

α-Glucosidase from *Saccharomyces cerevisiae*

α-D-Glucosidase; Maltase from yeast Hydrolysis of terminal, non-reducing 1→4-linked D-glucose residues with release of D-glucose.

For the determination of α-amylase and the synthesis of various 1'-O-sucrose and 1-O-fructose esters

- recombinant, expressed in unspecified host, lyophilized powder, activity: ≥125 units/mg protein

Lyophilized powder containing potassium phosphate buffer salt pH 7.15 and approx. 70% lactose

One unit will liberate 1.0 μmole of D-glucose from p-nitrophenyl α-D-glucoside per min at pH 6.8 at 37 °C.

G0660-750UN 750 units

- lyophilized powder, activity: ≥10 units/mg protein (using p-nitrophenyl α-D-glucoside as substrate.)

Sold on basis of p-nitrophenyl α-D-glucoside units.

One unit will liberate 1.0 μmole of D-glucose from p-nitrophenyl α-D-glucoside per min at pH 6.8 at 37 °C.

View more information on enzymes for complex carbohydrate analysis at sigma-aldrich.com/enzymeexplorer

G5003-100UN	100 units
G5003-1KU	1000 units

α-Glucosidase ▲

Ethanol Analysis Reagents

ASTM™ D5501 Denatured Fuel Ethanol Standards Kit

Use this quantitative calibration standard kit to determine if ethanol and gasoline fuel blends comply with federal and state laws. A certificate of analysis accompanies each kit.

suitable for D5501 per ASTM

Components

Sol.#1 Ethanol:Heptane: Methanol (92%, 7.40%, 0.60%)
 Sol.#2 Ethanol:Heptane: Methanol (93%, 6.50%, 0.50%)
 Sol.#3 Ethanol:Heptane: Methanol (94%, 5.60%, 0.40%)
 Sol.#4 Ethanol:Heptane: Methanol (95%, 4.70%, 0.30%)
 Sol.#5 Ethanol:Heptane: Methanol (96%, 3.80%, 0.20%)
 Sol.#6 Ethanol:Heptane: Methanol (97%, 2.90%, 0.10%)
 Sol.#7 Ethanol:Heptane: Methanol (98%, 1.95%, 0.05%)

40361-U 1 pkg

Fuel Ethanol Residual Saccharides Mix

Quantitative standard for monitoring unfermented sugars present during the process of transforming corn feedstock into ethanol.

Components

Glycerol 1 % (w/v)
 D-(+)-Glucose 2 % (w/v)
 Maltotriose (DP₃) 1 % (w/v)
 Maltose monohydrate 2 % (w/v)
 L-(+)-Lactic acid .3 % (w/v)
 Acetic acid .3 % (w/v)
 Dextrin 3.25 % (w/v)
 Ethanol 12 % (w/v)

48468-U 10 × 2 mL

Ethanol standards 10% (v/v)


Ethanol solution; Ethyl alcohol [64-17-5] C₂H₅OH FW 46.07

CH₃CH₂OH

- aqueous solution

Flame-sealed flint ampule

E2385-10AMP 10 amp



Sigma-Aldrich® and GeneGo

Sigma-Aldrich and GeneGo have integrated Sigma-Aldrich product listings into the MetaCore™ and MetaDrug™ software suites.

This partnership associates Sigma-Aldrich products to over 700 interactive maps of Canonical pathways located in MetaCore and MetaDrug.

Users of MetaCore and MetaDrug are able to link to Sigma-Aldrich products from molecules displayed in

the Drug Detail Pages. This functionality will accelerate productivity for researchers by showing molecule availability during the information gathering process of experiments.

Learn more about GeneGo and their MetaCore and MetaDrug software suites by visiting www.genego.com



MetaCore and MetaDrug are trademarks of GeneGo, Inc. Patent Chemistry Database: Copyright © 2007-2008, Elsevier Inc., licensed to Elsevier Information Systems GmbH. PharmaPendium: Copyright © 2007-2008, Elsevier Inc., licensed to Elsevier Information Systems GmbH.

sigma-aldrich.com

SIGMA-ALDRICH®



Enzymes and Reagents for BioDiesel Research

Lipases for Transesterification

Lipases for Transesterification

Amano Lipase PS

► (immobilized on diatomite)

708011-10G 10 g

Amano Lipase PS, from *Burkholderia cepacia*

► activity: ≥30,000 U/g, pH 7.0, 50 °C (Optimum pH and temperature)

534641-10G 10 g
534641-50G 50 g

Amano Lipase PS-C II (immobilized on ceramic)

One unit will produce 1.0 micromole of 1-Phenethyl alcohol to 1-Phenethyl acetate per min. at 25 °C in the presence of vinyl acetate

534889 Inquire

Amano lipase PS-IM (immobilized on diatomaceous earth)

709603-10G 10 g

Amano Lipase PS-C I (immobilized on ceramic)

One unit will produce 1.0 micromole of 1-Phenethyl alcohol to 1-Phenethyl acetate per min. at 25 °C in the presence of vinyl acetate

534897 Inquire

Lipase from *Candida* sp.

► Novozymes® CALB L recombinant, expressed in *Aspergillus niger*

aqueous solution

Minimum 5,000 LU/G of liquid

L3170-50ML 50 mL

Lipase acrylic resin from *Candida antarctica*

Novozym® 435
[9001-62-1]

Immobilized preparation of a thermostable lipase, particularly useful in the synthesis of esters and amides, and has a broad substrate specificity.

Lipase from (B lipase) *Candida antarctica* produced by submerged fermentation of a genetically modified *Aspergillus oryzae* microorganism and adsorbed on a macroporous resin.

► activity: ≥10,000 U/g, recombinant, expressed in *Aspergillus niger*

A product of Novozyme Corp.

L4777-3G 3 g
L4777-10G 10 g

Phospholipase C for Degumming Oils

▼ Phospholipase C

Phosphatidylcholine cholinophosphohydrolase; Lecithinase C; Lipophosphodiesterase I; PC-PLC [9001-86-9]

Hydrolyzes the phosphate bond on phosphatidylcholine and other glycerophospholipids yielding diacylglycerol; this enzyme will also hydrolyze the phosphate bonds of sphingomyelin, cardiolipin, choline plasmalogen and ceramide phospholipids.

Phospholipase C from *Clostridium perfringens* (*C. welchii*)

One unit will liberate 1.0 μmole of water soluble organic phosphorus from egg yolk L-α-phosphatidylcholine per min at pH 7.3 at 37 °C.

► lyophilized powder, activity: 10–50 units/mg protein

P7633-25UN	25 units
P7633-125UN	125 units
P7633-500UN	500 units

► lyophilized powder, activity: ≥150 units/mg protein

Lyophilized powder in buffered salts
Chromatographically purified

P4039-125UN	125 units
P4039-500UN	500 units

Phospholipase C ▲


Phospholipase C from *Bacillus cereus*

► activity: ≥200 units/mg protein


Lyophilized powder containing approx. 10% protein. Remainder: trehalose, zinc sulfate, and potassium phosphate

One unit will liberate 1.0 μmole of water soluble organic phosphorus from L-α-phosphatidylcholine per min at pH 7.3 at 37 °C.

P6621-250UN 250 units



Transesterification of triglycerides utilizing commercial lipase preparations is currently under investigation as a low energy input alternative to current "chemical" methods.



Degumming of feedstock oils by hydrolysis of phospholipids can be catalyzed by phospholipase C. Treatment of phosphatidylethanol, phosphatidylcholine and phosphatidylinositol with phospholipase C enables the hydrolysis and phase separation of the hydrophilic phosphate-containing contaminants.

BioDiesel Analysis

GC Columns for BioDiesel Analysis

MET-Biodiesel


This rugged metal column was designed specifically for the determination of free and total glycerin in B100 biodiesel samples. Features and benefits of this column include:

- Metal is used as the column material, virtually eliminating accidental column breakage during handling
- Metal columns do not require a coating to provide strength (NOTE: the protective polyimide coating on the outside of fused silica columns slowly burns off at 380 °C, exposing bare fused silica which is fragile and susceptible to breakage)
- Includes a guard that both protects the analytical column from excess reagent and non-volatile compounds, extending column life, and acts as a retention gap, minimizing peak broadening
- The guard is an integrated guard, thereby providing protection with a leak-free connection (the guard and analytical column are one continuous piece of tubing; there is no union between the guard and analytical column)
- Low bleed characteristic, even at 380 °C
- Provides good peak shape and resolution for all glyceride impurities of interest
- Able to separate glycerin in addition to mono-/diglycerides (as methyl esters) plus triglycerides from the FAMEs
- A maximum temperature of 380 °C (isothermal) and 430 °C (programmed) exceeds the temperature limitations specified in biodiesel methods such as ASTM D6584 and EN 14105

GC capillary columns

Temp. Limits:


- -60 °C to 380 °C (isothermal)
- -60 °C to 430 °C (programmable)



Derivatization

N-Methyl-N-(trimethylsilyl)trifluoroacetamide

MSTFA; N-Trimethylsilyl-N-methyl trifluoroacetamide
[24589-78-4] CF₃CON(CH₃)Si(CH₃)₃
FW 199.25



Silylating agent used to form volatile derivatives for GC-MS analysis.²

Improved silylation procedure for simultaneous detection of steroid hormones 17- α -ethynylestradiol and estrone.¹

Lit cited: 1. *J. Chromatogr.* **1108**, 121 (2006); 2. *Phosphorus, Sulfur Silicon Relat. Elem.* **88**, 53 (1994).

394866-10X1ML 10 × 1 mL

394866-5ML 5 mL

394866-25ML 25 mL

BioDiesel Calibration Standards

Supelco glycerin impurities in biodiesel reference solutions meet the specifications called out in ASTM D6584 and DIN EN 14105. Our reference solutions have been carefully formulated to save 4–6 hours in prep time, minimize the possibility of human error, and eliminate the need to prepare pyridine dilutions. Our standards offer:

- Raw materials tested for identity and purity
- Stock solutions for preparing your own multi-component solutions
- Multilevel, multi-component solutions ready for derivatization
- Monoglyceride commercial mixture for EN14105
- Internal standard solutions
- Certificate of composition shipped with each standard
- Instructions for derivatization with each kit

Our standards are manufactured under the Quality Management System which is registered under ISO 9001:2000.

Biodiesel Calibration Standards

Name	Suitability	Composition	Cat. No.
ASTM® D6584 Individual Standard Solution and Internal Standards Kit	D6584 ASTM	ASTM® D6584 Standard Solution 1 (Supelco 44899-U), 1 mL ASTM® D6584 Standard Solution 2 (Supelco 44914-U), 1 mL ASTM® D6584 Standard Solution 3 (Supelco 44915-U), 1 mL ASTM® D6584 Standard Solution 4 (Supelco 44916-U), 1 mL ASTM® D6584 Standard Solution 5 (Supelco 44917-U), 1 mL ASTM D6584 1,2,4-Butanetriol Solution, Internal Standard #1 (Supelco 44896-U), 5 mL ASTM D6584 Tricaprin Solution, Internal Standard #2 (Supelco 44897-U), 5 mL	44918-U
ASTM® D6584 Standard Solution 1	D6584 ASTM	1,3-Diolein, 50 µg/mL Glycerol, 5 µg/mL Monoolein, 100 µg/mL Glyceryl trioleate, 50 µg/mL	44899-U
ASTM® D6584 Standard Solution 2	D6584 ASTM	1,3-Diolein, 100 µg/mL Glycerol, 15 µg/mL Monoolein, 250 µg/mL Glyceryl trioleate, 100 µg/mL	44914-U
ASTM® D6584 Standard Solution 3	D6584 ASTM	1,3-Diolein, 200 µg/mL Glycerol, 25 µg/mL Monoolein, 500 µg/mL Glyceryl trioleate, 200 µg/mL	44915-U
ASTM® D6584 Standard Solution 4	D6584 ASTM	1,3-Diolein, 350 µg/mL Glycerol, 35 µg/mL Monoolein, 750 µg/mL Glyceryl trioleate, 350 µg/mL	44916-U
ASTM® D6584 Standard Solution 5	D6584 ASTM	1,3-Diolein, 500 µg/mL Glycerol, 50 µg/mL Monoolein, 1000 µg/mL Glyceryl trioleate, 500 µg/mL	44917-U
EN 14105:2003 Monoglyceride Stock Solution	EN 14105 DIN	Monoolein Monopalmitin Monostearin	49446-U
EN 14105:2003 Standard Solution 1	EN 14105 DIN	Butanetriol, 80 µg/mL 1,3-Diolein, 50 µg/mL Glycerol, 5 µg/mL Monoolein, 250 µg/mL Tricaprin, 800 µg/mL Triolein, 50 µg/mL	49441-U
EN 14105:2003 Standard Solution 2	EN 14105 DIN	Butanetriol, 80 µg/mL Diolein, 200 µg/mL Glycerol, 20 µg/mL Monoolein, 600 µg/mL Tricaprin, 800 µg/mL Triolein, 150 µg/mL	49442-U
EN 14105:2003 Standard Solution 3	EN 14105 DIN	Butanetriol, 80 µg/mL Diolein, 350 µg/mL Glycerol, 35 µg/mL Monoolein, 950 µg/mL Tricaprin, 800 µg/mL Triolein, 300 µg/mL	49443-U
EN 14105:2003 Standard Solution 4	EN 14105 DIN	Butanetriol, 80 µg/mL Diolein, 500 µg/mL Glycerol, 50 µg/mL Monoolein, 1250 µg/mL Tricaprin, 800 µg/mL Triolein, 400 µg/mL	49444-U
EN 14105:2003 Standard Solution Kit	EN 14105 DIN	EN 14105:2003 Standard Solution 1 (Supelco 49441-U), 1 mL EN 14105:2003 Standard Solution 2 (Supelco 49442-U), 1 mL EN 14105:2003 Standard Solution 3 (Supelco 49443-U), 1 mL EN 14105:2003 Standard Solution 4 (Supelco 49444-U), 1 mL	49445-U

Karl Fischer Titration—HYDRANAL® Reagents for the Determination of Water in Biodiesel

The presence of water in biodiesel poses problems for a number of reasons. Some of the water is residual from processing while some comes from condensation in

the storage tanks. Water content can be determined accurately and reproducibly by Karl Fischer Titration (volumetric, coulometric, and KF-oven methods) using reliable, superior, and proven HYDRANAL reagents (EN 14214).

One-Component Volumetric Reagents

Biodiesel contains components of different chain lengths. Therefore, a solubilizer is often required to dissolve or disperse it in order to extract the moisture. A number of HYDRANAL media are available for effective dispersion of biodiesel. HYDRANAL-Composite 2 is used as the titration agent.

One Component Reagents

Name	Grade	Cat. No.
HYDRANAL®-Composite 1	one-component reagent for volumetric Karl Fischer titration	34827-500ML-R 34827-6X500ML-R 34827-1L-R 34827-6X1L-R
HYDRANAL®-Composite 2	one-component reagent for volumetric Karl Fischer titration	34806-500ML 34806-6X500ML 34806-1L 34806-6X1L 34806-2.5L 34806-4X2.5L
HYDRANAL®-Composite 5	one-component reagent for volumetric Karl Fischer titration	34805-500ML-R 34805-6X500ML-R 34805-1L-R 34805-6X1L-R 34805-2.5L-R 34805-4X2.5L-R
HYDRANAL®-Composite 5 K	one-component reagent for volumetric Karl Fischer titration in ketones and aldehydes	34816-500ML-R 34816-6X500ML-R 34816-1L-R 34816-6X1L-R 34816-2.5L-R 34816-4X2.5L-R
HYDRANAL®-CompoSolver E	–	34734-1L-R 34734-6X1L-R 34734-2.5L-R 34734-4X2.5L-R
HYDRANAL®-KetoSolver	–	34738-500ML-R 34738-6X500ML-R 34738-1L-R 34738-6X1L-R
HYDRANAL®-LipoSolver CM	–	37855-1L 37855-6X1L
HYDRANAL®-LipoSolver MH	–	37856-1L 37856-6X1L
HYDRANAL®-Medium K	–	34698-1L-R 34698-6X1L-R
HYDRANAL®-Methanol dry	–	34741-1L-R 34741-6X1L-R 34741-2.5L-R 34741-4X2.5L-R
HYDRANAL®-Methanol Rapid	–	37817-1L-R 37817-6X1L-R 37817-2.5L-R 37817-4X2.5L-R
HYDRANAL®-Solver (Crude) oil	–	34697-1L-R 34697-6X1L-R 34697-2.5L-R 34697-4X2.5L-R
HYDRANAL®-Working Medium K	–	34817-1L 34817-6X1L
Karl-Fischer reagent	for titrimetric determination of water with one solution	36115-1L 36115-6X1L

Two-Component Volumetric Reagents

Name	Grade	Cat. No.
HYDRANAL®-Solvent	–	34800-1L-R 34800-6X1L-R 34800-2.5L-R 34800-4X2.5L-R
HYDRANAL®-Solvent E	–	34730-500ML-R 34730-6X500ML-R 34730-1L-R 34730-6X1L-R 34730-2.5L-R 34730-4X2.5L-R
HYDRANAL®-Solvent CM	–	34812-1L 34812-6X1L 34812-2.5L 34812-4X2.5L
HYDRANAL®-Solvent Oil	–	34749-1L 34749-6X1L
HYDRANAL®-Titrant 2	–	34811-500ML 34811-6X500ML 34811-1L 34811-6X1L
HYDRANAL®-Titrant 5	–	34801-500ML 34801-6X500ML 34801-1L 34801-6X1L 34801-2.5L 34801-4X2.5L
HYDRANAL®-Titrant 2 E	–	34723-1L-R 34723-6X1L-R
HYDRANAL®-Titrant 5 E	–	34732-100ML-R 34732-500ML-R 34732-6X500ML-R 34732-1L-R 34732-6X1L-R 34732-2.5L-R 34732-4X2.5L-R
Karl-Fischer reagent	for titrimetric determination of water with two separate solutions	36116-1L 36116-6X1L
Karl-Fischer reagent	for titrimetric determination of water with two separate solutions	36117-1L 36117-6X1L

Coulometric Titration Reagents

The coulometric procedure is more sensitive than the volumetric titration method. Determination of water is carried out in a coulometric cell (with or without diaphragm). The precision of the coulometric cell and the sample manipulation can be tested by means of HYDRANAL Water Standard 0.10.

Name	Grade	Cat. No.
HYDRANAL®-Coulomat A	anolyte for coulometric Karl Fischer titration (for cells with diaphragm)	34807-500ML 34807-6X500ML
HYDRANAL®-Coulomat E	–	34726-500ML-R 34726-6X500ML-R
HYDRANAL®-Coulomat AD	reagent for coulometric Karl Fischer titration for cells without diaphragm	34810-500ML-R 34810-6X500ML-R
HYDRANAL®-Coulomat AF 7	anolyte for coulometric Karl Fischer titration (suitable for coulometer AF7)	34829-1L-R 34829-6X1L-R
HYDRANAL®-Coulomat AG	anolyte for coulometric Karl Fischer titration. Suitable for cells with and without diaphragm	34836-500ML-R 34836-6X500ML-R 34836-1L-R 34836-6X1L-R
HYDRANAL®-Coulomat AG-H	anolyte for coulometric Karl Fischer titration in long-chained hydrocarbons (free of halogenated hydrocarbons)	34843-500ML 34843-6X500ML
HYDRANAL®-Coulomat AG-Oven	–	34739-500ML-R 34739-6X500ML-R
HYDRANAL®-Coulomat AK	anolyte for coulometric Karl Fischer titration in ketones (for cells with and without diaphragm)	34820-500ML-R 34820-6X500ML-R
HYDRANAL®-Coulomat CG	catholyte for coulometric Karl Fischer titration. Free of halogenated hydrocarbons	34840-50ML-R 34840-6X50ML-R
HYDRANAL®-Coulomat CG-K	–	34821-50ML-R 34821-6X50ML-R
HYDRANAL®-Coulomat Oil	–	34868-100ML-R 34868-6X100ML-R 34868-500ML-R 34868-6X500ML-R


Karl Fischer Oven products

This method is suitable for samples that release water at a higher temperature. In this method, HYDRANAL-Coulomat CG is placed in the cathode chamber and HYDRANAL Coulomat AG Oven is added to the anode chamber of a coulometric cell with diaphragm. A cell without a diaphragm requires only the Coulomat AG Oven. HYDRANAL Coulomat AG-H or HYDRANAL Coulomat AD may be used as a replacement for HYDRANAL Coulomat AG Oven. HYDRANAL Molecular Sieve 0.3 mm is a drying medium for the carrier gas.

Name	Grade	Cat. No.
HYDRANAL®-Coulomat AG-Oven	–	34739-500ML-R 34739-6X500ML-R
HYDRANAL®-Coulomat CG	catholyte for coulometric Karl Fischer titration. Free of halogenated hydrocarbons	34840-50ML-R 34840-6X50ML-R
HYDRANAL®-Molecular sieve 0.3 nm	–	34241-250G-R 34241-1KG-R
HYDRANAL®-Water Standard KF-Oven, 140–160 °C	–	34693-10G-R
HYDRANAL®-Water Standard KF-Oven 220 °C–230 °C	–	34748-10G-R 34748-6X10G-R




Metabolite Standards for Mevalonate and Isoprenoid Pathway Analysis



Name	Cat. No.
Acetoacetyl coenzyme A sodium salt hydrate	A1625-5MG A1625-10MG A1625-25MG
Acetyl coenzyme A sodium salt, powder	A2056-1MG A2056-5MG A2056-10MG
Acetyl coenzyme A trilithium salt	A2181-1MG A2181-5MG A2181-10MG
Butyryl coenzyme A lithium salt hydrate	B1508-5MG B1508-10MG B1508-25MG
1-Deoxy-D-xylulose-5-phosphate sodium salt	13368-1MG 13368-5MG
γ,γ -Dimethylallyl pyrophosphate triammonium salt	D4287-1VL D4287-5VL
Farnesyl monophosphate ammonium salt solution, methanol:ammonia solution	F1803-5VL
Farnesyl pyrophosphate ammonium salt solution, methanol:ammonia solution	F6892-1VL F6892-5VL
Geranylgeranyl monophosphate ammonium salt solution, methanol:ammonia solution	G2543-5VL
Geranylgeranyl pyrophosphate ammonium salt solution, methanol:ammonia solution	G6025-1VL G6025-5VL
Geranyl monophosphate ammonium salt solution, methanol:ammonia solution	G2293-1VL G2293-5VL
Geranyl pyrophosphate ammonium salt	G6772-1VL G6772-5VL
DL-Glyceraldehyde 3-phosphate solution	G5251-25MG G5251-100MG
DL-3-Hydroxy-3-methylglutaryl coenzyme A sodium salt	H6132-5MG H6132-10MG
Isopentenyl monophosphate ammonium salt solution	I1157
Isopentenyl pyrophosphate triammonium salt solution	I0503-1VL I0503-5VL
(R)-Mevalonic acid sodium salt	41288
(\pm)-Mevalonic acid 5-phosphate trilithium salt hydrate	79849-10MG 79849-50MG
(R)-Mevalonic acid 5-pyrophosphate tetralithium salt	77631
(\pm)-Mevalonic acid 5-pyrophosphate tetralithium salt	94259-10MG 94259-50MG
(R)-(–)-Mevalonolactone	68519-100MG
(\pm)-Mevalonolactone	M4667-1G
Phospho(enol)pyruvic acid monopotassium salt	P7127-100MG
Phospho(enol)pyruvic acid monosodium salt hydrate	P0564-100MG
Phospho(enol)pyruvic acid trisodium salt hydrate	P7002-100MG
n-Propionyl coenzyme A lithium salt	P5397-5MG
Pyruvic acid	107360-25G
Sodium pyruvate	15990-25G


Renewable and Alternative Energy Web Portal



The Portal offers a comprehensive set of links about environmentally benign production, storage and conversion of energy.

- Solar Energy
- Hydrogen Economy
- Batteries
- Supercapacitors
- Alternative Fuels
- Relevant events and conferences are also highlighted.

We focus on materials, so you can focus on results. For more information, visit sigma-aldrich.com/renewable



Have you seen these
Material Matters™ issues
on Alternative Energy?
Get your complimentary
subscription today:
sigma-aldrich.com/mm

Sigma-Aldrich® Worldwide Offices

Argentina

Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

Chile

Tel: (+56) 2 495 7395
Fax: (+56) 2 495 7396

China

Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

Tel: (+45) 43 56 59 00
Fax: (+45) 43 56 59 05

Finland

Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 88
Fax: (+33) 474 95 68 08

Germany

Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Hungary

Ingyenes telefonszám: 06 80 355 355
Ingyenes fax szám: 06 80 344 344
Tel: (+36) 1 235 9063
Fax: (+36) 1 269 6470

India

Telephone

Bangalore: (+91) 80 6621 9400
New Delhi: (+91) 11 4358 8000
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 4015 5488
Kolkata: (+91) 33 4013 8003

Fax

Bangalore: (+91) 80 6621 9550
New Delhi: (+91) 11 4358 8001
Mumbai: (+91) 22 4087 2364
Hyderabad: (+91) 40 4015 5488
Kolkata: (+91) 33 4013 8000

Ireland

Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 402 20370
Fax: (+353) 402 20375

Israel

Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

Free Tel: 800 827 018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: (+52) 722 276 1600
Fax: (+52) 722 276 1601

The Netherlands

Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

Tel: (+47) 23 17 60 00
Fax: (+47) 23 17 60 10

Poland

Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

Tel: (+7) 495 621 5828
Fax: (+7) 495 621 6037

Singapore

Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

Slovakia

Tel: (+421) 255 571 562
Fax: (+421) 255 571 564

South Africa

Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313

United States

Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Vietnam

Tel: (+84) 3516 2810
Fax: (+84) 6258 4238

Internet

sigma-aldrich.com



*Accelerating Customers'
Success through Innovation and
Leadership in Life Science,
High Technology and Service*

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052
Technical Service (800) 325-5832 • sigma-aldrich.com/techservice
Development/Custom Manufacturing Inquiries **SAFC® (800) 244-1173**
Safety-related Information sigma-aldrich.com/safetycenter

World Headquarters
3050 Spruce St.
St. Louis, MO 63103
(314) 771-5765
sigma-aldrich.com

©2010 Sigma-Aldrich Co. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH, FLUKA, and SUPELCO are trademarks belonging to Sigma-Aldrich Co. and its affiliate Sigma-Aldrich Biotechnology, L.P. Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip. ASTM is a registered trademark of American Society for Testing and Materials. Biomek is a registered trademark of Beckman Instruments, Inc. HypoThermosil is a registered trademark of BioLife Solutions, Inc. CryoStar is a trademark of BioLife Solutions, Inc. CellBIND, Corning and HYPER Flask are registered trademarks of Corning, Inc. TWEEN is a registered trademark of Croda International PLC. Triton is a registered trademark of Dow Chemical Co. Nafion is a registered trademark of E.I. du Pont de Nemours & Co., Inc. Eppendorf is a registered trademark of Eppendorf-Netheler-Hinz GmbH. Cy and Sepharose are registered trademarks of GE Healthcare. nanohA, Select-HA are trademarks of Hyalose LLC. Coomassie is a registered trademark of Imperial Chemical Industries Ltd. Legato is a trademark of KDScientific. Diselase is a registered trademark of Kyowa Hakko Kogyo Co. Ltd. Alexa, SYBR and Texas Red are registered trademarks of Novozymes Corp. BAN is a trademark of Novozymes Corp. Shimadzu is a registered trademark of Shimadzu Corporation. Aldrich, HYDRANAL and Prestige Antibodies are registered trademarks of Sigma-Aldrich Biotechnology LP and Sigma-Aldrich Co. and Ampliflu, EZBlue, GlycoProfile, HumanKine, Material Matters, ProteoMass, ProteoSilver and Sigma are trademarks of Sigma-Aldrich Biotechnology LP and Sigma-Aldrich Co.

MPC
73457-510295
1050