

Millistak+® HC Pod depth filter system

High performance modular depth filters for primary and secondary clarification operations in biopharmaceutical production

Millistak+® HC depth filter media is offered in a scalable, disposable format: the Pod filter system. Accommodating applications from lab to pilot to process-scale, the Pod format offers greater flexibility because of its modular and 100% disposable design.

The Millistak+® Pod system is ideal for a wide variety of primary and secondary clarification applications, including mammalian cell cultures, yeast and *E. coli* lysates post centrifuge, *E. coli* refolds, media, vaccines, plasma proteins and sera.

Millistak+® Pod filters are available in three distinct series of media grades in order to meet your specific application needs.

Millistak+® DE, CE and HC (high capacity) media deliver optimal performance through a gradient density matrix, as well as positive surface charge properties.



Benefits

- Low hold-up volume for greater product yield.
- Broad range of media types offered in single and multi-layer products.
- Millistak+® HC double layer media improves prefiltration and compresses clarification.
- Flexible, modular format offers scalability up to 20,000 L bioreactor volumes.
- Patented disposable design eliminates need for housing, CIP or cleaning validation.
- Self-contained Pod filters protect operators from exposure to biohazards.
- Robust construction is easy to use and set up.
- Smaller footprint facilitates use in tight spaces.
- 6 cm² and 20 cm² small surface area devices for screening and scalability studies for low feedstream volume requirements.
- Available in regionalized Optimized Packaging.

Easy to Use

With the compact, modular design of our Millistak+® Pod system, you can increase productivity and shorten cycle times from bench to manufacturing scale.

Installation and setup of the Pod system is simple and straightforward. The unique design of the disposable adapters and disposable manifolds makes it easy to connect the Pod filters to the rest of the unit operations in the process. The self-contained and disposable system eliminates maintenance as well as cleaning validation requirements.

Configurations

- NanoPod NP6 filters (X0 grades) – 6 cm².
- Micro 20 filters (C0/D0/X0 grades) – 20 cm².
- µPod® filters (CE/DE, A0/A1/B1/F0 grades) – 23 cm².
- Lab-scale Pod filters – 0.027 m² and 0.054 m².
- Process-scale Pod filters:
 - Millistak+® DE and CE media – 0.11 m², 0.77 m² and 1.4 m².
 - Millistak+® HC media – 0.11 m², 0.55 m² and 1.1 m².
- Pilot-scale Pod holder – accommodates up to two process-scale Pod filters for configurations from 0.11 m² to 2.8 m², depending on media type. An optional accessory kit expands capacity to five process-scale Pod filters.
- Process-scale Pod holder – accepts from five to ten process-scale Pod filters per rack. Up to three racks can be stacked for process flexibility.

- Disposable adapter – connects Pod filters to process piping, creating a disposable flow path.
- Disposable diverter plate – enables two media grades to be processed on a single rack.

Millistak+® Depth Filter Media


Millistak+® Pod filters incorporate multiple graded-density layers and adsorptive, positively charged filter media. Composed of select grade cellulose fiber and diatomaceous earth, the Millistak+® DE series not only improves the manufacturing process, but also increases contaminant holding. The Millistak+® CE series consists of single-layer media with cellulose fibers that are suitable for coarse filtration applications.

The Millistak+® HC series improves productivity by combining multiple media grades into one device, enabling compression of multiple filtration stages downstream of the bioreactor.

Millistak+® Pod Depth Filters are Supported by the Emprove® Program – The Smart Way to Master Compliance and Control

Complementing our product portfolio, the Emprove® Program provides convenient access to reliable technical and regulatory information in Emprove® Dossiers to support your risk assessment continuum. A subscription to our Emprove® Suite can help you stay current: In addition to accessing the Emprove® Dossiers, you can also receive notification updates to document changes, as well as generate metrics and reports.

For more information, please visit <https://sigmaaldrich.com/emprove>



Typical Extractables

	Millistak+® Single-Layer Pods (CE & DE Media)	Millistak+® HC Media (C0HC, D0HC, A1HC, A0HC, B1HC)	Millistak+® X0HC/F0HC Media
Conductivity	1.52–1.94 $\mu\text{S}/\text{cm}$ post autoclave (1 cycle of 60 minutes at 123 °C) and pure water flush of 5 L/ft ² (50 L/m ²) of media surface area	3.64–10.5 $\mu\text{S}/\text{cm}$ post autoclave (1 cycle of 60 minutes at 123 °C) and pure water flush of 10 L/ft ² (100 L/m ²) of media surface area	19.34–53.2 $\mu\text{S}/\text{cm}$ (X0HC) and 21.4–42.7 $\mu\text{S}/\text{cm}$ (F0HC) post autoclave (1 cycle of 60 minutes at 123 °C) and pure water flush of 10 L/ft ² (100 L/m ²) of media surface area
NVR Gravimetric Extractables	Not Tested	420–750 mg/ft ² (process-scale) and 630–1251 mg/m ² (lab-scale) per 24 hour static soak in pure water (type 1 DI water) post autoclave (1 cycle of 60 minutes at 123 °C) and pure water flush of 10 L/ft ² (100 L/m ²) of media surface area	Not Tested
TOC	910–1800 ppb post autoclave (1 cycle of 60 minutes at 123 °C) and pure water flush of 5 L/ft ² (50 L/m ²) of media surface area	720–4600 ppb post autoclave (1 cycle of 60 minutes at 123 °C) and pure water flush of 10 L/ft ² (100 L/m ²) of media surface area	1200–2800 ppb (X0HC) and 460–3200 ppb (F0HC) post autoclave (1 cycle of 60 minutes at 123 °C) and pure water flush of 10 L/ft ² (100 L/m ²) of media surface area
Metals	Per 24-hour static soak in 10 liters pure water post autoclave (1 cycle of 30 minutes at 123 °C) and pure water flush of 10 liters per ft ² (100 L/m ²) of surface area. Values based on worse case data from both process-scale and lab-scale pod devices.		



Nanopod NP6 filter:

HC (X0): 6 cm²




Micro 20:

HC (C0/D0/X0): 20 cm²



μPod® filter:

HC (A0/A1/B1/F0),
DE, CE: 23 cm²



Millistak+® Screening and Scaling devices specifications

	NanoPod NP6	Micro 20	μPod® Filter
Surface Area	6 cm ²	20 cm ²	23 cm ²
Media Grade	X0 media	C0, D0, X0 media	A0/A1/B1/F0, CE and DE media
Materials of Construction	Cellulose fibers with inorganic filter aid (CE Media contains cellulose only)		
Filter Media			
Filter Membrane	N/A	N/A	Mixed esters of cellulose (grades A1HC and B1HC only)
Pod Housings	Polypropylene		Glass-Filled Polypropylene
Inlet, Vent and Outlet Connections	Inlet: Female Luer-lock Outlet: Male Luer slip	Female Luer-lock	Female Luer-lock
Pod Dimensions			
Diameter	3.6 cm (1.4 in.)	6.4 cm (2.5 in.)	
Length			8.9 cm (3.5 in.)
Height			6.6 cm (2.6 in.)
Thickness	3.6 cm (1.4 in.)	4.3 cm (1.7 in.)	4.1 cm (1.6 in.)
Maximum Operating Pressure	30 psig (2.1 bar) at <40 °C	30 psig (2.1 bar) at <40 °C	50 psig (3.4 bar) at ≤40 °C
Maximum Differential Pressure			
Forward	30 psid (2.1 bar) at 25 °C		30 psid (2.1 bar) at 40 °C
Reverse	30 psid (2.1 bar) at 25 °C		15 psid (1.0 bar) at 40 °C
Sterilization	Integrity is maintained after 2 autoclave cycles of 60 minutes at 123 °C.		
Indirect Food Additive	All component materials meet the FDA Indirect Food Additive Requirements cited in 21 CFR 177-182.		
Biological Reactivity	All component materials meet the criteria for Biological Reactivity Testing. These tests can be any or a combination of the following test methods: USP <88> Class VI (<i>in vivo</i>), USP <87> (<i>in vitro</i>), ISO 10993-5 (<i>in vitro</i>).		
Bacterial Endotoxin	An aqueous extraction contained less than 0.25 EU/mL as determined using the Limulus Amoebocyte Lysate (LAL) kinetic turbidimetric test or the gel-clot test as a secondary method (on filter media only), aligned with USP <85>, Ph. Eur. 2.6.14, and JP 4.01.		



Millistak+® Lab devices and Process Scale Pod filters specifications

	Lab-scale Pod Filter		Multi-layer and Single-Layer Process-scale Pod Filter	Multi-layer Process-scale Pod Filter	Single-Layer Process-scale Pod Filter	Multi-layer Process-scale Pod Filter	Single-Layer Process-scale Pod Filter
Surface Area	0.027 m² (0.29 ft²)	0.054 m² (0.58 ft²)	0.11 m² (1.2 ft²)	0.55 m² (5.9 ft²)	0.77 m² (8.3 ft²)	1.1 m² (11.8 ft²)	1.4 m² (15.1 ft²)
Media Grade	HC, CE and DE media	HC, CE and DE media	HC, CE and DE media	HC media	CE and DE media	HC media	CE and DE media
Materials of Construction							
Filter Media	Cellulose fibers with inorganic filter aid (CE Media contains cellulose only)						
Filter Membrane	Mixed esters of cellulose (grades A1HC and B1HC only)						
Pod Housings	Glass-Filled Polypropylene						
Adapters	Glass-Filled Polypropylene*						
Gaskets and Plugs	Thermo Plastic Elastomer (TPE)*						
Inlet, Vent and Outlet Connections	Hose Barb		Flat seal				
Pod Dimensions							
Length	8.5 in. (22 cm)		24.2 in. (62 cm)				
Height	5.3 in. (14 cm)		12.5 in. (32 cm)				
Thickness	2.9 in. (7.4 cm)	3.7 in. (9.4 cm)	1.2 in. (3.0 cm)	2.8 in. (7.1 cm)	3.1 in. (7.9 cm)	4.8 in. (12.2 cm)	5.0 in. (12.7 cm)
Maximum Operating Pressure	30 psig (2.1 bar) at 25 °C		50 psig (3.5 bar) at 25 °C; 15 psig (1.0 bar) at 80 °C				
Maximum Differential Pressure							
Forward	30 psid (2.1 bar) at 4 °C; 30 psid (2.1 bar) at 37 °C		30 psid (2.1 bar) at 25 °C; 15 psid (1.0 bar) at 80 °C				
Reverse	30 psid (2.1 bar) at 37 °C		30 psid (2.1 bar) at 25 °C				
Sterilization	Integrity is maintained after 2 autoclave cycles of 60 minutes at 123 °C		Integrity is maintained after 1 autoclave cycle of 60 minutes at 123 °C				
Indirect Food Additive	All component materials meet the FDA Indirect Food Additive Requirements cited in 21 CFR 177-182.						
Biological Reactivity	All component materials meet the criteria for Biological Reactivity Testing. These tests can be any or a combination of the following test methods: USP <88> Class VI (<i>in vivo</i>), USP <87> (<i>in vitro</i>), ISO 10993-5 (<i>in vitro</i>).						
Bacterial Endotoxin	For Media Only: <0.25 EU/mL as determined by the Limulus Amoebocyte Lysate (LAL) test, aligned with USP <85>, Ph. Eur. 2.6.14 and JP 4.01.						
Pressure Equipment Directive 2014/68/EU	Devices and associated holders are designed and manufactured in accordance with the sound engineering practices (SEP) cited in Article 4(3) of 2014/68/EU.*						


* Process scale only

Select Millistak+® filtration media contain up to 30% pre-consumer fiber.

Conductivity/TOC



Water Permeability



Water Flow Rates

Millistak+® HC Filter Series		Millistak+® CE Filter Series		Millistak+® DE Filter Series	
Media Type/ Grade	Water Flow Rate L/min/m ² at 10 psid, 21 °C	Media Type/ Grade	Water Flow Rate L/min/m ² at 10 psid, 21 °C	Media Type/ Grade	Water Flow Rate L/min/m ² at 10 psid, 21 °C
D0HC	CE25 1425.8–2619.1	CE15	3054.9–5611.5	DE25	1425.8–2619.1
	DE40 454.6–835.1	CE20	2082.6–3825.5	DE30	974.1–1789.3
C0HC	DE30 974.1–1789.3	CE25	1425.8–2619.1	DE35	665.5–1222.4
	DE60 99.0–181.9	CE30	974.1–1789.3	DE40	454.6–835.1
B1HC	DE50 212.2–389.8	CE35	665.5–1222.4	DE45	310.6–570.5
	DE75 31.6–58.0	CE40	454.6–835.1	DE50	212.2–389.8
A0HC	DE60 99.0–181.9	CE45	310.6–570.5	DE55	145.0–266.3
	DE75 31.6–58.0	CE50	212.2–389.8	DE60	99.0–181.9
F0HC	DE60 99.0–181.9			DE65	67.7–124.3
	IM75 39.0–54.6			DE70	46.2–84.9
A1HC	DE60 99.0–181.9			DE75	31.6–58.0
	DE75 31.6–58.0				
X0HC	IM75 39.0–54.6				
	IM83 21.2–27.7				


Choose the Right Media

Media Grade	Application	Characteristics	Media Construction
Single-layer CE	Primary (coarse) clarification	Cellulose fibers	CE15 to 50
Single-layer DE	Primary or secondary clarification	Cellulose fibers + inorganic filter aid	DE25 to 75
Double-layer D0HC	Primary clarification directly out of the bioreactor	A more open CE layer and DE media combination	CE25 + DE40
Double-layer C0HC	Perfusion bioreactor fluid	Two layers of a more open DE media	DE30 + DE60
Double-layer A0HC	Centrate clarification	Two layers of a tighter DE media	DE60 + DE75
Double-layer F0HC	Secondary clarification of pretreated harvest by acid precipitation or flocculation, <i>E. coli</i> and yeast	Two DE layers. Provides sterile filter protection without an RW01 membrane	DE60 + IM75
Double-layer X0HC	Secondary clarification of bioreactor harvests, primarily for cell cultures	Two DE layers. Provides sterile filter protection without an RW01 membrane	IM75 + IM83
Triple-layer A1HC	Post-TFF (Prostak™ system) clarification fluids and primary or secondary clarification	Tight media combination with an additional membrane layer to protect downstream membrane filters	DE60 + DE75 + RW01
Triple-layer B1HC	Post-centrifuge or settled permeate containing cellular particulate and primary or secondary clarification	A more open first layer with an additional membrane layer to protect downstream membrane filters	DE50 + DE75 + RW01


Note: For clarification of serum, plasma, vaccines, cell culture or other fluids, choice of media grade should be based on small-scale trials. Reference our Clarification Portfolio Guide (PB1241EN00) for additional information.

Ordering Information


Nanopod NP6



Micro 20



μPod[®] Filter



μPod[®] Filter Accessories

μPod[®] Tubing Kit Catalogue No. MTUBEKITL1

Gauge 0–60 psi and Connection Fittings Catalogue No. XXPXLGAGE

Ordering Information (continued)

Lab-scale Pod Filter

M							H	1
Millistak+® Filter	Media Type and Grade					Size	Connection	Package Quantity
	A1HC	F0HC	CE30	CE50	DE40	DE60	H = ¼ in. (6 mm)	1 = 1/pack
	B1HC	X0HC	CE35	DE25	DE45	DE65	Hose Barb	
	A0HC	CE15	CE40	DE30	DE50	DE70		
	C0HC	CE20	CE45	DE35	DE55	DE75		
	D0HC	CE25						

Multi-layer Process-scale Pod Filter

M							F	S	1
Millistak+® Filter	Media Type and Grade					Size	Connection		Package Quantity
	A1HC	C0HC	X0HC			01 = 0.11 m ²	FS = Flat Seal		1 = 1/pack
	B1HC	D0HC	F0HC			05 = 0.55 m ²			
	A0HC					10 = 1.1 m ²			

Single-Layer Process-scale Pod Filter

M							F	S	1
Millistak+® Filter	Media Type and Grade					Size	Connection		Package Quantity
	CE15	CE35	DE25	DE45	DE65	01 = 0.11 m ²	FS = Flat Seal		1 = 1/pack
	CE20	CE40	DE30	DE50	DE70	07 = 0.77 m ²			
	CE25	CE45	DE35	DE55	DE75	13 = 1.4 m ²			
	CE30	CE50	DE40	DE60					

Optimized Packaging

To improve sustainability of the packaging and shipping of filter products, optimized pack solutions for several grades of Millistak+® Process Scale Pod filters have been developed to optimize transport and receiving processes as well as to reduce waste.

Optimized packaging configurations include different quantities of pod filters distributed in 3 boxes on one regional standardized pallet size.

The key benefits are:

- 42% average reduction in corrugated packaging material per product.
- 29% decrease in the number of pallets further reducing energy use and emissions.
- 75% reduction in operator time to open and manage the product and packaging.



Optimized Packaging – North American format

Description	Item Description	Qty/Pk	Item Number
US Pallet format HDPE pallet 40 in. × 48 in. × 5.9 in.	Millistak+® HC Pod A0HC 1.1 m ²	30	MA0HC10FS30
	Millistak+® HC Pod A1HC 1.1 m ²	30	MA1HC10FS30
	Millistak+® HC Pod B1HC 1.1 m ²	30	MB1HC10FS30
	Millistak+® HC Pod C0HC 1.1 m ²	30	MC0HC10FS30
	Millistak+® HC Pod D0HC 1.1 m ²	30	MD0HC10FS30
	Millistak+® HC Pod F0HC 1.1 m ²	30	MF0HC10FS30
	Millistak+® HC Pod X0HC 1.1 m ²	30	MX0HC10FS30

Optimized Packaging – European format

Description	Item Description	Qty/Pk	Item Number
Euro Pallet format EU Modified Heat treated (EU) Wooden CP2 1200 mm × 800 mm	Millistak+® HC Pod A0HC 1.1 m ²	27	MA0HC10FS27EU
	Millistak+® HC Pod A1HC 1.1 m ²	27	MA1HC10FS27EU
	Millistak+® HC Pod B1HC 1.1 m ²	27	MB1HC10FS27EU
	Millistak+® HC Pod C0HC 1.1 m ²	27	MC0HC10FS27EU
	Millistak+® HC Pod D0HC 1.1 m ²	27	MD0HC10FS27EU
	Millistak+® HC Pod F0HC 1.1 m ²	27	MF0HC10FS27EU
	Millistak+® HC Pod X0HC 1.1 m ²	27	MX0HC10FS27EU

Optimized Packaging – APAC format

Description	Item Description	Qty/Pk	Item Number
Asian Pallet format HDPE Pallet 1100 mm × 1100 mm	Millistak+® HC Pod A1HC 1.1 m ²	36	MA1HC10FS36AS
	Millistak+® HC Pod B1HC 1.1 m ²	36	MB1HC10FS36AS
	Millistak+® HC Pod C0HC 1.1 m ²	36	MC0HC10FS36AS
	Millistak+® HC Pod D0HC 1.1 m ²	36	MD0HC10FS36AS
	Millistak+® HC Pod F0HC 1.1 m ²	36	MF0HC10FS36AS
	Millistak+® HC Pod X0HC 1.1 m ²	36	MX0HC10FS36AS

The PODs provided in bulk packaging will all be from one catalog item, all from one lot. Millistak+® improved bulk packaging configuration was tested per the ISTA 3E and real-time ship testing. Devices were then visually inspected, and integrity tested according to the purchase specification. Testing was successful.

Catalog Numbering for Disposable Adapters

Connect Millistak+® Process-Scale Pods to process piping, creating a disposable flow path.

MP0DADAPT – disposable adapter kit with 3 through adapters and 3 blind adapters

MP0DADPTF – disposable adapter kit with 6 through adapters, required if using MP0DDIVERTR

Catalog Numbering for Disposable Diverter Plates

Enable more than one media grade on a single rack

MP0DDIVERTR – disposable diverter plates, 10/pk

Process-scale Pods require a pilot or process-scale Pod holder.

Lab-scale Pods and µPod® filters do not require a holder.

Please contact your local sales representative for more information.

For additional information, please visit SigmaAldrich.com

To place an order or receive technical assistance, please visit
SigmaAldrich.com/offices

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany

We have built a unique collection of life science brands with
unrivalled experience in supporting your scientific advancements.

Millipore® **Sigma-Aldrich®** **Supelco®** **Milli-Q®** **SAFC®** **BioReliance®**

