

## **Technical Bulletin**

# 9CFR AVA Testing of Bovine Sera

SAFC Biosciences is regulated by the Food and Drug Administration (FDA) and must follow the practices and procedures outlined in the United States Pharmacopeia (USP). The USP designates how raw materials should be tested before they are used in the manufacture of finished products; however, it does not provide specific rules for adventitious viral agents (AVA) testing of products of animal origin, namely sera products.

Title 9 of the Code of Federal Regulations (9CFR) includes the standard sampling and testing requirements for animal sera used in the manufacture of biologics. SAFC Biosciences performs the entire battery of AVA tests listed in the 9CFR, Part 113, Section 53(c), on every lot of serum, where applicable.

Specifically for bovine sera, 9CFR, Part 113, Section 53(c), requires that 2 prescreened cell lines, Vero and a line of bovine origin, typically Madin-Darby Bovine Kidney (MDBK), be subcultured in the presence of test and control sera. After at least 2 subcultures over 21 days, the monolayers are examined for the presence of AVA by 3 industry-accepted methods: visualization of cytopathic effect (CPE), a hemadsorption assay (HA) and fluorescent antibody (FA) techniques.

#### **CPE**

Many viruses cause a noticeable change in the physical appearance of the cells they infect. The cells may become abnormally large, grainy, have inclusion bodies or may lyse altogether. These cytopathic effects, commonly called CPE, can be observed when viewing an infected cell monolayer under a microscope and are dependable indicators of the presence of AVA. A cytological stain is used to enhance the visual differences between cells with and without CPE.

### HA

Some viruses do not cause a noticeable change in the cellular appearance, but can be identified by a hemadsorption assay (HA). Cell monolayers are washed and covered with a suspension of guinea pig and chicken red blood cells (RBCs). After an appropriate incubation period, uninfected cells will appear normal, while cells that contain hemadsorbing agents will take up the RBCs and appear dark.

#### FA

One of the most sensitive methods for detecting AVA is labeling with fluorescently tagged antibodies (FA). Antibodies against various viral agents are first tagged with a fluorescent marker. If the monolayers are infected, the antibodies will specifically adhere to the viral particles or infected cells and will glow under a microscope emitting ultraviolet light. The antibodies are easily washed off of uninfected cultures, which will not exhibit the same fluorescence under the UV light.

The combination of these 3 techniques provides a powerful tool that is used for the detection of the following adventitious viral agents:

Bovine Viral Diarrhea (BVD) Virus Parainfluenza Type 3 (PI<sub>3</sub>) Virus **Bovine Parvovirus Rabies Virus** Reovirus Infectious Bovine Rhinotracheitis (IBR) Virus Bovine Respiratory Syncytial Virus (BRSV) Blue Tonque Virus **Bovine Adenovirus** Vesticular Stomatitis Virus

For more information about this subject or other SAFC Biosciences' products and services, please contact our Technical Services department.

#### Warranty, Limitation of Remedies

WARRANTY OF HERGALTER SAFE BIOSCIENCE SAFE BIOSCIENCE SAFE BIOSCIENCES WHITENDED FOR PURPOSES DESCRIBED ONLY AND IS NOT INTENDED FOR PURPOSES DESCRIBED ONLY AND IS NOT INTENDED. ANY HUMAN OR THERAPEUTIC USE.

Additional Terms and Conditions are contained in the product Catalog, a copy of which is available upon request.

© 2006 SAFC Biosciences, Inc.

Issued April 2006 T020 0103 1105 0106

**United States** 

SAFC Biosciences, Inc. 13804 W. 107th Street Lenexa, Kansas 66215 USA

Phone Toll free-USA +1 913-469-5580 1 800-255-6032 +1 913-469-5584 Fax

Europe

SAFC Biosciences Ltd. Smeaton Road, West Portway Andover, Hampshire SP10 3LF UNITED KINGDOM

+44 (0)1264-333311 +44 (0)1264-332412 F-mail info-eu@sial.com

Asia Pacific

SAFC Biosciences Pty. Ltd. 18-20 Export Drive Brooklyn, Victoria 3025 AUSTRALIA

+61 (0)3-9362-4500 1 800-200-404 +61 (0)3-9315-1656 Phone Toll free-AUS Fax