

Analysis of Four Aromatic Carboxylic Acids using Different Phosphorylcholine Type HILIC Columns and UV Detection

Gisela Jung, R&D Scientist, Benjamin Peters, R&D Principal Scientist

Introduction

Aromatic carboxylic acids are characterized by the presence of an aromatic ring covalently bonded to a carboxylic acid functional group. The carboxylic acid moiety contributes significantly to the overall polarity of these compounds. Due to their inherent polarity, aromatic carboxylic acids can be effectively separated using hydrophilic interaction liquid chromatography (HILIC) mode in high-performance liquid chromatography (HPLC).

The objective of this study is to assess and compare the separation efficiency of four aromatic carboxylic acids

(**Figure 1**) using a fully porous particle (FPP) ZIC®-cHILIC column and a newly developed superficially porous particle (SPP) ZIC®-cHILIC column (**Figure 2**). The compounds investigated were benzoic acid, 3-hydroxybenzoic acid, 3,4-dihydroxybenzoic acid (protocatechuic acid), and gallic acid. The ZIC®-cHILIC columns are characterized by a grafted zwitterionic poly(phosphorylcholine) modification of the stationary phase, which contributes to the chromatographic retention and separation of the analytes under investigation.




Figure 1. Aromatic carboxylic acids in this study.




Figure 2. Illustration of FPP SeQuant® ZIC®-cHILIC (left) and SPP Ascentis® Express ZIC®-cHILIC column material (right).

Experimental

A standard mixture containing four aromatic carboxylic acids at concentrations ranging from 74 to 109 $\mu\text{g/mL}$ was prepared in the mobile phase (**Table 1**). The mixture was analyzed using a fully porous particle (FPP)

SeQuant® ZIC®-cHILIC column and a superficially porous particle (SPP) Ascentis® Express ZIC®-cHILIC column. The chromatographic conditions are summarized in **Table 1**.


Table 1. Chromatographic conditions used for the analysis of aromatic carboxylic acids

LC Conditions	
Columns:	SeQuant® ZIC®-cHILIC 100 Å, 3 µm, 150 x 2.1 mm I.D. (1.50658) Ascentis® Express 160 Å ZIC®-cHILIC, 2.7 µm, 150 x 2.1 mm I.D. (50504-U)
Mobile phase:	Acetonitrile/ 25 mM ammonium acetate buffer (adjusted to pH 4.5 with acetic acid), 80/20 (v:v)
Elution mode:	Isocratic
Flow rate:	0.1 mL/min
Pressure:	FPP: 42 bar (609 psi); SPP: 36 bar (522 psi)
Columns temp.:	23 °C
Detector:	UV, 220 nm
Injection:	0.1 µL
Sample:	Aromatic carboxylic acid standard mix - Benzoic acid (8.2 mg), 3-Hydroxybenzoic acid (10.5 mg), 3,4-Dihydroxybenzoic acid (7.4 mg), and Gallic acid (10.9 mg) in 100 mL mobile phase

Results & Discussion

A standard mixture of four aromatic carboxylic acid was analyzed under identical conditions using a SeQuant ZIC®-cHILIC column (**Figure 3, Table 2**) and an Ascentis® Express ZIC®-cHILIC column (**Figure 4, Table 3**). The compounds eluted in order of increasing number of hydroxyl groups in the molecules. The elution order of the analytes remained unchanged on both columns, as anticipated, reflecting the similar ZIC®-cHILIC surface chemistry.

Using the superficially porous particle (SPP) Ascentis® Express ZIC®-cHILIC column resulted in a significant reduction in retention time (approximately 30–55%) and improved chromatographic performance. This improvement was demonstrated by an increase in separation efficiency (plate count ~1.4 to 2-fold). The higher separation efficiency resulted in peaks with greater heights, thereby improving analyte detectability.

**Figure 3.** Separation of aromatic carboxylic acids on SeQuant® ZIC®-cHILIC 100 Å, 3 µm, 150 x 2.1 mm.**Table 2.** Chromatographic results on SeQuant® ZIC®-cHILIC 100 Å, 3 µm, 150 x 2.1 mm

Peak no.	Compound	Retention time (min)	Plates (USP)	Area (mAU*min)	Tailing Factor (USP)
1	Benzoic acid	5.51	78,527	9.503	1.27
2	3-Hydroxybenzoic acid	7.97	84,333	9.059	1.12
3	3,4-Dihydroxybenzoic acid	13.72	79,067	6.419	1.10
4	Gallic acid	23.04	75,020	7.734	1.13


**Figure 4.** Separation of aromatic carboxylic acids on Ascentis® Express 160 Å ZIC®-cHILIC, 2.7 µm, 150 x 2.1 mm.

Table 3. Chromatographic results on Ascentis® Express 160 Å ZIC®-cHILIC, 2.7 µm, 150 x 2.1 mm

Peak no.	Compound	Retention time (min)	Plates (USP)	Area (mAU*min)	Tailing Factor (USP)
1	Benzoic acid	3.93	143,020	10.818	1.10
2	3-Hydroxybenzoic acid	5.11	167,193	9.171	1.08
3	3,4-Dihydroxybenzoic acid	6.31	142,847	9.050	1.30
4	Gallic acid	10.18	106,973	12.191	1.60

Conclusion

The superficially porous particle Ascentis® Express ZIC®-cHILIC column was evaluated using a mixture of aromatic carboxylic acids and compared with a fully porous particle SeQuant® ZIC®-cHILIC column. Under the chromatographic conditions employed, the elution order remained consistent on both columns, as expected,

owing to their similar ZIC®-cHILIC surface chemistry. The effect of particle architecture was evident in a significant reduction in retention time (approximately 30–55%) and a substantial increase in separation efficiency (plates), resulting in higher peak heights and improved analyte detectability/method sensitivity.

Featured & Related Products

Description	Cat. No.
HPLC Columns	
Ascentis® Express 160 Å ZIC®-cHILIC, 2.7 µm, 150x2.1 mm I.D.	50504-U
SeQuant® ZIC®-cHILIC 100 Å, 3 µm, 150x2.1 mm I.D.	1.50658
Solvents, Reagents, and Reference Materials	
Acetonitrile gradient grade for liquid chromatography LiChrosolv®	1.00030
Water from a suitable Milli-Q® system	ZIQ7000T0C
Ammonium acetate, for HPLC LiChropur	5.43834
Acetic acid or HPLC LiChropur™	5.43808
Benzoic acid, Pharmaceutical Secondary Standard; Certified Reference Material, 1 g	PHR1050
3,4-Dihydroxybenzoic acid, analytical standard, 50 mg	08992
Gallic acid, certified reference material, TraceCERT®, 100 mg	91215
3-Hydroxybenzoic acid, ReagentPlus®, 99%	H20008

See more information on our HPLC column lines at
SigmaAldrich.com/HPLC



To place an order or receive technical assistance:
SigmaAldrich.com/support



For local contact information:
SigmaAldrich.com/offices

Merck KGaA
 Frankfurter Strasse 250
 64293 Darmstadt, Germany
SigmaAldrich.com

We have built a unique collection of life science brands with unrivalled experience in supporting your scientific advancements.

Millipore® Sigma-Aldrich® Supelco® Milli-Q® SAFC® BioReliance®

© 2026 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, BioReliance, Millipore, Milli-Q, SAFC, Sigma-Aldrich, Supelco, Ascentis, SeQuant, ZIC, LiChrosolv, LiChropur, TraceCERT, and ReagentPlus are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MK_AN14405EN Ver. 1.0 70764 06/2026

