

An Alternative to the USP Method for the Analysis of Riboflavin (Vitamin B₂) and Impurities

Using a Titan C18 Column with LC/MS/MS Detection

Nicolas J. Hauser and Carmen T. Santasania
carmen.santasania@sial.com

An alternative to the cumbersome USP and EP methods using LC/MS/MS detection and a Titan C18 UHPLC column to identify riboflavin and its impurities is described here. This new method is also compatible with UV-Vis detection.

Introduction

Riboflavin (vitamin B₂) is an essential micronutrient, synthesized by plants and many microorganisms, but not by higher animals. It is a precursor of coenzymes that are required for the enzymatic oxidation of carbohydrates, and is essential to basic metabolism. Riboflavin is a common food and feed additive, and a fortificant in infant formula and cereal. It is important to establish the purity of the riboflavin because of its widespread addition to food. The USP (US Pharmacopeia) and European Pharmacopeia (EP) have long-standing analytical methods for riboflavin impurity analysis, but the methods are time-consuming and labor intensive, and are not ideally suited to modern detection methods.^{1,2}

Problems with Current Regulated Methods for Riboflavin

The current USP method for riboflavin impurity analysis is non-quantitative, and calls for the use of ion-pair chromatography to separate the riboflavin from four major impurity peaks. Both EP and USP methods call for solubilization in strong base (0.1 M NaOH), which causes analyte degradation and formation of sodium adducts that compromise mass spectrometric (MS) detection.

Overview of Approach

The purpose of this work was to develop an MS-compatible alternative to the USP and EP methods. In addition to improving the specificity of the method with MS detection, reducing the analysis time was also a

desired outcome. Improvements were made on two levels: First, LC/MS compatibility was improved by replacing the sodium hydroxide solubilizing agent and eliminating ion-pair reagents. Second, fast gradients significantly reduced the analysis time without sacrificing resolution. Known impurities were confirmed, and other impurities were identified using reference standards.

Experimental


Riboflavin certified reference material (CRM) from Sigma-Aldrich/RTC was used in this study. Compound B (7,8-dimethylbenzo[g]pteridine-2,4(1*H*,3*H*)-dione) was not available as an individual standard; however, it was present in the riboflavin CRM. Other compounds were purchased from commercial sources. MRM (multiple reaction monitoring) transitions were determined and are reported in **Table 1**. Three transitions were used for each ion. **Figure 1** shows the structures, names, monoisotopic masses, and molecular formulas of the compounds studied. Structures for each of these major impurities were confirmed by MS. Because these compounds are not soluble in methanol or water, a number of experiments were performed to identify an appropriate diluent. Ultimately, DMSO was chosen because it completely solubilized riboflavin and impurities. A stock solution of riboflavin was prepared in DMSO at 1 mg/mL with 5 minutes of sonication. This stock solution was diluted as necessary and used for all future experiments. It was stored in the freezer when not in use.

LC and MS conditions are shown in **Figure 2**. The LC column was a Titan C18, 10 cm × 2.1 mm I.D., packed with 1.9 μm totally porous, monodisperse silica particles. The mobile phase was a simple gradient of acetonitrile–water in 0.1% formic acid. A tune solution of 1 μg/mL in 50:50 (0.1% formic acid in water: 0.1% formic acid in acetonitrile) was used for riboflavin and Compound A. Compounds C and D used the same solvent system at 0.6 μg/mL.

Table 1. MRM Transitions Used in the Study

Certified Reference Material	Transition 1	Transition 2	Transition 3
Riboflavin	377.10 > 43.23	377.10 > 172.24	377.10 > 243.14
Lumiflavine (7,8,10-trimethylbenzo[g]pteridine-2,4(3 <i>H</i> ,10 <i>H</i>)-dione), Compound A	257.12 > 156.28	257.12 > 171.16	257.12 > 186.24
7,8-Dimethylbenzo[g]pteridine-2,4(1 <i>H</i> ,3 <i>H</i>)-dione, Compound B	MS only at m/z=243.41		
6,7-Dimethyl-8-[(2 <i>S</i> ,3 <i>S</i> ,4 <i>R</i>)-2,3,4,5-tetrahydroxypentyl]pteridine-2,4(3 <i>H</i> ,8 <i>H</i>)-dione, Compound C	327.15 > 43.23	327.15 > 57.27	327.15 > 193.22
8-(Hydroxymethyl)-7-methyl-10-[(2 <i>S</i> ,3 <i>S</i> ,4 <i>R</i>)-2,3,4,5-tetrahydroxypentyl]benzo[g]pteridine-2,4(3 <i>H</i> ,10 <i>H</i>)-dione, Compound D	393.16 > 43.23	393.16 > 57.27	393.16 > 259.13

Figure 1. Structures of Riboflavin and Related Compounds




Results and Discussion

Literature methods for the analysis of riboflavin and impurities from the USP and EP report run times in excess of one hour and use MS-inhibiting sodium hydroxide and ion-pair reagents. The method reported here provides a viable alternative. **Figure 2** shows the MRM results for riboflavin and the related impurities. The combination of the Titan C18 column and rapid acetonitrile

gradient provided excellent resolution and short analysis time under conditions compatible with MS detection. The replacement of sodium hydroxide with DMSO and elimination of the ion-pair reagents permitted MS detection and the concurrent improvement in specificity and selectivity. The method is also compatible with UV detection.

(continued on next page)

Figure 2. MRM Results Chromatogram



column: Titan C18, 10 cm × 2.1 mm I.D., 1.9 μm particles (577124-U)
 mobile phase: 0.1% formic acid in water in (A) water; (B) acetonitrile
 gradient: 5 to 25% B in 6 min; held at 25% B for 0.1 min; to 5% B in 4 min
 flow rate: 0.5 mL/min
 pressure: 5,000 psi (345 bar)
 column temp.: 35 °C
 detector: ESI (+), MRM and TIC mode or UV, 276 nm
 injection: 2 μL
 sample: compounds in initial mobile phase

Other MS Conditions:

capillary (kV): 2.29
 cone (V): 54.95
 extractor (V): 5.13
 source temperature (°C): 151
 desolvation temperature (°C): 349
 cone gas flow (L/Hr): 2
 desolvation gas flow (L/hr): 646
 collision gas flow (mL/min): ON

References


1. USP Method Official Monograph of August 1, 2012.
2. EP Method 7.0 Official Monograph.

Featured Products

Description	Cat. No.
HPLC Column	
Titan C18, 10 cm × 2.1 mm I.D., 1.9 μm particles	577124-U
Mobile Phase Solvents and Additives	
Acetonitrile, tested for UHPLC-MS, 1 L, 2 L	14261
Water, tested for UHPLC-MS, 1 L, 2 L	14263
Formic acid, LC/MS Ultra eluent additive, 1 mL, 2 mL	14265
Standards	
Riboflavin, pharmaceutical secondary standard; traceable to USP and PhEur	PHR1054

Unfold a Variety of Trusted Options

Our newly updated Thin Layer Chromatography Plates Brochure offers a broad selection of standard, high performance and preparative TLC Plates. With the recent addition of twenty-three TLC plates from Merck KGaA, this brochure is a complete guide to all your thin layer chromatography needs.



Download our Thin Layer Chromatography Brochure (T412113A – OWF) at sigma-aldrich.com/tlc

TLC plates from Merck KGaA.

Supelco Molded Thermogreen LB-2 Septa



The superior quality and high performance you expect from Supelco septa are available in a molded design.

Don't be fooled by other septa that advertise a maximum temperature of 400 °C! To make septa with high thermal limits, they must be made stiffer, resulting in septa that are harder to pierce and easier to core! Our molded Thermogreen LB-2 septa offer the perfect combination of low bleed, thermal stability, and easy puncturability.

To find out more about molded Thermogreen LB-2 septa, visit sigma-aldrich.com/molded-septa