


SOAK ITUP

Chemizorb® absorbents – the fast, safe and easy way to clear away chemical spills



The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.



SPILLS HAPPEN

Chemizorb® absorbents for spilled liquids

Supelco[®] Chemizorb[®] absorbents, you can clear away aggressive or unpleasant liquids quickly and safely. Our fast-acting absorbents consist of porous mineral or synthetic copolymers that are chemically inert, and capable of absorbing up to 400% of their own weight in liquid.

Handling

All Chemizorb[®] absorbents are used in a similar way:

- 1. Cover spilled liquids with sufficient absorbent
- 2. Use a spatula, spoon, or small shovel to mix the absorbent and spilled liquid
- 3. Wait until neutralization and absorption processes are complete

Specific usage instructions are provided on each product's label.

Disposal

The used absorbent should be collected in a polyethylene bag, and sent for disposal in accordance with the company regulations and national guidelines for the hazardous products in question.

Features & Benefits

Easy dosing

Wide bottle neck and integrated handle facilitate

Prompt action

The red Chemizorb[®] absorbents bottle is quickly recognizable in emergencies

Rapid safety

Chemizorb[®] absorbents can be distributed directly over the spilled liquid, and requires no accessories

Personal protection

Quick adsorption reduces time exposed to the spilled chemical

Cost efficiency

Thanks to its high capacity, less absorbents is needed for each application

Simple disposal

4. Collect the used absorbent in a polyethylene (PE) bag

5. Thoroughly clean contaminated

surface with plenty of water

Used absorbed is placed in a PE bag and disposed of according to company and national guidelines

Clear process monitoring

The pH indicator allows you to easily keep track of the neutralization process

Packaging

The HPDE bottle for Chemizorb[®] absorbents is specially designed to facilitate cleanup work. Its wide neck allows easy dispensing straight from the bottle, while its integrated handle ensures comfortable transport and dosage. In emergencies, its bright red color stands out among other HPDE bottles.

To avoid errors, the packaging's label and shape clearly indicate the product type and name, even during the application process. The different product types of Chemizorb[®] absorbents can be easily and safely identified by their specific label colors.



Salunes Gramizora Granules

Learn more SigmaAldrich.com/che<u>mizorb</u>

PRODUCT INFORMATION

		lled sub:										Тур	e
Allenis College	(concertainty)	Acids (climeter)	Concie	4 Colored Colored	Hydronuo.	49ueous Solutious	Heating A	Ordenie	Providents	Processing and a second	olis Dolinger	Genuley.	Neutralies.
	•		•	•		•	•	•	•	•	•		
	•		•	•		•	•	•	•			•	
•	•					(●) ¹					•		•
		•	•	•		(●) ²						•	•
		•			•						•		•

¹ alkaline aqueous solution

² acetic aqueous solution

Neurophise	Color cry	Mineral Conge	Synthese Synthese Copolitients	rine de la comercia	b. Delice	Get No.	
		I			Allrounder		
		•		More suitable to use in calm areas or inside (Dust formation!). Bromine must first be treated with thiosulfate solution. Hydrofluoric acid with sodium hydroxide solution before applying.	Chemizorb® Powder	102051	Page 06
		•		More suitable to use e.g. in draughty rooms or outside. Bromine must first be treated with thiosulfate solution. Hydrofluoric acid with sodium hydroxide solution before applying.	Chemizorb® Granules	101568	Page 06
				Not suitable for viscous oils			
					Specialists		
					Specialists		
Acidic salt	Blue to pale blue	•			Chemizorb [®] OH ⁻	101596	Page 08
	pale	•		Not suitable for Hydrofluoric acid		101596 103874	Page 08 Page 09
salt Calcium	pale blue Bright yellow to pale	•	•	Not suitable for Hydrofluoric acid	Chemizorb [®] OH ⁻		

Additional information

Neutralization

	All-in-one kit for mercury			
Chemizorb [®] Hg is an all-inclusive set of reagents and auxiliaries for safe, complete removal of mercury drops, and traces of elementary mercury.	Chemizorb® HG	112576	Page 10	

ONE FOR ALL

Quick help for multiple spillages

Supelco[®] Chemizorb[®] powder and granules are insoluble in water and all other media that are liquid at room temperature. These "all-rounders" are suitable for removing almost any type of liquid spills, such as acids, alkalis, organic solvents and oils.

Dosage

The required amount of Chemizorb[®] absorbents powder or granules depends on the type and quantity of the spilled liquid. Chemizorb[®] absorbents should be added until the spill is completely absorbed.



Chemizorb[®] Powder

Absorbent for spilled liquids

Absorbance capacity

Aqueous solution 1-2 times its own weight

Organic solvent/oil 2-4 times its own weight

Ordering information

Pack size	Packaging	Ord. No.
500 g	PE bottle	1.02051.0500
25 kg	Fibre carton	1.02051.9025

General information

The powder form possess a slightly higher absorbance capacity than the granules.

Suitable for: concentrated acids and alkalis, acid chlorides, aqueous solutions, organic solvents, and paraffin oils.

Chemizorb® Granules

Absorbent for spilled liquids

Absorbance capacity

100% of its own weight

Ordering information

Pack size	Packaging	Ord. No.
1 kg	PE bottle	1.01568.1000
5 kg	Bucket, plastic	1.01568.5000
20 kg	Paper sack	1.01568.9020
20 kg	PE drum	1.01568.9021

General information

The granulates are more suitable for areas where the powder form cannot be used due to its finer particle size, e.g. draughty rooms or outdoors.

Suitable for: concentrated acids and alkalis, acid chlorides, aqueous solutions, organic solvents, and paraffin oils.

Not suitable for: viscous oils



BULK UP to meet the challenge!

Whether in the lab, production plant or warehouse, any chemical spill is a problem. And, the bigger the spill, the greater the risk of personal injury and other hazards.



With Supelco[®] Chemizorb[®] absorbents, you're always prepared. Our absorbents are available in different pack sizes and materials to help you overcome the challenges of

MAKE USE OF OUR EXPERIENCE

Use the specialists

We offer specific absorbents for alkalis, acids, and hydrofluoric acid. Each contains special carrier materials, water-soluble neutralizers, and pH indicators. Supelco® Chemizorb® "Specialists" consist of a mineral or synthetic copolymer as the absorbent, and an acidic or alkaline salt as the neutralizer. The admixed pH indicators allow you to visually monitor the neutralization process of neutralization of the spilled acid or alkalis.

Dosage

The amount of Chemizorb[®] absorbents needed for alkalis, acids, and hydrofluoric acid depends on the type, quantity and concentrations of the spilled chemical.

Recommended absorbent quantities are provided for different liquids. However, sufficient Chemizorb[®] absorbents should be added in each case until the spilled liquid is completely absorbed.

Please note that the reaction may generate heat and gas.



Chemizorb[®] OH⁻

Absorbent and neutralizer for spilled alkalis, with indicator

Absorbance capacity

100 g NaOH 25% requires ~120 g Chemizorb® OH-

Ordering information

Pack size	Packaging	Ord. No.
1 kg	PE bottle	1.01596.1000

General information Suitable for: diluted and concentrated alkalis

Color change Blue to pale blue



Chemizorb[®] H⁺

Absorbent and neutralizer for spilled acids, with indicator

Absorbance capacity

100 g HCl 25% requires ~100 g Chemizorb[®] H⁺

Ordering information

Pack size	Packaging	Ord. No.
2 kg	PE bottle	1.03874.2000
5 kg	Bucket	1.03874.5000

General information

Suitable for: diluted and concentrated acids Not suitable for: Hydrofluoric acid

Color change

Bright yellow to pale yellow



Chemizorb[®] HF

Absorbent and neutralizer for spilled hydrofluoric acid, with indicator

Absorbance capacity

100 g HNO₃ 25% requires ~160 g Chemizorb[®] HF

Ordering information

Pack size	Packaging	Ord. No.
1 kg	PE bottle	1.01591.1000

General information

Suitable for: hydrofluoric acid, and other diluted or semi-concentrated acids

Chemizorb[®] HF should be added until a dry powder mixture is obtained.

Color change

Red to rose

BE PREPARED

All-in-one kit for mercury

Supelco[®] Chemizorb[®] Hg is an all-inclusive set of reagents and auxiliaries for safe, complete removal of mercury drops, and traces of elementary mercury. The reagents included in the kit are sufficient for decontamination an area of roughly one square meter.

Handling

- 1. Suction off all mercury droplets using the pipette
- 2. Empty pipette contents into the mercury bottle
- 3. Strew a layer of **reagent 1** on all remaining mercury
- 4. Spray the reagent-covered area with reagent 2
- 5. Allow 15 to 30 minutes' reaction time

- Using a small shovel and spatula, remove the mercury-containing absorbent and place it in the small tub
- 7. Wipe away any remaining material
- 8. Place all tools and working materials into the large can

Disposal

The mercury-containing material should be disposed of as special waste in accordance with official regulations for hazardous products.

Features and Benefits

Easy-to-use

Kit includes clear, visual instructions for safe decontamination

Ready

All reagents and accessories needed are stored in one case

Portable

Case's light weight and integrated handle facilitate transportation

Safe

Accessories help prevent any contact with the hazardous product









Chemizorb[®] Hg

Reagents and accessories for absorption of mercury

Ordering information

Product	Pack size	Packaging	Ord. No.
Chemizorb® Hg Reagents and accessories for absorption of mercury	1 set	PE case	1.12576.0001
Complete kit consiting of			
500 g of reagent 1, 100 mL of reagent 2, one small tub, one large dis	nocal can inn	otoctivo alovor	

500 g of reagent 1, 100 mL of reagent 2, one small tub, one large disposal can, protective gloves, ...

Product	Pack size	Packaging	Ord. No.
Chemizorb [®] Hg Reagents refill pack for Cat. No. 1.12576.0001	1 set	PE can	1.01569.0001
Reagents refill set consisting of			
500 g of reagent 1, 100 mL of reagent 2			





MilliporeSigma 400 Summit Drive Burlington, MA 01803

SigmaAldrich.com

To place an Order or Receive Technical Assistance



Find contact information for your country at SigmaAldrich.com/Offices



For Technical Service, please visit SigmaAldrich.com/TechService

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

© 2020 Merck KGaA, Darmstadt, Germany and its affiliates. All Rights Reserved. MilliporeSigma, the vibrant M, Chemizorb and Supelco are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

Lit. No. MS_BR6299EN 07/2020