



Milli-Q[®] Services: Preventive Maintenance Visit

Your purification system is a sophisticated laboratory instrument that contains highly sensitive and complex technologies to accurately measure water quality. Critical mechanical and electrical components are in regular contact with water and must be routinely inspected, calibrated and serviced to prevent damage and guarantee your water meets manufacturer's specifications.





The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

Receive a complete review of your water purification system

During the Preventive Maintenance Visit, a Milli-O[®]-certified Field Service Engineer performs a complete review of your water purification system following auditable Standard Operating Procedures. The visit includes:

- Comprehensive check of system specifications
- Replacement of aging parts using Maintenance Kit
- · Hydraulics and mechanics inspection
- Control of electronic settings & monitoring devices
- Review and testing of events, alarms and alerts
- Software updates
- Operator training
- Application assistance
- Record of service history with traceability of system parameters (Diagnostic Report)

Optimize the performance of your water system

During an annual Preventive Maintenance Visit, the Milli-Q[®]-certified Field Service Engineer performs a full checkup of all critical components and replaces aging connections and tubing, including:

Verification of electrical components	Adjustments of power supplies & set points	
Display, touchscreen and/or keypad(s)	Motor voltage(s)	
Main and interface board(s)	EDI module power supply voltage(s)	
Power supply board	UV lamp ballast voltage	
UV ballast	Pressure set point(s)	
Solenoid valve(s)	Conductivity and resistivity set points	
Pressurization/recirculation/distribution pump(s)	Rejection set points	
Pressure sensor(s)		
Resistivity/Conductivity cell(s) and thermistor(s)	Measuring & recording	
Pack Tag Reader(s) or Pack detection microswitch(s)	Feed water conductivity	
Flow meter(s)	Product resistivity	
Level sensor(s)	Total Organic Carbon (TOC)	
Alarm relay	Export of system history	
RC/Link and POD connections	Motor(s) voltage(s)	
	EDI module power supply voltages	
Verification of mechanical/hydraulic components	UV ballast voltage	
Bodies of solenoid valve, 3-way valve and motorized valve		
Pressurization/recirculation/distribution pump(s)	Software updates	
Pressure regulator(s) and pressure gauge(s)	Software upgrade for qualified systems	
Sanitary overflow	Software upgrade for non-qualified systems	

Reverse osmosis housing and chevron seal

Verification & regulation of hydraulics

Product/distribution water flow rates

Rejection flow rates

Pump pressure

Software upgrade for non-qualified systems

Maintenance Kit to replace aging parts

Internal/external tubing

Internal/external fittings

Pump inlet/outlet fittings

Solenoid valve coils

Sanitary overflow solution

O-ring(s) for cleaning port and point-of-use

O-ring(s) and strainer(s) for conductivity cell(s) or inlet fittings

Check valve(s)

We also offer a set of tailor-made options for your lab's specific needs, including:

- Verification
- Calibration
- Pharmacopeia suitability tests
- Extended guarantees

• Customized user training

Scheduled shipments of consumables and parts

Sanitization

Additional Preventive Maintenance Visits

As part of our Qualification program, Preventive Maintenance Visits that are done within the framework of a Service Pharma[™] Plan, are performed and documented according to a Maintenance Procedure (MP) in compliance with GLP and cGMP.

Service history traceability

Via the provided Diagnostic Report and Service Visit Report, you benefit from full traceability and archiving of your system's service history and system parameters.

Extract of the Diagnostic Report

Milli-Q®		Systen	n Diagnosti	cs Report	
Lab Water Solution	าร				
nformation					
Work Order Number					
Work Order Type					
Business Field Name					
Product Name					
Serial Number					
Scheduled Date Time					
Work to be performed					
Preventive Maintenance V					
	en according to o or				
Work Performed					
				ensors, solenoid va	lves, wires and connectors, tubing & fittings connections. Replacement of
maintenance kit. Tank San	itization. TOC Cleaning and	RO Cleaning Soft	ware Upgrade		
Reverse Osmosis Diag	nosis and Maintenand	e	la i		
Parameters	System Specification	Pre Maintenance	Post	Unite	
Parameters	System Specification	Value	Maintenance Value	Units	Comments
Regulator Pressure	1 b <=P<= 2.5 b	2	2	bar	
RO Feed Temperature	5° C <=T<= 35 °C	18	18	°C	
RO Feed Conductivity	1 µS <=C<= 2000 µS	733	735	μS	
RO Membrane Pressure	3 b <=P<= 7.5 b	2.5	5		Performed a RO Cleaning
RO Permeate Conductivity		70	25	μS	
RO Reject Flowrate	12 l/h <=Q<= 22 l/h	11.5	17.5	1/h	
RO Rejection RO Pump Voltage	rej> 92 % 4 V <=U<= 40 V	90.45 23	96.60 23	% I	RO rejection is back to normal after cleaning
KO Pump voltage]4 V <=U<= 40 V	23	125	v	
UV Lamp Diagnosis a	nd Maintenance				
Parameters	System Specification	Pre Maintenance Value	e Maintenance Value	Units	Comments
254nm UV Lamp Current	I > 21 mA	22	23	mA	
185nm UV lamp Current		72	72	mA	
Point of Use Performa	nce & Quality				
		Maintenance	Post Maintenance Value	Units	Comments
	Q>4.25 l/h		5.2	l/h	
	R> 5 MOhm.cm@25°C		10	MOhm.cm@25°	c
Product temperature	5° C <=T<= 35 °C		20.4	°C	
Milli-O Product resistivity	R= 18.2 MOhm.cm@25°C	18.2	18.2	MOhm.cm@25°	c
	TOC<= 5 ppb*	7	2	ррь	Performed a TOC Cleaning
Total Organic Carbon (TOC)				1/min	
Total Organic Carbon (TOC) E-Pod Flowrate	1.3 l/mn <=Q<= 2.16 l/mn		2	l/min	
Total Organic Carbon (TOC)			2 2	l/min	

For more information, visit our website:

EMDMillipore.com/Milli-QServices

MilliporeSigma 400 Summit Drive Burlington, MA 01803

© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. MilliporeSigma, the vibrant M, and Milli-Q are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.