DNA Extraction from Agarose Gels with the Ultrafree[®]-DA Centrifugal Filter Device

The Ultrafree-DA device is designed to recover 100 to 10,000 bp DNA from agarose gel slices in one 10minute spin. It consists of a preassembled sample filter cup with agarose Gel Nebulizer, and a microcentrifuge vial. The device utilizes gel compression to extract DNA from the agarose. Centrifugal force collapses the gel structure, drives the agarose through a small orifice in the Gel Nebulizer and the resultant gel slurry is sprayed into the sample filter cup. As the agarose is compressed at 5,000xg, DNA is extruded from the gel's pores. The gel matrix is retained by the microporous membrane, and the DNA passes freely through the membrane. DNA can then be recovered in the filtrate vial.


DNA prepared with this device requires no further purification for most applications, including cloning and radioisotopic or fluorescent DNA sequencing. Since agarose gel electrophoresis has high resolving power, the small and large non-specific amplification products that frequently interfere with cloning and sequencing after PCR (polymerase chain reaction) are completely removed from the product.

Materials

- Microcentrifuge
- Pre-assembled Ultrafree-DA centrifugal filter device
- Modified TAE* electrophoresis buffer (40 mm Tris-acetate, pH 8.0, 0.1 mm Na₂EDTA)
- SeaKem[™] LE agarose (FMC BioProducts; Rockland, ME) or equivalent
- Long-wavelength UV lamp
- Scalpel or razor blade
- * Modified TAE rather than TBE is recommended for the following reasons: (1) TBE buffer strongly inhibits DNA sequencing reactions while modified TAE buffer does not. (2) Modified TAE has 0.1 mm Na₂EDTA while regular TAE has 0.0 mm Na₂EDTA. The EDTA level at 0.1 mm Na₂EDTA will not interfere with the magnesium concentration in sequencing reactions and other downstream enzymatic treatments, many of which are dependent on magnesium.

Procedure

- Electrophorese 30 µL of PCR product or other DNA through a <1.25% ordinary agarose gel, prepared in modified TAE buffer with ethidium bromide (0.5 µg/ml).
- Locate band of interest with a long wavelength UV lamp or transilluminator. With a razor blade or scalpel, cut out the slice of agarose (<100 µL or 100 mg) containing the band of interest. Trim any excess agarose away from band.



 Place gel slice into Gel Nebulizer sample cup assembly and seal device with the cap attached to vial.

MILLIPORE

- 4. Spin at 5,000xg for 10 minutes. Centrifugation forces the agarose through the Gel Nebulizer, converting it to a fine slurry that is captured by the sample filter cup. Extruded DNA in electrophoresis buffer passes through the microporous membrane in the sample filter cup and collects in the filtrate vial.
- 5. DNA in the filtrate is now ready for sequencing or cloning without further purification. After discarding the filter cup and Gel Nebulizer sample cup assembly, the DNA can be stored in the capped filtrate vial.


Typical DNA Recoveries

Gel Compression is a quick and easy technique for recovering DNA from an agarose gel slice. The Ultrafree-DA includes a Gel Nebulizer that allows the addition of an intact gel slice to the Nebulizer/Sample Cup stack. There is no need to manually cut up the gel slice since gel disintegration is automatic. With Ultrafree-DA the gel is completely macerated by the Nebulizer, thus recoveries of DNA are typically higher than when this component is not incorporated.

Typical DNA Recoveries

	Percent DNA Recovered	
DNA size (bp)	Intact Gel	Gel Nebulizer Disrupted Gel
100	74	78
400	39	ND
700	43	71
1000	55	77
2027	ND	47
4361	14	35
9416	ND	32
23130	ND	29

ND = Not detectable.



To Place an Order or Receive Technical Assistance

For additional information call your nearest Millipore office: In the U.S. and Canada, call toll-free 1-800-MILLIPORE (1-800-645-5476) In the U.S., Canada and Puerto Rico, fax orders to 1-800-MILLIFX (1-800-645-5439) On the Internet http://www.millipore.com/amicon E-mail: tech_service@millipore.com

MILLIPORE

Millipore is a registered trademark of Millipore Corporation or an affiliated company. Ultrafree is a trademark of Millipore Corporation or an affiliated company. SeaKem is a trademark of FMC BioProducts. Lit. No. AN024EN00

Printed in U.S.A. 8/99 99-182 Copyright 1999 Millipore Corporation or an affiliated company, Bedford, MA. All rights reserved.

Printed with soy inks on recycled paper.

Millipore Worldwide

AUSTRALIA Tel. 1 800 222 111 or (02) 9888 8999 Fax (02) 9878 0788

AUSTRIA, CENTRAL EUROPE, C.I.S., AFRICA, MIDDLE EAST AND GULF Tel. (43) 1 877-8926

Fax (43) 1 877-1654 **BALTIC COUNTRIES** Tel. +358 9 804 5110 Fax +358 9 256 5660

BELGIUM AND LUXEMBOURG Tel. (02) 726 88 40 Fax (02) 726 98 84

BRAZIL Tel. (011) 548-7011 Fax (011) 548-7923

CANADA Tel. 1-800-645-5476 Fax 1-800-645-5439

CHINA, PEOPLE'S REPUBLIC OF Beijing:

Beijing: Tel. (8610) 6500-8063 Fax (8610) 6500-7372 Guangzhou: Tel. (8620) 8755-4021 Fax (8620) 8755-4350 Hong Kong: Tel. (852) 2803-9111 Fax (852) 2513-0313 Shanghai: Tel. (8621) 5306-9100 Fax (8621) 5306-0838 Tel. 02-205 138 41 02-205 138 42 Fax 02-205 14294 **DENMARK** Tel. 70 10 00 23

CZECH REPUBLIC

Fax 70 10 13 14 FINLAND Tel. (09) 804 5110 Fax (09) 256 5660

FRANCE Tel. (01) 30.12.70.00 Fax (01) 30.12.71.80

GERMANY Tel. (06196) 494-0 Fax (06196) 43901

HUNGARY Tel. 1-205 9784 Fax 1-205 9792

INDIA Tel. (91) 80-839 46 57 Fax (91) 80-839 63 45

ITALY Vimodrone (Milano): Tel. (02) 25.07.81 Fax (02) 26.50.324 Roma: Tel. (06) 52.03.600 Fax (06) 52.95.735 JAPAN Tel. (03) 5442-9711 Fax (03) 5442-9711 Fax (03) 5442-9736 Analytical 9737 BioProcess

9734 Lab Water

KOREA Tel. (822) 551-0990 Fax (822) 551-0228

MALAYSIA Tel. 603-757 1322 Fax 603-757 1711

MEXICO Tel. (525) 576-9688 Fax (525) 576-8706

THE NETHERLANDS Tel. 076-5022000 Fax 076-5022436

NORWAY Tel. 22 67 82 53 Fax 22 66 04 60

POLAND Tel. 22-669 12 25 22-663 70 31 Fax 22-663 70 33

PUERTO RICO Tel. (787) 273-8495 Fax (787) 747-6553

SINGAPORE Tel. (65) 842 1822 Fax (65) 842 4988

SPAIN AND PORTUGAL Madrid:

Tel. 917 283 960 Fax 917 292 909 Barcelona: Tel. 934 525 530 Fax 934 516 048

SWEDEN Tel. 08-628 69 60

Fax 08-628 64 57

SWITZERLAND

Volketswil: Tel. (01) 908-30-60 Fax (01) 908-30-80 Lausanne: Tel. (021) 648 51 18 Fax (021) 648 51 37

TAIWAN

Taipei: Tel. (886-2) 700-1742 Fax (886-2) 755-742 Hsin Chu City: Tel. (886-3) 571-0178 Fax (886-3) 572-9520

U.K. AND IRELAND

Tel. (01923) 816375 Fax (01923) 818297

U.S.A. Tel. (781) 533-6000 Fax (781) 533-3110

IN ALL OTHER

COUNTRIES Millipore Intertech (U.S.A.) Tel. +1 (781) 533-8622 Fax +1 (781) 533-8630