

Data Sheet

GFP, FLAG-tagged SNRNP70 mRNA

Synthetic RNA

03-903**Pack Size 50 µg****Store at -80 °C****FOR RESEARCH USE ONLY****Not for use in diagnostic procedures. Not for Human or Animal Consumption.**

Background

Both mRNAs and non-coding RNAs interact with RNA-binding proteins and regulate critical cellular processes and their dysregulation contributes to multiple diseases. Thus, analyzing RNA-protein interactions are critical for understanding the mechanism of diseases and designing clinical strategies.

RNA-binding protein immunoprecipitation (RIP) has been used to interrogate RNA binding proteins (RBP) associated with specific RNA molecules. However, the experiments are not feasible because of the unavailability of antibodies for the particular RBP or low expression of the RBP. To overcome these issues epitope tags have been used for RIP experiments.

The GFP, FLAG-tagged SNRNP70 mRNA was developed to enable to perform control RIP experiments with RIPAb+ FLAG (03-901) or RIPAb+ GFP (03-902) RIP validated antibody and primer kits, and Magna RIP® Kit (17-700).

The high amount of GFP-FLAG-tagged SNRNP70 protein is expressed immediately after transfection to cells and the well-known specific interaction of SNRNP70 and U1snRNA can be used for control RIP experiments with anti-GFP or anti-FLAG antibodies.

GFP-tag and FLAG-tag sequences are fused to the 3' of the SNRNP70 mRNA. The GFP tag uses TagGFP2 an improved variant of the Aequorea macrodactyla GFP-like protein. TagGFP2 exhibits bright green fluorescence comparable to that of EGFP, with excitation/emission maxima at 483 and 506 nm, respectively. The anti-Green Fluorescent Protein Antibody in the RIPAb+ GFP (03-902) interact with both GFP and tagGFP2.

Source

In vitro synthesized RNA

Figure 1. Structure of GFP, FLAG-tagged SNRNP70 mRNA.

Storage and Handling

Store at: -80 °C

Presentation

GFP, FLAG-tagged SNRNP70 mRNA: (03-903) One (1) vial containing 50 μ L of RNA (1 μ g/ μ L)

Quality Control Testing




Figure 2. Quality of the RNA with Bioanalyzer and the Expression of TagGFP2 in HeLa Cells

(A) GFP, FLAG-tagged SNRNP70 mRNA was analyzed with Bioanalyzer. L: RNA ladder, 1: mRNA, (B, C) GFP, FLAG-tagged SNRNP70 mRNA was transfected into HeLa cells. One day after transfection, GFP expression was analyzed by flow cytometry and fluorescence microscopy.

Materials Required (Not Supplied)

- RiboJuice™ mRNA Transfection Kit (TR-1013)
- RIPAb+ FLAG (03-901) or RIPAb+ GFP (03-902) RIP validated antibody and primer kits
- Magna RIP® Kit (17-700)

Representative Data




Figure 3. RNA Binding Protein Immunoprecipitation

RIP Lysate was prepared from HeLa cells (1×10^7 cell equivalents per IP) expressing SNRNP70-GFP-FLAG (03-903) by RNA transfection. The lysate was subjected to immunoprecipitation using (A) 5 μ g of either a normal mouse IgG (CS200621) or Anti-FLAG M2 antibody (CS226432), and (B) 5 μ g of either a Rabbit IgG purified (PP64B) or Anti-Green Fluorescent Protein Antibody (AB3080) with the Magna RIP® RNA-Binding Protein Immunoprecipitation Kit (17-700). Successful immunoprecipitation of SNRNP70-GFP-FLAG associated RNA was verified by qPCR using RIP Primers, U1 snRNA, (CS203215)

Note: Please refer to the Magna RIP® (17-700) or EZ-Magna RIP® (17-701) protocol for experimental details.




Figure 4. Immunoprecipitation from RIP lysate

RIP Lysate was prepared from HeLa cells (1×10^6 cell equivalents per IP) expressing SNRNP70-GFP-FLAG (03-903) by RNA transfection. The lysate was subjected to immunoprecipitation using (A) 0.5 μ g of either normal rabbit IgG, (CS200621) or Anti-FLAG M2 Antibody (CS226432), (B) 0.5 μ g of either Rabbit IgG purified (PP64B) or Anti-Green Fluorescent Protein Antibody (AB3080). Precipitated proteins were resolved by electrophoresis, transferred to nitrocellulose, and probed with (A) anti-FLAG M2 (1.0 μ g/mL) or (B) Anti-Green Fluorescent Protein Antibody (1.0 μ g/mL). Proteins were visualized using (A) a goat anti-mouse secondary antibody conjugated to HRP (AP124P) or (B) a goat anti-Rabbit secondary antibody conjugated to HRP (AP187P), and a chemiluminescence detection system.

Transfection Protocol

Transfection of GFP, FLAG-tagged SNRNP70 mRNA has been validated using RiboJuice™ mRNA Transfection Kit (TR-1013). The amounts of RNAs and transfection reagents required may vary depending on the target cells. Here is an example to transfect HeLa cells on a 10 cm plate ($\sim 1.0 \times 10^7$ cells). 5 to 10 μ g of the mRNA per plate would be used.

1. Plate cells to reach 50-90% confluency at the time of transfection.
2. Wash cells once with EMEM (no serum, no antibiotics) and add 5 mL/plate of EMEM (no serum, no antibiotics). No serum condition increases the transfection efficiency. However, it is possible to use 1-10% serum depending upon cell types. Incubate cells in a 37 °C, 5% CO₂ incubator (10-20 minutes).
3. Set up transfection reactions in sterile Eppendorf tubes following the table. Follow the order of additions. Mix gently by pipetting during each addition of RNA and reagent. Do not vortex.

RiboJuice™ mRNA Transfection Kit

Component	5 μ g	10 μ g
Opti-MEM® I Reduced-Serum Medium (Life Technologies, 31985062)	500 μ L	1000 μ L
B18R-E3L RNA (1 μ g/ μ L)	5 μ L	10 μ L
RiboJuice™ mRNA Boost Reagent	10 μ L	20 μ L
RiboJuice™ mRNA Transfection Reagent	10 μ L	20 μ L
Total Volume	525 μL	1050 μL

4. Incubate the tubes at room temperature for 5 minutes and add the RNA-transfection reagent complex dropwise into plates containing cells.
5. Incubate the plate in a 37 °C, 5% CO₂ incubator for 2-4 hours. Maximum transfection efficiency is obtained with 4 hours of incubation using RiboJuice™ mRNA Transfection Kit with HeLa cells.
6. Aspirate the transfection medium and add 10 mL per plate of EMEM containing 10% FBS. Incubate in a 37 °C, 5% CO₂ incubator overnight.
7. Next day, Check the expression of GFP-FLAG tagged protein with a fluorescent microscope. Harvest the cells and perform RIP assay with Magna RIP® Kits and RIPAb+ antibodies.

Related Products

Description	Catalogue Number	Description	Catalogue Number
RIPAb+ EZH2, clone AC22	03-900	Anti-Ago1, clone 6D8.2	04-083
RIPAb+ FLAG	03-901	Anti-Ago 2	07-590
RIPAb+ GFP	03-902	Anti-Ago family	04-085
Magna RIP® RNA-Binding Protein Immunoprecipitation Kit	17-700	Anti-pan Ago, clone 2A8	MABE56
Magna RIP® Quad RNA-Binding Protein Immunoprecipitation Kit	17-704	Anti-Ago1, clone 4G7-E12	MABE143
EZ-Magna RIP® RNA-Binding Protein Immunoprecipitation Kit	17-701	Anti-Ago3, clone 4B1-F6	MABE144
		RIPAb+™ pan Ago	03-248

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck operates
as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates.
All other trademarks are the property of their respective owners. Detailed information on
trademarks is available via publicly accessible resources.

© 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

Document Template 20306518 Ver 6.0

23070403 Ver 1.0, Rev 23MAY2024, ab

