

Technical Bulletin

EX-CELL[™] CD CHO: A Chemically Defined, Animal-Component Free, Serum-Free Medium for CHO Cells

Introduction

SAFC Biosciences has developed EX-CELL[™] CD CHO, a chemically defined, animal-component free, serum-free medium optimized for the growth of CHO-derived cell lines producing recombinant proteins, including monoclonal antibodies (MAb). EX-CELL[™] CD CHO is hydrolysate-free and contains no animal-derived components.

The following study was undertaken to demonstrate the ability of EX-CELL[™] CD CHO to support the growth and MAb production in a CHO-derived (DXB11) cell line producing human MAb immunoglobulin G (IgG₄). Additionally, the ability of EX-CELL[™] CD CHO to support the growth of a non-producing CHO-K1 cell line was assessed. Both cell lines were adapted to EX-CELL[™] CD CHO using direct adaptation from serum-containing suspension cultures prior to the evaluation of growth and MAb production.

Materials

Cells

- DXB11, American Type Culture Collection, ATCC No. CRL-11397
- CHO-K1, American Type Culture Collection, ATCC No. CRL-61

Media and Supplements

- EX-CELL[™] CD Cho, with hypoxanthine and thymidine, SAFC Biosciences, Catalog No. 14361
- Iscove's Modified Dulbecco's Medium (IMDM), American Type Culture Collection, ATCC No. 30-2005
- Dulbecco's Modified Eagle's Medium/Ham's Nutrient Mixture F12 (DMEM/F12), SAFC Biosciences, Catalog No. 51448
- Fetal Bovine Serum Gamma Irradiated (FBS), SAFC Biosciences, Catalog No. 12107
- 200 mM L-glutamine, SAFC Biosciences, Catalog No. 59202

Antibody Assay Kit

• Immuno-Tek Quantitative IgG ELISA Kit, ZeptoMetrix Corporation, Catalog No. 0801182

Methods

Media/Supplement Preparation and Storage

EX-CELL[™] CD CHO was supplemented with 8 mM L-glutamine at point of use. All media were stored at 2 to 8 C protected from light. Cultures were maintained using aseptic technique with no antibiotic or fungicide supplementation.

Culture Techniques

Prior to adaptation, CHO-K1 and DXB11 cell lines were maintained as suspension cultures in 125 mL shaker flasks in DMEM/F12 (CHO-K1) and IMDM (DXB11) supplemented with 10% FBS and 4 mM L-glutamine. After adaptation to serum-free medium, the cells were transferred to 250 mL shaker flasks and routinely subcultured every 3 - 4 days at a density of 3 x 10⁵ cells/mL and 4 x 10⁵ cells/mL respectively (60 mL volume per 250 mL shaker flask). The flasks were shaken on an orbital shaker at 120 - 130 rpm and were maintained at 37 C in a humidified atmosphere with 10% CO₂. Cell densities and viabilities were determined by trypan blue exclusion.

Adaptation

The CHO-K1 and DXB11 cell lines were adapted to EX-CELLTM CD CHO by direct adaptation. During the first 2 - 3 passages, suspension cultures previously grown in basal medium + 10% FBS were seeded directly into pre-warmed EX-CELLTM CD CHO at a density of 5 x 10⁵ cells/mL. Cells were subcultured every 3 - 4 days and were considered fully adapted when cell densities achieved > 2.0 x 10⁶ cells/mL and viabilities of > 95%.

United States SAFC Biosciences, Inc. 13804 W. 107th Street Lenexa, Kansas 66215 USA Phone +1 913-469-5580 Toll free-USA 1 800-255-6032 Fax +1 913-469-5584 E-mail info-na@sial.com Europe SAFC Biosciences Ltd. Smeaton Road, West Portway Andover, Hampshire SP10 3LF UNITED KINGDOM Phone +44 (0)1264-333311 Fax +44 (0)1264-332412 E-mail info-eu@sia.com
 Asia Pacific

 SAFC Biosciences Pty. Ltd.

 18-20 Export Driv∈

 Brooklyn, Victoria 3025

 AUSTRALIA

 Phone
 +61 (0)3-9362-4500

 Toll free-AUS
 1 800-200-404

 Fax
 +61 (0)3-9315-1656

 E-mail
 info-ap@sial.com

www.safcbiosciences.com

Growth Studies and Antibody Production

Growth studies were initiated after adaptation and were monitored over an additional 6 passages. The cells were subcultured one additional time and daily cell counts were determined until culture viabilities dropped below 60%. Additionally, during the last passage, daily aliquots from the DXB11 cell suspension were taken for IgG determination. Each aliquot was micro-centrifuged (1000 rpm) for 2 minutes, followed by collection and storage (-20 C) of the supernatant. Antibody (human IgG) production was determined by ELISA (ZeptoMetrix Human IgG ELISA kit). The appropriate sample dilutions were prepared in dilution buffer supplied with the kit. The absorbance was read at 405 nm on a VersaMax[™] microplate reader and calculations were performed using SoftMax[®] Pro. 4.0 software (both from Molecular Devices Corporation).

Results


Growth Studies

Both the CHO-K1 and DXB11 cell lines typically reached densities of at least 3.5×10^6 and 2.5×10^6 respectively when subcultured every 4 days in EX-CELLTM CD CHO. Figure 1 depicts the typical growth of CHO-K1 and DXB11 cells over the course of 6 passages in EX-CELLTM CD CHO. Figure 2 illustrates the culture longevity of CHO-K1 and DXB11 cells over a period of 7 - 8 days in EX-CELLTM CD CHO medium.


CHO-K1 densities increase steadily over the first 4 days, reach a plateau by day 5 and steadily decline after day 5. Viabilities start to drop after day 5 and are below 60% on day 7. DXB11 densities increase steadily over the first 5 days, reach a plateau by day 6 and decline after day 6. Viabilities start to drop after day 5 and are below 60% on day 7.


Monoclonal Antibody Production


A human IgG ELISA assay kit was used to titer the production of IgG by the DXB11 cell line. Figure 3 depicts a typical viability and concurrent IgG production in the DXB11 cell line cultured in EX-CELL[™] CD CHO. Viabilities remain above 90% until day 5 and proceed to decline thereafter. IgG production steadily increases with peak productivity on day 8. Additional evaluations conducted by SAFC Biosciences demonstrate that IgG production in EX-CELL[™] CD CHO is comparable to levels typically experienced in a basal formulation + 10% FBS (Figure 4). These studies demonstrate that EX-CELL[™] CD CHO supports high density cell growth and MAb production in a CHO-derived cell line.



www.safcbiosciences.com







www.safcbiosciences.com

Warranty, Limitation of Remedies

Warranty, Limitation of Remeates SAFC Biosciences warrants to the purchaser for a period of one year from date of delivery that this product conforms to its specifications. Other terms and conditions of this warranty are contained in SAFC Biosciences' written warranty, a copy of which is available upon request. ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABLITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE EXCLUDED. In no case will SAFC Biosciences be liable for any special, incidental, or consequential damages arising out of this product or the use of this product by the customer or any third party based upon breach of warranty gainst claims by any third party by way of infringement or the like. THIS PRODUCT IS INTENDED FOR PURPOSES DESCRIBED ONLY AND IS NOT INTENDED FOR ANY HILMAN OR THERAPTITIC LISE. ANY HUMAN OR THERAPEUTIC USE.

Additional Terms and Conditions are contained in the product Catalog, a copy of which is available upon request.

EX-CELL[™] is a trademark of SAFC Biosciences, Inc. SoftMax[®] is a registered trademark of Molecular Devices Corporation. VersaMax[™] is a trademark of Molecular Devices Corporation.

© 2006 SAFC Biosciences, Inc.

Issued March 2006 T067 0404 1105

www.safcbiosciences.com

United States

SAFC Biosciences, Inc. 13804 W. 107th Street Lenexa, Kansas 66215 USA Phone Toll free-USA +1 913-469-5580 1 800-255-6032 Fax E-mail +1 913-469-5584 info-na@sial.com

Europe

SAFC Biosciences Ltd. Smeaton Road, West Portway Andover, Hampshire SP10 3LF UNITED KINGDOM Phone +44 (0)1264-333311 Fax +44 (0)1264-332412 E-mail info-eu@sial.com

Asia Pacific

Fax E-mail

SAFC Biosciences Pty. Ltd. 18-20 Export Drive Brooklyn, Victoria 3025 AUSTRALIA Phone Toll free-AUS +61 (0)3-9362-4500 1 800-200-404 +61 (0)3-9315-1656 info-ap@sial.com