Instructions for Use

Trichrome Stains (Gomori)

Procedure No. HT10

Intended Use

Trichrome Stains (Gomori) are intended for use in the histochemical study of connective tissue, muscle and collagen fibers. Trichrome Stain reagents are for "In Vitro Diagnostic Use". For professional use only. The data obtained from this manual, qualitative procedure demonstrates muscle and collagen fibers in tissue samples of human specimens. This data when reviewed in conjunction with other diagnostic tests and information may be used as an aid to diagnosis of increases in collagenous tissue associated with liver disfunction.

"Trichrome" stains are used primarily for distinguishing collagen from muscle tissue.¹ In general, they consist of nuclear, collagenous and cytoplasmic dyes in mordants such as phosphotungstic or phosphomolybdic acid. Historically, the first trichrome system was attributed to Mallory.^{2,3} Further modifications were introduced by Masson, Foot and Gomori.^{3,4} The procedure described here is based on the work of Gomori and is a one-step system combining the cytoplasmic and connective fiber stain in a phosphotungstic acid/acetic acid solution.⁴ Tissue sections are treated with Bouin's solution to intensify the final coloration. Nuclei are then stained with Weigert's iron hematoxylin, followed by chromotrope 2R (stains cytoplasm and muscle fibers) and fast green FCF or aniline blue (stains collagen fibers). Rinsing in acetic acid after staining makes the shades of color more delicate and transparent.⁴ Included is a Gomori trichrome stain technique for rapid staining in microwave ovens.⁵⁻⁸

Reagents

Trichrome Stain LG Solution (Cat. No. HT10316-500ML)

Chromotrope 2R, 0.6% (w/v), C.I. 16570, fast green FCF, (certified), 0.1% (w/v), C.I. 42053, phosphotungstic acid, 0.8% (w/v) and acetic acid, 1.0% (v/v). Reagent stains collagen green

Trichrome Stain AB Solution (Cat. No. HT10516-500ML)

Chromotrope 2R, 0.6% (w/v), C.I. 16570, aniline blue, (certified), 0.3% (w/v), C.I. 42780, phosphotungstic acid, 0.8% (w/v) and acetic acid, 1.0% (v/v). Reagent stains collagen blue.

Special Materials Required but Not Provided

- Positive control slides, such as Trichrome TISSUE-TROL™, (Cat. No. TTR012) should be included in each run
- 1N Acetic Acid Solution
- Bouin's Solution (Cat. No. HT10132-1L or HT101128-4L) picric acid, 1%, formaldehyde, 9% and acetic acid, 5%

For Standard Procedure Only

- Weigert's Iron Hematoxylin Set (Cat. No. HT1079)
- Weigert's Iron Hematoxylin Solution, Part A (Cat. No. HT107-500ML) 1% certified hematoxylin in ethanol
- Weigert's Iron Hematoxylin Solution, Part B (Cat. No. HT109-500ML), Ferric chloride 1.2% (w/v) and hydrochloric acid, 1% (v/v)

For Microwave Procedure Only

- Hematoxylin Solution, Gill No. 3 (Cat. No. GHS3-100ML)
- Scott's Tap Water Substitute Concentrate (Cat. No. S5134-6x100ML)
- Microwave Oven

Storage and Stability

Trichrome Stain solutions, Bouin's solution, Weigert's Iron Hematoxylin, Trichrome TISSUE-TROL™, Acetic Acid solution, Scott's Tap Water Substitute Concentrate, and Scott's Tap Water Substitute Working Solution should be stored at room temperature (18–26°C).

Reagents are stable until expiration date shown on the label.

Weigert's Iron Hematoxylin Working Solution should be prepared fresh for each use.

Preparation

Trichrome Stain LG Solution, Trichrome Stain AB and Bouin's solution are ready to use.

Weigert's Iron Hematoxylin Working Solution is prepared by mixing equal parts of Weigert's Iron Hematoxylin Solution, Part A and Weigert's Iron Hematoxylin Solution, Part B.

Scott's Tap Water Substitute working solution is prepared by mixing 1 part of Scott's Tap Water Substitute Concentrate with 9 parts deionized water.

Prepare 0.5% Acetic Acid by diluting 4.4 mL 1N acetic acid with 45.6 mL water.

Prepare 1% Acetic Acid by diluting 8.8 mL 1N acetic acid with 41.2 mL water.

Precautions

These IVDs are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Trichrome TISSUE-TROL™ control slides are paraffin embedded human tissue containing collagen and muscle tissue and should be considered potentially infectious.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Any well fixed paraffin section cut at 5 micron or less may be used.

Notes

If a microwave oven is used, please see the Owner's Manual for instructions.

Procedure

Standard Procedure

- Prepare Weigert's Iron Hematoxylin Working Solution according to instructions.
- Deparaffinize slides and hydrate to deionized water.
- Place slides in preheated Bouin's Solution at 56°C for **15 minutes** or at **room temperature overnight**.
- Cool slides in tap water (18–26°C) contained in a Coplin jar.
- Wash in running tap water to remove yellow color from sections.
- Stain in Weigert's Iron Hematoxylin Working Solution for **5 minutes**.
- Wash in running tap water for **5 minutes**.
- Stain in Trichrome Stain (LG or AB) Solution for **5 minutes**.
- Place in 0.5% Acetic Acid for **1 minute**.
- Rinse slides, dehydrate through alcohol, clear in xylene and mount.

Microwave Procedure

- Deparaffinize slides and hydrate to deionized water.
- Place slides in **40 mL** of Bouin's Solution in plastic Coplin jar. Loosely cover jar with lid or use lids with vented caps.
- Microwave on **600 watts** for **25 seconds**. Agitate solution with applicator stick or beral pipet then incubate for **5 minutes** in a fume hood or well-ventilated area.
- Rinse in running tap water until yellow color disappears.
- Place slides in **40 mL** of Hematoxylin Solution, Gill No. 3, contained in a plastic Coplin jar.
- Microwave on **800 watts** for **5 seconds**.
- Rinse well in tap water for **30 seconds** to **1 minute**.
- Blue in Scott's Tap Water Substitute working solution at **room temperature**.
- Rinse well in running tap water.
- Place slides in **40 mL** of Trichrome Stain (LG or AB) Solution contained in a plastic Coplin jar.
- Microwave on **800 watts** for **15 seconds**. Gently agitate solution using an applicator stick or beral pipet. Let incubate for **1 minute**.
- Rinse quickly in tap water.
- Place slides in 1% Acetic Acid for **30 seconds** to **1 minute** at room temperature.
- Rinse slides, dehydrate through alcohol, clear in xylene and mount.

Performance Characteristics

Target Structure	Staining Result
Nuclei	Blue to Black
Collagen	Green (LG) or Blue (AB)
Cytoplasm, Muscle	Red

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
HT10316	Trichrome Stain LG Solution	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
		Collagen	3 of 3	3 of 3	3 of 3	3 of 3
		Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
		Muscle	3 of 3	3 of 3	3 of 3	3 of 3
HT10516	Trichrome Stain AB Solution	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
		Collagen	3 of 3	3 of 3	3 of 3	3 of 3
		Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
		Muscle	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

HT10316, HT10516:

Not a hazardous substance or mixture.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

1. Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
2. Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
3. Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
4. Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 20:661, 1950
5. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
6. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
7. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
8. Histotechnology: A Self Instructional Text. Carson, ASCP Press, Chicago, 1990
9. Kok and Boon: Microwave Cookbook for Microscopists. Couloomb Press, Leyden. Leyden 1992

Contact Information

To place an order, please visit our web site at [SigmaAldrich.com](https://www.sigmaaldrich.com). For Technical Service, please visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Revision History

Rev. 5.0	2016
Rev. 6.0	2022
Rev. 7.0	2022 Transferred to new template with current branding. Specified for professional use in intended use and precautions. Moved aid to diagnosis statement to intended use. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Added a revision history table. Added Warnings and Hazards.
Rev. 8.0	2025 Updated naming conventions. Updated Intended Use statement. Added CH-REP and importer information.

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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The Initial M, TISSUE-TROL, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

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Referenzen

1. Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, S. 61
2. Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2. Auflage, Mosby, St. Louis, (MO), 1980
3. Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
4. Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 20:661, 1950
5. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
6. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
7. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
8. Histotechnology: A Self Instructional Text. Carson, ASCP Press, Chicago, 1990
9. Kok and Boon: Microwave Cookbook for Microscopists. Couloumb Press, Leyden. Leyden 1992

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Références

- Conn's Biological Stains, 10th ed., R.W. Horobin and J.A. Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
- Theory and Practice of Histotechnology, Edited by D.C. Sheehan and B.B. Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
- Lillie R.D.: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- Gomori G.: A rapid one-step trichrome stain. Am J. Clin Pathol 20:661, 1950
- Leong A.S.-Y., Milios J.: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J. Pathol 148:183, 1986
- Brinn N.T.: Rapid metallic histologic staining using the microwave oven. J. Histotechnol 6:125, 1983
- Valle S.: Special stains in the microwave oven. J. Histotechnol 9:237, 1986
- Histotechnology: A Self Instructional Text. Carson, ASCP Press, Chicago, 1990
- Kok and Boon: Microwave Cookbook for Microscopists. Coulomb Press, Leyden. Leyden 1992

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Bibliografia

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
- Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 20:661, 1950
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
- Histotechnology: A Self Instructional Text. Carson, ASCP Press, Chicago, 1990
- Kok and Boon: Microwave Cookbook for Microscopists. Couloomb Press, Leyden. Leyden 1992

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Referencias

1. Conn’s Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
2. Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
3. Lillie RD: Further experiments with the Masson trichrome modification of Mallory’s connective tissue stain. Stain Technol 15:17, 1940
4. Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 20:661, 1950
5. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
6. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
7. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
8. Histotechnology: A Self Instructional Text. Carson, ASCP Press, Chicago, 1990
9. Kok and Boon: Microwave Cookbook for Microscopists. Couloomb Press, Leyden. Leyden 1992

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- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
- Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 20:661, 1950
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
- Histotechnology: A Self Instructional Text. Carson, ASCP Press, Chicago, 1990
- Kok and Boon: Microwave Cookbook for Microscopists. Coulomb Press, Leyden. Leyden 1992

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- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
- Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 20:661, 1950.
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
- Histotechnology: A Self Instructional Text. Carson, ASCP Press, Chicago, 1990
- Kok and Boon: Microwave Cookbook for Microscopists. Coulomb Press, Leyden. Leyden 1992

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Reference

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, str. 61
- Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed., Mosby, St. Louis, (MO), 1980
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 20:661, 1950
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986
- Histotechnology: A Self Instructional Text. Carson, ASCP Press, Chicago, 1990
- Kok and Boon: Microwave Cookbook for Microscopists. Couloomb Press, Leyden. Leyden 1992

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Referanser

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
- Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 20:661, 1950.
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
- Histotechnology: A Self Instructional Text. Carson, ASCP Press, Chicago, 1990
- Kok and Boon: Microwave Cookbook for Microscopists. Couloumb Press, Leyden. Leyden 1992

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Referanslar

- Conn's Biological Stains, 10. baskı, RW Horobin and JA Kiernan, Editörler, Taylor & Francis, NY, 2002, s. 61
- Theory and Practice of Histotechnology, Editörler: DC Sheehan ve BB Hrapchak, 2. baskı, Mosby, St. Louis, (MO), 1980
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 20:661, 1950
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
- Histotechnology: A Self Instructional Text. Carson, ASCP Press, Chicago, 1990
- Kok and Boon: Microwave Cookbook for Microscopists. Couloomb Press, Leyden. Leyden 1992

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Referenties

- 1. Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
- 2. Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
- 3. Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- 4. Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 20:661, 1950.
- 5. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- 6. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- 7. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
- 8. Histotechnology: A Self Instructional Text. Carson, ASCP Press, Chicago, 1990
- 9. Kok and Boon: Microwave Cookbook for Microscopists. Couloomb Press, Leyden. Leyden 1992

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Reference

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
- Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 20:661, 1950
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
- Histotechnology: A Self Instructional Text. Carson, ASCP Press, Chicago, 1990
- Kok and Boon: Microwave Cookbook for Microscopists. Couloumb Press, Leyden. Leyden 1992

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Reference

- Conn's Biological Stains, 10. izd., RW Horobin i JA Kiernan, urednici, Taylor & Francis, NY, 2002., str. 61
- Theory and Practice of Histotechnology, uredili DC Sheehan i BB Hrapchak, 2. izd., Mosby, St. Louis, (MO), 1980.
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940.
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 20:661, 1950.
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986.
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983.
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986.
- Histotechnology: A Self Instructional Text. Carson, ASCP Press, Chicago, 1990.
- Kok i Boon: Microwave Cookbook for Microscopists. Couloumb Press, Leyden. Leyden 1992.

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Literatura

- 1. Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
- 2. Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
- 3. Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- 4. Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 20:661, 1950
- 5. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- 6. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- 7. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
- 8. Histotechnology: A Self Instructional Text. Carson, ASCP Press, Chicago, 1990
- 9. Kok and Boon: Microwave Cookbook for Microscopists. Coulomb Press, Leyden. Leyden 1992

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Atsauces

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, 61. lpp.
- Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 20:661, 1950
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
- Histotechnology: A Self Instructional Text. Carson, ASCP Press, Chicago, 1990
- Kok and Boon: Microwave Cookbook for Microscopists. Couloumb Press, Leyden. Leyden 1992

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MDSS GmbH
MDSS CH GmbH
Chemie GmbH