

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

ProductInformation

Lectin from *Phaseolus vulgaris* Phytohemagglutinin PHA-P

Product Number **L 9017** Storage Temperature 2-8 °C

Product Description

PHA-P is a mixture of PHA-E (MW = 128 kDa) and PHA-L (MW = 126 kDa). This product is affinity purified to remove additional protein content. Purity is determined by SDS-PAGE analysis.

Lectins are proteins or glycoproteins of non-immune origin that agglutinate cells and/or precipitate complex carbohydrates. Lectins are capable of binding glycoproteins even in presence of various detergents.¹ The agglutination activity of these highly specific carbohydrate-binding molecules is usually inhibited by a simple monosaccharide, but for some lectins, di, tri, and even polysaccharides are required.

Lectin PHA-P is not inhibited easily by monosaccharides, but may be inhibited by oligosaccharides.

Lectins are isolated from a wide variety of natural sources, including seeds, plant roots and bark, fungi, bacteria, seaweed and sponges, mollusks, fish eggs, body fluids of invertebrates and lower vertebrates, and from mammalian cell membranes. The precise physiological role of lectins in nature is still unknown, but they have proved to be very valuable in a wide variety of applications *in vitro*, including:

- 1. blood grouping and erythrocyte polyagglutination studies.
- 2. mitogenic stimulation of lymphocytes.
- 3. lymphocyte subpopulation studies.
- 4. fractionation of cells and other particles.
- 5. histochemical studies of normal and pathological conditions.

Sigma offers a range of lectins suitable for the above applications. Most Sigma lectins are highly purified by affinity chromatography, but some are offered as purified or partially purified lectins, suitable for specific applications. Many of the lectins are available conjugated to (conjugation does not alter the specificity of the lectin):

- 1. fluorochromes (for detection by fluorimetry).
- 2. enzymes (for enzyme-linked assays).
- 3. insoluble matrices (for use as affinity media).

Please refer to the table for general information on the most common lectins.

Procedure

A general agglutination procedure using this lectin with 96 well plates is as follows:

- 1. Prepare a lectin solution of 1 mg/ml in PBS buffer, pH 6.8.
- 2. Pipette 50 μ l of fresh PBS into each well and add 50 μ l of the lectin solution into the first well.
- 3. Serial dilutions are made by pipetting 50 μl from each successive well into the next well.
- 4. Blood type A with a 2% hematocrit is used as the substrate.
- 5. Pipette 50 µl of blood into each well.
- 6. Visually determine agglutination.

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This lectin is soluble in phosphate buffered saline, pH 7.2 (1 mg/ml).

Storage/Stability

Aggregation is thought to occur in the presence of high concentrations of 2-mercaptoethanol.

| Lectin | MW (kDa) | Subunits | Specificity Blood Group Sugar | | Mitogenic Activity |
|---|----------------------|------------------|----------------------------------|---------------------------------|-----------------------|
| Abrus precatorius | | | _ | | + |
| Agglutinin | 134 | 4 | | gal | |
| Abrin A (toxin) | 60 | 2 | | gal | |
| Abrin B (toxin) | 63.8 | 2(αβ) | | gal | |
| Agarius bisporus | 58.5 | _ | _ | β-gal(1→3)galNAc | |
| Anguilla anguilla | 40 | 2 | Н | α-L-Fuc | |
| Arachis hypogaea | 120 | 4 | Т | β-gal(1→3)galNAc | |
| Artocarpus integrifolia Bandeiraea simplicifolia | 42 | 4 | T | α-gal→OMe | + |
| BS-I | 114 | 4 | А, В | α-gal, α-galNAc | |
| BS-I-A ₄ | 114 | 4 | A | α -galNAc | |
| BS-I-B ₄ | 114 | 4 | В | α-gal | |
| BS-II | 113 | 4 | acq, B, Tk, T | glcNAc | |
| Bauhinia purpurea | 195 | 4 | ас ч , b, тк, т – | β-gal(1→3)galNAc | + |
| | 60; 120 ^a | 4 2/4 | - | | т |
| Caragana arborescens | 60; 120 44 | | - | galNAc fotuin | |
| Cicer arietinum | | 2 | _ | fetuin | |
| Codium fragile | 60 | 4 | - | galNAc | |
| Concanavalin A | 102 | 4 | - | α-man, α-glc | + |
| Succinyl-Concanavalin A | 51 | 2 | - | α -man, α -glc | + ^b |
| Cytisus scoparius | - | - | - | galNAc, gal | |
| Datura stramonium | 86 | 2(αβ) | - | (glcNAc) ₂ | |
| Dolichos biflorus | 140 | 4 | A ₁ | α-galNAc | |
| Erythrina corallodendron | 60 | 2 | - | β-gal(1→4)glcNAc | + |
| Erythrina cristagalli | 56.8 | 2(αβ) | - | β-gal(1→4)glcNAc | |
| Euonymus europaeus | 166 | 4(αβ) | В, Н | α-gal(1→3)gal | + |
| Galanthus nivalis | 52 | 4 | (h) | non-reduc. α-man | |
| Glycine max | 110 | 4 | _ | galNAc | + ^c |
| lelix aspersa | 79 | _ | А | galNAc | |
| lelix pomatia | 79 | 6 | А | galNAc | |
| .athyrus odoratus | 40-43 | 4(αβ) | _ | α-man | + |
| ens culinaris | 49 | 2 | _ | α-man | + |
| imulus polyphemus | 400 | _ 18 | _ | NeuNAc | |
| Bacterial agglutinin | - | _ | _ | galNAc, glcNAc | |
| ycopersicon esculentum | 71 | _ | _ | (glcNAc) ₃ | |
| laackia amurensis | 130 | 2(αβ) | 0 | sialic acid | |
| | 40-43 | | 0 | | Ŧ |
| Naclura pomifera | | $2(\alpha\beta)$ | - | α -gal, α -galNAc | |
| <i>Iomordica charantia</i> | 115-129 | 4(αβ) | - | gal, galNAc | |
| laja mocambique mocambique | | - | - | - | |
| laja naja kaouthia | - | _ | — | - | |
| larcissus pseudonarcissus | 26 | 2 | (h) | α-D-man | |
| Perseau americana | - | - | - | - | |
| Phaseolus coccineus | 112 | 4 | - | - | |
| Phaseolus limensis | 247(II) 124(III) | 8 4 | A | galNAc | + |
| Phaseolus vulgaris | | | | | |
| PHA-E | 128 | 4 | - | oligosaccharide | + |
| PHA-L PHA-P | 128 | 4 | - | oligosaccharide | + |
| PHA-M | | | | | |

| Lectin | MW (kDa) | Specificity Subunits Blood Group Sugar | | | Mitogenic Activity |
|-----------------------------|----------------------|---|----------------|--------------------------------|-----------------------|
| Phytolacca americana | 32 | _ | | (glcNAc) ₃ | + |
| Pisum sativum | 49 | 4(αβ) | _ | α-man | + |
| Pseudomonas aeruginosa PA-I | | - | _ | gal | + ^c |
| Psophocarpus tetragonolobus | 35 | 1 | _ | galNAc, gal | • |
| Ptilota plumosa | 65; 170 | _ | В | α -gal | |
| Ricinus communis | 00, 110 | | 2 | a gai | |
| Toxin, RCA ₆₀ | 60 | 2 | _ | galNAc, β-gal | |
| Toxin, RCA ₁₂₀ | 120 | 4 | _ | β-gal | |
| Sambucus nigra | 140 | 4(αβ) | _ | αNeuNAC(2→6)gal | + ^c |
| 3 | | (-1-) | | galNAc | |
| Solanum tuberosum | 50; 100 ^a | 1, 2 | _ | (glcNAc) ₃ | |
| Sophora japonica | 133 | 4 | А, В | β-galNAc | |
| Tetragonolobus purpureas | 120(A) | 4 | Н | α-L-fuc | |
| c | 58(BA) | 2 | Н | α-L-fuc | |
| | 117(C) | 4 | Н | α-L-fuc | |
| Triticum vulgaris | 36 | 2 | _ | (glcNAc) ₂ , NeuNAc | + |
| Ulex europaeus | | | | | |
| UEA I | 68 | _ | Н | α-L-fuc | |
| UEA II | 68 | _ | - | (glcNAc) ₂ | |
| Vicia faba | 50 | 4(αβ) | - | man, glc | + |
| Vicia sativa | 40 | 4(αβ) | - | glc, man | + |
| Vicia villosa | 139 | 4 | $A_{1+}T_n$ | galNAc | |
| A_4 | 134 | 4 | A ₁ | galNAc | |
| B ₄ | 143 | 4 | T _n | galNAc | |
| Vigna radiata | 160 | 4 | - | α-gal | |
| Viscum album | 115 | 4(αβ) | - | β-gal | |
| Wisteria floribunda | 68 | 2 | _ | galNAc | |

^a Concentration-dependent molecular weight

^b Non-agglutinating and mitogenic

[°] Mitogenic for neuraminidase-treated lymphocytes

References

 Rueben, L., et al., Activities of lectins and their immobilized derivatives in detergent solutions. Implications on the use of lectin affinity chromatography for the purification of membrane glycoproteins. Biochemistry, **16**, 1787-1794 (1977).

IRB/MWM/JRC/NSB/SAG 11/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.