

# Novabiochem® Peptide Delta Masses

## Amino-acid substitutions

-129.16	Trp to Gly	-42.09	Arg to Asn	-18.04	Met to Leu
-115.13	Trp to Ala	-42.08	Leu to Ala	-17.97	Asp to Pro
-106.12	Tyr to Gly	-42.08	Ile to Ala	-17.10	Met to Asn
-99.14	Arg to Gly	-42.08	Val to Gly	-16.98	Asn to Pro
-99.13	Trp to Ser	-42.04	Glu to Ser	-16.11	Met to Asp
-92.09	Tyr to Ala	-41.10	Arg to Asp	-16.06	Cys to Ser
-90.13	Phe to Gly	-41.09	Lys to Ser	-16.04	Ile to Pro
-89.09	Trp to Pro	-41.05	Gln to Ser	-16.04	Leu to Pro
-87.08	Trp to Val	-40.07	Pro to Gly	-16.00	Ser to Ala
-85.10	Trp to Thr	-40.02	His to Pro	-15.99	Tyr to Phe
-85.11	Arg to Ala	-39.03	Trp to Phe	-15.98	Phe to Met
-83.07	Trp to Cys	-38.01	His to Val	-15.96	Asp to Val
-80.09	His to Gly	-36.03	His to Thr	-15.96	Glu to Ile
-76.10	Phe to Ala	-35.04	Tyr to Gln	-15.96	Glu to Leu
-76.09	Tyr to Ser	-35.00	Tyr to Lys	-15.02	Glu to Asn
-74.15	Met to Gly	-34.08	Met to Pro	-15.01	Lys to Ile
-73.05	Trp to Ile	-34.05	Tyr to Ile	-15.01	Lys to Leu
-73.05	Trp to Leu	-34.02	Phe to Ile	-14.97	Gln to Ile
-72.11	Trp to Asn	-34.02	Phe to Leu	-14.97	Gln to Leu
-72.07	Glu to Gly	-34.00	His to Cys	-14.97	Asn to Val
-71.12	Trp to Asp	-33.08	Phe to Asn	-14.07	Lys to Asn
-71.12	Lys to Gly	-32.09	Phe to Asp	-14.03	Gln to Asn
-71.08	Gln to Gly	-32.07	Met to Val	-14.03	Ile to Val
-69.11	Arg to Ser	-32.06	Cys to Ala	-14.03	Leu to Val
-66.06	His to Ala	-32.00	Glu to Pro	-14.03	Ala to Gly
-66.05	Tyr to Pro	-31.97	Tyr to Met	-14.03	Glu to Asp
-64.04	Tyr to Val	-31.05	Lys to Pro	-14.03	Thr to Ser
-62.06	Tyr to Thr	-31.01	Gln to Pro	-13.98	Asp to Thr
-60.12	Met to Ala	-30.09	Met to Thr	-13.08	Lys to Asp
-60.10	Phe to Ser	-30.03	Ser to Gly	-13.04	Gln to Asp
-60.03	Tyr to Cys	-30.03	Thr to Ala	-12.99	Asn to Thr
-59.07	Arg to Pro	-30.02	Trp to Arg	-12.06	Ile to Thr
-58.08	Trp to Gln	-29.99	Glu to Val	-12.05	Leu to Thr
-58.04	Trp to Lys	-29.04	Lys to Val	-12.05	Val to Ser
-58.04	Asp to Gly	-29.00	Gln to Val	-11.95	Asp to Cys
-58.04	Glu to Ala	-28.06	Arg to Gln	-10.96	Asn to Cys
-57.09	Trp to Glu	-28.06	Met to Cys	-10.04	Phe to His
-57.09	Lys to Ala	-28.05	Val to Ala	-10.04	Pro to Ser
-57.06	Arg to Val	-28.02	Arg to Lys	-10.02	Ile to Cys
-57.05	Gln to Ala	-28.01	Glu to Thr	-10.02	Leu to Cys
-57.05	Asn to Gly	-28.01	Asp to Ser	-9.01	Arg to Phe
-56.11	Ile to Gly	-27.07	Arg to Glu	-9.01	His to Gln
-56.11	Leu to Gly	-27.06	Lys to Thr	-8.97	His to Lys
-55.08	Arg to Thr	-27.02	Gln to Thr	-8.02	His to Glu
-55.01	Trp to Met	-27.02	Asn to Ser	-6.98	Tyr to Arg
-53.05	Arg to Cys	-26.08	Ile to Ser	-6.02	Cys to Pro
-50.06	Phe to Pro	-26.08	Leu to Ser	-5.94	His to Met
-50.06	His to Ser	-26.04	Pro to Ala	-4.01	Cys to Val
-50.01	Tyr to Ile	-26.03	Tyr to His	-3.99	Gln to Pro
-50.01	Tyr to Leu	-25.98	Glu to Cys	-3.07	Met to Gln
-49.07	Trp to His	-25.03	Lys to Cys	-3.03	Met to Lys
-49.07	Tyr to Asn	-24.99	Arg to Met	-2.08	Met to Glu
-48.08	Tyr to Asp	-24.99	Gln to Cys	-2.03	Cys to Thr
-48.05	Phe to Val	-23.98	His to Ile	-2.01	Val to Pro
-46.09	Cys to Gly	-23.98	His to Leu	-1.98	Thr to Val
-46.07	Phe to Thr	-23.04	His to Asn	-1.93	Asp to Ile
-44.12	Met to Ser	-23.04	Trp to Tyr	-1.93	Asp to Leu
-44.06	Thr to Gly	-22.05	His to Asp	-0.99	Asp to Asn
-44.04	Phe to Cys	-19.05	Arg to His	-0.99	Glu to Lys
-44.01	Asp to Ala	-19.05	Arg to Gln	-0.95	Glu to Tyr
-43.03	Arg to Ile	-19.01	Phe to Lys	-0.94	Asn to Ile
-43.03	Arg to Leu	-18.06	Phe to Glu	-0.94	Asn to Leu
-43.02	Asn to Ala	-18.04	Met to Ile	-0.04	Lys to Gln

## Single amino-acid and dipeptide deletions

-372	TrpTrp	-248	PheThr	-206	CysCys
-349	TrpTyr	-246	AspMet	-204	CysThr
-342	ArgTrp	-246	PheVal	-204	GlyPhe
-333	PheTrp	-245	AsnMet	-202	AlaMet
-326	TyrTrp	-244	AspGlu	-202	AspSer
-323	HisTrp	-244	IleMet	-202	CysVal
-319	ArgTyr	-244	LeuMet	-202	ThrThr
-317	MetTrp	-244	PhePro	-201	AsnSer
-315	GluTrp	-243	pTyr	-200	AlaGlu
-314	GlnTrp	-243	Tyr(SO <sub>3</sub> H)	-200	CysPro
-314	LysTrp	-243	ArgSer	-200	IleSer
-312	ArgArg	-243	AsnGlu	-200	LeuSer
-310	PheTyr	-243	AspGln	-200	ThrVal
-303	ArgPhe	-243	AspLys	-199	AlaGln
-301	AspTrp	-243	GlyTrp	-199	AlaLys
-300	AsnTrp	-242	AsnGln	-198	ProThr
-300	HisTyr	-242	AsnLys	-198	ValVal
-299	IleTrp	-242	GluIle	-196	ProVal
-299	LeuTrp	-242	GluLeu	-194	GlyHis
-294	MetTyr	-241	GlnIle	-194	ProPro
-294	PhePhe	-241	GlnLeu	-190	CysSer
-293	ArgHis	-241	IleLys	-188	GlyMet
-292	GluTyr	-241	LeuLys	-188	SerThr
-291	GlnTyr	-240	CysHis	-186	Trp
-291	LysTyr	-238	HisThr	-186	AlaAsp
-289	CysTrp	-236	HisVal	-186	GluGly
-287	ArgMet	-234	AlaTyr	-186	SerVal
-287	ThrTrp	-234	CysMet	-185	AlaAsn
-285	ArgGlu	-234	HisPro	-185	GlnGly
-285	TrpVal	-234	PheSer	-185	GlyLys
-284	ArgGln	-232	CysGlu	-184	Arg(Me <sub>2</sub> )
-284	ArgLys	-232	MetThr	-184	Alalle
-284	HisPhe	-231	CysGln	-184	AlaLeu
-283	ProTrp	-231	CysLys	-184	ProSer
-278	AspTrp	-230	AspAsp	-181	pThr
-278	MetPhe	-230	GluThr	-174	AlaCys
-277	AsnTyr	-230	MetVal	-174	SerSer
-276	GluPhe	-229	AsnAsp	-172	AlaThr
-276	IleTyr	-229	GlnThr	-172	AspGly
-276	LeuTyr	-229	LysThr	-171	AsnGly
-275	GlnPhe	-228	AsnAsn	-171	Lys(Me <sub>2</sub> )
-275	LysPhe	-228	AspIle	-170	Arg(Me)
-274	HisHis	-228	AspLeu	-170	AlaVal
-273	SerTrp	-228	GluVal	-170	GlyIle
-271	ArgAsp	-228	MetPro	-170	GlyLeu
-270	ArgAsn	-227	AlaArg	-168	AlaPro
-269	ArgIle	-227	AsnIle	-167	pSer
-269	ArgLeu	-227	AsnLeu	-163	Tyr
-268	HisMet	-227	GlnVal	-160	CysGly
-266	CysTyr	-227	LysVal	-158	AlaSer
-266	GluHis	-226	GluPro	-158	GlyThr
-265	GlnHis	-226	IleIle	-157	Cit
-265	HisLys	-226	IleLeu	-156	Arg
-264	ThrTyr	-226	LeuLeu	-156	GlyVal
-262	AspPhe	-225	GlnPro	-156	Lys(Me <sub>2</sub> )
-262	MetMet	-225	LysPro	-154	GlyPro
-262	TyrVal	-224	HisSer	-153	Cha
-261	AsnPhe	-220	GlyTyr	-147	Phe
-260	GluMet	-218	AlaPhe	-144	GlySer
-260	IlePhe	-218	AspCys	-142	Lys(Me)
-260	LeuPhe	-218	MetSer	-142	AlaAla
-260	ProTyr	-217	AsnCys	-137	His
-259	ArgCys	-216	AspThr	-131	Met
-259	GlnMet	-216	CysIle	-129	Glu
-259	LysMet	-216	CysLeu	-128	Gln, Lys
-258	GluGlu	-216	GluSer	-128	AlaGly
-257	AlaTrp	-215	AsnThr	-115	Asp
-257	ArgThr	-215	GlnSer	-114	Asn, Orn
-257	GluLys	-215	LysSer	-114	GlnGly
-256	GlnGln	-214	AspVal	-113	Ile, Leu, Hyp, Nle, ε-Ahx
-256	GlnLys	-214	IleThr	-103	Cys
-256	LysLys	-214	LeuThr	-101	hSer, Thr
-255	ArgVal	-213	ArgGly	-99	Val
-253	ArgPro	-213	AsnVal	-97	Pro
-252	AspHis	-212	AspPro	-95	dPro
-251	AsnHis	-212	IleVal	-87	Ser
-250	CysPhe	-212	LeuVal	-85	Abu, γ-Abu, Aib
-250	HisIle	-211	AsnPro	-71	Ala, β-Ala, Sar
-250	HisLeu	-210	IlePro	-57	Gly
-250	SerTyr	-210	LeuPro		
		-208	AlaHis		

## Adducts, modifications, and protecting groups

-98	H <sub>3</sub> PO <sub>4</sub> ; artifact observed during MS of pSer and pThr peptides	82	Dmcp
-80	Dephosphorylation: avoid prolonged storage of phosphopeptides in acid aqueous media	84	Mpe
-42	Orn formation from Arg	86	Bum
-34	Dehydroalanine from Cys	88	StBu
-32	Lanthionine formation from cystine	90	Bzl; uncleaved phosphoamino acid protecting group or alkylation adduct of Cys, Trp or Tyr.
-18	Pyr formation from Glu; avoid using benzyl or allyl protection for N-terminal Glu residues	91	Bzl alkylation of Met.
-18	Dehydration (-H <sub>2</sub> O)	92	EDT disulfide adduct; reduce with TCEP
-18	Ser to dehydroalanine	96	Tfa; treat peptide with 0.1 M ammonium bicarbonate to cleave Tfa esters
-18	Cyanoalanine formation from Asn; use side-chain protected Asn derivatives	97	Pro insertion; avoid double coupling of Pro
-18	Aspartimide formation from Asp; use Fmoc-Asp(OtBu)-(Dmb)Gly-OH or Fmoc-Asp(OMpe)-OH	98	Tetramethylguanidium termination by-product on amine; use phosphonium reagents instead of uronium
-18	Pyrophosphate formation with peptide containing 2 adjacent phosphoresidues	100	Boc
-17	Pyr formed from Gln; avoid N-terminal Glu residues	104	MeBzl
-17	Aspartimide formation from Asn	104	Bz
-16	H-Phosphonate formation if a phosphopeptide; use benzyl protected phosphoramidites and anhydrous t-butyl hydroperoxide for oxidation	106	HMP adduct; use an HMPA resin, such as NovaSyn TGA instead of a Wang resin, or Rink Amide AM resin instead of Rink Amide resin
-13	N-Piperidinylalanine formation from pSer or pThr; avoid using microwave heating for Fmoc removal	111	Pyr insertion; use Pyr-OPcp for Pyr introduction
-2	Disulfide bond formation from cystine	117	Succinyl
-1	Amide formation (C terminus)	118	PhiPr
1	Deamidation of Asn and Gln to Asp and Glu; avoid leaving peptide in aqueous media for prolonged periods	120	MeOBzl
1	Cit formation from Arg	134	Z
2	Reduction of cystine	134	Ada
2	Reduction of indole double bond of Trp; use Trp(Boc) and avoid use of TES	136	Hmb
12	Cys to thioproline conversion, due to formaldehyde contamination of acetic acid	150	Dmb
12	Formaldehyde adduct of Trp	152	DDT disulfide adduct
14	Methyl ester; ensure peptide resin is free from methanol before cleavage	153	Nps
16	Hydrazide formation from free acid	154	Tos
16	Oxidation of Met to sulfoxide	154	Npys
16	Sulfenic Acid from Cys	163	Methylphenoxyacetamide - adduct of Knorr linker; use Rink Amide AM resin instead of Rink Amide MBHA resin
17	Ammonium adduct	166	Dnp
22	Sodium adduct	172	Ethanedithiol/TFA cyclic adduct; Use Trp(Boc) instead of Trp and avoid prolonged cleavage times when EDT is used as scavenger
28	Ethyl	178	AcHmb
28	Formyl; use good quality DMF and avoid heating peptides dissolved in formic acid	180	Tmob
32	Oxidation of Met to sulfone	180	Xan
38	Potassium adduct	182	Mts
40	Uncleaved pseudoproline; extend cleavage and add EDT	188	Alkylation of Trp by Pbf chroman
41	MeCN adduct of Lys	189	α-Cyanocinnamic acid adduct
42	Ac	201	HMP (hydroxymethylphenyl)/TFA adduct
44	Disodium adduct	204	S-Farnesyl
44	CO <sub>2</sub> adduct of Trp; dissolve in 1% AcOH aq. and freeze dry.	206	iVdde
45	Nitro	212	Mtr
48	Cysteic acid formation by oxidation of Cys	216	Mca
51	β-Piperidinylalanine formation from C-terminal Cys; use trityl resin such as 2-chlorotrityl chloride resin or NovaSyn TGT	222	Fmoc
56	tBu; alkylation adduct of Cys, Met, Trp, Tyr	226	Biotinyl
57	Gly insertion; avoid Gly symmetrical anhydride	226	Mbh
60	Sodium-potassium adduct	233	Dansyl
67	Asp piperidide formation; use Fmoc-Asp(OtBu)-(Dmb)Gly-OH or Fmoc-Asp(OMpe)-OH	238	Bpoc
71	Acm	238	Palmitoylation
76	Mercaptoethanol adduct	242	Trt
80	Sulfonation (SO <sub>3</sub> H) (Arg, Ser, Thr from Pbf or Pmc); recombine with TMSBr/thioanisole/EDT/TFA.	243	Rink amide linker adduct; use Rink Amide AM resin.
80	Sulfation of Tyr	251	Dabicyl
80	Phosphorylation	252	Pbf
82	cHex	256	Mtt
		264	EDANS
		266	Stearoylation
		266	Pmc
		268	Biotinyl-aminoethylamino (from Biotin NovaTag) addition to mass of peptide acid
		272	Mmt
		299	Knorr linker adduct; use Rink Amide AM resin instead of Rink Amide MBHA
		304	Glutathione disulfide adduct
		312	Dmab
		358	FAM
		412	TAMRA
		428	Biotinyl-PEG (from Biotin PEG NovaTag resin addition to mass of peptide acid