

## Technical Bulletin

# L-Glutamate Quick Test Strips

**Catalog Number MAS004**

## Product Description

Glutamate is an important chemical in general metabolism. It is also widely used as a flavor enhancer in the food industry. Consumption of food containing monosodium glutamate (MSG) has been known to cause headaches, flushing, sweating, nausea, and more.

The L-Glutamate Quick Test Strips are based on the glutamate dehydrogenase-catalyzed oxidation of L-glutamate in which the formed NADH reduces a chromogenic reagent. The intensity of the product color is directly proportional to the glutamate concentration in the sample. The semi-quantitative detection range of the kit is 0-300 mg/L (undiluted) L-glutamate.

This kit is suitable for the detection of L-glutamate in a variety of food and beverage samples.

## Components

The kit is sufficient for 10 tests

- L-Glutamate Test Strips (10 strips)  
Catalog Number MAS004A 1 Each
- Sample Development Tubes (400 µL of Development Reagent per tube)  
Catalog Number MAS004B 10 Each

## Equipment Required but Not Provided

- Pipetting devices and accessories

## Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

## Storage/Stability

The kit is shipped at room temperature. Store components at 2-8 °C. Keep strips dry and out of direct sunlight.

## Preparation Instructions

### Sample Preparation

For most instant soups and acidic samples, we strongly recommend diluting samples 21-fold dilution. Most soy sauces and fish sauces will require a 210-fold dilution. Milk and samples that are not expected to have very high levels of L-glutamate should be diluted 5-fold.

For samples requiring a 210-fold dilution, first make a 10-fold dilution by carefully transferring 100  $\mu\text{L}$  of Sample into 900  $\mu\text{L}$  of purified water. Mix dilution thoroughly and use for Procedure.

### Procedure

1. Unscrew the cap of one of the Sample Development Tubes and add the Sample volume required according to Table 1.

**Table 1.**  
Sample Dilutions

Final Dilution	Volume of Sample to Add to Tube
210-fold	20 $\mu\text{L}$ of 10-fold dilution of Sample
21-fold	20 $\mu\text{L}$
5-fold	100 $\mu\text{L}$

2. Replace cap on the Sample Development Tube, securely close the tube, and invert the tube 3-4 times to mix the diluted Sample.
3. Unscrew cap and dip in one of the L-Glutamate Test Strips, making sure to fully submerge the reaction pad at the end of the strip. Leave the strip submerged for 5 seconds and then remove the strip and shake gently to remove any excess liquid.
4. Allow the color to develop on the strip for 5 minutes.
5. Compare the color of the reaction pad of the strip with the provided L-Glutamate Chart. Multiply the concentration on the chart by the Sample dilution factor (i.e., 5, 21, or 210) to determine the concentration of L-Glutamate in the original Sample.



## Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

## Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://SigmaAldrich.com/offices).

## Technical Service

Visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice).

## Standard Warranty

The applicable warranty for the products listed in this publication may be found at [SigmaAldrich.com/terms](https://SigmaAldrich.com/terms).