Poster # T3604

Abstract

Purpose To develop cellular invasion assays, using cancer cells, in a 96 well format for application in cancer drug discovery and in *in vitro* cytotoxicity assays. Methods

drug discovery and in in three optionizity assys. Method
Tumor cell invasion assays were performed using polycarbonate hottomed 96 well plates
(MultiScene^{m-Au}RC). MDAME231, a highly invasion brassy conditions. The cells were
contract in PRIAC. MDAME231, a highly invasion brassy conditions. The cells were
contract in PRIAC MODAME251, a highly invasion brassy conditions. The cells were
contract in PRIAC MODAME251, a highly invasion brassy conditions. The cells were
to a set the descent of the transmission brassy conditions. The cells were
the prevent wells. Conditioned medium from NH1571 Biroblasts or serum containing medium was added to
the lower wells as chemonattractum. The ceffects of cell density, ECM concentration and invasion time
were evaluated on the MultiScenem^{Au}MIC plates with 3, 5 and 8 µm pore size polycarbonate
membrane. Set lists that ali winded through the ECM to the underside of the merbrane. Were stained
using a Hema-3 stain kit. Cells were enumerated microscopically over a portion (10-1594) of each
membrane as the basis for determining the total number of invading cells per vell. Fluorescent assay,
using DNA binding fluorescent probes, were also performed to quantify the invaded cells.
Result
Cell density at 1.65x 10¹/cm², ECM concentration at 30 µg/well and 72 hrs incubation at 37C were
optimal for invasion to occur. The MDAME323 cell mysion index, calculated relative to the noninvasive MCT² cells under similar conditions, was 5 fold or higher. Visual enumeration provided the noninvasive fuel cells were ensitive to only 1000 cells or more.
During the provide the cells were sensitive to only 1000 cells or more.
During the provide the fuel were sensitive to only 1000 cells or more.
During the provide the terming the of the outperformed to resensitive detection with a lower of 96 well formatin the transmiter fuel time of the outperformater
During the provide the terming the outperformater terming the outperformater
During the provide the terming the outperformateresent assay, alth

Conclusions Our results den cancer drugs. nonstrate the utility of a 96 well format in high throughput screening (HTS) of anti-

Introduction

Introduction Cell-based assays are gaining tremendous importance in pre-screening of compounds in order to generate target leads for the Absorption, Distribution, Metabolism, Excretion and Toxicity assays (ADMETox). Cancer drug discovery efforts have increasingly focused on incorporating functional cell-based assays in the pre-screening stages of lead compounds. Many drugs under development are directed at altering the migration and invasion properties of cancer cells. Development of HTS cell-based assays that are designed to be able to measure the migration, chemotaxis and invasion potential of cancer cells are therefore of interest. Chemotaxis is defined as the movement of cells in response to a concentration gradient set up by a chemoattractant. Invasion is defined as the chemotaxis of cancer cells across an extracellular matrix barrier. We have developed a 96 well plate to support these HTS migration, Invasion and Chemotaxis (MIC) assays. The data presented in this poster indicate that MultiScreen^m-MIC plates are useful for anti-cancer drug HTS and can serve as an alternative to existing, lower-throughput products in the market.

Methods

Invasion assays. Invasion assays were set up on Millipore MultiScreen⁷⁰⁴-MIC plates as described by Albini et al and Kannah et al. (1,2). MDAME231 and MCP7 (invasive and non-invasive adherent mammary admontory for an exist, regarding the set of the form ATCC and your by unlined in RPMI to 1480 were groups to 190% confluency and fastered overheading in serum. Free medium containing 0.1% BAS prior to initiation of invasion assays. The upper wells were coated with varying concentrations of Muriged⁴ - another derived extracted mutrix (ECK observed for two hours. Conditioned medium derived from MIH373 cells (passage 34 & 61 S2) (most for both sets. Obtained from ATCC) or serum-containing medium was added to the bottom wells. Invasion was allowed to occur evel 48 and 72 hr time periods. The non-invade cells were emoved from the upper wells by swiping with coton swabs followed by insing twice with PBS.

Table 1. Parameters Evaluated.

Invasior

edium

abeling

8, 5 and 3 µm MDAMB231 and MCF7 cells 25,000 to 100,000

cells Conditioned medium and serum containing

20, 30 and 40 µg/well 48 and 72 hrs

labeling and DNA binding fluorescent

copy imaging Whole cell fluore:

25

С

1

新闻是这个句子。

"ast as the

Assay evaluated

ore size Cell lines tested

emoattractants

cubation period

method 2


CM concentra

Cell number

ype of

Potection Assys Following the removal of cells from the upper wells, cells that had invaded to the membrane underside were emanerated by microscopy or fluorescent (FL) assay. For microscopic enumeration, membranes were stained with Hema-3 stain kit (Fisher Scientific). Stained cells were imaged using a Fuji Fingers ST to objetical carmer anomed on an Otympus BH-2 microscope and counted using the Optimas version 6.1 software. The cells in 16-15% of the stafface area of cash membrane well were counted and the number of (mvading cells in the entire for the objetical counter of the stafface area of cash member of the stafface area of cash for the provide the stafface area of cash were constrained and the number of (mvading cells in the entire forumber of cells invaded in the presence of matrigedrumber of cells imgrated in the absence of matricefs) through the DMAMB231 cells relative to the non-invasive MCF7 cells under similar conditions.

For fluorescent (FL) analysis, the cells were detuched from the membranes with PBS-EDTA or Trypsis-EDTA and quantified by whole cell fluorescent labeling (BCCF, Molecular Probes) or DNA binding (Inverseent probe (VO-RO-10 Idold, Molecular Probes). The number of invaded cells was calculated from standard curves generated for the FL probes with the cell lines.



Development of Cancer Cell Invasion Assays in a 96 well Format

Lakshmi Kamath, David Bell and Jeanne Phillips Millipore Corporation, 17 Cherry Hill Drive, Danvers, MA 01923 Phone: 978 762 5373, Email: lakshmi_kamath@millipore.com

4000 3500 2500 1500 500 0 0 and 5 µm MIC plates. 25,000 cells were seeded in . . and 5 µm MIC plates: 25,000 cells were seeded in wells coated with 30 µg of extracellular matrix. Conditioned medium was applied as a chemoattractant and plates were incubated for 72 hrs (n=2, r=2 or 4). Detection technique was whole cell labeling with BCECF. Representative experiment illustrated in figure. 40000 60000 80000 100000 12000 Cell number MDAMB231 YOPRO STANDARD CURVE MCF7 YOPRO STANDARD CURVE INVASION ON 8 AND 5 µM MIC PLATES 7 6 5 4 3 2 1 0 100000 V Ι 12000 80000 80000 80000 60000 400 2000 20000 40000 60000 80000 100000 120000 8 micror 40000 60000 80000 100000 12 0 20000 Cell number Cell number Figure 5. Percent basal invasion of MDAMB231 (invasive breast cancer cell line) in comparison to MCF7 (non-invasive breast cancer cell line) using 8 and 5 µm MIC plates. 50,000 cells were seeded in wells coated with 30 µg of extracellular matrix. Conditioned medium was applied as a chemoattractant and plates were incubated for 72 hrs. (n=2, r=2 or 4). Detection technique was DNA labeling with YOPRO. Representative experiment illustrated in figure. 8. Conclusions

MDAMB231 BCECF STANDARD CURVE (LOW)

Cell number

MCF7 BCECF STANDARD CURVE

8000

¹⁰

FL units

1. J. W.

24

and the second second second

" of " & S &

ting.

Figure 3. MDAMB231 (Panel A & B) and MCF7 cells (Panel C & D) migrated or invaded to the underside of an 8 µm MIC plate membrane.

Star Paula.

D

MDAMB231 BCECF STANDARD CURVE (HIGI

20000 40000 60000 80000 100000

Figure 4. Percent basal invasion of MDAMB231 (invasive breast cancer cell line) in comparison to MCF7 (non-invasive breast cancer cell line) using

Invasion assays were successfully optimized in a 96 well format using the MultiScreen¹⁰⁻MIC (Figure 1) plate. Several parameters were tested (Table 1) during protocol development.

The 8 µm pore size plates were most suitable for invasion assays for the cell lines tested compared to 5 µm or 3 µm [data not shown and unpublished work (3)].

Invasion index of MDAMB231 was typically two fold or more than the MCF7 cells Visual Microscopy and Imaging was the most sensitive method to detect invaded cells. The lower detection sensitivity of the fluorescent assays could be attributed to additional processing steps involved. Protocol development to optimize fluorescent detection method is in progress.

Our results demonstrate that the MIC plates are highly applicable to invasion assays and lend themselves to HTS of such assays with variety of other cancer cells.

INVASION ON 8 AND 5 µM MIC PLATES

Plate type

5 micron

MDAMB231

5 micron

MCF7

Plate type

MCF7

8 micron

References (1) Albini et al. A rapid in vitro assay for quantitating the invasive potential of tumor cells. Cancer Res. 47: 3239-3245, 1987. (2) Kamath et al. Signaling from Protease-activated Receptor-1 Inhibits Migration and invasion of Breast Cancer Cells. Cancer Res. 61, 5933-5940. 2001. (3) Kamath et al. Development of Cancer Cell Migration Assays in a 96 well Format. Unpublished work.

Acknowledgment Thanks to the Imaging facility at Millipore Corporation, Bedford for use of image analysis equipment.