

For life science research only.  
Not for use in diagnostic procedures.



# DIG-High Prime

 **Version: 10**

Content Version: November 2020

For the nonradioactive labeling of DNA with DIG-11-dUTP, alkali-labile using random oligonucleotides as primers.

**Cat. No. 11 585 606 910**    160 µl  
40 labeling assays

**Store the product at –15 to –25°C.**

<b>1.</b>	<b>General Information</b> .....	<b>3</b>
1.1.	Contents.....	3
1.2.	Storage and Stability.....	3
	Storage Conditions (Product).....	3
1.3.	Additional Equipment and Reagent required.....	3
1.4.	Application.....	3
<b>2.</b>	<b>How to Use this Product</b> .....	<b>4</b>
2.1.	Before you Begin.....	4
	Sample Materials.....	4
	General Considerations.....	4
	Template DNA requirements.....	4
	Safety Information.....	4
	For customers in the European Economic Area.....	4
	Working Solution.....	4
2.2.	Protocols.....	5
	Random primed DNA labeling.....	5
	Labeling DNA isolated from agarose.....	5
	Labeling reaction yield.....	5
	Detection of DIG-labeled DNA.....	6
	Semi-quantitative determination of labeling efficiency.....	6
2.3.	Parameters.....	7
	Chemical Name.....	7
	Structural formula.....	7
<b>3.</b>	<b>Additional Information on this Product</b> .....	<b>7</b>
3.1.	Test Principle.....	7
<b>4.</b>	<b>Supplementary Information</b> .....	<b>8</b>
4.1.	Conventions.....	8
4.2.	Changes to previous version.....	8
4.3.	Ordering Information.....	8
4.4.	Trademarks.....	9
4.5.	License Disclaimer.....	9
4.6.	Regulatory Disclaimer.....	9
4.7.	Safety Data Sheet.....	9
4.8.	Contact and Support.....	9

# 1. General Information

## 1.1. Contents

Vial / Bottle	Label	Function / Description	Content
1	DIG-High Prime, 5x conc.	Random prime labeling mixture: Premixed solution of 1 U/ $\mu$ l Klenow polymerase, labeling grade, 1 mM dATP, dCTP, dGTP each, 0.65 mM dTTP, 0.35 mM DIG-11-dUTP, alkali-labile, 5x stabilized reaction buffer in 50% (v/v) glycerol.	1 vial, 160 $\mu$ l

## 1.2. Storage and Stability

### Storage Conditions (Product)

When stored at  $-15$  to  $-25^{\circ}\text{C}$ , the product is stable through the expiry date printed on the label.

Vial / Bottle	Label	Storage
1	DIG-High Prime, 5x conc.	Store at $-15$ to $-25^{\circ}\text{C}$ . ⚠ <b>Avoid repeated freezing and thawing.</b> ⚠ <b>To avoid contamination, aliquot and store the solution in 2 to 3 vials.</b>

## 1.3. Additional Equipment and Reagent required

### For random primed DNA labeling

- 0.2 M EDTA, pH 8.0
- Autoclaved, double-distilled water
- Water bath
- Ice bath

### For labeling DNA isolated from agarose

- High Pure PCR Product Purification Kit\*

### For detection of DIG-labeled DNA

- Anti-Digoxigenin-AP, Fab fragments\*
- DIG Nucleic Acid Detection Kit\*
- DIG Luminescent Detection Kit\*

## 1.4. Application

DIG-High Prime-labeled probes are used in a variety of hybridization reactions:

- Southern blots
- Northern blots
- Dot/slot blots
- Screening of gene libraries
- *In situ* hybridizations

## 2. How to Use this Product

### 2.1. Before you Begin

#### Sample Materials

##### Templates for labeling reaction

- DNA fragments of at least 100 bp.
- Linearized plasmid, cosmid, or λDNA.
- Supercoiled DNA
- Minimal amounts of DNA (10 ng), such as DNA restriction fragments isolated from gels or in molten agarose.

#### General Considerations

##### Template DNA requirements

Feature	Detail
Purity	<ul style="list-style-type: none"> <li>▪ For plasmid DNA, use the High Pure Plasmid Isolation Kit* for purification.</li> <li>▪ When other commercially available purification kits are used, perform an additional phenol/chloroform extraction to remove residual protein.</li> </ul> <p><b>i</b> <i>This step is also necessary when templates have been treated with restriction or other modifying enzymes before labeling.</i></p>
Size	<ul style="list-style-type: none"> <li>▪ To obtain optimal results, template DNA should be linearized and have a size of <math>\geq 100</math> bp.</li> <li>▪ Template DNA <math>&gt;5</math> kb should be restriction digested using a 4 bp cutter prior to labeling.</li> </ul>
Amount	<p>For the Random primed DNA labeling protocol, 0.01 to 3 <math>\mu</math>g of template can be labeled.</p> <p><b>i</b> <i>Larger amounts can be labeled by scaling up of all components and volumes. If single-copy gene detection in complex genomes is performed, at least 300 ng of template DNA (probe concentration: 25 ng/ml hybridization solution) should be labeled.</i></p>

#### Safety Information

##### For customers in the European Economic Area

Contains SVHC: octyl/nonylphenol ethoxylates. For use in research and under controlled conditions only – acc. to Art. 56.3 and 3.23 REACH Regulation.

#### Working Solution

Solution	Composition	Use	Storage and Stability
Water	Autoclaved, double-distilled water	Dilution of DNA.	Store at +15 to +25°C.
EDTA	0.2 M ethylenediaminetetraacetic acid, pH 8.0	Stops the reaction.	Store at +15 to +25°C.

## 2.2. Protocols

### Random primed DNA labeling

Perform the standard random primed DNA labeling according to the following steps.

- 1 To a reaction vial, add 1 µg template DNA (linear or supercoiled) and autoclaved, double-distilled water to a final volume of 16 µl.

---

- 2 Denature the DNA by heating in a boiling water bath for 10 minutes.
  - Quickly chill in an ice water bath.

*i* Full denaturation is essential for efficient labeling.

---

- 3 Mix DIG-High Prime thoroughly and add 4 µl to the denatured DNA; mix and centrifuge briefly.
  - Incubate for 1 hour or overnight at +37°C.

*i* Longer incubations up to 20 hours increase the yield of DIG-labeled DNA, see Table, **Labeling reaction yield**.

---

- 4 Stop the reaction by adding 2 µl 0.2 M EDTA, pH 8.0, and/or by heating to +65°C for 10 minutes.
 

*i* The length of the DIG-labeled fragments obtained with DIG-High Prime ranges from 200 bp to ≥1,000 bp, depending on the length of the original template.

### Labeling DNA isolated from agarose

- 1 For hybridization of genomic Southern blots, separate the template insert DNA from the vector by agarose gel electrophoresis.

---

- 2 Isolate DNA from the gel using the High Pure PCR Product Purification Kit \* or an agarose gel DNA extraction kit for DNA fragments in the range of 400 bp to 5 kbp.
  - The kit can be used for standard agarose gels as well as low-melting point agarose gels.

*i* The DNA fragments are efficiently labeled with digoxigenin without further purification. However, labeled probes should be purified with the High Pure PCR Product Purification Kit\* to remove residual agarose particles.


### Labeling reaction yield

The labeling efficiency is shown in the following table and in Figure 1.

Template DNA [ng]	Template DNA [ng] and Labeling Time	
	1 Hour	20 Hours
10	45	600
30	130	1,050
100	270	1,500
300	450	2,000
1,000	850	2,300
3,000	1,350	2,650

- i* Reactions were performed with increasing amounts of different template DNA for 1 hour and 20 hours. The yield of DIG-labeled DNA was determined by incorporation of a radioactive tracer and confirmed by a dot blot. Numbers shown are the average of 10 independent labeling assays.

## 2. How to Use this Product



**Fig. 1:** Yield of DIG-labeled DNA from different amounts of template DNAs for 1 h and 20 h incubation of the DIG-High Prime reaction at +37°C.

### Detection of DIG-labeled DNA

Detect DIG-labeled DNA using Anti-DIG-AP, Fab fragments\* which catalyze a color or a chemiluminescent reaction. Special kits are available for color detection (DIG Nucleic Acid Detection Kit\*) or chemiluminescent detection (DIG Luminescent Detection Kit\*). Alternatively, especially for *in situ* applications, DIG-labeled hybrids can also be detected by anti-DIG antibodies conjugated to different fluorochromes.

### Semi-quantitative determination of labeling efficiency

Determination of the yield of DIG-labeled DNA is most important for optimal and reproducible hybridization results. Too high of a probe concentration in the hybridization step causes background, while too low of a concentration leads to weak signals.

The preferred method for quantification of labeled probes is the direct detection method.

- ① A series of dilutions of DIG-labeled DNA is applied to a small strip of Nylon Membrane, positively charged\*.
  - Part of the nylon membrane is preloaded with defined dilutions of DIG-labeled Control DNA\* which are used as standards.
- ② The nylon membrane is subjected to immunological detection with Anti-Digoxigenin-AP, Fab fragments\* and CSPD, ready-to-use\*.
  - The intensities of the dilution series of DIG-labeled DNA and control DNA are compared by exposure to an imaging device, X-ray film, or Lumi-Film\*.

### Probe quantification

To prepare the dilution series shown below, dilute labeled probes and the DIG-labeled Control DNA\* to 1 ng/μl, according to the expected yield of synthesized nucleic acid. The expected yield of DIG-labeled DNA in your probe can best be estimated by referring to the Table in, **Labeling reaction yield**. The yield depends on the starting amount of template and incubation time.

- i** The yields given in the Table, **Labeling reaction yield** were achieved under optimal conditions using highly purified template DNA.

## Dilution series

Prepare a dilution series of your labeled probe and your control DNA as described.

**i** See the Instructions for Use of the DIG-High Prime DNA Labeling and Detection Starter Kit II\* for the protocol and working solutions.

Tube	DNA [μl]	From Tube No.	DNA Dilution Buffer [μl]	Dilution	Final Concentration
1	–	original	–	–	1 ng/μl
2	5	1	495	1:100	10 pg/μl
3	15	2	35	1:3.3	3 pg/μl
4	5	2	45	1:10	1 pg/μl
5	5	3	45	1:10	0.3 pg/μl
6	5	4	45	1:10	0.1 pg/μl
7	5	5	45	1:10	0.03 pg/μl
8	5	6	45	1:10	0.01 pg/μl
9	0	–	50	–	0


## Analyzing the results

Compare the intensity of the spots from your labeling reaction to the control and calculate the amount of DIG-labeled DNA. If the 0.1 pg dilution spots of your probe and of the control are visible, then the labeled probe has reached the expected labeling efficiency, see Table, **Labeling reaction yield**, and can be used in the recommended concentration in the hybridization.

## 2.3. Parameters

### Chemical Name

### Structural formula



**Fig. 2:** Structure of alkali-labile DIG-11-dUTP.

## 3. Additional Information on this Product



### 3.1. Test Principle

DIG-labeled DNA probes are generated with DIG-High Prime according to the random primed labeling technique. DIG-High Prime is a specifically developed reaction mixture containing DIG-11-dUTP, and all reagents necessary for random-primed labeling. The premixed DIG-High-Prime reduces pipetting steps and increases yield, reproducibility, and convenience.

## 4. Supplementary Information

### 4.1. Conventions

To make information consistent and easier to read, the following text conventions and symbols are used in this document to highlight important information:

Text convention and symbols	
	<i>Information Note: Additional information about the current topic or procedure.</i>
	<b>Important Note: Information critical to the success of the current procedure or use of the product.</b>
① ② ③ etc.	Stages in a process that usually occur in the order listed.
1 2 3 etc.	Steps in a procedure that must be performed in the order listed.
* (Asterisk)	The Asterisk denotes a product available from Roche Diagnostics.

### 4.2. Changes to previous version

Layout changes.

Editorial changes.

New information added related to the REACH Annex XIV.

### 4.3. Ordering Information

Product	Pack Size	Cat. No.
Reagents, kits		
Agarose Gel DNA Extraction Kit	1 kit, up to 100 reactions	11 696 505 001
High Pure PCR Product Purification Kit	1 kit, up to 50 purifications	11 732 668 001
	1 kit, up to 250 purifications	11 732 676 001
DIG Nucleic Acid Detection Kit	1 kit, Detection of 40 blots of 10 cm x 10 cm	11 175 041 910
DIG Luminescent Detection Kit	1 kit, 50 blots with a size of 10 x 10 cm <sup>2</sup>	11 363 514 910
Nylon Membranes, positively charged	10 sheets, 20 x 30 cm	11 209 272 001
	20 sheets, 10 x 15 cm	11 209 299 001
	1 roll, 0.3 x 3 m	11 417 240 001
High Pure Plasmid Isolation Kit	1 kit, 50 purifications	11 754 777 001
	1 kit, 250 purifications	11 754 785 001
Lumi-Film Chemiluminescent Detection Film	100 films, 7.1 x 9.4 inches, 18 x 24 cm, <i>Not available in US</i>	11 666 916 001
CSPD, ready-to-use	2 x 50 ml	11 755 633 001
DIG-labeled Control DNA	50 µl, 5 µg/ml DIG-labeled plasmid DNA	11 585 738 910
DIG-High Prime DNA Labeling and Detection Starter Kit II	1 kit, 12 labeling reactions of 10 ng to 3 µg DNA and detection of 24 blots of 100 cm <sup>2</sup>	11 585 614 910



## 4.4. Trademarks

All product names and trademarks are the property of their respective owners.

## 4.5. License Disclaimer

For patent license limitations for individual products please refer to:

**List of biochemical reagent products.**

## 4.6. Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

## 4.7. Safety Data Sheet

Please follow the instructions in the Safety Data Sheet (SDS).

## 4.8. Contact and Support

To ask questions, solve problems, suggest enhancements or report new applications, please visit our **Online Technical Support Site.**

To call, write, fax, or email us, visit **sigma-aldrich.com**, and select your home country. Country-specific contact information will be displayed.

