



# Watercol™ Capillary GC Columns

Convenient Analysis of Water



The Life Science business  
of Merck operates as  
MilliporeSigma in the  
U.S. and Canada.

**Supelco®**  
Analytical Products

Watercol™ 1910 is a specialized ionic liquid column designed for gas chromatography, particularly effective for moisture analysis. Its chemical structure enhances water retention and detection sensitivity.

This column excels in separating water from various volatile organic compounds, making it suitable for analyzing complex matrices.



Compared to traditional methods like Karl Fischer titration, Watercol™ 1910 offers improved accuracy and speed, significantly reducing analysis time while maintaining high precision.

Its stability against water and thermal degradation further solidifies its role in analytical chemistry, including pharmaceutical, petroleum, and environmental applications.

### Water: Column Interactions

## Interaction Strength

The fundamental property driving the performance of Watercol™ series columns is the proper combination and strength of water, i.e. column interactions.



## Ionic Liquids: A Versatile Solution for Gas Chromatography

Ionic liquids (ILs) are a unique class of solvents characterized by their liquid state at room temperature, consisting entirely of ions. They have gained significant attention in the field of gas chromatography (GC) due to their distinct properties and advantages.

**Tunable Properties:** The physicochemical properties of ILs can be easily manipulated by varying the cation and anion, allowing for customization to suit specific analytical needs.

**High Thermal Stability:** ILs exhibit superior thermal stability compared to traditional GC stationary phases, enabling their use at elevated temperatures without degradation. Example, SLB®-IL60i has a higher maximum temperature than most PEG/WAX columns.

**Low Vapor Pressure:** Their negligible vapor pressure minimizes contamination risks and enhances safety during GC applications.

**Dual Polarity:** ILs can function as both polar and nonpolar solvents, providing excellent selectivity for a wide range of analytes, including polar and nonpolar compounds.

**Enhanced Separation Efficiency:** ILs improve peak shapes and resolution, facilitating the analysis of complex mixtures and trace components.

With these attributes, ionic liquids represent a promising alternative for GC stationary phases, offering scientists innovative solutions for advanced analytical challenges.

For more information visit: [SigmaAldrich.com/il-gc-inert](https://www.sigmaaldrich.com/il-gc-inert)

## Watercol™ Checklist

### Sharp Peak Shape:

Water can be integrated and quantified – see **Figure 1** and **Figure 3**.

### Linear Response:

$r^2$  value of 0.9923 achieved from a 5-point calibration curve – see **Figure 2**.

### Great Sensitivity:

- 100 ppm using a thermal conductive detector (TCD)
- <100 ppm using a mass spectrometer detector (MSD) or a vacuum ultraviolet spectroscopy detector (VUV)
- < 10 ppb using a barrier discharge ionization detector (BID)

### Small Sample Size:

GC results can be obtained using small sample sizes, such as <2 mL gasoline and <0.2 g ibuprofen.

### Selectivity Options:

Whether water is an analyte or the injection solvent – see **Figure 3** and **Figure 4**.

## Suitability

Narrow peak widths and optimal peak heights are also produced for other small polar analytes. Watercol™ columns are suitable for applications where:


- Water is an analyte: can be integrated and quantified
- Water is the injection solvent: because it does not tail, it does not interfere chromatographically with other analytes

## Column Selection

Quantify water as an analyte or quantify small polar analytes in water. The chromatograms in this brochure cover a variety of analytes and compound classes.


**Figure 1. Water Standard (0.05% in Ethanol) on Watercol™ 1910**

<b>Column</b>	Watercol™ 1910, 30 m × 0.25 mm I.D., 0.20 μm ( <b>29711-U</b> )
<b>Oven</b>	80 °C (10 min)
<b>Inj. temp.</b>	250 °C
<b>Detector</b>	TCD, 200 °C
<b>Carrier gas</b>	helium, 26 cm/sec
<b>Injection</b>	0.5 μL, 100:1 split
<b>Liner</b>	4 mm I.D., split type, cup design ( <b>2048225</b> )
<b>Sample</b>	water at 0.05% (v/v) in ethanol



**Figure 2. Water Calibration Curve (0.01–0.5%) on Watercol™ 1910**

<b>Column</b>	Watercol™ 1910, 30 m × 0.25 mm I.D., 0.20 μm ( <b>29711-U</b> )
<b>Oven</b>	80 °C (10 min)
<b>Inj. temp.</b>	250 °C
<b>Detector</b>	TCD, 200 °C
<b>Carrier gas</b>	helium, 26 cm/sec
<b>Injection</b>	0.5 μL, 100:1 split
<b>Liner</b>	4 mm I.D., split type, cup design ( <b>2048225</b> )
<b>Sample</b>	5 standards, water at 0.01, 0.05, 0.1, 0.3, and 0.5% (v/v), each in ethanol




## Water Measurement Determination

Karl Fischer Titration (KFT) has long been the standard for water determination. Gas Chromatography (GC) with Watercol™ columns represents a modern alternative, offering enhanced sensitivity and accuracy for water analysis, paving the way for more reliable and efficient future applications.


Feature/Aspect	Gas Chromatography with Watercol™ 1910	Karl Fischer Titration (KFT)
Sensitivity	Higher sensitivity (1-100 ppm)	Moderate sensitivity (5-500 ppm)
Sample Size	Smaller sample volume required	Larger sample volume needed
Hands-On Time	Less hands-on time	More hands-on time
Interference from Reagents	Minimal unwanted side reactions	Prone to interferences from organic compounds and acids
Chemical Waste	Lower volume of chemical waste	Higher volume of chemical waste
Solubility Issues	No solubility issues	May encounter solubility issues
Simultaneous Analysis	Can analyze water and volatiles	Typically analyzes only water
Worker Exposure	Reduced exposure to harmful chemicals	Higher exposure risk
Legacy vs. Modern	Modern, innovative method	Legacy method

**Figure 3. Water as an Analyte**


<b>Column</b>	Watercol™ 1910, 30 m x 0.25 mm I.D., 0.20 µm (29711-U)
<b>Oven</b>	96 °C
<b>Inj. temp.</b>	250 °C
<b>Detector</b>	TCD, 200 °C
<b>Carrier gas</b>	helium, 26 cm/sec
<b>Injection</b>	1 µL, 100:1 split
<b>Liner</b>	4 mm I.D., split type, cup design
<b>Sample</b>	6 analytes in acetone

**Figure 4. Water as the Injection Solvent**

<b>Column</b>	Watercol™ 1910, 30 m x 0.25 mm I.D., 0.20 µm (29711-U)
<b>Oven</b>	35 °C, 4 °C/min to 125 °C (2 min)
<b>Inj. temp.</b>	250 °C
<b>Detector</b>	TCD, 300 °C
<b>Carrier gas</b>	helium, 25 cm/sec
<b>Injection</b>	1 µL, 100:1 split
<b>Liner</b>	4 mm I.D., split type, cup design (2048225)
<b>Sample</b>	8-component solvent mix in water

**Figure 5. Water in a 60-Component Solvent Mix (Chemical Application)**


<b>Column</b>	Watercol™ 1910, 30 m x 0.25 mm I.D., 0.20 µm (29711-U); two identical columns installed in the same inlet using a 2-hole ferrule, each going to a separate detector
<b>Oven</b>	50 °C (1 min), 5 °C/min to 180 °C (5 min)
<b>Inj. temp.</b>	250 °C
<b>Detector</b>	FID, 250 °C
<b>Detector</b>	TCD, 200 °C
<b>Carrier gas</b>	helium, 25 cm/sec (measured in each column at an oven temperature of 125 °C)
<b>Injection</b>	1 µL, 100:1 split
<b>Liner</b>	4 mm I.D., split/splitless type, wool-packed single-taper FocusLiner™ design (2879925-U)
<b>Sample</b>	mixture of 60 solvents, each at 0.2% (v/v), and water at 0.3% (v/v), in pentane




Additional chromatograms, product information, real-time availability, and ordering information is available 24 hours a day at [SigmaAldrich.com/watercol](https://www.sigmaaldrich.com/watercol)

**Figure 6. Water in Ibuprofen (Pharmaceutical Application)**


<b>Column</b>	Watercol™ 1910, 30 m x 0.25 mm I.D., 0.20 µm ( <b>29711-U</b> )
<b>Oven</b>	35 °C (4 min), 8 °C/min to 130 °C (1 min)
<b>Inj. temp.</b>	250 °C
<b>Detector</b>	TCD, 200 °C
<b>Carrier gas</b>	helium, 11.40 psi constant pressure (equal to 25 cm/sec), measured with oven at 35 °C
<b>Injection</b>	1 µL, 100:1 split
<b>Liner</b>	4 mm I.D., split/splitless type, wool-packed single-taper FocusLiner™ design ( <b>2879925-U</b> )
<b>Sample</b>	0.198 g ibuprofen weighed into a 1 mL volumetric flask, 0.5 mL acetone added, vortexed slowly for 30 seconds, acetone added to the 1 mL line, vortexed, transferred to an autosampler vial

**Figure 7. Water Impurity in C1-C6 Alcohol Mix**

<b>Column</b>	Watercol™ 1910, 30 m x 0.25 mm I.D., 0.20 µm ( <b>29711-U</b> )
<b>Oven</b>	40 °C (2 min), 4 °C/min to 125 °C (1 min)
<b>Inj. temp.</b>	300 °C
<b>Detector</b>	MSD, m/z = 18–500
<b>MSD interface</b>	260 °C
<b>Carrier gas</b>	helium, 1.1 mL/min
<b>Injection</b>	1 µL, 100:1 split
<b>Liner</b>	3.4 mm I.D., split/splitless type, wool-packed straight FocusLiner™ design ( <b>2877605-U</b> )
<b>Sample</b>	10-component alcohol mix in methylene chloride, 500 µg/mL total concentration (mixture adsorbed some water during storage)

**Figure 8. Aroma Fusel Alcohols in Tequila (Food and Beverage Application)**

<b>Column</b>	Watercol™ 1910, 30 m x 0.25 mm I.D., 0.20 µm ( <b>29711-U</b> ); two identical columns installed in the same inlet using a 2-hole ferrule, each going to a separate detector
<b>Oven</b>	50 °C (1 min), 5 °C/min to 180 °C (3 min)
<b>Inj. temp.</b>	250 °C
<b>Detector</b>	FID, 200 °C
<b>Detector</b>	TCD, 200 °C
<b>Carrier gas</b>	helium, 25 cm/sec (measured in each column at an oven temperature of 125 °C)
<b>Injection</b>	1 µL, 100:1 split
<b>Liner</b>	4 mm I.D., split/splitless type, wool packed single taper FocusLiner™ design ( <b>2879905-U</b> )
<b>Sample</b>	neat tequila



Watercol™ capillary GC columns contain innovative ionic liquid stationary phases that produce a sharp peak shape for water, allowing the convenient measurement of water by GC.

**Table 1. Watercol™ Series Column Specifications**

Watercol™ 1910	
USP Code	None
Phase	Non-bonded; 1,11-Di(3-hydroxyethylimidazolium) 3,6,9-trioxaundecane trifluoromethanesulfonate
Temp. Limits	30 °C to 180 °C (isothermal or programmed)

### Ordering Information

Description	Cat. No.
30 m × 0.25 mm I.D., 0.20 µm	29711-U
30 m × 0.32 mm I.D., 0.26 µm	29714-U

### Related Information

Additional chromatograms, product information, real-time availability, and ordering information is available 24 hours a day at [SigmaAldrich.com/watercol](http://SigmaAldrich.com/watercol)

## Did you know...?

The descriptive numbers (1910) indicate the Kovats Retention Index (KRI) value of water at 100 °C isothermal oven temperature. KRI values indicate the elution location of the analyte relative to n-alkanes, in which the n-alkane carbon number is multiplied by 100. For example, a KRI value of 1910 indicates that water elutes after nonadecane (nC19) and 10% of the interval before eicosane (nC20) elutes in 100 °C isothermal oven temperature.

Supelco®, a brand under the Life Science business of Merck KGaA, Darmstadt, Germany, has a legacy of providing high-quality chromatography products and solutions, particularly in gas chromatography. Tracing its roots to 1966, Supelco® Analytical Products has been recognized for its innovative stationary phases, columns, and accessories that have significantly advanced analytical chemistry with a focus on quality, reliability, and supporting you, our customers.

# Supelco®

Analytical Products

Merck KGaA  
Frankfurter Strasse 250  
64293 Darmstadt, Germany

**SigmaAldrich.com**



To place an order or receive  
technical assistance:  
**SigmaAldrich.com/support**



For local contact information:  
**SigmaAldrich.com/offices**

We have built a unique collection of life science brands with  
unrivalled experience in supporting your scientific advancements.

**Millipore® Sigma-Aldrich® Supelco® Milli-Q® SAFC® BioReliance®**

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, BioReliance, Millipore, Milli-Q, SAFC, Sigma-Aldrich, SLB, Watercol, and Supelco are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MK\_BR14203EN Ver. 1.0  
60243  
01/2025