User Guide

Amicon® Ultra-4 10K Centrifugal Filter Devices for volumes up to 4 mL

UFC801008D
UFC801024D
UFC801096D

IVD For *in vitro* diagnostic use

Languages available

English	2	Lietuviškai	68
Français	8	Magyar	74
Italiano	14	Čeština	80
Deutsch	20	Polski	86
Español	26	日本語	92
Português	32	中文	98
Ελληνικά	38	한국어	104
Nederlands	44	Norsk	110
Dansk	50	Slovenčina	116
Svenska	56	Türkçe	122
Latviski	62		

Introduction

Amicon® Ultra-4 10K centrifugal filter devices provide fast ultrafiltration, with the capability for high concentration factors and easy concentrate recovery from dilute and complex sample matrices. The vertical design and available membrane surface area provide fast sample processing, high sample recovery (typically greater than 90% of dilute starting solution), and the capability for 80-fold concentration. Typical processing time is 10 to 20 minutes. Solute polarization and subsequent fouling of the membrane are minimized by the vertical design, and a physical deadstop in the filter device prevents spinning to dryness and potential sample loss. The concentrate is collected from the filter device sample reservoir using a pipettor, while the ultrafiltrate is collected in the provided centrifuge tube. The device can be spun in a swinging-bucket (for optimal performance) or fixed-angle rotor. Amicon® Ultra-4 10K devices are supplied non-sterile and are for single use only.

Intended Use

Amicon® Ultra-4 centrifugal filter devices are disposable, single-use ultrafiltration filter devices for *in vitro* diagnostic use for processing biological solutions in the 1-4 mL volume range. They are intended for concentrating biological samples such as serum, urine, cerebrospinal fluid and other body fluids prior to *in vitro* diagnostic analysis. Device intended to be single-use and used by a laboratory professional.

Applications

- Concentration of biological samples containing antigens, antibodies, enzymes, nucleic acids (DNA/RNA samples, either single- or double-stranded), microorganisms, column eluates, and purified samples
- Purification of macromolecular components found in tissue culture extracts and cell lysates, removal of primer, linkers, or molecular labels from a reaction mix, and protein removal prior to HPLC
- Desalting, buffer exchange, or diafiltration

Materials Supplied


The Amicon® Ultra-4 10K device is supplied with a cap, a filter device, and a centrifuge tube.

Required Equipment

- Centrifuge with swinging-bucket (preferred) or fixed-angle rotor with wells/carriers that can accommodate 17 mm × 124 mm 15 mL conical-bottomed tubes

CAUTION: To avoid damage to the device during centrifugation, check clearance before spinning.

- Pipettor with 200 microliter (µL) tip for concentrate recovery

Suitability

Preliminary recovery and retention studies are suggested to ensure suitability for intended use. See the "How to Quantify Recoveries" section.

Storage and Stability

Refer to product label for storage conditions and shelf life.

Rinsing Before Use

The ultrafiltration membranes in Amicon® Ultra-4 10K devices contain trace amounts of glycerine. If this material interferes with analysis, rinse the device with buffer or Milli-Q® water before use. If interference continues, rinse with 0.1 N NaOH followed by a second spin of buffer or Milli-Q® water.

CAUTION: Do not allow the membrane in Amicon® Ultra filter devices to dry out once wet. If you are not using the device immediately after rinsing, leave fluid on the membrane until the device is used.

How to Use Amicon® Ultra-4 Centrifugal Filter Devices

1. Add up to 4 mL of sample (3.5 mL if using a 23° fixed-angle rotor) to the Amicon® Ultra filter device.
2. Place capped filter device into centrifuge rotor (swinging-bucket preferred); counterbalance with a similar device. Check Centrifuge clearance before use.
3. **When using a swinging-bucket rotor**, spin the device at $4,000 \times g$ maximum for approximately 10–20 minutes.
4. **When using a fixed-angle rotor**, orient the device with the membrane panel facing up and spin at $7,500 \times g$ maximum for approximately 10–20 minutes.

NOTE: Refer to Figure 1 and Table 1 for typical spin times.


5. To recover the concentrated solute, insert a pipettor into the bottom of the filter device and withdraw the sample using a side-to-side sweeping motion to ensure total recovery. The ultrafiltrate can be stored in the centrifuge tube.

NOTE: For optimal recovery, remove concentrated sample immediately after centrifugation.

NOTE: Filter may not function correctly if allowed to dry out after wetting.

Desalting or Diafiltration

Desalting, buffer exchange, or diafiltration are important methods for removing salts or solvents in solutions containing biomolecules. The removal of salts or the exchange of buffers can be accomplished in the Amicon® Ultra-4 device by concentrating the sample, then reconstituting the concentrate to the original sample volume with any desired solvent. The process of “washing out” can be repeated until the concentration of the contaminating microsolite has been sufficiently reduced. See example below.


Performance

Flow Rate

Factors affecting flow rate include sample concentration, starting volume, chemical nature of solute, relative centrifugal force, centrifuge rotor angle, membrane type, and temperature. Figure 1 and Table 1 can be used to estimate the time required to achieve a given volume of filtrate or concentrate. A typical spin time for a 4 mL sample is approximately 10 to 20 minutes. While most of the sample is filtered in the first 10 to 15 minutes of centrifugation, the lowest concentrate volume (30–75 μL) is reached after spinning for 15 to 20 minutes.

Figure 1. Typical Filtrate Volume vs. Spin Time

Spin Conditions: Swinging-bucket rotor, $4,000 \times g$, or fixed-angle rotor, $7,500 \times g$, room temperature, 4 mL starting volume. Protein marker used: Cytochrome c, $n=6$.

Table 1. Typical Concentrate Volume vs. Spin Time

Spin time (min)	Concentrate volume (μL)	
	Swinging-bucket rotor $4,000 \times g$	35° Fixed-angle rotor $7,500 \times g$
10	176	97
15	76	54
20	58	35

Spin Conditions: Room temperature, 4 mL starting volume. Protein marker used: Cytochrome c, $n=6$ (mean value of 3 device lots). Shaded volumes were used for the calculation of protein recovery in Table 3.

Protein Retention and Concentrate Recovery

The membranes used in Amicon® Ultra devices are characterized by a molecular weight cutoff (MWCO); that is, their ability to retain molecules above a specified molecular weight. Solutes with molecular weights close to the MWCO may be only partially retained. Membrane retention depends on the solute's molecular size and shape. For most applications, molecular weight is a convenient parameter to use in assessing retention characteristics. For best results, use a membrane with a MWCO at least two times smaller than the molecular weight of the protein solute that one intends to concentrate. Refer to Table 2.

Table 2. Typical Retention of Protein Markers

Marker/Concentration	Molecular Weight	Device MWCO	% Retention Swinging-bucket	% Retention Fixed-angle	Spin Time (min)
α -Chymotrypsinogen (1 mg/mL)	25,000	10K	> 95	> 95	15
Cytochrome c (0.25 mg/mL)	12,400		> 95	> 95	15
Vitamin B-12 (0.2 mg/mL)	1,350		< 15	< 15	15

Spin Conditions: Swinging-bucket rotor, $4,000 \times g$, or 35° fixed-angle rotor, $7,500 \times g$, 4 mL starting volume, room temperature, $n=6$ (mean value of 3 device lots).

Factors that determine sample recovery include the nature of the protein solute relative to the device MWCO chosen, starting concentration, and concentration factor. Table 3 provides typical recoveries for the Amicon® Ultra-4 10K device.

Table 3. Typical Concentrate Recovery

Marker/Concentration	Device MWCO	Spin Time (min)	Concentrate Volume (μL)		Concentration Factor (x)		Concentrate Recovery (%)	
			Swinging-bucket	Fixed-angle	Swinging-bucket	Fixed-angle	Swinging-bucket	Fixed-angle
Cytochrome c (0.25 mg/mL)	10K	15	76	54	52.3	76.6	97.3	98.5

Spin Conditions: Swinging-bucket rotor, $4,000 \times g$, or 35° fixed-angle rotor, $7,500 \times g$, 4 mL starting volume, room temperature, $n=6$ (mean value of 3 device lots). The shaded volumes were taken from Table 1.

Maximizing Sample Recovery

Low sample recovery in the concentrate may be due to adsorptive losses, over-concentration, or passage of sample through the membrane.

- Adsorptive losses depend upon solute concentration, its hydrophobic nature, temperature and time of contact with filter device surfaces, sample composition, and pH. To minimize losses, remove concentrated samples immediately after centrifugal spin.
- If the starting sample concentration is high, monitor the centrifugation process in order to avoid over-concentration of the sample. Over-concentration can lead to precipitation and potential sample loss.
- If the sample appears to be passing through the membrane, choose a lower MWCO Amicon® Ultra-4 device.

How to Quantify Recoveries

Calculate total recovery, percent concentrate recovery, and percent filtrate recovery using the method below. The procedure provides a close approximation of recoveries for solutions having concentrations up to roughly 20 mg/mL.

NOTE: Appropriate assay techniques include absorption spectrophotometry, radioimmunoassay, refractive index, and conductivity.

Direct Weighing Procedure

The density of most dilute proteins is nearly equal to the density of water (i.e., 1 g/mL). Using this property, the concentrate and filtrate volumes can be quantified by weighing them and converting the units from grams to milliliters. This technique is valid only for solutions with concentrations of approximately 20 mg/mL or less.

- Before use, separately weigh the empty filter device, the centrifuge tube, and an empty tube for concentrate collection.
- Fill filter device with solution and reweigh.
- Assemble device and centrifuge per instructions.
- Collect the concentrate with a pipettor and dispense it into the preweighed concentrate collection tube.
- Remove the device from the centrifuge tube and weigh the centrifuge tube and concentrate collection tube.
- Subtract weight of empty device/tubes to calculate weights of starting material, filtrate, and concentrate.
- Assay the starting material, filtrate, and concentrate to determine solute concentration.
- Calculate recoveries using the weight/volume data and the measured concentrations as follows:

$$\% \text{ concentrate recovery} = 100 \times \frac{W_c \times C_c}{W_o \times C_o}$$

$$\% \text{ filtrate recovery} = 100 \times \frac{W_f \times C_f}{W_o \times C_o}$$

$$\% \text{ total recovery} = \% \text{ concentrate recovery} + \% \text{ filtrate recovery}$$

W_c = total weight of concentrate before assay
 W_o = weight of original starting material
 W_f = weight of filtrate

C_c = concentrate concentration
 C_o = original starting material concentration
 C_f = filtrate concentration

Specifications

Maximum initial sample volume

Swinging-bucket and fixed-angle rotors (45° and 35°) 4.0 mL

Fixed-angle rotor (23°) 3.5 mL

Typical final concentrate volume 50–100 µL

Maximum relative centrifugal force

Swinging-bucket rotor 4,000 × g

Fixed-angle rotor 7,500 × g

Active membrane area 3 cm²

Dimensions

Filter device in tube (capped)

Length: 124 mm (4.9 in.) Diameter: 17.3 mm (0.7 in.)

Filter device

Length: 73.4 mm (2.9 in.) Diameter: 17.2 mm (0.7 in.)

Materials of construction

Filter device Copolymer styrene/butadiene

Membrane Ultracel® low binding regenerated cellulose

Filtrate tube Polypropylene

Filtrate cap Polypropylene

Filtrate liner Polyethylene

Filter MWCO 10K

Chemical Compatibility

Amicon® Ultra centrifugal devices are intended for use with biological fluids and aqueous solutions. Before use, check the sample for chemical compatibility with the device. Go to SigmaAldrich.com/FilterChemicalCompatibility for more information.

Symbol Definitions

Symbol	Definition	Symbol	Definition
	In vitro diagnostic medical device		Do not use if package is damaged and consult instructions for use
	Catalogue number		Date of manufacture
	Batch code		Manufacturer
	Consult electronic instructions for use		Importer
	Download product documentation online		CE conformity marking
	Non-sterile		UK Conformity Assessed
	Do not reuse		Swiss authorized representative
	Use-by date		Unique device identifier
	Temperature limit		

Product Ordering

Purchase products online at SigmaAldrich.com.

Description	Qty/Pk	Catalogue Number
	8	UFC801008D
Amicon® Ultra-4	24	UFC801024D
	96	UFC801096D

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

Collection and Disposal

All samples must be clearly labelled. Suitable instruments must be used for obtaining and preparing samples.

NOTE: Follow precautions for disposal of items contaminated with potentially infectious or hazardous biological material according to all applicable international, federal, state, and local regulations.

Technical Assistance

Visit the tech service page on our web site at SigmaAldrich.com/techservice.

Any serious incident of this device should be reported to manufacturer and competent authority of country where user is established.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

Revision History

2021-OCT	• IFU PR05779 Date of Issue OCT 2021 - Replaced PR05139. • IFU, Packaging Damage and Temperature Limit symbols added. • Chemical Compatibility and Ordering Information linked to website. • Disposal and Complaints information added. • UK Responsible Person and UKCA symbol information added
2024-MAR	• Figure 1 corrected • "Materials of construction": Filtrate cap corrected text to polypropelene
2024-OCT	• Added to Symbol Definitions table: CH-REP, Importer, UDI • Added UKCA symbol on title page

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Made in Ireland

Merck Millipore Ltd.
Tullagreen, Carrigtwohill,
Co. Cork, Ireland
T45 KD29

an affiliate of Merck KGaA,
Darmstadt, Germany

UK Responsible Person:

Millipore (UK) Ltd.
Kirkton Campus
Fleming Road,
Livingston, UK
EH54 7BN

Merck, Millipore, Amicon, Milli-Q, Ultracel, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.
© 2001-2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck operates
as MilliporeSigma in the U.S. and Canada.