

Data Sheet

HPDE6c7 Human Pancreatic Ductal Epithelial Cell Line

Immortalized Cell Line

SCC442**Pack Size ≥ 1x10⁶ viable cells/vial****Store at: Liquid nitrogen**FOR RESEARCH USE ONLY**Not for use in diagnostic procedures. Not for Human or Animal Consumption.**

Background

Pancreatic cancer is the fourth commonest cause of cancer death in North America and has one of the worst prognoses. Greater than 90% of these tumors arise from the pancreatic duct epithelium.¹ Even when the primary tumor is small and localized, the prognosis remains poor and chemotherapy or radiotherapy has demonstrated limited effectiveness.² The availability of dynamic models remains crucial to the study and understanding of the biological significance of the genetic changes underlying pancreatic cancer, especially in the context of pancreatic duct epithelial cell carcinogenesis.

The immortalized HPDE6c7 cell line demonstrates a near normal genotype and phenotype of pancreatic duct epithelial cells.³ HPDE6c7 cells retain normal Ki-ras, p53, C-myc, and p16INK4A genotypes but lack the p53 functional pathway.³ The HPDE6c7 cell line is useful for studies on the molecular basis of pancreatic duct cell carcinogenesis and islet cell differentiation. HPDE6c7 cells are a versatile model for the pancreatic ductal epithelium, facilitating the development of strategies for the chemoprevention of human pancreatic cancers.

Source

HPDE6c7 cells originated from normal human pancreatic duct of a 63-year-old female, having been genotypically altered by infection with retrovirus vector expressing E6E7 genes of human papilloma virus (HPV)-16.³

Short Tandem Repeat

D3S1358: 16	D18S51: 13, 18	TPOX: 8, 10
D7S820: 9, 10	D5S818: 10, 12	CSF1PO: 11, 12
vWA: 14, 17	D13S317: 11, 12	Amel: X
FGA: 21, 22	D16S539: 11, 13	Penta D: 9, 12
D8S1179: 14, 16	TH01: 9.3	Penta E: 5, 10
D21S11: 30, 33		

Cancer cell lines are inherently genetically unstable. Instability may arise in the form of loss of heterozygosity of alleles at one or more genetic sites with increased passages.

Quality Control Testing

- Each vial contains $\geq 1 \times 10^6$ viable cells.
- Cells are tested negative for infectious diseases by a Human Essential CLEAR Panel by Charles River Animal Diagnostic Services.
- Cells are verified to be of human origin and negative for inter-species contamination from rat, mouse, Chinese hamster, Golden Syrian hamster, and Non-human Primate (NHP) as assessed by a Contamination Clear panel by Charles River Animal Diagnostic Services.
- Cells are negative for mycoplasma contamination.

Storage and Handling

HPDE6c7 cells should be stored in liquid nitrogen. The cells can be cultured for at least 10 passages after initial thawing without significantly affecting functionality.

Representative Data




Figure 1. Bright-field images of HPDE6c7 cells one (A) and three (B) days after thaw in a T75 flask.

Protocols

Thawing Cells

1. Do not thaw the cells until the recommended medium is on hand. Cells can grow on normal tissue cultureware surfaces without any additional coating.
Cells are thawed and expanded in HPDE6c7 Expansion Medium comprised of Keratinocyte Serum-Free Growth Medium (Cat. No. 131-500A) with 50 $\mu\text{g}/\text{mL}$ bovine pituitary extract (Cat. No. P1476), 5 ng/mL EGF (Cat. No. GF316) and 1X penicillin/streptomycin (Cat. No. TMS-AB2-C, optional).
2. Remove the vial of frozen HPDE6c7 cells from liquid nitrogen and incubate in a 37 °C water bath. Closely monitor until the cells are completely thawed. Maximum cell viability is dependent on the rapid and complete thawing of frozen cells.
IMPORTANT: Do not vortex the cells.
3. As soon as the cells are completely thawed, disinfect the outside of the vial with 70% ethanol. Proceed immediately to the next step.
4. In a laminar flow hood, use a 1- or 2-mL pipette to transfer the cells to a sterile 50 mL conical tube. Be careful not to introduce any bubbles during the transfer process.

5. Using a 10 mL pipette, slowly add dropwise 9 mL of HPDE6c7 Expansion Medium (Step 1 above) to the 15 mL conical tube.
IMPORTANT: Do not add the entire volume of media all at once to the cells. This may result in decreased cell viability due to osmotic shock.
6. Gently mix the cell suspension by slowly pipetting up and down twice. Be careful not to introduce any bubbles.
IMPORTANT: Do not vortex the cells.
7. Centrifuge the tube at $300 \times g$ for 2-3 minutes to pellet the cells.
8. Decant as much of the supernatant as possible. Steps 5-8 are necessary to remove residual cryopreservative (DMSO).
9. Resuspend the cells in 15 mL of HPDE6c7 Expansion Medium.
10. Transfer the cell mixture to a T75 tissue culture flask.
11. Incubate the cells at 37 °C in a humidified incubator with 5% CO₂.

Subculturing Cells

Cells are strongly adherent. HPDE6c7 cells are sensitive to trypsin treatment and thus it is essential to neutralize the trypsin with the same volume of 1% soybean trypsin inhibitor. Cell should not be split too thinly. We recommend split ratio of 1:3-1:4.

1. HPDE6c7 cells should be passaged at ~80-85% confluency. Do not allow the cells to grow over 85% confluency.
2. Carefully remove the medium from the T75 tissue culture flask containing the 80% confluent layer of HPDE6c7 cells.
3. Rinse the flask with 10 mL 1X PBS. Aspirate after the rinse.
4. Apply 6 mL of 0.25% Trypsin/EDTA and incubate in a 37 °C incubator for 7-10 minutes.
Note: If cells do not detach, add a bit more trypsin and wait a little longer.
5. Inspect the flask and ensure the complete detachment of cells by gently tapping the side of the flask with the palm of your hand.
6. Neutralize the trypsin by adding the same volume of 1% soybean trypsin inhibitor (Cat. No. 10109886001).
7. Add 10-30 mL of HPDE6c7 Expansion Medium to the flask.
8. Gently rotate the flask to mix the cell suspension. Transfer the dissociated cells to a 50 mL conical tube.
9. Centrifuge the tube at $300 \times g$ for 3-5 minutes to pellet the cells.
10. Discard the supernatant, then loosen the cell pellet by tapping the tip of the tube with a finger.
11. Apply 2-5 mL of HPDE6c7 Expansion Medium to the conical tube and resuspend the cells thoroughly.
IMPORTANT: Do not vortex the cells.
12. Count the number of cells using a hemocytometer.
13. Plate the cells to the desired density. Typical split ratio should not exceed 1:3-1:4. The medium should be replaced every 2-3 days. HPDE6c7 cells should not be split too thinly.

Cryopreservation of Cells

HPDE6c7 cells may be frozen in HPDE6c7 Expansion Medium and 10% DMSO using a Nalgene® slow freeze Mr. Frosty® container.

References

1. Am J Pathol 1996, 148(6): 1763-1770.
2. Am J Pathol 1998, 153(1): 263-269.
3. Am J Pathol 2000, 157(5): 1623-1631

Academic Use Agreement

Subject to local law

THIS PRODUCT MAY ONLY BE USED BY INDIVIDUALS EMPLOYED BY AN ACADEMIC INSTITUTION AND IS INTENDED SOLELY TO BE USED FOR ACADEMIC RESEARCH, WHICH IS FURTHER DEFINED BELOW. BY OPENING THIS PRODUCT, YOU ("PURCHASER") HEREBY REPRESENT THAT YOU HAVE THE RIGHT AND AUTHORITY TO LEGALLY BIND YOURSELF AND/OR YOUR EMPLOYER INSTITUTION, AS APPLICABLE, AND CONSENT TO BE LEGALLY BOUND BY THE TERMS OF THIS ACADEMIC USE AGREEMENT. IF YOU DO NOT AGREE TO COMPLY WITH THESE TERMS, YOU MAY NOT OPEN OR USE THE PRODUCT AND YOU MUST CALL MILLIPORESIGMA ("SELLER") CUSTOMER SERVICE (1-800-645-5476) TO ARRANGE TO RETURN THE PRODUCT FOR A REFUND.

"Product" means HPDE6c7 Human Pancreatic Ductal Epithelial Cell Line (SCC442).

"Academic Research" means any internal in vitro research use by individuals employed by an academic institution. Academic Research specifically excludes the following uses of whatever kind or nature:

- Re-engineering or copying the Product
- Making derivatives, modifications, or functional equivalents of the Product
- Obtaining patents or other intellectual property rights claiming use of the Product
- Using the Product in the development, testing, or manufacture of a Commercial Product
- Using the Product as a component of a Commercial Product
- Reselling or licensing the Product
- Using the Product in clinical or therapeutic applications including producing materials for clinical trials
- Administering the Product to humans
- Using the Product in collaboration with a commercial or non-academic entity

"Commercial Product" means any product intended for: (i) current or future sale; (ii) use in a fee-for-service; or (iii) any diagnostic, clinical, or therapeutic use.

Access to the Product is limited solely to those officers, employees, and students of PURCHASER's academic institution who need access to the Product to perform Academic Research. PURCHASER shall comply with all applicable laws in its use and handling of the Product and shall keep it under reasonably safe and secure conditions to prevent unauthorized use or access.

These use restrictions will remain in effect for as long as PURCHASER possesses the Product.

COMMERCIAL OR NON-ACADEMIC ENTITIES INTERESTED IN PURCHASING OR USING THE PRODUCT MUST CONTACT licensing@emdmillipore.com AND AGREE TO SEPARATE TERMS OF USE PRIOR TO USE OR PURCHASE.

Genetically Modified Organisms (GMO)

This product contains genetically modified organisms.
Este producto contiene organismos genéticamente modificados.
Questo prodotto contiene degli organismi geneticamente modificati.
Dieses Produkt enthält genetisch modifizierte Organismen.
Ce produit contient des organismes génétiquement modifiés.
Dit product bevat genetisch gewijzigde organismen.
Tämä tuote sisältää geneettisesti muutettuja organismeja.
Denna produkt innehåller genetiskt ändrade organismer.

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck operates
as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates.
All other trademarks are the property of their respective owners. Detailed information on
trademarks is available via publicly accessible resources.

© 2023-2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

Document Template 20306518 Ver 6.0

20634641 Ver 3.0, Rev 10May2024, CJ

