

## Product Information

## Extract-N-Amp™ Tissue PCR Kit

## XNAT2, XNAT2R

## Product Description

The Extract-N-Amp™ Tissue PCR Kit for direct PCR contains the reagents needed to rapidly extract and amplify genomic DNA from mouse tails and other animal tissues, buccal swabs, hair shafts, and saliva. Briefly, the DNA is released from the starting material by incubating the sample with a mixture of the Extraction Solution and the Tissue Preparation Solution at room temperature for 10 minutes. There is no need for mechanical disruption, organic extraction, column purification, or precipitation of the DNA.

After adding Neutralization Solution B, the extract is ready for PCR. An aliquot of the neutralized extract is then combined with the Extract-N-Amp™ PCR Reaction Mix and user-provided PCR primers to amplify target DNA. The Extract-N-Amp™ PCR Reaction Mix is a 2X ready mix containing buffer, salts, dNTPs, and *Taq* polymerase. It is optimized specifically for use with the extraction reagents. It also contains the JumpStart™ *Taq* antibody for hot start PCR to enhance specificity but does not contain the inert red dye found in the REExtract-N-Amp™ PCR Reaction Mix.

Reagents Provided	Cat. No.	XNAT2	XNAT2R
		100 Preps, 100 PCRs	1000 Preps, 1000 PCRs
Extraction Solution	E7526	24 mL	240 mL
Tissue Preparation Solution	T3073	3 mL	30 mL
Neutralization Solution B	N3910	24 mL	240 mL
Extract-N-Amp™ PCR Reaction Mix			
This is a 2X PCR reaction mix containing buffer, salts, dNTPs, <i>Taq</i> polymerase, and JumpStart™ <i>Taq</i> antibody.	E3004	1.2 mL	12 mL

## Reagents and Equipment Required

(Not Provided)

- Microcentrifuge tubes (1.5 or 2 mL) or multi-well plate for extractions (200 µL minimal well volume)
- Small dissecting scissors
- Forceps (small to medium in size)
- Buccal swab - Sterile foam tipped applicator (Cat. No. WHAWB100032)
- Sample collection card - Bloodstain card (Cat. No. WHAWB100014)
- Tubes or plate for PCR
- Heat block or thermal cycler at 95 °C
- PCR Primers (Cat. No. OLIGO)
- Thermal cycler
- Water, PCR Reagent (Cat. No. W1754)

## Precautions and Disclaimer

This product is for R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

## Storage

The Extract-N-Amp™ Tissue PCR Kit can be stored at 2 to 8 °C for up to 3 weeks. For long-term storage, greater than 3 weeks, -20 °C is recommended. Do not store in a "frost-free" freezer.

## Procedure

All steps are carried out at room temperature unless otherwise noted.

### DNA Extraction from Mouse Tails, Animal Tissues, Hair, or Saliva

1. Pipette 100 µL of Extraction Solution into a microcentrifuge tube or well of a multi-well plate. Add 25 µL of Tissue Preparation Solution to the tube or well and pipette up and down to mix.

**Note:** If several extractions will be performed, sufficient volumes of Extraction and Tissue Preparation Solutions may be pre-mixed in a ratio of 4:1 up to 2 hours before use.

2. **For fresh or frozen mouse tails:** Rinse the scissors and forceps in 70% ethanol prior to use and between different samples. Place a 0.5–1 cm piece of mouse tail tip (cut end down) into the solution. Mix thoroughly by vortexing or pipetting. Ensure the mouse tail is in solution.

**Note:** For fresh mouse tails, perform extractions within 30 minutes of snipping the tail.

**For animal tissues:** Rinse the scissors or scalpel and forceps in 70% ethanol prior to use and between different samples. Place a 2–10 mg piece of tissue into the solution. Mix thoroughly by vortexing or pipetting. Ensure the tissue is in the solution.

**For hair shafts:** Rinse the scissors and forceps in 70% ethanol prior to use and between different samples. Trim excess off of the hair shaft leaving the root and place sample (root end down) into solution. Only one hair shaft, with root, is required per extraction.

**For Saliva:** Pipette 10 µL of saliva into the solution. Mix thoroughly by vortexing or pipetting.

**For saliva dried on card:** Pipette 50 µL of saliva onto collection card and allow the card to dry. Rinse the punch in 70% ethanol prior to use and between different samples. Punch a disk (preferably 1/8 inch or 3 mm) out of the card from the area with the dried saliva sample. Place disk into the solution. Tap tube or plate on hard surface to ensure disk is in solution for incubation period.

3. Incubate sample at room temperature for 10 minutes.
4. Incubate sample at 95 °C for 3 minutes.  
**Note:** Tissues will not be completely digested at the end of the incubations. This is normal and will not affect performance.
5. Add 100 µL of Neutralization Solution B to sample and mix by vortexing.
6. Store the neutralized tissue extract at 4 °C or use immediately in [PCR amplification](#).

**Note:** For long term storage, remove the undigested tissue or transfer the extracts to new tubes or wells. Extracts may now be stored at 4 °C for at least 6 months without notable loss in most cases.

### DNA Extraction for Buccal Swabs

1. Collect buccal cells on swab and allow the swab to dry. Drying time is approximately 10 to 15 minutes.  
**Note:** Due to the low volume of solution used for DNA extraction, a foam tipped swab should be used. Swabs with fibrous tips, such as cotton or Dacron®, should be avoided because the solution cannot be recovered efficiently.
2. Pipette 200 µL of Extraction Solution into a microcentrifuge tube. Add 25 µL of Tissue Preparation Solution to the tube and pipette up and down to mix.  
**Note:** If several extractions will be performed, sufficient volumes of Extraction and Tissue Preparation Solutions may be pre-mixed in a ratio of 8:1 up to 2 hours before use.
3. Place dried buccal swab into solution and incubate at room temperature for 1 minute.
4. Twirl swab in solution 10 times and then remove excess solution from the swab into the tube by twirling swab firmly against the side of the tube. Discard the swab. Close the tube and vortex briefly.
5. Incubate sample at room temperature for 10 minutes.
6. Incubate sample at 95 °C for 3 minutes.
7. Add 200 µL of Neutralization Solution B to sample and mix by vortexing.

- Store the neutralized extract at 4 °C or use immediately in PCR. Continue to PCR amplification.

**Note:** Extracts may be stored at 4 °C for at least 6 months without notable loss in most cases.

## PCR Amplification

The Extract-N-Amp™ PCR Reaction Mix contains JumpStart™ *Taq* antibody for specific hot start amplification. Therefore, PCR mixtures can be assembled at room temperature without premature *Taq* DNA polymerase activity.

Typical final primer concentrations are approximately 0.4 μM each. The optimal primer concentration and cycling parameters will depend on the system being used.

- Add the following reagents to a thin-walled PCR microcentrifuge tube or plate:

Reagent	Volume
Water, PCR grade	Variable
Extract-N-Amp™ PCR reaction mix	10 μL
Forward primer	Variable
Reverse primer	Variable
Tissue extract	4 μL*
Total volume	20 μL

\*The Extract-N-Amp™ PCR Reaction Mix is formulated to compensate for components in the Extraction, Tissue Preparation, and Neutralization Solutions. If less than 4 μL of tissue extract is added to the PCR reaction volume, use a 50:50 mixture of Extraction and Neutralization B

Solutions to bring the volume of tissue extract up to 4 μL.

- Mix gently.
- For thermal cyclers without a heated lid, add 20 μL of mineral oil on top of the mixture in each tube to prevent evaporation.
- Perform thermal cycling. The amplification parameters should be optimized for individual primers, template, and thermal cycler.

### Common cycling parameters:

Step	Temperature	Time	Cycles
Initial Denaturation	94 °C	3 minutes	1
Denaturation	94 °C	30 seconds	
Annealing	45 to 68 °C	30 seconds	30-35
Extension	72 °C	1-2 minutes (1 min/kb)	
Final Extension	72 °C	10 minutes	1
Hold	4 °C	Indefinitely	

- The amplified DNA can be loaded onto an agarose gel after the PCR is completed with the addition of a separate loading buffer/tracking dye such as Gel Loading Solution, Cat. No. G2526.

**Note:** PCR products can be purified, if desired, for downstream applications such as sequencing with the GenElute™ PCR Clean-Up Kit, Cat. No. NA1020.

## Troubleshooting Guide

Problem	Cause	Solution
Little or no PCR product is detected.	PCR reaction may be inhibited due to contaminants in the tissue extract.	Dilute the tissue extract with a 50:50 mix of Extraction and Neutralization Solutions. To test for inhibition, include a DNA control and/or spike a known amount of template (100-500 copies) into the PCR along with the tissue extract.
	Extraction is insufficient.	Incubate samples at 55 °C for 10 minutes instead of room temperature.
	A PCR component may be missing or degraded.	Run a positive control to ensure that components are functioning. A checklist is also recommended when assembling reactions.
	There may be too few cycles performed.	Increase the number of cycles (5-10 additional cycles at a time).
	The annealing temperature may be too high.	Decrease the annealing temperature in 2-4 °C increments.
	The primers may not be designed optimally.	Confirm the accuracy of the sequence information. If the primers are less than 22 nucleotides long, try to lengthen the primer to 25-30 nucleotides. If the primer has a GC content of less than 45%, try to redesign the primer with a GC content of 45-60%.
	The extension time may be too short.	Increase the extension time in 1-minute increments, especially for long templates.
Multiple products	Target template is difficult.	In most cases, inherently difficult targets are due to unusually high GC content and/or secondary structure. Betaine, Cat. No. B0300, has been reported to help amplification of high GC content templates at a concentration of 1.0-1.7 M.
	JumpStart™ <i>Taq</i> antibody is not working correctly.	Do not use DMSO or formamide with Extract-N-Amp™ PCR Reaction Mix. It can interfere with the enzyme-antibody complex. Other cosolvents, solutes (e.g., salts), and extremes in pH or other reaction conditions may reduce the affinity of the JumpStart™ <i>Taq</i> antibody for <i>Taq</i> polymerase and thereby compromise its effectiveness.
	Touchdown PCR may be needed.	"Touchdown" PCR significantly improves the specificity of many PCR reactions in various applications. Touchdown PCR involves using an annealing/extension temperature that is higher than the TM of the primers during the initial PCR cycles. The annealing/extension temperature is then reduced to the primer TM for the remaining PCR cycles. The change can be performed in a single step or in increments over several cycles.
Negative control shows a PCR product or "false positive" result.	Reagents are contaminated.	Include a reagent blank without DNA template be included as a control in every PCR run to determine if the reagents used in extraction or PCR are contaminated with a template from a previous reaction.
Tissue is not digested after incubations.	Tissue is not expected to be completely digested.	The REExtract-N-Amp™ Tissue PCR Kit does not require the tissue to be completely digested. Sufficient DNA is released for PCR without completely digesting the tissue.
Buccal swab absorbed all the solution.	The recommended type of swab was not used.	Due to the low volume of solution used for DNA extraction, a foam tipped swab should be used. Swabs with fibrous tips, such as cotton or Dacron®, should be avoided because the solution cannot be recovered efficiently.

## References

1. Dieffenbach, C.W., and Dveksler, G.S. (Eds.), PCR Primer: A Laboratory Manual, 2<sup>nd</sup> ed., Cold Spring Harbor Laboratory Press, New York (1995).
2. Don, R.H. et al., 'Touchdown' PCR to circumvent spurious priming during gene amplification. *Nucleic Acids Res.*, **19**, 4008 (1991).
3. Erlich, H.A. (Ed.), PCR Technology: Principles and Applications for DNA Amplification, Stockton Press, New York (1989).
4. Griffin, H.G., and Griffin, A.M. (Eds.), PCR Technology: Current Innovations, CRC Press, Boca Raton, FL (1994).
5. Innis, M.A., et al., (Eds.), PCR Strategies, Academic Press, New York (1995).
6. Innis, M., et al., (Eds.), PCR Protocols: A Guide to Methods and Applications, Academic Press, San Diego, California (1990).
7. McPherson, M.J. et al., (Eds.), PCR 2: A Practical Approach, IRL Press, New York (1995).
8. Newton, C.R. (Ed.), PCR: Essential Data, John Wiley & Sons, New York (1995).
9. Roux, K.H. Optimization and troubleshooting in PCR. *PCR Methods Appl.*, **4**, 5185-5194 (1995).
10. Saiki, R., PCR Technology: Principles and Applications for DNA Amplification, Stockton, New York (1989).

## Product Ordering

Order products online at [SigmaAldrich.com](https://www.sigmaaldrich.com)

Related Products	Cat. No.
Ethanol	E7148; E7023; 459836
Forceps, micro-dissecting	F4267
PCR Marker	P9577
PCR microtubes	Z374873; Z374962; Z374881
PCR multi-well plates	Z374903
Precast Agarose Gels	P6097
Sealing mats & tapes	Z374938; A2350
TBE Buffer	T4415, T6400, T9525

---

## Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

### Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice).

### Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at [SigmaAldrich.com/terms](https://SigmaAldrich.com/terms).

### Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://SigmaAldrich.com/offices).

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Merck, Extract-N-Amp, REExtract-N-Amp, JumpStart, GenElute and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

XNAT2pis Rev 09/22

