SIGMA-ALDRICH®

sigma-aldrich.com

3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

Product Information

StableCell™ Trypsin Solution

10×, 5.0 g porcine trypsin and 2 g EDTA • 4Na per liter of 0.9% sodium chloride BioReagent, sterile-filtered, suitable for cell culture

Catalog Number **T2610** Storage Temperature 2-8 °C

EC 3.4.21.4 Synonyms: Tryptase, Tryptar, Cocoonase, Parenzyme, Parenzymol

Product Description

Molecular mass:¹ 23.4 kDa Extinction Coefficient:² $E^{1\%}$ = 15.0 (280 nm) pl:^{1,2} 10.2–10.8

StableCell[™] Trypsin 10× Solution (5.0 g/L of porcine trypsin and 2.0 g/L of EDTA • 4Na in 0.9% sodium chloride, plus a proprietary protein stabilizing agent) is sterile-filtered and cell culture tested.

Trypsin consists of a single chain polypeptide of 223 amino acid residues. Trypsin is produced by the removal of the N-terminal hexapeptide from trypsinogen which is cleaved at the Lys⁶-Ile⁷ peptide bond. The amino acid sequence of trypsin is crosslinked by 6 disulfide bridges. This native form of trypsin is referred to as β -trypsin. Autolysis of β -trypsin (which is cleaved at Lys¹³¹-Ser¹³²) results in α -trypsin which is held together by disulfide bridges. Trypsin is a member of the serine protease family. The active site amino acid residues of trypsin include His⁴⁶ and Ser¹⁸³.^{1,3}

Trypsin will cleave peptides on the C-terminal side of lysine and arginine amino acid residues. The rate of hydrolysis is slower if an acidic residue is on either side of the cleavage site and no cleavage occurs if a proline residue is on the carboxyl side of the cleavage site. The pH optimum of trypsin is 7–9.² Trypsin will also hydrolyze ester and amide linkages of synthetic derivatives of amino acids such as: benzoyl L-arginine ethyl ester (BAEE), p-toluenesulfonyl-L-arginine methyl ester (TAME), tosyl-L-arginine methyl ester, N α -benzoyl-L-arginine p-nitroanilide (BAPNA), L-lysyl-p-nitroanilide, and benzoyl-L-arginamide.^{2–5} Assuming the pH and temperature are the same and using a molar extinction coefficient of 808 at 254 nm for BAEE, the following conversions are valid:

1 BAEE μ M Unit = 200 BAEE Units 1 TAME μ M Unit = 0.27 BAEE μ M Units 1 BAEE μ M Unit = 3.64 TAME Units 1 TAME μ M Unit = 55 BAEE A₂₅₃ Units 1 BAEE A₂₅₃ Unit = 0.018 TAME μ M Unit 1 TAME μ M Unit = 180 TAME A₂₄₇ Units 1 TAME A₂₄₇ Unit = 0.33 BAEE Units A USP Unit = Δ A₂₅₃ of 0.003 per minute 1 NF Unit = 3.3 A₂₅₃ BAEE Units.⁹

<u>Note</u>: These activity conversions were determined using bovine trypsin; however, they are thought to be similar for porcine trypsin.

Serine protease inhibitors that will inhibit trypsin include DFP (diisopropyl fluorophosphate), TLCK (N α -p-tosyl-L-lysine chloromethyl ketone), PMSF (phenylmethanesulfonyl fluoride), APMSF(4-amidinophenylmethanesulfonyl fluoride), AEBSEF (4-(2-aminoethyl)benzenesulfonyl fluoride), aprotinin, leupeptin, α_2 -macroglobulin, α_1 -antitrypsin, *p*-aminobenzamidine, benzamidine (reversible), soybean trypsin inhibitor, lima bean inhibitor, bovine pancreas trypsin inhibitor, chicken egg white inhibitor, and turkey egg white inhibitor.^{1,7}

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

Recommended storage is 2-8 °C upon arrival. During stability studies, data showed this product retains $\geq 90\%$ of its activity when stored at 37 °C for up to 8 weeks.

Procedure

StableCell[™] Trypsin Solution may be used to remove adherent cells from a culture surface. Cells are most commonly removed from the culture substrate by treatment with trypsin or trypsin • EDTA. StableCell[™] Trypsin Solutions can range from 0.05% to 0.5%. The reasons for the range of concentrations are as follows:
(1) Differences in trypsin activity or potency,
(2) Different incubation times, and (3) Different cell lines.

Incubating cells with too high a trypsin concentration for too long a time period will damage cell membranes and kill the cells. If unsure about the concentration of trypsin to use, use a low concentration. There can be lot-to-lot variation in dissociation times which is to be expected since the enzymatic activity of each lot will differ. If trypsin is being solubilized or diluted from a concentrated solution, it is important to use a buffered salt solution that contains no Ca²⁺ or Mg²⁺, such as Hank's Balanced Salt Solution, Modified (Catalog No. H9394). Adjust the pH of trypsin solution to 7.4–7.6.

- Remove medium from culture vessel by aspiration and wash the monolayer with Ca⁺² and Mg⁺²-free salt solution to remove all traces of serum. Remove salt solution by aspiration.
- Dispense enough StableCell[™] Trypsin Solution into culture vessel(s) to completely cover the monolayer of cells and place in 37 °C incubator for ~2 minutes.
- Remove the StableCell[™] Trypsin Solution by aspiration and return closed culture vessel(s) to incubator. The coated cells are allowed to incubate until cells detach from the surface. Progress can be checked by examination with an inverted microscope.

<u>Note</u>: The time required to remove cells from the culture surface is dependent on cell type, population density, serum concentration in the growth medium, potency of trypsin, and time since last subculture. Trypsin can cause cellular damage, thus time of exposure should be kept to a minimum.

4. When trypsinization process is complete the cells will be in suspension and appear rounded.

- 5. It is advisable to add serum or medium containing serum to the cell suspension as soon as possible to inhibit further tryptic activity which may damage cells. Soybean trypsin inhibitor (Catalog No. T6414) can also be added at an equimolar concentration to inhibit the trypsin that is present. Soybean trypsin inhibitor is used when culturing in serum-free conditions.
- Cells can be resuspended by gently pipetting the cell suspension to break up the clumps. Further dilution can be made, if required, for cell counts and/or subculturing.

References

- 1. Walsh, K.A., Trypsinogens and trypsins of various species. Meth. Enzymol., **19**, 41-63 (1970).
- Buck, F.F. et al., On the mechanism of enzyme action. LXXII. Studies on trypsins from beef, sheep, and pig pancreas. Arch. Biochem. Biophys., 97, 417-424 (1962).
- Keil, B., in The Enzymes, 3rd ed., Vol. III, Boyer, P. D., Academic Press (New York, NY: 1971), pp. 250-275.
- Burdon, R.H. et al., in Laboratory Techniques in Biochemistry and Molecular Biology, Vol. 9, 2nd ed., Allen, G., ed., Elsevier/North (New York, NY: 1989), pp. 73-104.
- Enzyme Handbook, Vol. II, Barman, T. E., Springer-Verlag (New York, NY: 1969), pp. 618-619.
- 6. USP, Vol. 23, p. 1611.
- Proteolytic Enzymes, A Practical Approach, Beynon, R. J., ed., IRL Press (New York, NY: 1989), p. 240.
- 8. Methods of Molecular Biology, Vol. 3, Smith, B. J., Humana Press, (New Jersey, 1988), pp. 57-69.
- Porter, W.H., and Preston, J.L., Retention of trypsin and chymotrypsin proteolytic activity in sodium dodecyl sulfate solutions. Anal. Biochem., 66, 69-77 (1975).

StableCell is a trademark of Sigma-Aldrich Co. LLC.

BG,MAM 10/17-1

©2017 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA-ALDRICH is a trademark of Sigma-Aldrich Co. LLC, registered in the US and other countries. Sigma brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.