





SIGMA-ALDRICH

3050 Spruce Street
 Saint Louis, Missouri 63103 USA
 Telephone 800-325-5832 • (314) 771-5765
 Fax (314) 286-7828
 email: techserv@sial.com
 sigma-aldrich.com

Product Information

SEPHAROSE AND SEPHAROSE CL GEL FILTRATION MEDIA Exact replacement for Product Code 84965

The information below is abstracted, for customer convenience, from a supplier technical manual and a supplier information sheet. Complete copies of supplier data are available upon request from Sigma Research Technical Service.

SEPHAROSE	SEPHAROSE CL
 <p>D-Galactose 3,6-Anhydro-L-Galactose Partial structure of Sepharose CL</p>	 <p>D-Galactose 3,6-Anhydro-L-Galactose Partial structure of Sepharose CL</p>
CHARACTERISTICS	
<ul style="list-style-type: none"> - Broad fractionation range - High exclusion limits - Negligible non-specific adsorption - Appearance: white suspension 	<ul style="list-style-type: none"> - Broad fractionation range - High exclusion limits - Negligible non-specific adsorption - Excellent chemical and physical stability - Appearance: white suspension
PRODUCT DESCRIPTION	
<p>Sepharose is a beaded agarose gel filtration medium with a broad fractionation range. Three different agarose contents are available: 2%,</p>	<p>Sepharose CL is a cross-linked derivative of Sepharose, prepared by reacting Sepharose with 2,3-dibromopropanol under strongly alkaline</p>

4% and 6%, designated 2B, 4B and 6B, respectively. As agarose concentration increases porosity decreases, thus increasing rigidity and altering the fractionation range; nucleic acids and polysaccharides with molecular weights up to $\sim 4 \times 10^7$ can be separated on Sepharose 2B.

conditions. After cross-linking, the gel is desulfated by alkaline hydrolysis under reducing conditions. The resulting crosslinked polysaccharide chains have a very low content of ionizable groups and better chemical and physical resistance than Sepharose, as well as improved flow properties. Porosity, however, is comparable to that of Sepharose. Three different agarose contents are available: 2%, 4% and 6%, designated CL-2B, CL-4B and CL-6B, respectively.

Property	2B	4B	6B	Property	CL-2B	CL-4B	CL-6B
MW Range:				MW Range:			
Globular Proteins	7X10 ⁴ -4X10 ⁷	6X10 ⁴ -2X10 ⁷	1X10 ⁴ -4X10 ⁶	Globular Proteins	7X10 ⁴ -4X10 ⁷	6X10 ⁴ -2X10 ⁷	1X10 ⁴ -4X10 ⁶
Dextrans	1X10 ⁵ -2X10 ⁷	3X10 ⁴ -5X10 ⁶	1X10 ⁴ -1X10 ⁶	Dextrans	1X10 ⁵ -2X10 ⁷	3X10 ⁴ -5X10 ⁶	1X10 ⁴ -1X10 ⁶
DNA exclusion limit	~ 1352 base pairs	~ 872 base pairs	$\sim 45-165$ base pairs	DNA exclusion limit	~ 1353 base pairs	~ 872 base pairs	$\sim 45-165$ base pairs
Diameter (wet bead)	$\sim 60-200 \mu\text{m}$	$\sim 45-165 \mu\text{m}$	$\sim 45-154 \mu\text{m}$	Diameter (wet bead)	$\sim 60-200 \mu\text{m}$	$\sim 45-165 \mu\text{m}$	$\sim 45-165 \mu\text{m}$
pH Range	4-9	4-9	4-9	pH Range	3-14	3-14	3-14
Max. Pressure	40 cm H ₂ O	80 cm H ₂ O	200 cm H ₂ O	Max. Pressure	50 cm H ₂ O	120 cm H ₂ O	>200 cm H ₂ O
Max. Volumetric Flow Rate	0.83 mL/min	0.96 mL/min	1.16 mL/min	Max. Volumetric Flow Rate	1.25 mL/min	2.17 mL/min	2.5 mL/min
Max. Linear Flow Rate	10 mL/cm ² h	11.5 mL/cm ² h	14 mL/cm ² h	Max. Linear Flow Rate	15 mL/cm ² h	26 mL/cm ² h	30 mL/cm ² h

STABILITY

Sepharose melts upon heating to 40°C, cannot be autoclaved, and the bead structure may be damaged upon freezing. Due to the presence of 3,6-anhydro-L-galactose, the matrix is resistant to biological degradation. Sepharose is stable in aqueous (including saline) solutions at pH 4-9. Use of dissociation media such as guanidine hydrochloride and urea, chaotropic salts such as KSCN, and oxidizing agents is not advisable because these reagents may disrupt the hydrogen bonds which stabilize the matrix. Sterilization can be done chemically (e.g., by treatment with DEPC).

Sepharose CL is stable in aqueous media at pH 3-14, and in organic solvents. Dissociating media and chaotropic salts may be used, but oxidizing conditions should be avoided in order to preserve the identity of the constituent sugar residues. Sepharose CL may be autoclaved repeatedly at 120°C and pH 7. The matrix is resistant to biological degradation. Flow rates may be up to 50% higher than on Sepharose with an equivalent fractionation range.

USAGE AND REGENERATION

Sepharose and Sepharose CL are supplied pre-swollen as suspensions in distilled water. Before packing a column, dilute the required amount of gel with starting buffer to form a thick slurry, about 75% of which is settled resin, then degas the slurry. The maximum operating pressures (shown above) should never be exceeded. Pass 2-3 column volumes (CV) of eluent through the gel to equilibrate the bed. Elution may be achieved by gravity feed or through

the use of a peristaltic pump; more consistent flow rates and more reproducible separations are obtained with a pump. Sepharose contains a small number of ionic sulfate and carboxyl groups which may cause adsorption of basic proteins at low ionic strengths. Therefore, eluents with ionic strengths exceeding 0.02 M are sometimes necessary (for example, tRNA species have been resolved in high concentrations of ammonium sulfate, and DNA has been separated from RNA in 1.5 N NaCl). The gels can be cleaned as indicated below and stored at 4-8°C in a suitable antimicrobial agent (e.g., 20% ethanol) for indefinite time periods.

Sepharose should be cleaned in the column or batchwise with a non-ionic detergent solution.

Sepharose CL can be cleaned in the column or batchwise. Wash with at least 1 CV of 0.5 N NaCl in 0.1 N NaOH (a step which also sterilizes the resin), then with 10 CV of water (or until the eluent is at neutral pH).