



# THE SELECTIVITY PROFILE OF 35 SMALL MOLECULE PROTEIN KINASE INHIBITORS

## INTRODUCTION

Protein kinases regulate most aspects of normal cellular life, such as differentiation and proliferation, whereas inappropriate phosphorylation is a cause of many diseases. As the regulation and roles of protein kinases have become better understood, the interest in developing orally active protein kinase inhibitors has increased, and has recently been validated by the success of several small molecule and antibody therapies (Gleevec®, Iressa®, Tarceva®, Herceptin®, Erbitux®, etc). With around 500 distinct protein kinases identified through the human genome project and around 30 kinase inhibitors in clinical trials or approved for human use, this class of protein continues to be a rich area for therapeutic drug development.

## KINASE INHIBITOR PROFILING

While kinase inhibitors hold great promise for the treatment of many diseases, target selectivity remains a common problem in the development process. Inhibition of off-target kinases could produce undesired side-effects or toxicity, and as the development of such compounds continues, questions are arising around how to design multi-targeted kinase inhibitors.

To facilitate this by allowing a better understanding of compound inhibition patterns and structure-activity relationships (SAR), Upstate has used its unparalleled experience in the cloning, manufacturing and assaying of protein kinases to assemble the **KinaseProfiler™** panel. This panel currently contains more than 240

kinases and is increasing in breadth continually. It already represents the largest and most disease-relevant panel in the industry, and we have augmented this with the introduction of clinically-relevant mutant versions of standard kinases. All assays are conducted in the gold-standard radiometric filter-binding format, so there is no deconvolution, interpretation or compound interference to worry about.

Customers can use the service to screen large numbers of compounds against a single kinase or screen that vital compound against the whole panel. They can observe and examine inhibition patterns to enable them to fully understand the SAR of their compound library, learn more about the specificities of their lead compounds, see if their leads hit any crucial, 'survival' kinases or determine rank orders of potency using our IC<sub>50</sub> estimation service. The overall aim is to give them the flexibility to more quickly design potent and targeted compounds to enable the production of more effective drugs.

In order to demonstrate the type of information which can be obtained using **KinaseProfiler™**, we challenged our entire panel with 35 commercially available compounds which are reported to be relatively selective inhibitors of various protein kinases, and this data is presented on pages 2 and 3.



KINASE INHIBITION DATA GENERATED BY KINASEPROFILER™

|                 | PP2<br>[1µM] | PP3<br>[1µM] | LY294002<br>[50µM] | SB203580<br>[10µM] | HA1077<br>[10µM] | Roscovitine<br>[10µM] | Rotterlin<br>[10µM] | EGCG<br>[10µM] | K252c<br>[10µM] | Staurosporine<br>[1µM] | SU6656<br>[1µM] | AG1296<br>[10µM] | Olomoucine<br>[50µM] | Ro31-8220<br>[1µM] | AG183<br>[10µM] | AG18<br>[10µM] | AG1024<br>[10µM] | AG538<br>[10µM] | AG1478<br>[10µM] | H89<br>[10µM] | Y27632<br>[10µM] | Alsterpaullone<br>[1µM] | Indinbin<br>monocine<br>[10µM] | CDK2/5<br>inhibitor<br>[10µM] | JNK<br>inhibitor II<br>[10µM] | KT5720<br>[1µM] | ST638<br>[10µM] | ZM336372<br>[10µM] | GW583340<br>[10µM] | GW441756<br>[10µM] | API-2<br>[10µM] | PD153035<br>[1µM] | JAK<br>inhibitor I *<br>[1µM] | UO126<br>[10µM] |    |
|-----------------|--------------|--------------|--------------------|--------------------|------------------|-----------------------|---------------------|----------------|-----------------|------------------------|-----------------|------------------|----------------------|--------------------|-----------------|----------------|------------------|-----------------|------------------|---------------|------------------|-------------------------|--------------------------------|-------------------------------|-------------------------------|-----------------|-----------------|--------------------|--------------------|--------------------|-----------------|-------------------|-------------------------------|-----------------|----|
| AMH(h)          | 96           | 0            | 0                  | 24                 | 0                | 0                     | 49                  | 0              | 1               | 93                     | 81              | 0                | 0                    | 0                  | 4               | 0              | 99               | 73              | 95               | 0             | 0                | 67                      | 92                             | 7                             | 64                            | 0               | 22              | 43                 | 0                  | 90                 | 0               | 78                | 0                             | 0               |    |
| AMPK(ε)         | 1            | 7            | 8                  | 2                  | 27               | 3                     | 0                   | 38             | 90              | 98                     | 86              | 7                | 10                   | 80                 | 51              | 13             | 65               | 46              | 31               | 99            | 22               | 56                      | 98                             | 25                            | 84                            | 23              | 12              | 0                  | 0                  | 75                 | 0               | 0                 | 75                            | 0               |    |
| Aurora-A(h)     | 7            | 0            | 3                  | 3                  | 2                | 0                     | 93                  | 10             | 94              | 95                     | 99              | 3                | 26                   | 3                  | 40              | 31             | 56               | 62              | 69               | 51            | 0                | 6                       | 100                            | 9                             | 94                            | 98              | 90              | 2                  | 17                 | 100                | 0               | 37                | 77                            | 0               |    |
| Axl(h)          | 0            | 1            | 21                 | 0                  | 3                | 11                    | 68                  | 47             | 78              | 93                     | 89              | 11               | 8                    | 57                 | 22              | 6              | 40               | 34              | 56               | 25            | 0                | 22                      | 93                             | 0                             | 78                            | 64              | 36              | 0                  | 14                 | 92                 | 0               | 34                | 46                            | 0               |    |
| Bms(h)          | 89           | 11           | 22                 | 39                 | 0                | 0                     | 0                   | 1              | 27              | 93                     | 0               | 0                | 0                    | 44                 | 16              | 0              | 0                | 0               | 51               | 3             | 0                | 0                       | 73                             | 0                             | 53                            | 47              | 0               | 0                  | 0                  | 22                 | 0               | 13                | 49                            | 0               |    |
| CDK1/cyclinB(h) | 0            | 0            | 16                 | 0                  | 0                | 77                    | 0                   | 2              | 65              | 98                     | 55              | 0                | 53                   | 96                 | 11              | 0              | 7                | 17              | 2                | 7             | 3                | 83                      | 92                             | 47                            | 80                            | 3               | 0               | 0                  | 0                  | 44                 | 0               | 0                 | 9                             | 0               |    |
| CDK2/cyclinA(h) | 2            | 0            | 18                 | 0                  | 3                | 93                    | 0                   | 0              | 81              | 98                     | 32              | 0                | 84                   | 97                 | 15              | 0              | 3                | 8               | 0                | 20            | 24               | 93                      | 94                             | 78                            | 81                            | 21              | 0               | 0                  | 0                  | 32                 | 0               | 0                 | 16                            | 2               |    |
| CDK2/cyclinE(h) | 12           | 0            | 28                 | 0                  | 10               | 96                    | 0                   | 23             | 77              | 98                     | 26              | 0                | 92                   | 98                 | 36              | 2              | 0                | 17              | 0                | 50            | 44               | 96                      | 99                             | 88                            | 86                            | 14              | 8               | 0                  | 0                  | 47                 | 0               | 0                 | 20                            | 0               |    |
| CDK5p35(h)      | 4            | 0            | 33                 | 8                  | 20               | 97                    | 7                   | 24             | 93              | 99                     | 60              | 0                | 91                   | 96                 | 17              | 1              | 2                | 21              | 0                | 48            | 32               | 95                      | 98                             | 89                            | 89                            | 19              | 0               | 0                  | 0                  | 90                 | 11              | 0                 | 23                            | 0               |    |
| CHK1(h)         | 0            | 0            | 0                  | 2                  | 6                | 0                     | 0                   | 0              | 88              | 93                     | 45              | 0                | 0                    | 82                 | 15              | 5              | 0                | 28              | 16               | 79            | 0                | 77                      | 89                             | 25                            | 86                            | 82              | 0               | 0                  | 0                  | 63                 | 0               | 6                 | 83                            | 0               |    |
| CHK2(h)         | 11           | 0            | 3                  | 5                  | 51               | 29                    | 82                  | 19             | 94              | 93                     | 80              | 0                | 73                   | 57                 | 39              | 4              | 11               | 51              | 70               | 74            | 4                | 90                      | 95                             | 7                             | 97                            | 57              | 50              | 48                 | 30                 | 90                 | 0               | 42                | 94                            | 0               |    |
| CK1_α(h)        | 92           | 45           | 79                 | 97                 | 13               | 85                    | 38                  | 26             | 47              | 21                     | 33              | 65               | 50                   | 8                  | 14              | 9              | 80               | 35              | 64               | 48            | 4                | 11                      | 96                             | 9                             | 96                            | 0               | 38              | 26                 | 2                  | 96                 | 17              | 39                | 23                            | 27              |    |
| CK2(h)          | 0            | 8            | 93                 | 0                  | 0                | 0                     | 0                   | 0              | 49              | 0                      | 0               | 0                | 10                   | 0                  | 71              | 85             | 0                | 71              | 1                | 0             | 0                | 0                       | 20                             | 0                             | 72                            | 0               | 0               | 0                  | 0                  | 37                 | 3               | 0                 | 17                            | 0               |    |
| ckr(h)          | 63           | 9            | 58                 | 0                  | 0                | 16                    | 0                   | 93             | 67              | 95                     | 79              | 84               | 12                   | 24                 | 82              | 0              | 6                | 39              | 25               | 0             | 0                | 5                       | 97                             | 5                             | 81                            | 81              | 0               | 33                 | 96                 | 0                  | 89              | 0                 | 4                             | 0               | 9  |
| CSK(h)          | 38           | 15           | 0                  | 11                 | 3                | 13                    | 10                  | 101            | 38              | 74                     | 13              | 0                | 6                    | 11                 | 79              | 0              | 0                | 51              | 15               | 7             | 5                | 1                       | 25                             | 4                             | 21                            | 0               | 3               | 6                  | 0                  | 2                  | 8               | 9                 | 0                             | 0               |    |
| cSRC(h)         | 97           | 40           | 26                 | 74                 | 2                | 9                     | 1                   | 21             | 47              | 96                     | 97              | 12               | 0                    | 12                 | 24              | 3              | 90               | 53              | 96               | 9             | 1                | 82                      | 97                             | 29                            | 62                            | 13              | 75              | 18                 | 11                 | 97                 | 0               | 81                | 21                            | 3               |    |
| DAPK1(h)        | 0            | 0            | 45                 | 6                  | 22               | 0                     | 0                   | 1              | 40              | 98                     | 0               | 0                | 15                   | 46                 | 9               | 1              | 59               | 18              | 10               | 46            | 0                | 0                       | 81                             | 0                             | 71                            | 0               | 4               | 0                  | 0                  | 11                 | 0               | 0                 | 95                            | 0               |    |
| DYRK2(h)        | 2            | 9            | 26                 | 1                  | 0                | 6                     | 0                   | 69             | 91              | 50                     | 0               | 0                | 64                   | 6                  | 15              | 49             | 12               | 74              | 53               | 7             | 0                | 41                      | 78                             | 14                            | 96                            | 7               | 0               | 0                  | 40                 | 33                 | 27              | 25                | 8                             |                 |    |
| EGFR(h)         | 51           | 5            | 17                 | 54                 | 13               | 0                     | 0                   | 97             | 24              | 37                     | 0               | 0                | 15                   | 6                  | 80              | 0              | 0                | 19              | 98               | 15            | 3                | 0                       | 19                             | 0                             | 0                             | 0               | 68              | 3                  | 97                 | 0                  | 0               | 98                | 0                             | 0               |    |
| EphB4(h)        | 95           | 9            | 16                 | 56                 | 0                | 1                     | 0                   | 98             | 3               | 36                     | 9               | 0                | 0                    | 5                  | 73              | 0              | 0                | 0               | 89               | 2             | 0                | 0                       | 85                             | 0                             | 23                            | 0               | 11              | 53                 | 0                  | 43                 | 0               | 75                | 0                             | 0               |    |
| Flt1(h)         | 74           | 0            | 13                 | 75                 | 30               | 0                     | 27                  | 41             | 94              | 91                     | 62              | 33               | 15                   | 83                 | 20              | 4              | 25               | 67              | 87               | 84            | 0                | 0                       | 50                             | 91                            | 18                            | 94              | 64              | 30                 | 7                  | 0                  | 94              | 0                 | 68                            | 59              | 10 |
| Flt3(h)         | 30           | 5            | 44                 | 26                 | 59               | 21                    | 40                  | 34             | 83              | 94                     | 75              | 71               | 33                   | 84                 | 26              | 6              | 51               | 66              | 71               | 68            | 6                | 69                      | 93                             | 46                            | 84                            | 83              | 32              | 0                  | 3                  | 97                 | 8               | 53                | 76                            | 10              |    |
| Flt4(h)         | 88           | 0            | 69                 | 62                 | 46               | 28                    | 62                  | 37             | 98              | 98                     | 82              | 79               | 70                   | 96                 | 35              | 13             | 77               | 76              | 96               | 88            | 0                | 71                      | 100                            | 11                            | 98                            | 86              | 62              | 78                 | 0                  | 0                  | 99              | 0                 | 84                            | 76              | 4  |
| GSK3_β(h)       | 0            | 0            | 65                 | 37                 | 0                | 0                     | 23                  | 0              | 32              | 91                     | 0               | 0                | 26                   | 95                 | 68              | 3              | 41               | 22              | 0                | 0             | 0                | 95                      | 92                             | 0                             | 50                            | 19              | 17              | 0                  | 0                  | 42                 | 0               | 0                 | 2                             | 0               |    |
| GSK3β(h)        | 0            | 0            | 56                 | 64                 | 0                | 0                     | 20                  | 0              | 21              | 89                     | 0               | 0                | 12                   | 89                 | 39              | 2              | 38               | 10              | 0                | 0             | 0                | 90                      | 86                             | 0                             | 29                            | 0               | 0               | 0                  | 0                  | 27                 | 0               | 0                 | 0                             | 0               |    |
| Hsk(h)          | 90           | 21           | 32                 | 51                 | 0                | 0                     | 76                  | 0              | 28              | 87                     | 80              | 0                | 0                    | 0                  | 11              | 0              | 83               | 46              | 81               | 0             | 0                | 36                      | 91                             | 16                            | 58                            | 24              | 40              | 80                 | 0                  | 84                 | 4               | 70                | 0                             | 2               |    |
| IGF-1R(h)       | 0            | 0            | 10                 | 0                  | 0                | 0                     | 0                   | 53             | 32              | 63                     | 0               | 0                | 0                    | 0                  | 15              | 0              | 0                | 0               | 0                | 0             | 0                | 85                      | 0                              | 30                            | 0                             | 16              | 0               | 0                  | 0                  | 55                 | 1               | 0                 | 7                             | 3               |    |
| IR(h)           | 0            | 0            | 5                  | 0                  | 0                | 7                     | 0                   | 0              | 29              | 41                     | 0               | 0                | 0                    | 0                  | 0               | 0              | 0                | 0               | 0                | 0             | 0                | 89                      | 0                              | 40                            | 0                             | 6               | 0               | 0                  | 0                  | 74                 | 0               | 0                 | 14                            | 0               |    |
| KDR(h)          | 59           | 8            | 36                 | 54                 | 14               | 9                     | 0                   | 36             | 93              | 93                     | 61              | 67               | 26                   | 68                 | 25              | 18             | 42               | 42              | 78               | 63            | 9                | 42                      | 91                             | 9                             | 90                            | 50              | 38              | 29                 | 3                  | 92                 | 16              | 60                | 34                            | 0               |    |
| Lsk(h)          | 96           | 0            | 1                  | 35                 | 0                | 0                     | 0                   | 77             | 0               | 98                     | 99              | 0                | 0                    | 1                  | 0               | 0              | 76               | 77              | 85               | 0             | 0                | 25                      | 94                             | 0                             | 7                             | 0               | 29              | 91                 | 0                  | 90                 | 0               | 52                | 32                            | 0               | 0  |
| Lyn(h)          | 99           | 55           | 62                 | 80                 | 12               | 11                    | 44                  | 95             | 95              | 99                     | 99              | 28               | 0                    | 32                 | 14              | 3              | 93               | 64              | 99               | 39            | 16               | 97                      | 99                             | 31                            | 96                            | 39              | 73              | 94                 | 0                  | 98                 | 0               | 97                | 31                            | 0               |    |
| MAPK1(h)        | 0            | 4            | 14                 | 4                  | 0                | 19                    | 32                  | 62             | 56              | 35                     | 0               | 0                | 41                   | 13                 | 1               | 0              | 0                | 3               | 18               | 26            | 0                | 8                       | 11                             | 0                             | 20                            | 19              | 0               | 0                  | 0                  | 6                  | 0               | 0                 | 56                            | 8               |    |
| MAPKAP-K2(h)    | 0            | 0            | 0                  | 0                  | 0                | 0                     | 56                  | 2              | 49              | 83                     | 0               | 0                | 0                    | 15                 | 0               | 0              | 0                | 0               | 0                | 26            | 0                | 45                      | 32                             | 0                             | 31                            | 44              | 0               | 0                  | 0                  | 0                  | 0               | 0                 | 44                            | 0               |    |
| MAPKAP-K3(h)    | 0            | 0            | 0                  | 0                  | 0                | 2                     | 0                   | 14             | 21              | 0                      | 0               | 0                | 0                    | 0                  | 0               | 0              | 0                | 0               | 0                | 0             | 8                | 0                       | 0                              | 0                             | 2                             | 0               | 0               | 0                  | 0                  | 0                  | 0               | 6                 | 0                             | 0               |    |
| Mei(h)          | 0            | 0            | 52                 | 42                 | 0                | 26                    | 0                   | 0              | 83              | 92                     | 81              | 2                | 0                    | 41                 | 26              | 0              | 28               | 19              | 42               | 0             | 0                | 11                      | 86                             | 0                             | 68                            | 38              | 0               | 0                  | 0                  | 53                 | 0               | 14                | 0                             | 0               |    |
| MSK1(h)         | 7            | 0            | 5                  | 15                 | 70               | 0                     | 23                  | 14             | 46              | 91                     | 0               | 0                | 19                   | 0                  | 91              | 0              | 0                | 9               | 8                | 0             | 0                | 89                      | 61                             | 3                             | 74                            | 0               | 45              | 40                 | 0                  | 0                  | 16              | 0                 | 0                             | 38              | 0  |
| MSK2(h)         | 26           | 0            | 17                 | 17                 | 81               | 0                     | 48                  | 25             | 84              | 88                     | 6               | 3                | 0                    | 90                 | 1               | 6              | 39               | 22              | 0                | 68            | 59               | 3                       | 93                             | 0                             | 68                            | 79              | 3               | 0                  | 0                  | 69                 | 0               | 0                 | 79                            | 0               |    |
| NEK2(h)         | 0            | 5            | 13                 | 0                  | 0                | 8                     | 2                   | 0              | 41              | 68                     | 54              | 10               | 1                    | 6                  | 21              | 0              | 4                | 51              | 18               | 0             | 0                | 0                       | 31                             | 3                             | 42                            | 0               | 34              | 0                  | 0                  | 74                 | 0               | 4                 | 4                             | 0               |    |
| p70S6K(h)       | 14           | 3            | 10                 | 24                 | 88               | 0                     | 14                  | 17             | 93              | 101                    | 10              | 0                | 0                    | 99                 | 4               | 31             | 75               | 23              | 14               | 100           | 37               | 0                       | 94                             | 0                             | 81                            | 40              | 1               | 0                  | 0                  | 78                 | 0               | 0                 | 45                            | 2               |    |
| PAK2(h)         | 8            | 6            | 14                 | 11                 | 3                | 5                     | 1                   | 15             | 37              | 98                     | 4               | 0                | 14                   | 21                 | 6               | 0              | 2                | 17              | 2                | 13            | 0                | 79                      | 67                             | 0                             | 32                            | 29              | 2               | 0                  | 2                  | 9                  | 0               | 0                 | 19                            | 0               |    |
| PDGFR_β(h)      | 0            | 0            | 2                  | 0                  | 9                | 55                    | 0                   | 0              | 80              | 98                     | 92              | 0                | 48                   | 31                 | 0               | 0              | 0                | 0               | 18               | 87            | 0                | 33                      | 97                             | 0                             | 91                            | 5               | 43              | 0                  | 0                  | 72                 | 0               | 3                 | 0                             | 0               |    |
| PDGFRβ(h)       | 75           | 13           | 6                  | 16                 | 5                | 5                     | 0                   | 103            | 34              | 105                    | 30              | 40               | 0                    | 26                 | 77              | 0              | 0                | 23              | 3                | 13            | 0                | 2                       | 80                             | 0                             | 17                            | 6               | 0               | 84                 | 8                  | 55                 | 0               | 0                 | 0                             | 0               |    |
| PDK1(h)         | 2            | 4            | 1                  | 2                  | 0                | 0                     | 14                  | 0              | 97              | 99                     | 0               | 0                | 0                    | 55                 | 0               | 0              | 0                | 0               | 0                | 4             | 0                | 16                      | 83                             | 0                             | 62                            | 94              | 0               | 0                  | 0                  | 7                  | 0               | 0                 | 91                            | 0               |    |
| Pim-1(h)        | 6            | 8            | 102                | 8                  | 13               | 3                     | 50                  | 84             | 98              | 99                     | 2               | 6                | 13                   | 101                | 79              | 60             | 84               | 89              | 82               | 36            | 7                | 22                      | 94                             | 0                             | 93                            | 62              | 38              | 9                  | 10                 | 84                 | 3               | 65                | 81                            | 10              |    |
| PKA(h)          | 27           | 11           | 0                  | 13                 | 75               | 6                     | 0                   | 6              | 27              | 99                     | 0               | 0                | 5                    | 72                 | 0               | 13             | 0                | 0               | 5                | 96            | 11               | 10                      | 23                             | 4                             | 9                             | 18              | 0               | 0                  | 0                  | 0                  | 0               | 0                 | 0                             | 0               |    |
| PKB_β(h)        | 0            | 0            | 37                 | 14                 | 38               | 0                     | 5                   | 2              | 21              | 95                     | 0               | 0                | 0                    | 81                 | 12              | 0              | 0                | 14              | 1                | 95            | 24               | 0                       | 42                             | 0                             | 21                            | 12              | 0               | 0                  | 0                  | 3                  | 0               | 0                 | 8                             | 0               |    |
| PKBβ(h)         | 0            | 0            | 19                 | 24                 | 0                | 0                     | 53                  | 22             | 25              | 96                     | 0               | 0                | 0                    | 59                 | 2               | 0              | 0                | 9               | 0                | 72            | 0                | 0                       | 15                             | 0                             | 0                             | 4               | 0               | 0                  | 0                  | 10                 | 0               | 0                 | 16                            | 0               |    |
| PKB_α(h)        | 19           | 1            | 25                 | 22                 | 65               | 0                     | 20                  | 4              | 38              | 98                     | 0               | 0                | 0                    | 100                | 0               | 0              | 6                | 0               | 0                | 97            | 9                | 0                       | 55                             | 0                             | 29                            | 13              | 0               | 0                  | 0                  | 4                  | 0               | 0                 | 53                            | 0               |    |
| PKCβ1(h)        | 3            | 0            | 13                 | 6                  | 27               | 0                     | 0                   | 11             | 88              | 96                     | 0               | 0                | 1                    | 100                | 4               | 1              | 0                | 0               | 0                | 64            | 14               | 0                       | 11                             | 0                             | 12                            | 0               | 0               | 0                  | 0                  | 0                  | 0               | 0                 | 33                            | 0               |    |
| PRAK(h)         | 0            | 3            | 36                 | 0                  | 0                | 0                     | 61                  | 9              | 71              | 79                     | 3               | 0                | 0                    | 0                  | 0               | 0              | 32               | 2               | 0                | 18            | 0                | 11                      |                                |                               |                               |                 |                 |                    |                    |                    |                 |                   |                               |                 |    |



the Serologicals® Group

|              |               |                  |
|--------------|---------------|------------------|
| <b>USA</b>   | <b>UK</b>     | <b>Australia</b> |
| 800 437 7500 | 0800 0190 333 | 03 9839 2000     |

|                       |                     |                     |
|-----------------------|---------------------|---------------------|
| <b>Europe</b>         | Austria 01 585 4116 | France 0805 300 576 |
| Germany 0800 759 8965 | Sweden 040 22 42 00 |                     |

To find your local distributor, please visit our website  
**[www.upstate.com](http://www.upstate.com)**