For life science research only. Not for use in diagnostic procedures.



Expand High Fidelity
 PCR System

Version: 10

Content Version: October 2020

Thermostable DNA polymerase blend.

Deoxynucleoside-triphosphate: DNA deoxynucleotidyltransferase

Cat. No. 03 300 226 001 500 U

2 x 250 U, 5 U/µl

200 reactions in a final volume of 50 µl

Cat. No. 03 300 234 001 2,500 U

10 x 250 U, 5 U/µl

1,000 reactions in a final volume of 50 µl

Store the kit at -15 to -25°C.

1.	General Information	3
1.1.	Contents	3
1.2.	Storage and Stability	
	Storage Conditions (Product)	
1.3.	Additional Equipment and Reagent required	4
1.4.	Application	4
2.	How to Use this Product	5
2.1.	Before you Begin	5
	Sample Materials	5
	Primers	
	Mg2+ Concentration	
	General ConsiderationsPrevention of Carryover Contamination	
	Safety Information	
	For customers in the European Economic Area	
2.2.	Protocols	
2.2.	Preparation of PCR master mix	
	PCR protocol	
2.3.	Parameters	9
	EC-Number	9
	Error Rate	
	Incorporation of Modified Nucleotides	
	Maximum Fragment Size	
	PCR Cloning Proofreading Activity	
	Temperature Optimum	
	Volume Activity	
	Working Concentration	9
3.	Troubleshooting	10
4.	Additional Information on this Product	11
4.1.	Quality Control	11
5.	Supplementary Information	12
5.1.	Conventions	12
5.2.	Changes to previous version	12
5.3.	Ordering Information	
5.4.	Trademarks	
5.5.	License Disclaimer	
5.6.	Regulatory Disclaimer	
5.7.	Safety Data Sheet	
5.7.	Contact and Support	13
.10	CONTROL AND ANDROLL	1 3

1. General Information

1.1. Contents

Vial / Bottle	Label	Function / Description	Catalog Number	Content
1	PCR System, 20 mM Tris-HCl, pH 8 (+4°C),	03 300 226 001	2 vials, 50 µl each	
		03 300 234 001	10 vials, 50 μl each	
2	Expand High Fidelity ^{PLUS} PCR System, Reaction Buffer, 5x conc. with MgCl ₂	Reaction buffer with 7.5 mM MgCl ₂ .	03 300 226 001	2 vials, 1.25 ml each
			03 300 234 001	10 vials, 1.25 ml each
3	Expand High Fidelity ^{PLUS} Reaction buffer without MgCl ₂ . PCR System,		03 300 226 001	2 vials, 1.25 ml each
	Reaction Buffer, 5x conc. without MgCl ₂		03 300 234 001	10 vials, 1.25 ml each
4	Expand High Fidelity ^{PLUS} PCR System, MgCl ₂ 25 mM Stock Solution	To adjust final Mg ²⁺ concentration.	03 300 226 001	2 vials, 1 ml each
			03 300 234 001	8 vials, 1 ml each

1.2. Storage and Stability

Storage Conditions (Product)

When stored at -15 to -25°C, the kit is stable through the expiry date printed on the label.

Vial / Bottle	Label	Storage
1	Enzyme Mix	Store at -15 to -25 °C.
2	Reaction Buffer, 5x conc. with MgCl ₂	Store at −15 to −25°C.
3	Reaction Buffer, 5x conc. without MgCl ₂	⚠ Thaw and equilibrate at +37 to +56°C before use; vortex thoroughly. If crystals have formed, incubate at +37 to +56°C until dissolved.
4	MgCl ₂ 25 mM Stock Solution	Store at −15 to −25°C.

1.3. Additional Equipment and Reagent required

Standard laboratory equipment

- Nuclease-free, aerosol-resistant pipette tips
- Pipettes with disposable, positive-displacement tips
- Sterile reaction tubes for preparing PCR mixes and dilutions
- PCR reaction vessels, such as thin-walled PCR tubes or plates
- Standard benchtop microcentrifuge
- Thermal block cycler

For PCR

- PCR primers
- Template DNA
- PCR Nucleotide Mix, PCR Grade*
- Water, PCR Grade*
- Mineral oil (optional)

For incorporation of nonradioactive modified nucleotides

- Biotin-16-dUTP*
- Fluorescein-12-dUTP*
- Digoxigenin-11-dUTP, alkali-stable*
- Digoxigenin-11-dUTP, alkali-labile*
- Deoxynucleoside Triphosphate Set, PCR Grade* (optional)

For carryover prevention (optional)

- Uracil-DNA Glycosylase, heat-labile*
- PCR Nucleotide Mix^{PLUS}, PCR Grade*

1.4. Application

The Expand High Fidelity PLUS PCR System is designed to:

- Amplify fragments up to 5 kb from all types of DNA with outstanding yield and fidelity.
 - The system is a blend of Taq DNA Polymerase and a thermostable proofreading protein that lacks polymerase activity.
- Amplify DNA with twofold greater replicational accuracy (fidelity) than the Expand High Fidelity PCR System, and sixfold greater fidelity than Taq DNA Polymerase alone.
 - *The synergy between this proofreading protein and the highly processive Taq DNA Polymerase is the key to the outstanding yield, specificity, sensitivity, and accuracy of the Expand High Fidelity*^{PLUS} System.
- Incorporate dUTP and, in combination with Uracil-DNA Glycosylase, can be used to safeguard PCR reactions from cross-contamination. Therefore, it is suited for simultaneous amplification of a large number of different targets in the same run.

These characteristics make the Expand High Fidelity^{PLUS} PCR System blend the product of choice if a larger number of targets need to be amplified at the same time with high yield in combination with improved accuracy.

• Efficiently label DNA fragments with radioactive or nonradioactive modified nucleotides.

2. How to Use this Product

2.1. Before you Begin

Sample Materials

Use any template DNA such as genomic or plasmid DNA, cDNA suitable for PCR in terms of purity, concentration, and absence of inhibitors. For reproducible isolation of nucleic acids, use:

- Either a MagNA Pure System together with a dedicated reagent kit (for automated isolation),
- or a High Pure Nucleic Acid Isolation Kit (for manual isolation).
- Use 5 to 500 ng human genomic DNA or 0.1 to 10 ng plasmid DNA/cDNA. Use a starting concentration of 250 ng genomic DNA or 1 ng plasmid DNA.

⚠ Store the template DNA in either Water, PCR Grade* or 5 to 10 mM Tris-HCl, pH 7 to 8. Avoid dissolving the template in TE buffer since EDTA chelates Mg²⁺.

Primers

Forward and reverse primer

- Start with 0.4 µM (final concentration) of each primer.
- For optimization, the concentration can vary between 0.2 and 0.6 μM.

Mg2+ Concentration

1.5 to 4 mM (as MgCl₂) (optimal)
1.5 mM (as MgCl₂) when used with 200 μM of each dNTP (standard)
i) The Mg²⁺ concentration must be optimized when dUTP is used.

General Considerations

The optimal conditions, including incubation times and temperatures, concentration of enzyme, template DNA, Mg²⁺ vary from system to system and must be determined for each individual experimental system. At the very least, titrate the Mg²⁺ concentration and the amount of enzyme mix used per assay to ensure optimal efficiency of DNA synthesis. As a starting point, use the following guidelines:

- Optimal enzyme concentration: 2.5 U/50 μl.
- Optimal Mg²⁺ concentration can vary between 1.5 to 4 mM. In most cases, a Mg²⁺ concentration of 1.5 mM will produce satisfactory results if you use 200 µM of each dNTP.
- dNTP concentration: Always use equal concentrations of all four dNTPs. The final concentration of each dNTP should be between 50 and 500 μM; the most commonly used concentration is 200 μM. If you increase the dNTP concentration, you must also increase the Mg²+ concentration.

Prevention of Carryover Contamination

Yes

Uracil-DNA Glycosylase (UNG) is suitable for preventing carryover contamination in PCR. This carryover prevention technique involves incorporating deoxyuridine triphosphate (dUTP instead of dTTP) into amplification products, then pretreating later PCR mixtures with UNG. If a dUTP-containing contaminant is present in the later PCRs, it will be cleaved by a combination of the UNG and the high temperatures of the initial denaturation step; it will not serve as a PCR template. Since your target DNA template contains thymidine rather than uridine, it is not affected by this procedure.

Adapt the protocol in section, Protocols, Preparation of PCR master mix as follows:

- Use the Expand High Fidelity^{PLUS} Reaction Buffer, 5x conc. without MgCl₂ (Vial 3) and MgCl₂, 25 mM stock solution (Vial 4) in concentrations ranging from 1.5 mM up to 4 mM final.
- Use the PCR Nucleotide Mix^{PLUS} which contains 200 mM dATP, dCTP, dGTP and 600 mM dUTP (final conc.) To prevent carryover contamination using heat-labile UNG:
- 1 Per 50 μl final reaction volume, add 2 μl heat-labile UNG to the RT-PCR mix.
- 2 Add DNA template and incubate the reaction mixture for 10 minutes at +15 to +25°C to destroy any contaminating template.
- 3 Inactivate heat-labile UNG by performing the Pre-Incubation step at +94°C for 2 minutes.

Safety Information

For customers in the European Economic Area

Contains SVHC: octyl/nonylphenol ethoxylates. For use in research and under controlled conditions only – acc. to Art. 56.3 and 3.23 REACH Regulation.

2.2. Protocols

Preparation of PCR master mix

- *See sections,* **Prevention of Carryover Contamination** *and* **Incorporation of Modified Nucleotides** *for additional information.*
- Thaw the reagents and store on ice.
 - Briefly vortex and centrifuge all reagents before setting up the reactions.
- 2 Prepare a 10x-concentrated solution of each respective primer.
 - i If you are using, for example, the final concentration of 0.4 μM for each primer, the 10x-concentrated solution would contain a 4 μM concentration of the respective primer.
- 3 To a sterile reaction tube on ice, add the components in the order listed for each 50 µl reaction:
 - *(i)* To prepare the PCR mix for more than one reaction, multiply the amounts in the Volume column by z, where z equals the number of reactions to be run plus one additional reaction.

Reagent	Volume [µl]	Final conc.
Water, PCR Grade*	add up to a final volume of 50	-
Reaction Buffer, 5x conc. with MgCl ₂	10	1x
PCR Grade Nucleotide Mix* (10 mM of each dNTP)	1	200 μM of each dNTP
Forward primer 1, 10x conc.	5	0.4 μM
Reverse primer 2, 10x conc.	5	0.4 μM
Template DNA	variable	5 – 500 ng genomic DNA, 100 pg – 10 ng plasmid DNA
Enzyme Mix (Vial 1)	0.5	2.5 U/reaction
Final Volume	50	

4 Gently vortex the mixture to produce a homogeneous reaction, then centrifuge briefly to collect the solution at the bottom of the tube.

PCR protocol

- 1 The following thermal profiles are an example. Different thermal cyclers may require different profiles.
- 1 Place your samples in a thermal block cycler and use either of the thermal profiles below to perform PCR.

 1 Thermal Profile A has a fixed extension time.

Step	Temperature [°C]	Time	Number of Cycles
Pre-Incubation	94 ⁽¹⁾	2 min	1
Denaturation	94 ⁽¹⁾	10 - 30 sec	25 – 35
Annealing	55 – 68 ⁽²⁾	30 sec	
Elongation	68 - 72 ⁽³⁾	30 sec – 4 min ⁽⁴⁾	
Final Elongation	68 – 72	7 min	1
Cooling	4	indefinitely	

Thermal Profile B has a gradually increasing extension time, ensuring a higher yield of amplification products.

Step	Temperature [°C]	Time	Number of Cycles
Pre-Incubation	94 ⁽¹⁾	2 min	1
Denaturation Annealing Elongation	94 ⁽¹⁾ 55 - 68 ⁽²⁾ 72	10 – 30 sec 30 sec 30 sec – 4 min ⁽⁴⁾	10
Denaturation Annealing Elongation	92 - 94 ⁽¹⁾ 55 - 68 ⁽²⁾ 68 - 72 ⁽³⁾	15 – 30 sec 30 – 60 sec 30 sec – 4 min ⁽⁴⁾ + 10 sec cycle elongation for each successive cycle ⁽⁵⁾	15 – 20
Final Elongation	68 - 72 ⁽³⁾	7 min	1
Cooling	4	indefinitely	

- 2 After cycling, use samples immediately or store them frozen for later use.
 - **Proposition** The PCR product on an agarose gel for size and specificity. Use an appropriate size marker. In addition, purify the PCR product with the High Pure PCR Product Purification Kit, for example, before performing nested PCR.
- (1) Optimal denaturation temperature and time depends upon the GC content of the template.
- ⁽²⁾ Optimal annealing temperature depends on the melting temperature of the primers and on the experimental system.
- (3) For PCR products up to 3 kb, elongation temperature should be approximately +72°C; for PCR products >3 kb, elongation temperature should be approximately +68°C.
- (4) Elongation time depends upon the length of the product to be amplified, approximately 1 minute per kb.
- (5) For example, cycle number 11 is 10 seconds longer than cycle 10. Cycle number 12 is 20 seconds longer than cycle 10. Cycle number 13 is 30 seconds longer than cycle 10, etc.

2.3. Parameters

EC-Number

EC 2.7.7.7

Error Rate

- Twofold more accurate compared to the Expand High Fidelity PCR System.
- Sixfold more accurate compared to Taq DNA polymerase.
- ? Relative fidelity determined by the lacl assay.

Incorporation of Modified Nucleotides

Accepts modified nucleotides, such as Digoxigenin-dUTP*, Biotin-dUTP*, or Fluorescein-dUTP*. To use the Expand High Fidelity^{PLUS} PCR system to incorporate nonradioactive modified nucleotides, exchange the PCR Nucleotide Mix with one of the following modified nucleotide mixes:

Nucleotide Mix, 10x conc.	Description/Preparation
Digoxigenin-11-dUTP	 0.7 mM Digoxigenin-11-dUTP*, 1.3 mM dTTP, 2 mM dCTP, 2 mM dATP, and 2 mM dGTP. Use a final MgCl₂ concentration of 3 mM, however it is best to optimize the concentration for each new reaction.
Biotin-16-dUTP	 0.7 mM Biotin-16-dUTP*, 1.3 mM dTTP, 2 mM dCTP, 2 mM dATP, and 2 mM dGTP. Use a final MgCl₂ concentration of 1.5 mM.
Fluorescein-12-dUTP	 0.5 mM Fluorescein-12-dUTP*, 1.5 mM dTTP, 2 mM dCTP, 2 mM dATP, and 2 mM dGTP. Use a final MgCl₂ concentration of 1.5 mM.

For convenience, use the Deoxynucleoside Triphosphate Set, PCR Grade which contains individual vials of dATP, dCTP, dGTP, dTTP at a concentration of 100 mM.*

Maximum Fragment Size

Up to 5 kb.

PCR Cloning

TA cloning

Enzyme adds a single, overhanging adenine (A).

Proofreading Activity

Yes

Temperature Optimum

+72°C (elongation) for amplicons up to 3 kb. For PCR products >3 kb, the optimal elongation temperature is +68°C.

Volume Activity

5 U/µl

Working Concentration

2.5 U per 50 µl reaction (standard).

3. Troubleshooting

Observation	Possible cause	Recommendation
Little or	Pipetting errors	Check all concentrations and storage conditions of reagents.
no PCR product.	Difficult templates, such as GC-rich templates.	Add DMSO (final concentration 8%), and reduce enzyme concentration down to 0.5 U per reaction.
	•	Perform PCR with the GC-RICH PCR System*.
	Poor DNA template quality	 Check quality and concentration of template: Analyze an aliquot on an agarose gel to check for possible degradation. Perform a control reaction on template using an established primer pair. Check or repeat template purification.
	Enzyme concentration too low.	Increase the amount of enzyme mix in 0.5 U steps to 2 U per 50 μ l reaction.
	MgCl ₂ concentration too low	Increase the MgCl ₂ concentration in 0.25 mM steps.
	Cycle conditions not	Decrease annealing temperature.
	optimal.	Increase cycle number.
		Make sure that the final elongation step is included in the program.
	Primer design not optimal.	Design alternative primers.
	Primer concentration not optimal.	Both primers must have the same concentration.
		Titrate primer concentration.
	Annealing temperature	Reduce annealing temperature.
	too high.	Determine the optimal annealing temperature by touchdown PCR.
	Primer specificity not optimal.	Perform nested PCR.
	Primer quality or storage problems.	If you use an established primer pair, check performance with a control template.
		Make sure that the primers are not degraded.
		Always store primers at −15 to −25°C.
Multiple bands or background	Annealing temperature too low.	Increase annealing temperature according to primer length.
smear.	Primer design or concentration not optimal.	Review primer design.
		Titrate primer concentration (0.1 to 0.6 μM).
		Both primers must have the same concentration.
		Perform nested PCR.
	Difficult templates, such	Perform PCR with the GC-RICH PCR System*.
	as GC-rich templates.	Add DMSO (final concentration 8%), and reduce enzyme concentration down to 0.5 U per reaction.
	DNA template problems.	Use serial dilution of template.
PCR products in	Carryover contamination	Replace all reagents, especially water.
negative control	present.	Use aerosol-resistant pipette tips.
experiments.		Set up PCR reactions in an area separate from that used for PCR product analysis.
		To eliminate carryover contaminants: Use dUTP (600 µM) instead of dTTP, and Uracil-DNA Glycosylase* (1 U/50 µl reaction). Increase Mg²+ concentration to a maximum of 4 mM.

Problems specific to RT-PCR.	No product, additional bands, background smear	The volume of cDNA template (RT reaction) should not exceed 10% of the final concentration of the PCR reaction.
	observed.	Follow all troubleshooting tips.
		Increase MgCl ₂ in 0.25 mM steps.

4. Additional Information on this Product

4.1. Quality Control

For lot-specific certificates of analysis, see section, **Contact and Support**.

5. Supplementary Information

5.1. Conventions

To make information consistent and easier to read, the following text conventions and symbols are used in this document to highlight important information:

Text convention and symbols				
information Note: Add	information Note: Additional information about the current topic or procedure.			
⚠ Important Note: Information critical to the success of the current procedure or use of the product.				
1 2 3 etc.	Stages in a process that usually occur in the order listed.			
1 2 3 etc. Steps in a procedure that must be performed in the order listed.				
* (Asterisk)	The Asterisk denotes a product available from Roche Diagnostics.			

5.2. Changes to previous version

Layout changes.

Editorial changes.

New information added related to the REACH Annex XIV.

5.3. Ordering Information

Product	Pack Size	Cat. No.
Reagents, kits		
PCR Nucleotide Mix	200 μ l, 500 reactions of 20 μ l final reaction volume	11 581 295 001
	$5 \ x \ 200 \ \mu l, \ 2,500$ reactions of 20 μl final reaction volume.	04 638 956 001
	10 x 200 μ l, 5,000 reactions of 20 μ l final reaction volume.	11 814 362 001
PCR Nucleotide MixPLUS	2 x 100 μl, 200 PCR reactions in 50 μl	11 888 412 001
Water, PCR Grade	25 ml, 25 x 1 ml	03 315 932 001
	25 ml, 1 x 25 ml	03 315 959 001
	100 ml, 4 x 25 ml	03 315 843 001
Uracil-DNA Glycosylase,	100 U, 1 U/μl	11 775 367 001
heat-labile	500 U, 1 U/μl	11 775 375 001
GC-RICH PCR System	100 U, 50 reactions in a final volume of 50 μl	12 140 306 001
Biotin-16-dUTP	custom fill	11 093 711 103
Fluorescein-12-dUTP	custom fill	11 375 601 103
Digoxigenin-11-dUTP,	25 nmol, 25 μl, 1 mM	11 093 088 910
alkali-stable	125 nmol, 125 μl, 1 mM	11 558 706 910
	5 x 125 nmol, 5x 125 μl, 1 mM	11 570 013 910
Digoxigenin-11-dUTP,	25 nmol, 25 μl, 1 mM	11 573 152 910
alkali-labile	125 nmol, 125 μl, 1 mM	11 573 179 910
Deoxynucleoside	4 x 250 μl, 4 x 25 μmol, 100 mM	11 969 064 001
Triphosphate Set	4 x 1,250 μl, 4 x 125 μmol, 100 mM	03 622 614 001

5.4. Trademarks

MAGNA PURE and EXPAND are trademarks of Roche.
All other product names and trademarks are the property of their respective owners.

5.5. License Disclaimer

For patent license limitations for individual products please refer to: **List of biochemical reagent products**.

5.6. Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

5.7. Safety Data Sheet

Please follow the instructions in the Safety Data Sheet (SDS).

5.8. Contact and Support

To ask questions, solve problems, suggest enhancements or report new applications, please visit our **Online Technical Support Site**.

To call, write, fax, or email us, visit **sigma-aldrich.com**, and select your home country. Country-specific contact information will be displayed.