

Data Sheet

BioTracker™ Dual-Color Mitochondria-Lysosome Live Cell Dye

Live Cell Dye

SCT244**Pack Size: 1 mg****Store at -20 °C**FOR RESEARCH USE ONLY**Not for use in diagnostic procedures. Not for human or animal consumption.**

Background

The BioTracker™ Dual-Color Mitochondria-Lysosome Live Cell Dye is a fluorescent probe used to detect mitochondria-lysosome interactions.

Mitochondrial-lysosomal synergy is critically important in the processes of cellular homeostasis. Mitochondria and lysosomes have unique interactions which involve the transfer of metabolites. Certain disease pathogenesis can be linked to this important cellular interaction. The BioTracker™ Dual-Color Mitochondria-Lysosome Live Cell Dye is a unique cell-permeable and viscosity-responsive molecular probe that is capable of monitoring mitochondrial and lysosomal intercommunication in living cells. Through the use of fluorescence conversion, this live cell probe can simultaneously tag mitochondria with blue fluorescence and lysosomes with red fluorescence.

Source

SCT244 does not contain genetically modified organisms.

Spectral Properties

Fluorescence images obtained by $\lambda_{ex} = 405$ nm and emission at 450 nm (Mitochondria) and $\lambda_{ex} = 561$ nm and emission at 590 nm (Lysosome).

Quality Control Testing

Purity: $\geq 98\%$ confirmed by HPLC, HNMR, LC-MS and elemental analysis. Molar Mass: 586.46 g/mol.

Storage and Handling

Store BioTracker™ Dual-Color Mitochondria-Lysosome Live Cell Dye at -20 °C, desiccated and protected from light.

Note: Centrifuge vial briefly to collect contents at bottom of vial before opening.

Presentation

Lyophilized. Brown solid.

Representative Data




Figure 1: Confocal microscopy images of Biotracker™ Dual-Color Mitochondria-Lysosome Live Cell Dye staining. Raw264.7 mouse macrophage cells were incubated with 10 μ M SCT244 solution (Red and blue). **A.** Mitochondrial staining is shown in the blue channel. **B.** Lysosomal staining is shown in red. **C.** A merged image is shown to display the interactions between mitochondrial and lysosomal channels.




Figure 2: Probe excitation and emission data. 7 μL of probe at stock concentration (10 mM) was diluted in 1 mL of DMSO before undergoing excitation and emission scans. Spectral scans were conducted using a Fluorescence Spectrophotometer (PerkinElmer FL8500). **A.** Mitochondrial probe spectral data. **B.** Lysosomal probe spectral characteristics.

References

1. Chen, Qixin, et al. "A dual-labeling probe to track functional mitochondria-lysosome interactions in live cells." *Nature communications* 11.1 (2020): 1-10.

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

Merck, Biotracker and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2008-2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

Document Template 20306518 Ver 6.0

00153745 Ver 1.0, Rev 17NOV2023, RC, AB

