



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

Chymopapain from papaya latex

Product Number **C 8526**
Storage Temperature -0 °C

Product Description

Molecular weight: 27 kDa¹
CAS Number: 9001-09-6
Enzyme Commission (EC) Number: 3.4.22.6
 λ_{max} = 280 nm
Extinction coefficient: $E^{1\%} = 18.7^1$

This enzyme cleaves peptide (amide) linkages. The carboxyl side of arginine residues are most susceptible,⁷ but it will also cleave on the carboxyl side of leucine, glutamic acid, and alanine residues, as well as arginine esters.¹

This enzyme has been used to detach cultured CD34⁺ cells from beads (130 units/ml for 2 minutes, repeated three times).² Other examples of the use of chymopapain for cell dissociation have been published. For digestion of human intestinal tissue, chymopapain and pronase were found to be most effective.³ For rabbit lung cells, chymopapain has been used at a concentration of 0.05% in calcium-magnesium-free Krebs serum substitute, along with pronase, collagenase, elastase, DNase, and catalase.⁴ Methods for the release of neuroblastoma cells from marrow bound to antibody-coated microspheres have been published.⁵

The method of enzyme assay used by Sigma for this product has been published.⁶ The pH dependence of the activity of this enzyme using several substrates has also been published.⁷ For example, the pH optimum for casein digestion is broad (pH 7-9); hemoglobin digestion is optimal at pH 4.0, but urea-denatured hemoglobin digestion is optimal at pH 7.0.⁷

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This product is soluble in water (5 mg/ml), yielding a clear, colorless solution.

Storage/Stability

The enzyme can be diluted in water if a concentrated solution is made (at least several mg/ml). However, to insure stability for a more dilute solution (1 mg/ml), dissolve the protein in 0.05 M sodium acetate buffer, pH 6.2, containing 1 – 5 mM cysteine and 0.1 – 1 mM EDTA.

References

1. Ebata, M., and Yasunobu, K. T., Isolation, crystallization, and preliminary characterization. *J. Biol. Chem.*, **237**(4), 1086-1094 (1962).
2. Issaad, C., and Vainchenker, W., Growth of erythroid colonies in chronic myelogenous leukemia is independent of erythropoietin only in the presence of steel factor. *Blood*, **84**(10), 3447-3456 (1994).
3. Fox, C. C., et al., Isolation and characterization of human intestinal mucosal mast cells. *J. Immunol.*, **135**(1), 483-491 (1985).
4. Gould, K. G. Jr., et al., Dispersal of rabbit lung into individual viable cells: a new model for the study of lung metabolism. *Science*, **178**(66), 1209-1210 (1972).
5. Gazitt, Y., et al., Expression of N-myc, c-myc, and MDR-1 proteins in newly established neuroblastoma cell lines: a study by immunofluorescence staining and flow cytometry. *Cancer Res.*, **52**(10), 2957-2965 (1992).

6. Arnon, R., Papain. Methods Enzymol., **19**, 226-244 (1970).

7. Kunimitsu, D.K. and Yasunobu, K.T., Chymopapain B. Methods Enzymol., **19**, 244-252 (1970).

IRB/RXR 10/03

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.