

BioTracker™ 560 Orange Lysosome Dye

Live Cell Dye

Cat. # SCT019

FOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSTIC PROCEDURES.
NOT FOR HUMAN OR ANIMAL CONSUMPTION.

pack size: 10x10µg

Store at -20°C



Data Sheet

page 1 of 2

Background

Lysosomes are membrane-enclosed organelles that contain an array of enzymes capable of breaking down all types of biological material including proteins, nucleic acids, carbohydrates, and lipids. Lysosomes function as the digestive system of the cell, serving both to degrade material taken up from outside the cell and to digest obsolete components of the cell itself.

The BioTracker™ 560 Orange Lysosome Dye is a live cell fluorescence imaging probe which enhances fluorescence dramatically in acidic environments. This probe can stain acidic organelles such as lysosome, late endosome and granule selectively. Namely, in the condition of pH 5.0, identical to the environment of acidic organelles, fluorescence intensity is 50-fold or more compared to physical pH 7.4. The BioTracker™ 560 Orange Lysosome Dye shows great photostability against irradiating excitation light. As it emits orange color fluorescent by excitation at 532nm or 514 nm, multicolor imaging is possible by combining with blue fluorescence (CFP, Hoechst, etc.), green fluorescence (GFP, fluorescein, etc.) and near-infrared fluorescence.

Storage

Store BioTracker 560 Orange Lysosome Dye at -20°C, desiccate and protect from light

Note: Centrifuge vial briefly to collect contents at bottom of vial before opening.

Spectral Properties

Absorbance: 535nm

Emission: 560nm




Figure 1. Fluorescence spectra of BioTracker 560 Orange Lysosome Dye was measured in phosphate buffer pH 5.0 or pH 7.4, respectively. Fluorescent intensity of BioTracker 560 Orange Lysosome Dye at pH 5.0 increased 50 fold to that in the pH 7.4 buffer. (λ_{exc} 532 nm / λ_{em} 568 nm)




Figure 2. HeLa cells expressing Mitochondria-GFP were stained with BioTracker 560 Orange Lysosome Dye and Hoechst33342.

Please visit www.millipore.com for additional product information and references.

Submit your published journal article, and earn credit toward future purchases. Visit www.millipore.com/publicationrewards to learn more!




Figure 3. Photostability of BioTracker 560 Orange Lysosome Dye. HeLa cells were continuously illuminated for 180 seconds and viewed by confocal microscopy. Although LysoTracker® Green DND-26 and LysoTracker® Red DND-99 discolored dramatically, fluorescence of BioTracker 560 Orange Lysosome Dye remained strong. LysoSensor® Green DND-189 was not suitable for continuous observation because fluorescence leaked to cytoplasm by excitation light irradiation.




Figure 4. Photostability of BioTracker 560 Orange Lysosome Dye vs LysoTracker Green and LysoTracker Red over 180 seconds of continuous light exposure using HeLa cells.

BioTracker™ is a trademark of Merck KGaA

■ antibodies ■ Multiplex products ■ biotools ■ cell culture ■ enzymes ■ kits ■ proteins/peptides ■ siRNA/cDNA products

Please visit www.millipore.com for additional product information, test data and references

EMD Millipore Corporation, 28820 Single Oak Drive, Temecula, CA 92590, USA 1-800-437-7500

Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures. Not for human or animal consumption. Purchase of this Product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

EMD Millipore®, the M mark, Upstate®, Chemicon®, Linco® and all other registered trademarks, unless specifically identified above in the text as belonging to a third party, are owned by Merck KGaA, Darmstadt, Germany. Copyright ©2008-2018 Merck KGaA, Darmstadt, Germany. All rights reserved.



We Buy 100% Certified Renewable Energy