

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

## **ProductInformation**

# Xanthine Oxidase from bovine milk

Product Number **X 4376** Storage Temperature 2-8 °C

### **Product Description**

Enzyme Commission (EC) Number: 1.1.3.22

CAS Number: 9002-17-9
Molecular Weight: 283 kDa<sup>1</sup>
Extinction Coefficient: E<sup>1%</sup> = 11.7

Synonyms: Xanthine:oxygen oxidoreductase, XOD

Xanthine oxidase from buttermilk is a homodimer consisting of two equal subunits of 140 kDa. Each subunit contains one mole of FAD, one atom of Mo, and four iron atoms. The enzyme catalyzes the following reaction:

Xanthine +  $H_2O + O_2 \rightarrow Uric Acid + H_2O_2$ 

Hypoxanthine, purine, acetaldehyde, salicylaldehyde, and benzaldehyde may also be utilized as substrates. Reported  $K_{\rm M}$  values are xanthine (1.7 mM), hypoxanthine (1.3 mM), and salicylaldehyde (1.1 mM).

Xanthine oxidase is inhibited by acacetin, allopurinol, ellagic acid, 4-hydroxybenzaldehyde, 3-methyl-quercitin, persicarin and luteolin.<sup>3</sup>

Under some reaction conditions, the product of the reaction catalyzed by xanthine oxidase is the superoxide radical  $O_2$ - rather than hydrogen peroxide. As such, xanthine oxidase can be used to generate superoxide radicals for use in the enzymatic assay of superoxide dismutase. The addition of xanthine oxidase as a superoxide donor to hepatic stellate cells (HSC) cell cultures strongly increased procollagen I synthesis.  $^5$ 

#### **Precautions and Disclaimer**

For Laboratory Use Only. Not for drug, household or other uses.

### **Preparation Instructions**

This enzyme is soluble in 50 mM potassium phosphate buffer (1 mg/ml), yielding a clear solution.

#### References

- Bray, R. C., in The Enzymes, 3rd ed., vol. XII, pt. B, Academic Press (New York, NY: 1975), pp. 303-388.
- 2. Methods of Enzymatic Analysis, 2nd ed., vol. 1, Bergmeyer, H. U., ed., Academic Press (New York, NY: 1974), pp. 521-522.
- 3. Zollner, H., Handbook of Enzyme Inhibitors, 2nd ed., pt. A, VCH (Weinheim, Federal republic of Germany: 1993), pp. 517-519.
- 4. McCord, J. M., and Fridovich, I., Superoxide dismutase. An enzymic function for erythrocuprein (hemocuprein). J. Biol. Chem., **244(22)**, 6049-6055 (1969).
- Casini, A, et al., Neutophil derived superoxide anion induces lipid peroxidation and stimulates collagen synthesis in human hepatic stellate cells: role of niotric oxide. Hepatology, 25(2), 361-367 (1997).

TMG/RXR 2/03