Instructions for Use

Iron Stain

Procedure No. HT20

Intended Use

Iron Stain is for use in the histologic staining of iron in tissue or blood and bone marrow films. Iron Stain reagents are for "In Vitro Diagnostic Use". For professional use only. The data obtained from this manual, qualitative procedure demonstrates iron deposits in blood and bone marrow films and tissue samples of human specimens. This data when reviewed in conjunction with other diagnostic tests and information may be used as an aid to diagnosis of hemochromatosis and hemosiderosis.

Iron Stain is based on the well known Prussian blue reaction¹ in which ionic iron reacts with acid ferrocyanide producing a blue color. Both the standard procedures and a procedure for rapid staining in the microwave²⁻⁴ are described.

Reagents

Potassium Ferrocyanide Solution (Cat. No. HT201-250ML)
Potassium ferrocyanide, 4% w/v.

Hydrochloric Acid Solution (Cat. No. HT202-250ML)
Hydrochloric acid, 1.2 mmol/L. Molarity can be verified by titration. Danger. Causes severe skin burns and eye damage. Wear protective gloves/protective clothing/eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician.

Pararosaniline Solution (Cat. No. HT203-30ML)
Pararosaniline hydrochloride, 1% w/v, in methanol with a stabilizer. Danger. Highly flammable liquid and vapour. Toxic if swallowed or in contact with skin. Causes serious eye irritation. Toxic if inhaled. May cause cancer. Causes damage to organs.

Special Materials Required but Not Provided

- Positive control slides, such as Iron TISSUE-TROL™ (Cat. No. TTR007) should be included in each run
- Ethanol, absolute, or Reagent Alcohol
- Acetone-free Methanol (NOTE: required for blood/bone marrow films only)

For Microwave Procedure Only

- Microwave Oven
- Staining Jars (Cat. No. S5641)

Storage and Stability

Store reagents at room temperature (18–26°C). Reagents are stable until expiration date shown on the labels.

Preparation

Working Iron Stain Solution is prepared by mixing equal volumes of Potassium Ferrocyanide Solution and Hydrochloric Acid Solution. Use once and then discard.

Working Pararosaniline Solution is prepared by adding 1 mL of Pararosaniline Solution to 50 mL water. Prepare fresh daily. Use once and discard.

Precautions

The IVDs included in this kit are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Iron TISSUE-TROL™ control slides are paraffin embedded human tissue containing iron and should be considered potentially infectious.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Any well prepared tissue section cut 5–6 microns and fixed in neutral buffered formalin is satisfactory for iron staining. Blood or bone marrow films are prepared in the usual fashion, air dried at least 30 minutes and fixed in absolute methanol for 7 minutes.⁵

Notes

Iron deposits may be removed by the use of acid containing fixatives.⁷

Procedure

Procedure for Tissue Sections

- Deparaffinize and hydrate tissue to deionized water.
- Place slides in Working Iron Stain Solution for **10 minutes**.
- Rinse in deionized water.
- Stain in Working Pararosaniline Solution for **3–5 minutes**.
- Rinse in deionized water.
- Rapidly dehydrate through alcohol and xylene and mount.

Procedure for Blood or Bone Marrow Films

- Place fixed films in Working Iron Staining Solution for **10 minutes**.
- Rinse in deionized water.
- Counterstain **5 minutes** in Working Pararosaniline Solution.
- Rinse in deionized water and air dry.

Microwave Procedure

- Deparaffinize slides and hydrate to deionized water. If staining peripheral blood or bone marrow smears fix in methanol for **7 minutes**, then hydrate to deionized water.
- Place slides in **40 mL** of Working Iron Staining Solution contained in plastic Coplin jar. Loosely cover with lid before placing in oven, or use Coplin jar lids with holes drilled into them.
- Microwave on **400 watts** for **30 seconds**. Discard solution after use.
- Rinse in deionized water.
- Place slides in **40 mL** of Working Pararosaniline Solution contained in plastic Coplin jar.
- Microwave on **800 watts** for **10 seconds**. Let incubate for **2 minutes**. Use once and discard.
- Rinse tissue sections in deionized water, dehydrate quickly through alcohols, clear in xylene and mount. Blood or bone marrow smears should be rinsed in deionized water and air dried.

Performance Characteristics

Tissue Sections

Target Structure	Staining Result
Iron Pigment	Bright Blue
Nuclei	Red
Cytoplasm	Light Pink

Note: Heavy iron deposits appear dark blue.

Blood or Bone Marrow Films

- Sideroblasts:** These are nucleated erythrocytes containing at least one small blue granule. If the blue granules surround the nucleus, the cell is a ringed sideroblast.
- Siderocytes:** These are non-nucleated erythrocytes containing at least one blue granule.
- Reticuloendothelial Iron:** Usually seen as blue particles on the marrow film or as blue particles in the cytoplasm or phagocytic cells.

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
HT201	Potassium Ferrocyanide Solution	Iron	3 of 3	3 of 3	3 of 3	3 of 3
		Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
HT202	Hydrochloric Acid Solution	Normality	3 of 3	3 of 3	3 of 3	3 of 3
		Purity	3 of 3	3 of 3	3 of 3	3 of 3
HT203	Pararosaniline Solution	Iron	3 of 3	3 of 3	3 of 3	3 of 3
		Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
		Nuclei	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

HT20:

H225: Highly flammable liquid and vapor.

H290: May be corrosive to metals.

H301 + H311 + H331: Toxic if swallowed, in contact with skin or if inhaled.

H350: May cause cancer.

H370: Causes damage to organs (Eyes, Central nervous system).

P202: Do not handle until all safety precautions have been read and understood.

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P280: Wear protective gloves/ protective clothing/ eye protection/ face protection.

P301 + P310: IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P304 + P340 + P311: IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer

References

1. Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, CV Mosby, St. Louis (MO) 1980 p 217
2. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
3. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
4. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
5. Carson: A Self Instructional Text. ASCP Press, Chicago (IL), 1990
6. Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. T Sun, C-Y Li, LT Yam, Editors, ASCP, Chicago (IL), 1985, p 205
7. Theory and Practice of Histological Techniques, 5th Edition JD Bancroft, M Gamble, Editors Churchill Livingstone, New York 2002 p 244

Contact Information

To place an order, please visit our web site at SigmaAldrich.com. For Technical Service, please visit the tech service page on our web site at SigmaAldrich.com/techservice.

Revision History

Rev. 4.0	2016
Rev. 5.0	2022
Rev. 6.0	2022

Transferred to new template with current branding. Specified for professional use in intended use and precautions. Moved aid to diagnosis statement to intended use. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Added Warnings and Hazards.

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271 Darmstadt,
Germany

The Initial M, TISSUE-TROL, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

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Referenzen

1. Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, CV Mosby, St. Louis (MO) 1980 p 217
2. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
3. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
4. Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986
5. Carson: A Self Instructional Text. ASCP Press, Chicago (IL), 1990
6. Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. T Sun, C-Y Li, LT Yam, Redakteure, ASCP, Chicago (IL), 1985, S. 205
7. Theory and Practice of Histological Techniques, 5th Edition JD Bancroft, M Gamble, Editors Churchill Livingstone, New York 2002 p 244

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Références

1. Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, CV Mosby, St. Louis (MO) 1980 p 217
2. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
3. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
4. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
5. Carson: A Self Instructional Text. ASCP Press, Chicago (IL), 1990
6. Atlas of Cytochemistry and Immunochemistry of Hematologic Neoplasms. T Sun, C-Y Li, LT Yam, Editors, ASCP, Chicago (IL), 1985, p 205
7. Theory and Practice of Histological Techniques, 5th Edition JD Bancroft, M Gamble, Editors Churchill Livingstone, New York 2002 p 244

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Referências

1. Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, CV Mosby, St. Louis (MO) 1980 p 217
2. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
3. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
4. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
5. Carson: A Self Instructional Text. ASCP Press, Chicago (IL), 1990
6. Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. T Sun, C-Y Li, LT Yam, Editors, ASCP, Chicago (IL), 1985, p 205
7. Theory and Practice of Histological Techniques, 5th Edition JD Bancroft, M Gamble, Editors Churchill Livingstone, New York 2002, p 244

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1. Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, CV Mosby, St. Louis (MO) 1980 p 217
2. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
3. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
4. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
5. Carson: A Self Instructional Text. ASCP Press, Chicago (IL), 1990
6. Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. T Sun, C-Y Li, LT Yam, Editors, ASCP, Chicago (IL), 1985, 205. o.
7. Theory and Practice of Histological Techniques, 5th Edition JD Bancroft, M Gamble, Editors Churchill Livingstone, New York 2002 p 244

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Reference

1. Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, CV Mosby, St. Louis (MO) 1980 p 217
2. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
3. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
4. Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986
5. Carson: A Self Instructional Text. ASCP Press, Chicago (IL), 1990
6. Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. T Sun, C-Y Li, LT Yam, Editors, ASCP, Chicago (IL), 1985, str. 205
7. Theory and Practice of Histological Techniques, 5th Edition JD Bancroft, M Gamble, Editors Churchill Livingstone, New York 2002 p 244

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Referanslar

1. Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, CV Mosby, St. Louis (MO) 1980 s. 217
2. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
3. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
4. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
5. Carson: A Self Instructional Text. ASCP Press, Chicago (IL), 1990
6. Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. T Sun, C-Y Li, LT Yam, Editors, ASCP, Chicago (IL), 1985, s. 205
7. Theory and Practice of Histological Techniques, 5th Edition JD Bancroft, M Gamble, Editors Churchill Livingstone, New York 2002 s. 244

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MDSS GmbH