


**MERCK**

# Millipore Express<sup>®</sup> Family Guide



The Life Science business  
of Merck operates as  
MilliporeSigma in the  
U.S. and Canada.

**Millipore<sup>®</sup>**

Expert Pharm/BioPharm  
Products & CTDMO Services

# Overview

Choosing the right membrane filter is critical to the success of your process. A trusted name in the industry, Millipore Express® filters contain hydrophilic polyethersulfone (PES) membranes that offer sterile filtration with high flux and capacity. Filters containing Millipore Express® membranes provide broad chemical compatibility across a wide pH range, are easy to wet and integrity test, and are available in a range of formats to meet the specific needs of different biomanufacturing operations.

Our Millipore Express® filters can be used in a wide range of applications including filtration of cell culture media and feeds and for sterile filtration at different steps in monoclonal antibodies (mAbs), vaccine, plasma, viral vector and large and small volume parenteral production.

## Our Millipore Express® Filter Family

### Millipore Express® SHR Filters

#### Sterile, High Retention

Sterilizing-grade filters for efficient filtration of cell culture media and feeds.

An optional, integrated PES membrane prefilter protects the high flux sterilizing 0.1 µm membrane from premature plugging. These filters offer mycoplasma reduction.

### Millipore Express® SHC Filters

#### Sterile, High Capacity

High capacity, sterilizing-grade filters for plugging streams. Contain two PES membrane layers; 0.5/0.2 µm.

### Millipore Express® SHF Filters

#### Sterile, High Flux

High flow, sterilizing-grade filters for critical process steps. Contain one layer of 0.2 µm PES membrane.

### Millipore Express® PHF Filters

#### Process Protection, High Flux

High flow filters for cost-effective sterile filtration of buffers. Contain one layer of sterilizing-grade 0.2 µm PES membrane.

# Selection & Applications Guide

We offer a full portfolio of filters to meet the needs of different bioprocess applications.

The table below provides a high-level overview of key applications and our filtration solutions as a starting point for development or optimization.

	mAb & r-protein intermediates	Plasma	Vaccines & viral vectors	Ophthalmics	SVPs	LVPs	Cell Culture Media Filtration	Buffer Filtration	Final Filtration	Gas	Colloids	Lipid Removal
<b>Particulate removal and sterile filter protection</b>												
Milligard PES® filters	✓	✓	✓	✓	✓	✓	✓	✓			✓	
Milligard® filters	✓	✓	✓		✓		✓				✓	
Polysep™ II filters	✓	✓	✓				✓				✓	✓
Lifegard™ filters		✓					✓				✓	✓
<b>Bioburden reduction</b>												
Milligard PES® filters	✓	✓	✓		✓	✓		✓				
Durapore® 0.45 µm filters		✓	✓		✓	✓						
<b>Sterile filtration</b>												
Millipore Express® SHC filters	✓	✓	✓		✓	✓	✓	✓	✓		✓	✓
Millipore Express® SHF filters			✓	✓	✓	✓			✓			
Millipore Express® PHF filters	✓	✓				✓		✓				
Durapore® 0.22 µm filters		✓	✓	✓	✓	✓			✓			
Durapore® Multilayer filters			✓								✓	✓
Aervent® filters										✓		
Aerex® filters										✓		
Millipore Express® SPG filters										✓		
<b>Mycoplasma removal and sterile filtration</b>												
Millipore Express® SHR filters							✓					

# Formats

Multiple formats offer flexibility and scalable solutions for both single-use and stainless steel operations.

	Format Size	Millipore Express® SHR Filters	Millipore Express® SHC Filters	Millipore Express® SHF Filters	Millipore Express® PHF Filters
<b>Single-use capsule filters</b>					
OptiScale® Capsules	25	<b>A</b>	<b>A</b>	<b>A</b>	<b>A</b>
	47	<b>A</b>	<b>A</b>	<b>A</b>	
Opticap® XL Capsules	150/300/600	<b>G, S</b>	<b>G, S</b>	<b>G, S</b>	<b>G, S</b>
	3/5/10	<b>G, S, A</b>	<b>G, S, A</b>	<b>G, S, A</b>	<b>G, S, A</b>
Opticap® XLT Capsules	10/20/30 Standard Area	<b>G, S, A</b>	<b>G, S, A</b>	<b>G, S, A</b>	<b>G, S, A</b>
	10/20/30 High Area	<b>G*</b>	<b>G</b>		
<b>Cartridge filters for stainless steel operations</b>					
Cartridge	5	<b>A*</b>	<b>A*</b>	<b>A*</b>	<b>A*</b>
	10/20/30 Standard Area	<b>A*</b>	<b>A*</b>	<b>A*</b>	<b>A*</b>
	10/20/30 High Area	<b>A**</b>	<b>A*</b>		

## Key

**G** = Gamma and X-ray compatible

**G\*** = Gamma and X-ray compatible: product is gamma and X-ray compatible and only available in Millipore Express® SHR with Prefilter capsules

**A** = Autoclavable: product can be autoclaved

**A\*** = Autoclavable: product can be autoclaved and is compatible with steam-in-place (SIP) sterilization methods

**A\*\*** = Autoclavable, compatible with SIP, only available in Millipore Express® SHR with Prefilter capsules

**S** = Sterile: product has been presterilized by gamma irradiation

## Quality Information

Filters with Millipore Express® membranes are designed, developed, and manufactured in accordance with a Quality Management System approved by an accredited registering body to an ISO 9001 Quality Systems Standard. Each Millipore Express® filter has a Certificate of Quality, which is available online. Millipore Express® SHR, Millipore Express® SHC and Millipore Express® SHF filters are integrity tested during manufacturing and are supported with an Emprove® documentation package or Validation Guide. For traceability and easy identification, each device is marked with the product name and identifying characteristics.

# The Emprove<sup>®</sup> Program

## Your fast track through regulatory challenges

Complementing our product portfolio, the Emprove<sup>®</sup> Program provides three types of dossiers to support different stages of development and manufacturing operations such as qualification, risk assessment and process optimization. The dossiers consolidate comprehensive product-specific testing data, quality statements and regulatory information in a readily-available format to simplify your compliance needs.

For more information, please visit:  
[SigmaAldrich.com/Emprove](https://SigmaAldrich.com/Emprove)




# Mobius<sup>®</sup> Single-use Solutions

## Millipore Express<sup>®</sup> filters are part of the Mobius<sup>®</sup> component library.

Whether you are looking to introduce single-use manufacturing components into your current process or investigating how you can implement a single-use process train, Mobius<sup>®</sup> products and solutions help meet your evolving process needs.

For more information, please visit:  
[SigmaAldrich.com/single-use-assemblies](https://SigmaAldrich.com/single-use-assemblies)



# OptiScale® Capsules

## For Filter Screening and Scaling

Our OptiScale® disposable capsule filters provide a convenient small-volume option for process development screening and scaling. They are ideal for quickly evaluating performance of different filters with various process streams.



OptiScale® Capsules

# Cartridge Filters

## For Pilot and Production-Scale Processing


Our cartridge filters are designed for pilot and production-scale processing in stainless steel housings. These filters provide high flow rates and throughput with minimal differential pressure and are designed to withstand multiple steam-in-place cycles. A full range of filtration areas are available for maximum flexibility.

**Standard area cartridge filters** contain membrane with a conventional pleat pattern.

**High area cartridge filters** contain membrane with an M-pleat pattern which doubles the membrane area compared to standard area filters. High area filters are designed to maximize filtration area while minimizing filter footprint.



Cartridge Filters



# Capsule Filters

For Pilot and Production-Scale Processing

## Opticap® XL and XLT Capsules

Opticap® XL and XLT capsule filters allow unparalleled hydraulic stress resistance in a single-use filter and are available from small-scale to 30-inch formats, enabling easy scale-up to pilot and production-scale processing. They are available with a range of inlet/outlet connections, and are offered in autoclavable, presterilized, gamma and X-ray compatible formats. These capsules minimize cleaning, assembly and validation requirements which translates to increased flexibility, more rapid turnaround and less downtime than maintaining stainless steel operations.

Opticap® XL 150, 300 and 600 capsule filters have the option to add a filling bell to protect an open container from airborne particles.

The T-line design of Opticap® XLT capsule filters accommodates series or parallel filtration, and a specially designed stand enables quick and easy integration into existing operations.

**Standard area capsule filters** contain membrane with a conventional pleat pattern.

**High area capsule filters** contain membrane with an M-pleat pattern which doubles the membrane area compared to standard area filters. High area filters are designed to maximize filtration area while minimizing filter footprint.



Opticap® XL and XLT Capsules shown as gamma and X-ray compatible



Opticap® XL 150, 300, and 600 capsules shown with optional filling bell attachment



Opticap® XLT Filters shown as autoclavable

# Millipore Express® SHR Filters

Sterilizing-grade filters for efficient filtration of cell culture media and feeds

Filters containing Millipore Express® SHR (Sterile, High Retention) membrane provide sterilizing-grade performance and mycoplasma removal from cell culture media, media additives and other biological solutions. These filters contain a 0.1 µm polyethersulfone (PES) membrane that provides high sterility assurance, broad chemical compatibility, high flow rates and extended throughput.

Millipore Express® SHR with Prefilter (SHRP) filters contain an integrated PES 0.5 µm membrane prefilter that protects the sterilizing-grade 0.1 µm membrane from premature plugging and extends filtration capacity in fouling streams.



## Benefits

- Reliable mycoplasma removal and sterilizing-grade performance
- Available with an on-board PES membrane prefilter for extended throughput in high-fouling solutions (SHRP)
- Broad chemical compatibility across a wide pH range
- 100% integrity tested during the manufacturing process

## Filter Formats

- Cartridge filters: standard and high area (SHRP only)
- Opticap® XL and XLT capsule filters
  - Sterile, gamma and X-ray compatible, or autoclavable
  - Standard and high area (XLT only)

# Millipore Express® SHR Filters

## Enhanced Sterility Assurance and Filter Capacity

Millipore Express® SHR filters are designed for the removal of mycoplasma and small microorganisms that could pass through 0.2 µm rated sterilizing-grade filters. These filters typically demonstrate Log Reduction Value (LRV) >6 with *Acholeplasma laidlawii* ATCC® 23206 using our validated test method.

Millipore Express® SHRP filters are designed to efficiently process cell culture media and feeds. The high permeability of these membranes minimizes filtration area requirements and their high capacity for particulates makes them an efficient option for plugging streams.

## Mobius® Single-use Solutions

Millipore Express® SHR filters are part of the Mobius® library providing you with the flexibility to design single-use assemblies that meet your specific processing requirements.


For more information, please visit:  
[SigmaAldrich.com/single-use-assemblies](http://SigmaAldrich.com/single-use-assemblies)

## The Emprove® Program – Your Fast Track through Regulatory Challenges

Complementing our product portfolio, the Emprove® Program provides three types of dossiers to support different stages of development and manufacturing operations such as qualification, risk assessment and process optimization.


For more information, please visit:  
[SigmaAldrich.com/Emprove](http://SigmaAldrich.com/Emprove)

## High Retention




**Figure 1.** Retention of *Acholeplasma laidlawii* by Millipore Express® SHR and Millipore Express® SHR with prefilter filters (n=3) at a challenge level of 10<sup>7</sup> CFU/cm<sup>2</sup> of filter area.

## High Permeability



**Figure 2.** Permeability per 10-inch cartridge. Tested in duplicate 47 mm discs at 10 psi and scaled to 10-inch cartridge. PVDF: polyvinylidene fluoride; PES: polyethersulfone.

## Superior Efficiency Reduces Filter Footprint



**Figure 3.** Number of standard 10-inch cartridges needed to filter 10,000 L of NSO Serum-free, Protein-free, or CHO Serum-free Cell Culture Medium in 2 hours at 10 psi (n=2).

# Millipore Express® SHR Filters

## OptiScale® Capsule and Cartridge Filter Specifications

Description	OptiScale® 25 Capsules	OptiScale® 47 Capsules	5-inch Standard Area Cartridge	Per Standard 10-inch Cartridge	Per High Area 10-inch Cartridge*
<b>Dimensions</b>					
Diameter	31 mm (1.21 in.)	70 mm (2.75 in.)	6.9 cm (2.7 in.)	6.9 cm (2.7 in.)	7.1 cm (2.8 in.)
Maximum Length	39 mm (1.52 in.)	94 mm (3.70 in.)	12.7 cm (5 in.)	25.4 cm (10 in.)	25.4 cm (10 in.)
<b>Filtration Area</b>					
Millipore Express® SHR	3.5 cm <sup>2</sup>	17.7 cm <sup>2</sup>	0.29 m <sup>2</sup> (3.1 ft <sup>2</sup> )	0.60 m <sup>2</sup> (6.5 ft <sup>2</sup> )	-
Millipore Express® SHR with Prefilter	3.5 cm <sup>2</sup>	17.7 cm <sup>2</sup>	0.23 m <sup>2</sup> (2.5 ft <sup>2</sup> )	0.49 m <sup>2</sup> (5.3 ft <sup>2</sup> )	1.0 m <sup>2</sup> (10.8 ft <sup>2</sup> )
<b>Materials of Construction</b>					
Filter membrane	Hydrophilic polyethersulfone (PES)		Hydrophilic PES		Hydrophilic PES
Film Edge	—	—	Polypropylene		Polypropylene
Supports	Polypropylene		Polypropylene		Polypropylene
Vent Cap	Polypropylene	Polyvinylidene fluoride	—		—
Structural Components	Polypropylene	Polycarbonate	Polypropylene		Polypropylene
Core	—	—	Polysulfone		Polyethersulfone
O-Rings	—	Fluoroelastomer	Silicone, EPDM, or Fluoroelastomer		Silicone, EPDM, or Fluoroelastomer
<b>Maximum Inlet Pressure</b>					
	4.1 bar (60 psi) at 25 °C	5.1 bar (80 psi) at 25 °C			
<b>Maximum Differential Pressure</b>					
Forward:	4.1 bar (60 psi) at 25 °C	5.5 bar (80 psi) at 25 °C	6.9 bar (100 psi) at 25 °C 1.7 bar (25 psi) at 80 °C 1 bar (15 psi) at 135 °C		6.9 bar (100 psi) at 25 °C 1.7 bar (25 psi) at 80 °C 300 mbar (5 psi) at 135 °C
Reverse:	0 bar (0 psi)	690 mbar (10 psi) at 25 °C	2.1 bar (30 psi) at 25 °C 69 mbar (1 psi) at 135 °C		2.1 bar (30 psi) at 25 °C 69 mbar (1 psi) at 135 °C
<b>70/30 IPA/Water Bubble Point</b>					
	—	—	≥ 2590 mbar (37.5 psi) with nitrogen		
<b>Air Diffusion</b>					
			<b>Through a water wet membrane at 3450 mbar (50 psi):</b>		
Millipore Express® SHR	—	—	≤ 15.9 cc/min.	≤ 33.3 cc/min.	-
Millipore Express® SHR with Prefilter	—	—	≤ 12.8 cc/min.	≤ 27.1 cc/min.	≤ 54.2 cc/min.
<b>Bacterial Retention</b>					
			Quantitative retention of 10 <sup>7</sup> CFU/cm <sup>2</sup> <i>Brevundimonas diminuta</i> ATCC® 19146 per ASTM F838 methodology.		
<b>Mycoplasma Removal</b>					
			Typical Log Reduction Value (LRV) > 6 using <i>Acholeplasma laidlawii</i> ATCC® 23206 and our validated test method.		
<b>Bacterial Endotoxin</b>					
	Aqueous extraction contains < 0.25 EU/mL as determined by the Limulus Ameboyte Lysate (LAL) Test as described in USP <85>. Specification not applicable to OptiScale® 47 capsules.				

# Millipore Express® SHR Filters

## OptiScale® Capsule and Cartridge Filter Specifications (cont.)

Description	OptiScale® 25 Capsules	OptiScale® 47 Capsules	5-inch Standard Area Cartridge	Per Standard 10-inch Cartridge	Per High Area 10-inch Cartridge*
<b>Total Organic Carbon (TOC)/ Conductivity</b>	Autoclaved filter effluent meets the WFI requirement of USP <643> for Total Organic Carbon and USP <645> for Water Conductivity at 25 °C after a WFI flush of:				
Millipore Express® SHR	15 mL	—	5.5 L	10 L	—
Millipore Express® SHR with Prefilter	15 mL	—	9.5 L	20 L	20 L
<b>Oxidizable Substances</b>					
Millipore Express® SHR	—	100 mL	—	—	—
Millipore Express® SHR with Prefilter	—	100 mL	—	—	—
<b>Sterilization</b>					
Millipore Express® SHR					
Autoclave	1 cycle at 123 °C for 60 min.	3 cycles at 126 °C for 60 min.	25x, 60 min. cycles at 126 °C		—
In-line Steam	—	—	25 forward cycles, 30 min., 135 °C at ≤ 340 mbar (5 psi) or, 22 forward cycles, 30 min., 135 °C at ≤ 340 mbar (5 psi) and 3 reverse cycles, 30 min. at ≤ 69 mbar (1 psi)		—
Millipore Express® SHR with Prefilter					
Autoclave	1 cycle at 123 °C for 60 min.	3 cycle at 126 °C for 60 min.	25x, 60 min. cycles at 126 °C		5x, 60 min. cycles at 126 °C
In-line Steam	—	—	25 forward cycles, 30 min., 135 °C at ≤ 340 mbar (5 psi) or, 22 forward cycles, 30 min., 135 °C at ≤ 340 mbar (5 psi) and 3 reverse cycles, 30 min. at ≤ 69 mbar (1 psi) or, 3 forward cycles, 30 min., 135 °C at ≤ 1.0 bar (15 psi) or, 3 forward cycles, 30 min, 145 °C at ≤ 69 mbar (1 psi)		5x forward cycles, 30 min., 135 °C at ≤ 340 mbar (5 psi) or 2x forward cycles, 30 min., 145 °C at ≤ 69 mbar (1 psi)
<b>Toxicity</b>	Component materials meet the criteria for Biological Reactivity Testing. These tests can be any one or a combination of the following test methods USP <88> Class VI ( <i>in vivo</i> ), USP <87> ( <i>in vitro</i> ), ISO 10993-5 ( <i>in vitro</i> ).				
<b>Particle Shedding</b>	Effluent meets the acceptance criteria set forth in USP <788> for large volume parenterals. Specification does not apply to OptiScale® 47 capsules.				
<b>Non-fiber Releasing</b>	Meets the criteria for a 'non-fiber releasing' filter as defined in 21 CFR 210.3 (b)(6).				
<b>Indirect Food Additive</b>	All component materials meet the FDA Indirect food additives requirement cited in 21 CFR 177-182 based on information provided by raw material suppliers.				

\* High Area available in Millipore Express® SHR with Prefilter only

# Millipore Express® SHR Filters

## Opticap® XL and XLT Autoclavable Capsule Filter Specifications

Description	Opticap® XL 3 Standard Area Capsules	Opticap® XL 5 Standard Area Capsules	Opticap® XL 10 Standard Area Capsules	Opticap® XLT 10 Standard Area Capsules	Opticap® XLT 20 Standard Area Capsules	Opticap® XLT 30 Standard Area Capsules
<b>Dimensions</b>						
Body Diameter:	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)
Maximum Width:	—	—	—	19.8 cm (7.8 in.)	19.8 cm (7.8 in.)	19.8 cm (7.8 in.)
Maximum Length:	17.3 cm (6.8 in.)	21.6 cm (8.5 in.)	33.5 cm (13.2 in.)	38.1 cm (15.0 in.)	62.5 cm (24.6 in.)	87.1 cm (34.3 in.)
<b>Filtration Area</b>						
Millipore Express® SHR	0.16 m <sup>2</sup> (1.7 ft <sup>2</sup> )	0.29 m <sup>2</sup> (3.1 ft <sup>2</sup> )	0.60 m <sup>2</sup> (6.5 ft <sup>2</sup> )	0.60 m <sup>2</sup> (6.5 ft <sup>2</sup> )	1.21 m <sup>2</sup> (13.0 ft <sup>2</sup> )	1.81 m <sup>2</sup> (19.5 ft <sup>2</sup> )
Millipore Express® SHR with Prefilter	0.13 m <sup>2</sup> (1.4 ft <sup>2</sup> )	0.23 m <sup>2</sup> (2.5 ft <sup>2</sup> )	0.49 m <sup>2</sup> (5.3 ft <sup>2</sup> )	0.49 m <sup>2</sup> (5.3 ft <sup>2</sup> )	0.98 m <sup>2</sup> (10.6 ft <sup>2</sup> )	1.48 m <sup>2</sup> (15.9 ft <sup>2</sup> )
<b>Materials of Construction</b>						
Filter Membrane	Hydrophilic polyethersulfone (PES)					
Film Edge	Polypropylene					
Supports	Polypropylene					
Core	Polysulfone					
Housing and Cage	Polypropylene					
O-Rings	Silicone					
<b>Vent/Drain</b>	6 mm (¼ in.) hose barb with double O-ring seal					
<b>Maximum Inlet Pressure</b>	6.9 bar (100 psi) intermittent at 25 °C 5.5 bar (80 psi) at 25 °C 2.8 bar (40 psi) at 60 °C 1 bar (15 psi) at 80 °C					
<b>Maximum Differential Pressure</b>						
Forward:	6.9 bar (100 psi) intermittent at 25 °C 5.5 bar (80 psi) at 25 °C 1 bar (15 psi) at 80 °C					
Reverse:	2.1 bar (30 psi) intermittent at 25°C					
<b>70/30 IPA/Water Bubble Point</b>	≥ 2590 mbar (37.5 psi) with nitrogen					
<b>Air Diffusion</b>						
	<b>Through a water wet membrane at 3450 mbar (50 psi):</b>					
Millipore Express® SHR	≤ 8.7 cc/min.	≤ 15.9 cc/min.	≤ 33.3 cc/min.	≤ 33.3 cc/min.	≤ 66.6 cc/min.	≤ 99.9 cc/min.
Millipore Express® SHR with Prefilter	≤ 7.2 cc/min.	≤ 12.8 cc/min.	≤ 27.1 cc/min.	≤ 27.1 cc/min.	≤ 54.2 cc/min.	≤ 81.3 cc/min.
<b>Bacterial Retention</b>	Quantitative retention of 10 <sup>7</sup> CFU/cm <sup>2</sup> <i>Brevundimonas diminuta</i> ATCC® 19146 per ASTM F838 methodology.					
<b>Mycoplasma Removal</b>	Typical Log Reduction Value (LRV) > 6 using <i>Acholeplasma laidlawii</i> ATCC® 23206 and our validated test method.					
<b>Bacterial Endotoxin</b>	Aqueous extraction contains < 0.25 EU/mL as determined by the Limulus Ameboyte Lysate (LAL) Test as described in USP <85>.					

## Millipore Express® SHR Filters


### Opticap® XL and XLT Autoclavable Capsule Filter Specifications (cont.)

Description	Opticap® XL 3 Standard Area Capsules	Opticap® XL 5 Standard Area Capsules	Opticap® XL 10 Standard Area Capsules	Opticap® XLT 10 Standard Area Capsules	Opticap® XLT 20 Standard Area Capsules	Opticap® XLT 30 Standard Area Capsules
<b>Total Organic Carbon (TOC)/Conductivity</b>	Autoclaved filter effluent meets the WFI requirement of USP <643> for Total Organic Carbon and USP <645> for Water Conductivity at 25 °C after a WFI flush of:					
Millipore Express® SHR	3.0 L	5.5 L	10 L	10 L	20 L	30 L
Millipore Express® SHR with Prefilter	5.5 L	9.5 L	20 L	20 L	40 L	60 L
<b>Sterilization</b>	May be autoclaved for 3 cycles for 60 minutes at 126 °C. Cannot be steam sterilized in-line.					
<b>Toxicity</b>	Component materials meet the criteria for Biological Reactivity Testing. These tests can be any one or a combination of the following test methods USP <88> Class VI ( <i>in vivo</i> ), USP <87> ( <i>in vitro</i> ), ISO 10993-5 ( <i>in vitro</i> ).					
<b>Particle Shedding</b>	Effluent meets the acceptance criteria set forth in USP <788> for large volume parenterals.					
<b>Non-fiber Releasing</b>	Meets the criteria for a 'non-fiber releasing' filter as defined in 21 CFR 210.3 (b)(6).					
<b>Indirect Food Additive</b>	All component materials meet the FDA Indirect food additives requirement cited in 21 CFR 177-182 based on information provided by raw material suppliers.					


# Millipore Express® SHR Filters

## Typical Clean Water Flow Rates – Cartridge Filters


Standard Area Cartridge Filters with Millipore Express® SHR 0.1 µm Membrane



Standard Area Cartridge Filters with Millipore Express® SHR 0.5/0.1 µm Membrane with Prefilter




High Area Cartridge Filters with Millipore Express® SHR 0.5/0.1 µm Membrane with Prefilter




# Millipore Express® SHR Filters

## Typical Clean Water Flow Rates – Opticap® XL and XLT Autoclavable Capsules Containing Millipore Express® SHR Membrane


**Opticap® XL 3 Capsule Filters with 0.1 µm Millipore Express® SHR Membranes**




**Opticap® XL 5 Capsule Filters with 0.1 µm Millipore Express® SHR Membranes**




**Opticap® XL 10 Capsule Filters with 0.1 µm Millipore Express® SHR Membranes**




**Opticap® XLT 10 Capsule Filters with 0.1 µm Millipore Express® SHR Membranes**



**Opticap® XLT 20 Capsule Filters with 0.1 µm Millipore Express® SHR Membranes**



**Opticap® XLT 30 Capsule Filters with 0.1 µm Millipore Express® SHR Membranes**



### Opticap® XL Capsule Legends Refer to Connection Type

- TT = 38 mm (1½ in.) sanitary flange inlet and outlet
- HH = 14 mm (½ in.) hose barb inlet and outlet
- TH = 38 mm (1½ in.) sanitary flange inlet and 14 mm (½ in.) hose barb outlet
- TB = 38 mm (1½ in.) sanitary flange inlet and 25 mm (1 in.) hose barb outlet


### Opticap® XLT Capsule Legends Refer to Connection Type

- TT = 38 mm (1½ in.) sanitary flange inlet and outlet
- TH = 38 mm (1½ in.) sanitary flange Inlet and 16 mm (⅝ in.) hose barb outlet
- HH = 16 mm (⅝ in.) hose barb inlet and outlet
- TB = 38 mm (1½ in.) sanitary flange inlet and 25 mm (1 in.) hose barb outlet
- BB = 25 mm (1 in.) hose barb inlet and outlet


# Millipore Express® SHR Filters

## Typical Clean Water Flow Rates – Opticap® XL and XLT Autoclavable Capsules Containing Millipore Express® SHR Membrane with Prefilter


**Opticap® XL 3 Capsule Filters with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**




**Opticap® XL 5 Capsule Filters with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**




**Opticap® XL 10 Capsule Filters with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**




**Opticap® XLT 10 Capsule Filters with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**



**Opticap® XLT 20 Capsule Filters with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**



**Opticap® XLT 30 Capsule Filters with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**



**Opticap® XL Capsule Legends Refer to Connection Type**

- TT = 38 mm (1½ in.) sanitary flange inlet and outlet
- HH = 14 mm (½ in.) hose barb inlet and outlet
- TH = 38 mm (1½ in.) sanitary flange inlet and 14 mm (½ in.) hose barb outlet
- TB = 38 mm (1½ in.) sanitary flange inlet and 25 mm (1 in.) hose barb outlet

**Opticap® XLT Capsule Legends Refer to Connection Type**

- TT = 38 mm (1½ in.) sanitary flange inlet and outlet
- TH = 38 mm (1½ in.) sanitary flange Inlet and 16 mm (⅝ in.) hose barb outlet
- HH = 16 mm (⅝ in.) hose barb inlet and outlet
- TB = 38 mm (1½ in.) sanitary flange inlet and 25 mm (1 in.) hose barb outlet
- BB = 25 mm (1 in.) hose barb inlet and outlet

# Millipore Express® SHR Filters

## Opticap® XL 150, 300 and 600 Sterile, Gamma and X-ray Compatible Capsule Filter Specifications


Description	Opticap® XL 150 Standard Area Capsules	Opticap® XL 300 Standard Area Capsules	Opticap® XL 600 Standard Area Capsules
<b>Dimensions</b>			
Body Diameter	5.6 cm (2.2 in.)	5.6 cm (2.2 in.)	5.6 cm (2.2 in.)
Maximum Length	9.7 cm (3.8 in.)	11.9 cm (4.7 in.)	16.5 cm (6.5 in.)
Maximum Length with Optional Bell	11.4 cm (4.5 in.)	13.7 cm (5.4 in.)	18.3 cm (7.2 in.)
<b>Filtration Area</b>			
Millipore Express® SHR	0.025 m <sup>2</sup> (0.268 ft <sup>2</sup> )	0.048 m <sup>2</sup> (0.517 ft <sup>2</sup> )	0.101 m <sup>2</sup> (1.090 ft <sup>2</sup> )
Millipore Express® SHR with Prefilter	0.015 m <sup>2</sup> (0.163 ft <sup>2</sup> )	0.028 m <sup>2</sup> (0.308 ft <sup>2</sup> )	0.062 m <sup>2</sup> (0.664 ft <sup>2</sup> )
<b>Materials of Construction</b>			
Filter Membrane	Hydrophilic polyethersulfone (PES)		
Supports	Polyethylene		
Core	Polysulfone		
Housing and Cage	Gamma and X-ray compatible polypropylene		
O-Rings	Silicone		
Optional Bell	Polycarbonate		
<b>Vent/Drain</b>	6 mm (¼ in.) hose barb with double O-ring seal		
<b>Maximum Inlet Pressure</b>	6.9 bar (100 psi) intermittent at 25 °C 5.5 bar (80 psi) at 25 °C 2.8 bar (40 psi) at 60 °C 1 bar (15 psi) at 80 °C		
<b>Maximum Differential Pressure</b>	Forward: 6.9 bar (100 psi) intermittent at 25 °C 5.5 bar (80 psi) at 25 °C 1 bar (15 psi) at 80 °C Reverse: 2.1 bar (30 psi) intermittent at 25 °C		
<b>70/30 IPA/Water Bubble Point</b>	≥ 2590 mbar (37.5 psi) with nitrogen		
<b>Air Diffusion</b>	<b>Through a water wet membrane at 3450 mbar (50 psi):</b>		
Millipore Express® SHR	≤ 1.4 cc/min.	≤ 2.8 cc/min.	≤ 5.8 cc/min.
Millipore Express® SHR with Prefilter	≤ 1.0 cc/min.	≤ 1.9 cc/min.	≤ 3.7 cc/min.
<b>Bacterial Retention</b>	Quantitative retention of 10 <sup>7</sup> CFU/cm <sup>2</sup> <i>Brevundimonas diminuta</i> ATCC® 19146 per ASTM F838 methodology.		
<b>Mycoplasma Removal</b>	Typical Log Reduction Value (LRV) > 6 using <i>Acholeplasma laidlawii</i> ATCC® 23206 and our validated test method.		
<b>Bacterial Endotoxin</b>	Aqueous extraction contains < 0.25 EU/mL as determined by the Limulus Ameboyte Lysate (LAL) Test as described in USP <85>.		
<b>Total Organic Carbon (TOC)/ Conductivity</b>	Sterilized filter effluent meets the WFI requirement of USP <643> for Total Organic Carbon and for USP <645> for Water Conductivity at 25 °C after a WFI flush of:		
Millipore Express® SHR	2 L	2.5 L	3 L
Millipore Express® SHR with Prefilter	1 L	2 L	3 L
<b>Sterilization</b>			
Gamma and X-ray Compatible Capsules	Compatible to 45 kGy and may be autoclaved for 3 cycles of 60 minutes at 123 °C. Cannot be steam sterilized in-line.		
Sterile Capsules	May be autoclaved for 3 cycles of 60 minutes at 123 °C. Cannot be steam sterilized in-line.		
<b>Sterility (Sterile Capsules)</b>	These capsules meet current USP and AAMI guidelines for sterility utilizing a validated sterilization cycle.		
<b>Toxicity</b>	Component materials meet the criteria for Biological Reactivity Testing. These tests can be any one or a combination of the following test methods USP<88> Class VI (in vivo), USP <87> (in vitro), ISO 10993-5 (in vitro).		
<b>Particle Shedding</b>	Effluent meets the acceptance criteria set forth in USP <788> for large volume parenterals.		
<b>Non-fiber Releasing</b>	Meets the criteria for a 'non-fiber releasing' filter as defined in 21 CFR 210.3 (b)(6).		
<b>Indirect Food Additive</b>	All component materials meet the FDA Indirect food additives requirement cited in 21 CFR 177-182 based on information provided by raw material suppliers.		

# Millipore Express® SHR Filters


## Typical Clean Water Flow Rates – Opticap® XL 150, 300 and 600 Sterile, Gamma and X-ray Compatible Capsules Containing Millipore Express® SHR Membrane

Filters tested post irradiation at 25–45 kGy and autoclaved at 123 °C for 60 minutes


**Opticap® XL 150 Capsule with 0.1 µm Millipore Express® SHR Membrane**



**Opticap® XL 300 Capsule with 0.1 µm Millipore Express® SHR Membrane**



**Opticap® XL 600 Capsule with 0.1 µm Millipore Express® SHR Membrane**



### Opticap® XL 150, 300 and 600 Capsule Connection Type


- FF = 19 mm (¾ in.) sanitary flange inlet and outlet
- FH = 19 mm (¾ in.) sanitary flange inlet and 14 mm (½ in.) hose barb outlet
- HH = 14 mm (½ in.) hose barb inlet and outlet

# Millipore Express® SHR Filters


## Typical Clean Water Flow Rates – Opticap® XL 150, 300 and 600 Sterile, Gamma and X-ray Compatible Capsules Containing Millipore Express® SHR Membrane with Prefilter

Filters tested post irradiation at 25–45 kGy and autoclaved at 123 °C for 60 minutes.


**Opticap® XL 150 Capsule with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**



**Opticap® XL 300 Capsule with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**



**Opticap® XL 600 Capsule with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**



### Opticap® XL 150, 300 and 600 Capsule Connection Type

FF = 19 mm (¾ in.) sanitary flange inlet and outlet

FH = 19 mm (¾ in.) sanitary flange inlet and 14 mm (⅜ in.) hose barb outlet

HH = 14 mm (⅜ in.) hose barb inlet and outlet

# Millipore Express® SHR Filters

## Opticap® XL and XLT Sterile, Gamma and X-ray Compatible Capsule Filter Specifications

Description	Opticap® XL 3 Standard Area Capsules	Opticap® XL 5 Standard Area Capsules	Opticap® XL 10 Standard Area Capsules	Per 10" Standard Area Opticap® XLT Capsules	Per 10" High Area Opticap® XLT Capsule*
<b>Dimensions</b>					
Body Diameter	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)
Maximum Width:	—	—	—	19.8 cm (7.8 in)	19.8 cm (7.8 in)
Maximum Length:	17.3 cm (6.8 in.)	21.6 cm (8.5 in.)	33.5 cm (13.2 in.)	—	—
10 inch	—	—	—	37.6 cm (14.8 in.)	37.6 cm (14.8 in.)
20 inch	—	—	—	62.5 cm (24.6 in.)	62.5 cm (24.6 in.)
30 inch	—	—	—	87.1 cm (34.3 in.)	87.1 cm (34.3 in.)
<b>Filtration Area</b>					
Millipore Express® SHR	0.17 m <sup>2</sup> (1.8 ft <sup>2</sup> )	0.31 m <sup>2</sup> (3.3 ft <sup>2</sup> )	0.69 m <sup>2</sup> (7.4 ft <sup>2</sup> )	0.69 m <sup>2</sup> (7.4 ft <sup>2</sup> )	—
Millipore Express® SHR with Prefilter Per 10 inch	—	—	—	0.54m <sup>2</sup> (5.8 ft <sup>2</sup> )	1.0 m <sup>2</sup> (10.8 ft <sup>2</sup> )
<b>Materials of Construction</b>					
Filter Membrane:	Hydrophilic polyethersulfone (PES)				Hydrophilic PES
Film Edge:	Polyethylene				Polyethylene
Supports:	Polyester				Polyester
Core:	Polysulfone				Polyethersulfone
Housing and Cage:	Gamma and X-ray compatible Polypropylene				Gamma and X-ray compatible Polypropylene
O-Rings:	Silicone				Silicone
<b>Vent/Drain</b>	6 mm (¼ in.) hose barb with double O-ring seal				
<b>Maximum Inlet Pressure</b>	6.9 bar (100 psi) intermittent at 25 °C 5.5 bar (80 psi) at 25 °C 2.8 bar (40 psi) at 60 °C 1 bar (15 psi) at 80 °C				
<b>Maximum Differential Pressure</b>	Forward: 6.9 bar (100 psi) intermittent at 25 °C 5.5 bar (80 psi) at 25 °C 1 bar (15 psi) at 80 °C Reverse: 2.1 bar (30 psi) intermittent at 25 °C				
<b>70/30 IPA/Water Bubble Point</b>	≥ 2590 mbar (37.5 psi) with nitrogen				
<b>Air Diffusion</b>	<b>Through a water wet membrane at 3450 mbar (50 psi):</b>				
Millipore Express® SHR	≤ 9.4 cc/min.	≤ 17.3 cc/min.	≤ 38.8 cc/min.	≤ 38.8 cc/min.	—
Millipore Express® SHR with Prefilter	≤ 7.3 cc/min.	≤ 13.6 cc/min.	≤ 30.4 cc/min.	≤ 30.4 cc/min.	≤ 54.2 cc/min.
<b>Bacterial Retention</b>	Quantitative retention of 10 <sup>7</sup> CFU/cm <sup>2</sup> <i>Brevundimonas diminuta</i> ATCC® 19146 per ASTM F838 methodology.				
<b>Mycoplasma Removal</b>	Typical Log Reduction Value (LRV) >6 using <i>Acholeplasma laidlawii</i> ATCC® 23206 and our validated test method.				
<b>Bacterial Endotoxin</b>	Aqueous extraction contains < 0.25 EU/mL as determined by the Limulus Ameboyte Lysate (LAL) Test as described in USP <85>.				

## Millipore Express® SHR Filters

### Opticap® XL and XLT Sterile, Gamma and X-ray Compatible Capsule Filter Specifications (cont.)

Description	Opticap® XL 3 Standard Area Capsules	Opticap® XL 5 Standard Area Capsules	Opticap® XL 10 Standard Area Capsules	Per 10" Standard Area Opticap® XLT Capsules	Per 10" High Area Opticap® XLT Capsule*
<b>Total Organic Carbon (TOC)/ Conductivity</b>	Sterilized filter effluent meets the WFI requirement of USP <643> for Total Organic Carbon and for USP <645> for Water Conductivity at 25 °C after a WFI flush of:				
Millipore Express® SHR	3.5 L	6.0 L	11 L	11 L	—
Millipore Express® SHR with Prefilter	5.5 L	9.5 L	21 L	21 L	21 L
<b>Sterilization</b>					
Gamma and X-ray Compatible Capsules	Compatible to 45 kGy and may be autoclaved for 3 cycles of 60 minutes at 123 °C. Cannot be steam sterilized in-line.				
Sterile:	May be autoclaved for 3 cycles of 60 minutes at 123 °C. Cannot be steam sterilized in-line.				
<b>Sterility (Sterile Capsules)</b>	These capsules meet current USP and AAMI guidelines for sterility utilizing a validated sterilization cycle.				
<b>Toxicity</b>	Component materials meet the criteria for Biological Reactivity Testing. These tests can be any one or a combination of the following test methods USP <88> Class VI ( <i>in vivo</i> ), USP <87> ( <i>in vitro</i> ), ISO 10993-5 ( <i>in vitro</i> ).				
<b>Particle Shedding</b>	Effluent meets the acceptance criteria set forth in USP <788> for large volume parenterals.				
<b>Non-fiber Releasing</b>	Meets the criteria for a 'non-fiber releasing' filter as defined in 21 CFR 210.3 (b)(6).				
<b>Indirect Food Additive</b>	All component materials meet the FDA Indirect food additives requirement cited in 21 CFR 177-182 based on information provided by raw material suppliers.				


\* Only available in Millipore Express® SHR with Prefilter gamma and X-ray compatible capsules

# Millipore Express® SHR Filters


## Typical Clean Water Flow Rates – Opticap® XL and XLT Sterile, Gamma and X-ray Compatible Capsules Containing Millipore Express® SHR Membrane

Filters tested post irradiation at 25–45 kGy and autoclaved at 123 °C for 60 minutes


Gamma and X-ray compatible Opticap® XL 3 Capsule Filters with 0.1 µm Millipore Express® SHR Membrane




Gamma and X-ray compatible Opticap® XL 5 Capsule Filters with 0.1 µm Millipore Express® SHR Membrane




Gamma and X-ray compatible Opticap® XL 10 Capsule Filters with 0.1 µm Millipore Express® SHR Membrane




Gamma and X-ray compatible Opticap® XLT 10 Capsule Filters with 0.1 µm Millipore Express® SHR Membrane



Gamma and X-ray compatible Opticap® XLT 20 Capsule Filters with 0.1 µm Millipore Express® SHR Membrane



Gamma and X-ray compatible Opticap® XLT 30 Capsule Filters with 0.1 µm Millipore Express® SHR Membrane



### Opticap® XL Capsule Legends Refer to Connection Type

- TT = 38 mm (1½ in.) sanitary flange inlet and outlet
- HH = 14 mm (⅝ in.) hose barb inlet and outlet
- TH = 38 mm (1½ in.) sanitary flange inlet and 14 mm (⅝ in.) hose barb outlet
- TB = 38 mm (1½ in.) sanitary flange inlet and 25 mm (1 in.) hose barb outlet

### Opticap® XLT Capsule Legends Refer to Connection Type


- TT = 38 mm (1½ in.) sanitary flange inlet and outlet
- TH = 38 mm (1½ in.) sanitary flange Inlet and 16 mm (⅝ in.) hose barb outlet
- HH = 16 mm (⅝ in.) hose barb inlet and outlet
- TB = 38 mm (1½ in.) sanitary flange inlet and 25 mm (1 in.) hose barb outlet
- BB = 25 mm (1 in.) hose barb inlet and outlet

# Millipore Express® SHR Filters


## Opticap® XL and XLT Sterile, Gamma and X-ray Compatible Capsules Containing Millipore Express® SHR Membrane with Prefilter

Filters tested post irradiation at 25–45 kGy and autoclaved at 123 °C for 60 minutes


**Gamma and X-ray compatible Opticap® XL 3 with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**



**Gamma and X-ray compatible Opticap® XL 5 with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**



**Gamma and X-ray compatible Opticap® XL 10 with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**




**Opticap® XL Capsule Legends Refer to Connection Type**

- TT = 38 mm (1½ in.) sanitary flange inlet and outlet
- HH = 14 mm (½ in.) hose barb inlet and outlet
- TH = 38 mm (1½ in.) sanitary flange inlet and 14 mm (½ in.) hose barb outlet
- TB = 38 mm (1½ in.) sanitary flange inlet and 25 mm (1 in.) hose barb outlet


**Opticap® XLT Capsule Legends Refer to Connection Type**

- TT = 38 mm (1½ in.) sanitary flange inlet and outlet
- TH = 38 mm (1½ in.) sanitary flange Inlet and 16 mm (⅝ in.) hose barb outlet
- HH = 16 mm (⅝ in.) hose barb inlet and outlet
- TB = 38 mm (1½ in.) sanitary flange inlet and 25 mm (1 in.) hose barb outlet
- BB = 25 mm (1 in.) hose barb inlet and outlet


**Gamma and X-ray compatible Standard Area Opticap® XLT 10 with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**




**Gamma and X-ray compatible Standard Area Opticap® XLT 20 with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**




**Gamma and X-ray compatible Standard Area Opticap® XLT 30 with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**




**Gamma and X-ray compatible High Area Opticap® XLT 10 with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**



**Gamma and X-ray compatible High Area Opticap® XLT 20 with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**




**Gamma and X-ray compatible High Area Opticap® XLT 30 with 0.5/0.1 µm Millipore Express® SHR Membrane with Prefilter**




# Millipore Express® SHR Filters

## Ordering Information

### OptiScale® Capsules




### Cartridge Filters




# Millipore Express® SHR Filters

## Ordering Information

### Opticap® XL 150, 300, 600 Capsule Filters




### Opticap® XL Capsule Filters



# Millipore Express® SHR Filters

## Ordering Information

### Opticap® XLT Capsule Filters



# Millipore Express® SHC Filters

High capacity, sterilizing-grade filters for plugging streams

Filters containing Millipore Express® SHC (Sterile, High Capacity) sterilizing-grade membrane provide superior throughput and capacity in applications where premature filter plugging is a concern. These filters contain two layers of polyethersulfone (PES) membrane (0.5  $\mu\text{m}$  and 0.2  $\mu\text{m}$ ) that provide sterility assurance, broad chemical compatibility, high flow rates and capacity.



## Benefits

- Sterilizing-grade membranes that are easy to wet and integrity test
- High flux, high capacity PES membrane that provides superior throughput in high fouling streams
- Broad chemical compatibility across a wide pH range
- 100% integrity tested during the manufacturing process


## Filter Formats

- OptiScale® capsules
- Cartridge filters: standard and high area
- Opticap® XL and XLT capsule filters:
  - Sterile, gamma and X-ray compatible, or autoclavable
  - Standard or high area (XLT only)


# Millipore Express® SHC Filters

## Extended Capacity and Fewer Filter Change

Millipore Express® SHC filters are designed to maximize the capacity of constrained filtration systems. With their high flux and superior capacity, these filters can double your output without adding to your filter footprint, Figures 1 and 2.



**Figure 1.** CHO Cell Culture Growth Media: No Prefiltration



**Figure 2.** Protein Concentrate Post Ultrafiltration

## Mobius® Single-use Solutions

Millipore Express® SHC filters are part of the Mobius® library providing you with the flexibility to design single-use assemblies that meet your specific processing requirements.

**For more information, please visit:**  
[SigmaAldrich.com/single-use-assemblies](https://SigmaAldrich.com/single-use-assemblies)

## The Emprove® Program – Your Fast Track through Regulatory Challenges

Complementing our product portfolio, the Emprove® Program provides three types of dossiers to support different stages of development and manufacturing operations such as qualification, risk assessment and process optimization.

**For more information, please visit:**  
[SigmaAldrich.com/Emprove](https://SigmaAldrich.com/Emprove)

# Millipore Express® SHC Filters

## OptiScale® Capsule and Cartridge Filter Specifications

Description	OptiScale® 25 Capsules	OptiScale® 47 Capsules	5-inch Standard Area Cartridges	Per 10-inch Standard Area Cartridge	Per 10-inch High Area Cartridge
<b>Dimensions</b>					
Diameter:	31 mm (1.21 in.)	70 mm (2.75 in.)	6.9 cm (2.7 in.)	6.9 cm (2.7 in.)	7.1 cm (2.8 in.)
Maximum Length:	39 mm (1.52 in.)	94 mm (3.70 in.)	12.7 cm (5 in.)	25.4 cm (10 in.)	25.4 cm (10 in.)
<b>Filtration Area</b>	3.5 cm <sup>2</sup>	17.7 cm <sup>2</sup>	0.23 m <sup>2</sup> (2.5 ft <sup>2</sup> )	0.49 m <sup>2</sup> (5.3 ft <sup>2</sup> )	1.0 m <sup>2</sup> (10.8 ft <sup>2</sup> )
<b>Materials of Construction</b>					
Filter membrane	Hydrophilic polyethersulfone (PES)		Hydrophilic PES		Hydrophilic PES
Film Edge	—	—	Polypropylene		Polypropylene
Supports	Polypropylene		Polypropylene		Polypropylene
Vent Cap	Polypropylene	Polyvinylidene fluoride	—		—
Structural Components	Polypropylene	Polycarbonate	Polypropylene		Polypropylene
Core	—	—	Polysulfone		Polyethersulfone
O-Rings	—	Fluoroelastomer	Silicone, EPDM, or Fluoroelastomer		Silicone, EPDM, or Fluoroelastomer
<b>Maximum Inlet Pressure</b>	4.1 bar (60 psi) at 25 °C	5.5 bar (80 psi) at 25 °C			
<b>Maximum Differential Pressure</b>					
Forward:	4.1 bar (60 psi) at 25 °C	5.5 bar (80 psi) at 25 °C	6.9 bar (100 psi) at 25 °C 1.7 bar (25 psi) at 80 °C 1 bar (15 psi) at 135 °C 0.34 bar (5 psi) at 135 °C	6.9 bar (100 psi) at 25 °C and 4 °C	1.7 bar (25 psi) at 80 °C  340 mbar (5 psi) at 135 °C
Reverse:	0 bar (0 psi)	690 mbar (10 psi) at 25 °C	2.1 bar (30 psi) at 25 °C 69 mbar (1 psi) at 135 °C	2.1 bar (30 psi) at 25 °C	69 mbar (1 psi) at 135 °C
<b>Bubble Point</b>	-	-	≥ 4000 mbar (58 psi) air with water ≥ 1280 mbar (18.5 psi) nitrogen with 70/30% IPA/water mixture		
<b>Air Diffusion</b>	-	-	Through a water wet membrane at 2800 mbar (40 psi): ≤ 13.3 cc/min.	≤ 28.2 cc/min.	≤ 56.4 cc/min.
<b>Bacterial Retention</b>	-	-	Quantitative retention of 10 <sup>7</sup> CFU/cm <sup>2</sup> <i>Brevundimonas diminuta</i> ATCC® 19146 per ASTM F838 methodology.		
<b>Bacterial Endotoxin</b>	Aqueous extraction contains < 0.25 EU/mL as determined by the Limulus Ameboyte Lysate (LAL) Test as described in USP <85>. Specification not applicable to OptiScale® 47 capsules.				
<b>Total Organic Carbon (TOC)/Conductivity</b>	Autoclaved filter effluent meets the WFI criteria for USP <643> Total Organic Carbon, and USP <645> Conductivity at 25 °C after a WFI flush of: 15 mL	-	9.5 L	20 L	20 L

# Millipore Express® SHC Filters

## OptiScale® Capsule and Cartridge Filter Specifications (cont.)

Description	OptiScale® 25 Capsules	OptiScale® 47 Capsules	5-inch Standard Area Cartridges	Per 10-inch Standard Area Cartridge	Per 10-inch High Area Cartridge
<b>Oxidizable Substances</b>	-	Meets the USP Oxidizable Substance Test requirements for sterile purified water after a water flush of:  100 mL	-	-	-
<b>Sterilization</b>					
Autoclave	1 cycle at 123 °C for 60 min.	3 cycles at 126 °C for 60 min.	25 cycles at 126 °C for 60 min.		5 cycles at 126 °C, for 60 min.
In-line Steam			25x forward cycles, 30 min., 135 °C at ≤ 340 mbar (5 psi) or 22x forward cycles, 30 min., 135 °C at ≤ 340 mbar (5 psi) and 3x reverse cycles, 30 min., 135 °C at ≤ 69 mbar (1 psi) or, 3x forward cycles, 30 min., 135 °C at ≤ 1.0 bar (15 psi), or 3x forward cycles, 30 min., 145 °C at ≤ 69 mbar (1 psi)		5x forward cycles, 30 min., 135 °C at ≤ 340 mbar (5 psi), or 5x forward cycles, 30 min., 145 °C at ≤ 69 mbar (1 psi), or 1x reverse cycle, 30 min., 125 °C at ≤ 69 mbar (1 psi)
<b>Toxicity</b>	Component materials meet the criteria for Biological Reactivity Testing. These tests can be any one or a combination of the following test methods USP<88> Class VI ( <i>in vivo</i> ), USP <87> ( <i>in vitro</i> ), ISO 10993-5 ( <i>in vitro</i> ).				
<b>Particle Shedding</b>	Effluent meets the acceptance criteria set forth in USP <788> for large volume parenterals.				
<b>Non-fiber Releasing</b>	Meets the criteria for a 'non-fiber releasing' filter as defined in 21 CFR 210.3 (b)(6).				
<b>Indirect Food Additive</b>	All component materials meet the FDA Indirect food additives requirement cited in 21 CFR 177-182 based on information provided by raw material suppliers.				

# Millipore Express® SHC Filters


## Opticap® XL and XLT Autoclavable Capsule Filter Specifications

Description	Opticap® XL 3 Standard Area Capsules	Opticap® XL 5 Standard Area Capsules	Opticap® XL 10 Standard Area Capsules	Opticap® XLT 10 Standard Area Capsules	Opticap® XLT 20 Standard Area Capsules	Opticap® XLT 30 Standard Area Capsules
<b>Dimensions</b>						
Body Diameter	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)
Maximum Width	-	-	-	19.8 cm (7.8 in)	19.8 cm (7.8 in)	19.8 cm (7.8 in)
Maximum Length	17.3 cm (6.8 in.)	21.6 cm (8.5 in.)	33.5 cm (13.2 in.)	38.1 cm (15.0 in.)	62.5 cm ( 24.6 in.)	87.1 cm (34.3 in.)
<b>Filtration Area</b>	0.13 m <sup>2</sup> (1.4 ft <sup>2</sup> )	0.23 m <sup>2</sup> (2.5 ft <sup>2</sup> )	0.49 m <sup>2</sup> (5.3 ft <sup>2</sup> )	0.49 m <sup>2</sup> (5.3 ft <sup>2</sup> )	0.98 m <sup>2</sup> ( 10.6 ft <sup>2</sup> )	1.48 m <sup>2</sup> (15.9 ft <sup>2</sup> )
<b>Materials of Construction</b>						
Filter Membrane	Hydrophilic polyethersulfone (PES)					
Film Edge	Polypropylene					
Supports	Polypropylene					
Core	Polysulfone					
Housing and Cage	Polypropylene					
O-Rings	Silicone					
<b>Vent/Drain</b>	6 mm (¼ in.) hose barb with double O-ring seal					
<b>Maximum Inlet Pressure</b>	6.9 bar (100 psi) intermittent at 25 °C 5.5 bar (80 psi) at 25 °C 2.8 bar (40 psi) at 60 °C 1 bar (15 psi) at 80 °C					
<b>Maximum Differential Pressure</b>						
Forward:	6.9 bar (100 psi) intermittent at 25 °C 5.5 bar (80 psi) at 25 °C 1 bar (15 psi) at 80 °C					
Reverse:	2.1 bar (30 psi) intermittent at 25 °C					
<b>Bubble Point</b>	≥ 4000 mbar (58 psi) air with water ≥ 1280 mbar (18.5 psi) nitrogen with 70/30% IPA/water mixture					
<b>Air Diffusion</b>	Through a water wet membrane at 2800 mbar (40 psi). ≤ 7.3 cc/min.      ≤ 13.3 cc/min.      ≤ 28.2 cc/min.      ≤ 28.2 cc/min.      ≤ 56.3 cc/min.      ≤ 84.5 cc/min.					
<b>Bacterial Retention</b>	Quantitative retention of 10 <sup>7</sup> CFU/cm <sup>2</sup> <i>Brevundimonas diminuta</i> ATCC® 19146 per ASTM F838 methodology.					
<b>Bacterial Endotoxin</b>	Aqueous extraction contains < 0.25 EU/mL as determined by the Limulus Ameboyte Lysate (LAL) Test as described in USP <85>.					
<b>Total Organic Carbon (TOC)/ Conductivity</b>	Autoclaved filter effluent meets the WFI requirement of USP <643> for Total Organic Carbon, and USP <645> for Water Conductivity at 25 °C after a WFI flush of: 5.5 L                      9.5 L                      20 L                      20 L                      40 L                      60 L					
<b>Sterilization</b>	May be autoclaved for 3 cycles for 60 minutes at 126 °C. Cannot be steam sterilized in-line.					
<b>Toxicity</b>	Component materials meet the criteria for Biological Reactivity Testing. These tests can be any one or a combination of the following test methods USP<88> Class VI ( <i>in vivo</i> ), USP <87> ( <i>in vitro</i> ), ISO 10993-5 ( <i>in vitro</i> )					
<b>Particle Shedding</b>	Effluent meets the acceptance criteria set forth in USP <788> for large volume parenterals.					
<b>Non-fiber Releasing</b>	Meets the criteria for a 'non-fiber releasing' filter as defined in 21 CFR 210.3 (b)(6).					
<b>Indirect Food Additive</b>	All component materials meet the FDA Indirect food additives requirement cited in 21 CFR 177-182 based on information provided by raw material suppliers.					


# Millipore Express® SHC Filters

## Typical Clean Water Flow Rates – Cartridge Filters

Standard Area Cartridge Filters with 0.5/0.2 µm Millipore Express® SHC Membranes




High Area Cartridge Filters with 0.5/0.2 µm Millipore Express® SHC Membranes




# Millipore Express® SHC Filters

## Typical Clean Water Flow Rates – Opticap® XL and XLT Autoclavable Capsules


Opticap® XL 3 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane




Opticap® XL 5 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane




Opticap® XL 10 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane




Opticap® XLT 10 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane



Opticap® XLT 20 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane



Opticap® XLT 30 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane



### Opticap® XL Capsule Legends Refer to Connection Type

- TT = 38 mm (1½ in.) sanitary flange inlet and outlet
- HH = 14 mm (9/16 in.) hose barb inlet and outlet
- TH = 38 mm (1½ in.) sanitary flange inlet and 14 mm (9/16 in.) hose barb outlet
- TB = 38 mm (1½ in.) sanitary flange inlet and 25 mm (1 in.) hose barb outlet

### Opticap® XLT Capsule Legends Refer to Connection Type

- TT = 38 mm (1½ in.) sanitary flange inlet and outlet
- TH = 38 mm (1½ in.) sanitary flange Inlet and 16 mm (5/8 in.) hose barb outlet
- HH = 16 mm (5/8 in.) hose barb inlet and outlet
- TB = 38 mm (1½ in.) sanitary flange inlet and 25 mm (1 in.) hose barb outlet
- BB = 25 mm (1 in.) hose barb inlet and outlet



## Millipore Express® SHC Filters


### Opticap® XL 150, 300 and 600 Sterile, Gamma and X-ray Compatible Capsule Filter Specifications (cont.)

Description	Opticap® XL 150 Standard Area Capsules	Opticap® XL 300 Standard Area Capsules	Opticap® XL 600 Standard Area Capsules
<b>Non-fiber Releasing</b>	Meets the criteria for a 'non-fiber releasing' filter as defined in 21 CFR 210.3 (b)(6).		
<b>Indirect Food Additive</b>	All component materials meet the FDA Indirect food additives requirement cited in 21 CFR 177-182 based on information provided by raw material suppliers.		


# Millipore Express® SHC Filters

## Typical Clean Water Flow Rates – Opticap® XL 150, 300, 600 Sterile, Gamma and X-ray Compatible Capsules


Opticap® XL 150 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane



Opticap® XL 300 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane



Opticap® XL 600 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane



### Opticap® XL 150, 300 and 600 Capsule Connection Type

- FF:** 19 mm (3/4 in.) sanitary flange inlet and outlet
- FH:** 19 mm (3/4 in.) sanitary flange inlet and 14 mm (9/16 in.) hose barb outlet
- HH:** 14 mm (9/16 in.) hose barb inlet and outlet

# Millipore Express® SHC Filters

## Opticap® XL and XLT Sterile, Gamma and X-ray Compatible Capsule Filter Specifications

Description	Opticap® XL 3 Standard Area Capsules	Opticap® XL 5 Standard Area Capsules	Opticap® XL 10 Standard Area Capsules	Standard Area Opticap® XLT Capsule	Per High Area Opticap® XLT Capsule
<b>Dimensions</b>					
Body Diameter	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)
Maximum Width	-	-	-	19.8 cm (7.8 in)	19.8 cm (7.8 in)
Maximum Length	17.3 cm (6.8 in.)	21.6 cm (8.5 in.)	33.5 cm (13.2 in.)		
10 inch				38.1 cm (15.0 in.)	38.1 cm (15.0 in.)
20 inch				62.5 cm (24.6 in.)	62.5 cm (24.6 in.)
30 inch				87.1 cm (34.3 in.)	87.1 cm (34.3 in.)
<b>Filtration Area</b>	0.13 m <sup>2</sup> (1.4 ft <sup>2</sup> )	0.24 m <sup>2</sup> (2.6 ft <sup>2</sup> )	0.54 m <sup>2</sup> (5.8 ft <sup>2</sup> )	0.54 m <sup>2</sup> (5.8 ft <sup>2</sup> )	1.0 m <sup>2</sup> ( 10.8 ft <sup>2</sup> ) per 10-inch
<b>Materials of Construction</b>					
Filter Membrane	Hydrophilic polyethersulfone (PES)				Hydrophilic PES Polyethylene Polyester Polyethersulfone Gamma and X-ray compatible Polypropylene Silicone
Film Edge	Polyethylene				
Supports	Polyester				
Core	Polysulfone				
Housing and Cage	Gamma and X-ray compatible Polypropylene				
O-Rings	Silicone				
<b>Vent/Drain</b>	6 mm (¼ in.) hose barb with double O-ring seal				
<b>Maximum Inlet Pressure</b>	6.9 bar (100 psi) intermittent at 25 °C 5.5 bar (80 psi) at 25 °C 2.8 bar (40 psi) at 60 °C 1 bar (15 psi) at 80 °C				
<b>Maximum Differential Pressure</b>					
Forward:	6.9 bar (100 psi) intermittent at 25 °C 5.5 bar (80 psi) at 25 °C 1 bar (15 psi) at 80 °C				
Reverse:	2.1 bar (30 psi) intermittent at 25 °C				
<b>Bubble Point</b>	≥ 4000 mbar (58 psi) air with water ≥ 1280 mbar (18.5 psi) nitrogen with 70/30% IPA/water mixture				
<b>Air Diffusion</b>	Through a water wet membrane at 2800 mbar (40 psi): ≤ 7.6 cc/min.                      ≤ 14.0 cc/min.                      ≤ 31.2 cc/min.                      ≤ 31.2 cc/min.                      ≤ 56.4 cc/min.				
<b>Bacterial Retention</b>	Quantitative retention of 10 <sup>7</sup> CFU/cm <sup>2</sup> <i>Brevundimonas diminuta</i> ATCC® 19146 per ASTM F838 methodology.				
<b>Bacterial Endotoxin</b>	Aqueous extraction contains < 0.25 EU/mL as determined by the Limulus Ameboyte Lysate (LAL) Test as described in USP <85>.				
<b>Total Organic Carbon (TOC)/ Conductivity</b>	Sterilized filter effluent meets the WFI requirement of USP <643> for Total Organic Carbon, and for USP <645> for Water Conductivity at 25 °C after a WFI flush of: 5.0 L                      9.5 L                      21 L                      21 L                      21 L				
<b>Sterilization</b>					
Gamma and X-ray Compatible Capsules	Compatible to 45 kGy and may be autoclaved for 3 cycles of 60 minutes at 123 °C. Cannot be steam sterilized in-line.				
Sterile Capsules	May be autoclaved for 3 cycles of 60 minutes at 123 °C. Cannot be steam sterilized in-line.				

## Millipore Express® SHC Filters

### Opticap® XL and XLT Sterile, Gamma and X-ray Compatible Capsule Filter Specifications (cont.)


Description	Opticap® XL 3 Standard Area Capsules	Opticap® XL 5 Standard Area Capsules	Opticap® XL 10 Standard Area Capsules	Standard Area Opticap® XLT Capsule	Per High Area Opticap® XLT Capsule
<b>Sterility (Sterile Capsules)</b>	These capsules meet current USP and AAMI guidelines for sterility utilizing a validated sterilization cycle.				
<b>Toxicity</b>	Component materials meet the criteria for Biological Reactivity Testing. These tests can be any one or a combination of the following test methods USP<88> Class VI ( <i>in vivo</i> ), USP <87> ( <i>in vitro</i> ), ISO 10993-5 ( <i>in vitro</i> ).				
<b>Particle Shedding</b>	Effluent meets the acceptance criteria set forth in USP <788> for large volume parenterals.				
<b>Non-fiber Releasing</b>	Meets the criteria for a 'non-fiber releasing' filter as defined in 21 CFR 210.3 (b)(6).				
<b>Indirect Food Additive</b>	All component materials meet the FDA Indirect food additives requirement cited in 21 CFR 177-182 based on information provided by raw material suppliers.				

\* Only available in gamma and X-ray compatible capsules


# Millipore Express® SHC Filters

## Typical Clean Water Flow Rates – Opticap® XL and XLT Sterile, Gamma and X-ray Compatible Capsules


**Opticap® XL 3 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane**



**Opticap® XL 5 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane**



**Opticap® XL 10 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane**



**Opticap® XL Capsule Legends Refer to Connection Type**

FF = 19 mm (¾ in.) sanitary flange inlet and outlet

FF = 19 mm (¾ in.) sanitary flange inlet and 14 mm (9/16 in.) hose barb outlet


TT = 38 mm (1 ½ in.) sanitary flange inlet and outlet

HH = 14 mm (9/16 in.) hose barb inlet and outlet


TH = 38 mm (1 ½ in.) sanitary flange inlet and 14 mm (9/16 in.) hose barb outlet

TB = 38 mm (1 ½ in.) sanitary flange inlet and 25 mm (1 in.) hose barb outlet


**Standard Area Opticap® XLT 10 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane**



**Standard Area Opticap® XLT 20 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane**



**Standard Area Opticap® XLT 30 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane**



**Opticap® XLT Capsule Legends Refer to Connection Type**

TT = 38 mm (1 ½ in.) sanitary flange inlet and outlet


TH = 38 mm (1 ½ in.) sanitary flange inlet and 16 mm (5/8 in.) hose barb outlet

HH = 16 mm (5/8 in.) hose barb inlet and outlet


TB = 38 mm (1 ½ in.) sanitary flange inlet and 25 mm (1 in.) hose barb outlet

BB = 25 mm (1 in.) hose barb inlet and outlet


**High Area Opticap® XLT 10 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane**



**High Area Opticap® XLT 20 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane**




**High Area Opticap® XLT 30 Capsules with 0.5/0.2 µm Millipore Express® SHC Membrane**




# Millipore Express® SHC Filters

## Ordering Information

### OptiScale® Capsules




### Cartridge Filters




# Millipore Express® SHC Filters

## Ordering Information

### Opticap® XL 150, 300, 600 Capsule Filters




### Opticap® XL Capsule Filters



# Millipore Express® SHC Filters

## Ordering Information

### Opticap® XLT Capsule Filters



# Millipore Express® SHF Filters

High flow, sterilizing-grade filters for critical process steps

Filters containing Millipore Express® SHF (Sterile, High Flux) sterilizing-grade membrane provide exceptionally high flow rates for low plugging streams. These filters contain a 0.2 µm polyethersulfone (PES) membrane that provides broad chemical compatibility for critical process steps that require validated sterility assurance.



## Benefits

- Sterilizing-grade membranes that are easy to wet and integrity test
- High flux PES membrane that provides superior processing efficiency
- Broad chemical compatibility across a wide pH range
- 100% integrity tested during the manufacturing process


## Filter Formats

- OptiScale® capsules
- Cartridge filters
- Opticap® XL and XLT capsule filters: sterile, gamma and X-ray compatible or autoclavable

## Millipore Express® SHF Filters

Millipore Express® SHF filters provide faster flow rates allowing you to process equivalent volumes with reduced filtration area, thus delivering savings in filtration costs.

### Millipore Express® SHF Membrane Flux Relative to Competitors



### Mobius® Single-use Solutions

Millipore Express® SHF filters are part of the Mobius® library providing you with the flexibility to design single-use assemblies that meet your specific processing requirements.

For more information, please visit: [SigmaAldrich.com/single-use-assemblies](http://SigmaAldrich.com/single-use-assemblies)



### The Emprove® Program – Your Fast Track through Regulatory Challenges

Complementing our product portfolio, the Emprove® Program provides three types of dossiers to support different stages of development and manufacturing operations such as qualification, risk assessment and process optimization.

For more information, please visit: [SigmaAldrich.com/Emprove](http://SigmaAldrich.com/Emprove)

### Filter Sizing

Testing with Millipore Express® SHF filters across a range of commonly used buffers and cleaning solutions (pH 1-14) showed an average permeability of 450 liters per hour per psi per 10-inch cartridge. Filters were tested after irradiation at 25-37 kGy and autoclaved at 123° C for 60 minutes.



# Millipore Express® SHF Filters

## OptiScale® Capsule and Cartridge Filter Specifications

Description	OptiScale® 25 Capsules	OptiScale® 47 Capsules	5-inch Cartridge	Per 10-inch Cartridge
<b>Dimensions</b>				
Diameter:	31 mm (1.21 in.)	70 mm (2.75 in.)	6.9 cm (2.7 in.)	6.9 cm (2.7 in.)
Maximum Length:	39 mm (1.52 in.)	94 mm (3.70 in.)	12.7 cm (5 in.)	25.4 cm (10 in.)
<b>Filtration Area</b>	3.5 cm <sup>2</sup>	17.7 cm <sup>2</sup> *	0.29 m <sup>2</sup> (3.1 ft <sup>2</sup> )	0.54 m <sup>2</sup> (5.8 ft <sup>2</sup> )
<b>Materials of Construction</b>				
Filter Membrane	Hydrophilic polyethersulfone (PES)	-	-	-
Film Edge	-	-	Polypropylene	-
Supports	Polypropylene	-	Polypropylene	-
Structural Components	Polypropylene	Polycarbonate	Polypropylene	-
Core	-	-	Polysulfone	-
Vent cap	Polypropylene	Polyvinylidene fluoride	-	-
O-Rings	-	Fluoroelastomer	Silicone, EPDM, or Fluoroelastomer	-
<b>Maximum Inlet Pressure</b>	4.1 bar (60 psi) at 25 °C	5.5 bar (80 psi) at 25 °C	-	-
<b>Maximum Differential Pressure</b>				
Forward:	4.1 bar (60 psi) at 25 °C	5.5 bar (80 psi) at 25 °C	6.9 bar (100 psi) at 25 °C 1.7 bar (25 psi) at 80 °C 340 mbar (5 psi) at 135 °C	-
Reverse:	0 bar (0 psi)	690 mbar (10 psi) at 25 °C	2.1 bar (30 psi) at 25 °C 69 mbar (1 psi) at 135 °C	-
<b>Bubble Point</b>	-	-	≥ 4000 mbar (58 psi) air with water ≥ 1280 mbar (18.5 psi) nitrogen with 70/30% IPA/ water mixture	-
<b>Air Diffusion</b>				
	-	-	Through a water wet membrane at 2800 mbar (40 psi):	-
	-	-	≤ 16.4 cc/min.	≤ 30 cc/min.
<b>Bacterial Retention</b>	-	-	Quantitative retention of 10 <sup>7</sup> CFU/cm <sup>2</sup> <i>Brevundimonas diminuta</i> ATCC® 19146 per ASTM F838 methodology.	-
<b>Bacterial Endotoxin</b>	Aqueous extraction contains < 0.25 EU/mL as determined by the Limulus Ameboyte Lysate (LAL) Test as described in USP <85>. Specification not applicable to OptiScale® 47 capsules.			
<b>Total Organic Carbon (TOC)/ Conductivity</b>	Autoclaved filter effluent meets the WFI requirement of USP <643>, for Total Organic Carbon and USP <645> for Water Conductivity at 25 °C after a WFI flush of: 15 mL	-	5.5 L	10 L
<b>Oxidizable Substances</b>	-	Meets the USP Oxidizable Substances Test requirements for sterile purified water after a water flush of: 100 mL	-	-
<b>Sterilization</b>				
Autoclave	1 cycle at 123 °C for 60 min.	3 cycles at 126 °C for 60 min.	15 cycles of 60 minutes at 126 °C	
In-line Steam	-	-	Forward for 25x , 30 min cycles at 135 °C, or 22x (forward) and 3x (reverse), 30 min cycles at 135 °C	
<b>Toxicity</b>	Component materials meet the criteria for Biological Reactivity Testing. These tests can be any one or a combination of the following test methods USP<88> Class VI ( <i>in vivo</i> ), USP <87> ( <i>in vitro</i> ), ISO 10993-5 ( <i>in vitro</i> ).			

## Millipore Express® SHF Filters

### OptiScale® Capsule and Cartridge Filter Specifications (cont.)

Description	OptiScale® 25 Capsules	OptiScale® 47 Capsules	5-inch Cartridge	Per 10-inch Cartridge
<b>Particle Shedding</b>	Effluent meets the acceptance criteria set forth in USP <788> for large volume parenterals. Specification does not apply to OptiScale® 47 capsules.			
<b>Non-fiber Releasing</b>	Meets the criteria for a 'non-fiber releasing' filter as defined in 21 CFR 210.3 (b) (6).			
<b>Indirect Food Additive</b>	All component materials meet the FDA Indirect food additives requirement cited in 21 CFR 177-182 based on information provided by raw material suppliers.			

\*OptiScale® 47 is recommended for screening purposes only, please refer to OptiScale® 25 capsule for scaling

# Millipore Express® SHF Filters


## Opticap® XL and XLT Autoclavable Capsule Filter Specifications

Description	Opticap® XL 3 Capsules	Opticap® XL 5 Capsules	Opticap® XL 10 Capsules	Opticap® XLT 10 Capsules	Opticap® XLT 20 Capsules	Opticap® XLT 30 Capsules
<b>Dimensions</b>						
Body Diameter	10.7 (4.2 in)	10.7 (4.2 in)	10.7 (4.2 in)	10.7 (4.2 in)	10.7 (4.2 in)	10.7 (4.2 in)
Maximum Width	-	-	-	19.8 cm (7.8 in)	19.8 cm (7.8 in)	19.8 cm (7.8 in)
Maximum Length	17.3 cm (6.8 in.)	21.6 cm (8.5 in.)	33.5 cm (13.2 in.)	38.1 cm (15.0 in.)	62.5 cm (24.6 in.)	87.1 cm (34.3 in.)
<b>Filtration Area</b>	0.16 m <sup>2</sup> (1.7 ft <sup>2</sup> )	0.29 m <sup>2</sup> (3.1 ft <sup>2</sup> )	0.54 m <sup>2</sup> (5.8 ft <sup>2</sup> )	0.54 m <sup>2</sup> (5.8 ft <sup>2</sup> )	1.08 m <sup>2</sup> (11.6 ft <sup>2</sup> )	1.62 m <sup>2</sup> (17.4 ft <sup>2</sup> )
<b>Materials of Construction</b>	Hydrophilic polyethersulfone (PES)					
Filter Membrane	Hydrophilic polyethersulfone (PES)					
Film Edge	Polypropylene					
Supports	Polypropylene					
Core	Polysulfone					
Housing and Cage	Polypropylene					
O-Rings	Silicone					
<b>Vent/Drain</b>	6 mm (¼ in.) hose barb with double O-ring seal					
<b>Maximum Inlet Pressure</b>	6.9 bar (100 psi) intermittent 25 °C					
<b>Maximum Differential Pressure</b>						
Forward:	5.5 bar (80 psi) at 25 °C 6.9 bar (100 psi) intermittent 25 °C 1 bar (15 psi) at 80 °C					
Reverse:	2.1 bar (30 psi) intermittent at 25 °C					
<b>Bubble Point</b>	≥ 4000 mbar (58 psi) air with water ≥ 1280 mbar (18.5 psi) nitrogen with 70/30% IPA/water mixture					
<b>Air Diffusion</b>	Through a water wet membrane at 2800 mbar (40 psi) at: ≤ 9.1 cc/min.    ≤ 16.4 cc/min.    ≤ 30 cc/min.    ≤ 30 cc/min.    ≤ 60 cc/min.    ≤ 90 cc/min.					
<b>Bacterial Retention</b>	Quantitative retention of 10 <sup>7</sup> CFU/cm <sup>2</sup> <i>Brevundimonas diminuta</i> ATCC® 19146 per ASTM F838 methodology.					
<b>Bacterial Endotoxin</b>	Aqueous extraction contains < 0.25 EU/mL as determined by the Limulus Ameboyte Lysate (LAL) Test as described in USP <85>.					
<b>Total Organic Carbon (TOC)/ Conductivity</b>	Autoclaved filter effluent meets the WFI requirement of USP <643> for Total Organic Carbon, and USP <645> for Water Conductivity at 25 °C after a WFI flush of: 3.0 L                      5.5 L                      10 L                      10 L                      20 L                      30 L					
<b>Sterilization</b>	May be autoclaved for 3 cycles for 60 minutes at 126 °C. Cannot be steam sterilized in-line.					
<b>Non-fiber Releasing</b>	Meets the criteria for a 'non-fiber releasing' filter as defined in 21 CFR 210.3 (b) (6).					
<b>Toxicity</b>	Component materials meet the criteria for Biological Reactivity Testing. These tests can be any one or a combination of the following test methods USP<88> Class VI ( <i>in vivo</i> ), USP <87> ( <i>in vitro</i> ), ISO 10993-5 ( <i>in vitro</i> )					
<b>Indirect Food Additive</b>	All component materials meet the FDA Indirect food additives requirement cited in 21 CFR 177-182 based on information provided by raw material suppliers.					

# Millipore Express® SHF Filters


## Typical Clean Water Flow Rates – Cartridge Filters

Cartridge Filters with Millipore Express® SHF Hydrophilic Membrane




## Typical Clean Water Flow Rates – Opticap® XL and XLT Autoclavable Capsules


Opticap® XL 3 Capsule Filters with 0.2 µm Millipore Express® SHF Membrane



Opticap® XL 5 Capsule Filters with 0.2 µm Millipore Express® SHF Membrane




Opticap® XL 10 Capsule Filters with 0.2 µm Millipore Express® SHF Membrane




Opticap® XL Capsule Legends Refer to Connection Type

- TT = 38 mm (1 1/2 in.) Sanitary Flange Inlet and Outlet
- HH = 14 mm (9/16 in.) Hose Barb Inlet and Outlet
- TH = 38 mm (1 1/2 in.) Sanitary Flange Inlet and 14 mm (9/16 in.) Hose Barb Outlet
- TB = 38 mm (1 1/2 in.) Sanitary Flange Inlet and 25 mm (1 in.) Hose Barb Outlet


Opticap® XLT 10 Capsule Filters with 0.2 µm Millipore Express® SHF Membrane



Opticap® XLT 20 Capsule Filters with 0.2 µm Millipore Express® SHF Membrane



Opticap® XLT 30 Capsule Filters with 0.2 µm Millipore Express® SHF Membrane



Opticap® XLT Capsule Legends Refer to Connection Type

- TT = 38 mm (1 1/2 in.) Sanitary Flange Inlet and Outlet
- TH = 38 mm (1 1/2 in.) Sanitary Flange Inlet and 16 mm (5/8 in.) Hose Barb Outlet
- HH = 16 mm (5/8 in.) Hose Barb Inlet and Outlet
- BB = 25 mm (1 in.) Hose Barb Inlet and Outlet
- TB = 38 mm (1 1/2 in.) Sanitary Flange Inlet and 25 mm (1 in.) Hose Barb Outlet



## Millipore Express® SHF Filters

### Opticap® XL 150, 300 and 600 Sterile, Gamma and X-ray Compatible Capsule Filter Specifications (cont.)


Description	Opticap® XL 150 Capsules	Opticap® XL 300 Capsules	Opticap® XL 600 Capsules
<b>Particle Shedding</b>	Effluent meets the acceptance criteria set forth in USP <788> for large volume parenterals.		
<b>Non-fiber Releasing</b>	Meets the criteria for a 'non-fiber releasing' filter as defined in 21 CFR 210.3 (b)(6).		
<b>Indirect Food Additive</b>	All component materials meet the FDA Indirect food additives requirement cited in 21 CFR 177-182 based on information provided by raw material suppliers.		

# Millipore Express® SHF Filters


## Typical Clean Water Flow Rates – Opticap® XL 150, 300, 600 Sterile, Gamma and X-ray Compatible Capsules

Filters tested post irradiation at 45–65 kGy and autoclaved at 123 °C for 60 minutes


### Opticap® XL 150 Capsule Filters with 0.2 µm Millipore Express® SHF Membrane



### Opticap® XL 300 Capsule Filters with 0.2 µm Millipore Express® SHF Membrane



### Opticap® XL 600 Capsule Filters with 0.2 µm Millipore Express® SHF Membrane



### Opticap® XL 150, 300 & 600 Capsule Connection Type

FF = 19 mm (¾ in.) sanitary flange inlet and outlet

FH = 19 mm (¾ in.) sanitary flange inlet and  
14mm (9/16 in.) hose barb outlet

HH = 14 mm (9/16 in.) hose barb inlet and outlet

## Millipore Express® SHF Filters

### Opticap® XL and XLT Sterile, Gamma and X-ray Compatible Capsule Filter Specifications


Description	Opticap® XL 3 Capsules	Opticap® XL 5 Capsules	Opticap® XL 10 Capsules	Opticap® XLT 10 Capsules	Opticap® XLT 20 Capsules	Opticap® XLT 30 Capsules
<b>Dimensions</b>						
Body Diameter	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)
Maximum Width	-	-	-	19.8 cm (7.8 in.)	19.8 cm (7.8 in.)	19.8 cm (7.8 in.)
Maximum Length	17.3 cm (6.8 in.)	21.6 cm (8.5 in.)	33.5 cm (13.2 in.)	38.1 cm (15.0 in.)	62.5 cm (24.6 in.)	87.1 cm (34.3 in.)
<b>Filtration Area</b>	0.17 m <sup>2</sup> (1.8 ft <sup>2</sup> )	0.31 m <sup>2</sup> (3.3 ft <sup>2</sup> )	0.57 m <sup>2</sup> (6.1 ft <sup>2</sup> )	0.57 m <sup>2</sup> (6.1 ft <sup>2</sup> )	1.14 m <sup>2</sup> (12.3 ft <sup>2</sup> )	1.71 m <sup>2</sup> (18.4 ft <sup>2</sup> )
<b>Materials of Construction</b>						
Filter Membrane	Hydrophilic polyethersulfone (PES)					
Film Edge	Polyethylene					
Supports	Polyester					
Core	Polysulfone					
Housing and Cage	Gamma and X-ray compatible polypropylene					
O-Rings	Silicone					
<b>Vent/Drain</b>	6 mm (¼ in.) hose barb with double O-ring seal					
<b>Maximum Differential Pressure</b>						
Forward:	5.5 bar (80 psi) at 25 °C 6.9 bar (100 psi) intermittent 25 °C 1 bar (15 psi) at 80 °C					
Reverse:	2.1 bar (30 psi) intermittent at 25 °C					
<b>Maximum Inlet pressure</b>	6.9 bar (100psi) intermittent at 25 °C 5.5 bar (80 psi) at 25 °C 2.8 bar (40 psi) at 60 °C 1 bar (15 psi) at 80 °C					
<b>Bubble Point</b>	≥ 4000 mbar (58 psi) air with water ≥ 1280 mbar (18.5 psi) nitrogen with 70/30% IPA/water mixture					
<b>Air Diffusion</b>	Through a water wet membrane at 2800 mbar (40 psi): ≤ 9.5 cc/min.      ≤ 17.4 cc/min.      ≤ 32.7 cc/min.      ≤ 32.7 cc/min.      ≤ 65.5 cc/min.      ≤ 98.2 cc/min.					
<b>Bacterial Retention</b>	Quantitative retention of 10 <sup>7</sup> CFU/cm <sup>2</sup> <i>Brevundimonas diminuta</i> ATCC® 19146 per ASTM F838 methodology.					
<b>Bacterial Endotoxin</b>	Aqueous extraction contains < 0.25 EU/mL as determined by the Limulus Ameboyte Lysate (LAL) Test as described in USP <85>.					
<b>Total Organic Carbon (TOC)/ Conductivity</b>	Sterilized filter effluent meets the WFI requirement of USP <643>, for Total Organic Carbon and for USP <645> for Water Conductivity at 25 °C after a WFI flush of: 3.5 L                  6.0 L                  11 L                  11 L                  22 L                  33 L					
<b>Sterilization</b>						
Gamma and X-ray compatible Capsules	Compatible to 45 kGy and may be autoclaved for 3 cycles of 60 minutes at 123 °C. Cannot be steam sterilized in-line.					
Sterile Capsules	May be autoclaved for 3 cycles of 60 minutes at 123 °C. Cannot be steam sterilized in-line.					
<b>Sterility (Sterile Capsules)</b>	These capsules meet current USP and AAMI guidelines for sterility utilizing a validated sterilization cycle.					
<b>Toxicity</b>	Component materials meet the criteria for Biological Reactivity Testing. These tests can be any one or a combination of the following test methods USP<88> Class VI ( <i>in vivo</i> ), USP <87> ( <i>in vitro</i> ), ISO 10993-5 ( <i>in vitro</i> ).					
<b>Non-fiber Releasing</b>	Meets the criteria for a 'non-fiber releasing' filter as defined in 21 CFR 210.3 (b)(6).					
<b>Indirect Food Additive</b>	All component materials meet the FDA Indirect food additives requirement cited in 21 CFR 177-182 based on information provided by raw material suppliers.					

# Millipore Express® SHF Filters


## Typical Clean Water Flow Rates – Opticap® XL and XLT Sterile, Gamma and X-ray Compatible Capsules

Filters tested post irradiation at 25–45 kGy


**Opticap® XL 3 Capsule Filters with 0.2 µm Millipore Express® SHF Membrane**




**Opticap® XL 5 Capsule Filters with 0.2 µm Millipore Express® SHF Membrane**




**Opticap® XL 10 Capsule Filters with 0.2 µm Millipore Express® SHF Membrane**




**Opticap® XLT 10 Capsule Filters with 0.2 µm Millipore Express® SHF Membrane**



**Opticap® XLT 20 Capsule Filters with 0.2 µm Millipore Express® SHF Membrane**




**Opticap® XLT 30 Capsule Filters with 0.2 µm Millipore Express® SHF Membrane**




# Millipore Express® SHF Filters

## Ordering Information

### OptiScale® Capsules




### Cartridge Filters




# Millipore Express® SHF Filters

## Ordering Information

### Opticap® XL 150, 300, 600 Capsule Filters




### Opticap® XL Capsule Filters



# Millipore Express® SHF Filters

## Ordering Information

### Opticap® XLT Capsule Filters



# Millipore Express® PHF Filters

High flow filters for cost-effective filtration of buffers

Filters containing Millipore Express® PHF (Process protection, High Flux) membrane are a cost-effective option for sterile filtration of buffers. These filters contain a 0.2 µm polyethersulfone (PES) membrane and provide high flow rates and extended throughput for superior process efficiency.

## Benefits


- High flux sterilizing-grade PES membrane that reduces process footprint
- Broad chemical compatibility across a wide pH range
- Fully scalable product offerings

## Available Pore Sizes

- OptiScale® capsules
- Cartridge filters
- Opticap® XL and XLT capsule filters: sterile, gamma and X-ray compatible or autoclavable



## Relative Water Permeability of 10-inch Cartridges



## Mobius® Single-use Solutions

Millipore Express® PHF filters are part of the Mobius® library providing you with the flexibility to design single-use assemblies that meet your specific processing requirements.

**For more information, please visit:**  
[SigmaAldrich.com/single-use-assemblies](https://www.sigmaaldrich.com/single-use-assemblies)

# Millipore Express® PHF Filters

## OptiScale® Capsule and Cartridge Filter Specifications

Description	OptiScale® 25 Capsules	5-inch Cartridge	Per Standard Area 10-inch Cartridge
<b>Dimensions</b>			
Diameter:	31 mm (1.21 in.)	6.9 cm (2.7 in.)	6.9 cm (2.7 in.)
Length:	39 mm (1.52 in.)	12.5 cm (5 in.)	25.4 cm (10 in.)
<b>Filtration Area</b>	3.5 cm <sup>2</sup>	0.29 m <sup>2</sup> (3.1 ft <sup>2</sup> )	0.54 m <sup>2</sup> (5.8 ft <sup>2</sup> )
<b>Materials of Construction</b>			
Filter Membrane	Hydrophilic polyethersulfone (PES)		
Film Edge	—	Polypropylene	
Supports	Polypropylene	Polypropylene	
Structural components	Polypropylene	Polypropylene	
Core	—	Polysulfone	
Vent Cap	Polypropylene	—	
O-Ring	—	Silicone, EPDM, or Fluoroelastomer	
<b>Maximum Inlet Pressure</b>	4.1 bar (60 psi) at 25 °C	5.5 bar (80 psi) at 25 °C	—
<b>Maximum Differential Pressure</b>			
Forward:	4.1 bar (60 psi) at 25 °C	5.5 bar (80 psi) at 25 °C 1.7 bar (25 psi) at 80 °C 300 mbar (5 psi) at 135 °C	
Reverse:	0 bar (0 psi)	1.4 bar (20 psi) at 25 °C 69 mbar (1 psi) at 135 °C	
<b>Air Diffusion</b>			
	—	Through a water wet membrane at 2800 mbar (40 psi).	
	—	≤ 16.4 cc/min.	≤ 30 cc/min.
<b>Bacterial Retention</b>			
	—	Quantitative retention of 10 <sup>7</sup> CFU/cm <sup>2</sup> <i>Brevundimonas diminuta</i> ATCC® 19146 per ASTM F838 methodology.	
<b>Bacterial Endotoxin</b>			
	Aqueous extraction contains < 0.25 EU/mL as determined by the Limulus Ameboyte Lysate (LAL) Test as described in USP <85>.		
<b>Total Organic Carbon (TOC)/ Conductivity</b>			
	Autoclaved filter effluent meets the WFI requirement of USP <643> for Total Organic Carbon, and USP <645> for Water Conductivity at 25 °C after a WFI water flush of: 15 mL	5.5 L	10 L
<b>Sterilization</b>			
Autoclave	1 cycle at 123 °C for 60 min.	15 cycles of 60 minutes at 126 °C	
In-line Steam	—	Forward for 15x, 30 min. cycles at 135 °C, or 12x (forward) and 3x (reverse), 30 min. cycles at 135 °C	
<b>Toxicity</b>			
	Component materials meet the criteria for Biological Reactivity Testing. These tests can be any one or a combination of the following test methods USP<88> Class VI ( <i>in vivo</i> ), USP <87> ( <i>in vitro</i> ), ISO 10993-5 ( <i>in vitro</i> ).		
<b>Particle Shedding</b>			
	Effluent meets the acceptance criteria set forth in USP <788> for large volume parenterals.		
<b>Non-fiber Releasing</b>			
	Meets the criteria for a 'non-fiber releasing' filter as defined in 21 CFR 210.3 (b)(6).		
<b>Indirect Food Additive</b>			
	All component materials meet the FDA Indirect food additives requirement cited in 21 CFR 177-182 based on information provided by raw material suppliers.		

# Millipore Express® PHF Filters


## Opticap® XL and XLT Autoclavable Capsule Filter Specifications

Description	Opticap® XL 3 Capsules	Opticap® XL 5 Capsules	Opticap® XL 10 Capsules	Opticap® XLT 10 Capsules	Opticap® XLT 20 Capsules	Opticap® XLT 30 Capsules
<b>Dimensions</b>						
Body Diameter	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)
Maximum Width	-	-	-	19.8 cm (7.8 in.)	19.8 cm (7.8 in.)	19.8 cm (7.8 in.)
Maximum Length	17.3 cm (6.8 in.)	21.6 cm (8.5 in.)	33.5 cm (13.2 in.)	38.1 cm (15.0 in.)	62.5 cm (24.6 in.)	87.1 cm (34.3 in.)
<b>Filtration Area</b>	0.16 m <sup>2</sup> (1.7 ft <sup>2</sup> )	0.29 m <sup>2</sup> (3.1 ft <sup>2</sup> )	0.54 m <sup>2</sup> (5.8 ft <sup>2</sup> )	0.54 m <sup>2</sup> (5.8 ft <sup>2</sup> )	1.08 m <sup>2</sup> (11.6 ft <sup>2</sup> )	1.62 m <sup>2</sup> (17.4 ft <sup>2</sup> )
<b>Materials of Construction</b>						
Filter Membrane	Hydrophilic polyethersulfone (PES)					
Film Edge	Polypropylene					
Supports	Polypropylene					
Core	Polysulfone					
Housing and Cage	Polypropylene					
O-rings	Silicone					
<b>Vent/Drain</b>	6 mm (¼ in.) hose barb with double O-ring seal					
<b>Maximum Inlet Pressure</b>	5.5 bar (80 psi) at 25 °C					
<b>Maximum Differential Pressure</b>						
Forward:	5.5 bar (80 psi) at 25 °C; 1.0 bar 15 psi at 80 °C					
Reverse:	1.4 bar (20 psi) at 25 °C					
<b>Air Diffusion</b>						
	Through a water wet membrane at 2800 mbar (40 psi):					
	≤ 9.1 cc/min.	≤ 16.4 cc/min.	≤ 30 cc/min.	≤ 30 cc/min.	≤ 60 cc/min.	≤ 90 cc/min.
<b>Bacterial Retention</b>	Quantitative retention of 10 <sup>7</sup> CFU/cm <sup>2</sup> <i>Brevundimonas diminuta</i> ATCC® 19146 per ASTM F838 methodology.					
<b>Bacterial Endotoxin</b>	Aqueous extraction contains < 0.25 EU/mL as determined by the Limulus Ameboyte Lysate (LAL) Test as described in USP <85>.					
<b>Total Organic Carbon (TOC)/ Conductivity</b>						
	Autoclaved filter effluent meets the WFI requirement of USP <643> for Total Organic Carbon, and USP <645> for Water Conductivity at 25 °C after a WFI flush of:					
	3 L	5.5 L	10 L	10 L	20 L	30 L
<b>Sterilization</b>	May be autoclaved for 3 cycles for 60 minutes at 126 °C. Cannot be steam sterilized in-line.					
<b>Toxicity</b>						
	Component materials meet the criteria for Biological Reactivity Testing. These tests can be any one or a combination of the following test methods USP<88> Class VI ( <i>in vivo</i> ), USP <87> ( <i>in vitro</i> ), ISO 10993-5 ( <i>in vitro</i> )					
<b>Non-fiber Releasing</b>	Meets the criteria for a 'non-fiber releasing' filter as defined in 21 CFR 210.3 (b)(6).					
<b>Indirect Food Additive</b>	All component materials meet the FDA indirect food additives requirement cited in 21 CFR 177-182 based on information provided by raw material suppliers.					

# Millipore Express® PHF Filters


## Typical Clean Water Flow Rates – Cartridge Filters

Cartridge Filters with 0.2 µm Millipore Express® PHF Hydrophilic Membrane




## Typical Clean Water Flow Rates – Opticap® XL and XLT Autoclavable Capsules


Opticap® XL 3 Capsule Filters with 0.2 µm Millipore Express® PHF Membrane



Opticap® XL 5 Capsule Filters with 0.2 µm Millipore Express® PHF Membrane




Opticap® XL 10 Capsule Filters with 0.2 µm Millipore Express® PHF Membrane




Opticap® XL Capsule Legends Refer to Connection Type

- TT = 38 mm (1 in.) sanitary flange inlet and outlet
- HH = 14mm (<sup>9</sup>/<sub>16</sub> in.) hose barb inlet and outlet
- TH = 38mm (1 ½ in.) sanitary flange inlet and 14mm (<sup>9</sup>/<sub>16</sub> in.) hose barb outlet
- TB = 38mm (1 ½ in.) sanitary flange inlet and 25mm (1 in.) hose barb outlet


Opticap® XLT 10 Capsule Filters with 0.2 µm Millipore Express® PHF Membrane



Opticap® XLT 20 Capsule Filters with 0.2 µm Millipore Express® PHF Membrane



Opticap® XLT 30 Capsule Filters with 0.2 µm Millipore Express® PHF Membrane



Opticap® XLT Capsule Legends Refer to Connection Type

- TT = 38 mm (1 ½ in.) sanitary flange inlet and outlet
- TH = 38mm (1 ½ in.) sanitary flange inlet and 16mm (<sup>5</sup>/<sub>8</sub> in.) hose barb outlet
- HH = 16mm (<sup>5</sup>/<sub>8</sub> in.) hose barb inlet and outlet
- BB = 25mm (1 in.) hose barb inlet and outlet
- TB = 38mm (1 ½ in.) sanitary flange inlet and 25mm (1 in.) hose barb outlet

# Millipore Express® PHF Filters

## Opticap® XL 150, 300 and 600 Sterile, Gamma and X-ray Compatible Capsule Filter Specifications


Description	Opticap® XL 150 Capsules	Opticap® XL 300 Capsules	Opticap® XL 600 Capsules
<b>Dimensions</b>			
Body Diameter	5.6 cm (2.2 in.)	5.6 cm (2.2 in.)	5.6 cm (2.2 in.)
Maximum Length	9.7 cm (3.8 in.)	11.9 cm (4.7 in.)	16.5 cm (6.5 in.)
<b>Filtration Area</b>	0.022 m <sup>2</sup> (0.240 ft <sup>2</sup> )	0.048 m <sup>2</sup> (0.514 ft <sup>2</sup> )	0.097 m <sup>2</sup> (1.046 ft <sup>2</sup> )
<b>Materials of Construction</b>			
Filter Membrane	Hydrophilic polyethersulfone (PES)		
Supports	Polyethylene		
Core	Polysulfone		
Housing and Cage	Gamma and X-ray compatible polypropylene		
O-Rings	Silicone		
<b>Vent/Drain</b>	6 mm (¼ in.) hose barb with double O-ring seal		
<b>Maximum Inlet Pressure</b>	5.5 bar (80 psi) at 25 °C		
<b>Maximum Differential Pressure</b>			
Forward:	5.5 bar (80 psi) at 25 °C; 1.0 bar (15 psi) at 80 °C		
Reverse:	1.4 bar (20 psi) at 25 °C		
<b>Air Diffusion</b>	Through a water wet membrane at 2800 mbar (40 psi):		
	≤ 1.4 cc/min.	≤ 2.8 cc/min.	≤ 5.8 cc/min.
<b>Bacterial Retention</b>	Quantitative retention of 10 <sup>7</sup> CFU/cm <sup>2</sup> <i>Brevundimonas diminuta</i> ATCC® 19146 per ASTM F838 methodology.		
<b>Bacterial Endotoxin</b>	Aqueous extraction contains < 0.25 EU/mL as determined by the Limulus Ameboyte Lysate (LAL) Test as described in USP <85>.		
<b>Total Organic Carbon (TOC)/ Conductivity</b>	Sterilized filter effluent meets the WFI requirement of USP <643> for Total Organic Carbon, and for USP <645> for Water Conductivity at 25 °C after a WFI flush of:		
	2.0 L	2.5 L	3.0 L
<b>Sterilization</b>			
Gamma and X-ray Compatible Capsules	Compatible to 40 kGy and may be autoclaved for 3 cycles of 60 minutes at 123 °C. Cannot be steam sterilized in-line.		
Sterile Capsules	May be autoclaved for 3 cycles of 60 minutes at 123 °C. Cannot be steam sterilized in-line.		
<b>Sterility (Sterile Capsules)</b>	These capsules meet current USP and AAMI guidelines for sterility utilizing a validated sterilization cycle.		
<b>Toxicity</b>	Component materials meet the criteria for Biological Reactivity Testing. These tests can be any one or a combination of the following test methods USP<88> Class VI ( <i>in vivo</i> ), USP <87> ( <i>in vitro</i> ), ISO 10993-5 ( <i>in vitro</i> )		
<b>Particle Shedding</b>	Effluent meets the acceptance criteria set forth in USP <788> for large volume parenterals.		
<b>Non-fiber Releasing</b>	Meets the criteria for a 'non-fiber releasing' filter as defined in 21 CFR 210.3 (b)(6).		
<b>Indirect Food Additive</b>	All component materials meet the FDA Indirect food additives requirement cited in 21 CFR 177-182 based on information provided by raw material suppliers.		

# Millipore Express® PHF Filters


## Typical Clean Water Flow Rates – Opticap® XL 150, 300, 600 Sterile, Gamma and X-ray Compatible Capsules

Filters tested post irradiation at 25–40 kGy and autoclaved at 123 °C for 60 minutes


### Opticap® XL 150 Capsule Filters with 0.2 µm Millipore Express® PHF Membrane



### Opticap® XL 300 Capsule Filters with 0.2 µm Millipore Express® PHF Membrane



### Opticap® XL 600 Capsule Filters with 0.2 µm Millipore Express® PHF Membrane



### Opticap® XL 150, 300 & 600 Capsule Connection Type

FF = 19 mm (¾ in.) sanitary flange inlet and outlet

FH = 19 mm (¾ in.) sanitary flange inlet and  
14mm (9/16 in.) hose barb outlet

HH = 14 mm (9/16 in.) hose barb inlet and outlet

## Millipore Express® PHF Filters

### Opticap® XL and XLT Sterile, Gamma and X-ray Compatible Capsule Filter Specifications


Description	Opticap® XL 3 Capsules	Opticap® XL 5 Capsules	Opticap® XL 10 Capsules	Opticap® XLT 10 Capsules	Opticap® XLT 20 Capsules	Opticap® XLT 30 Capsules
<b>Dimensions</b>						
Body Diameter	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)
Maximum Width	-	-	-	19.8 cm (7.8 in)	19.8 cm (7.8 in)	19.8 cm (7.8 in)
Maximum Length	17.3 cm (6.8 in.)	21.6 cm (8.5 in.)	33.5 cm (13.2 in.)	38.1 cm (15.0 in.)	62.5 cm (24.6 in.)	87.1 cm (34.3 in.)
<b>Filtration Area</b>	0.17 m <sup>2</sup> (1.8 ft <sup>2</sup> )	0.31 m <sup>2</sup> (3.3 ft <sup>2</sup> )	0.57 m <sup>2</sup> (6.1 ft <sup>2</sup> )	0.57 m <sup>2</sup> (6.1 ft <sup>2</sup> )	1.14 m <sup>2</sup> (12.3 ft <sup>2</sup> )	1.71 m <sup>2</sup> (18.4 ft <sup>2</sup> )
<b>Materials of Construction</b>						
Filter Membrane	Hydrophilic polyethersulfone (PES)					
Film Edge	Polyethylene					
Supports	Polyester					
Core	Polysulfone					
Housing and Cage	Gamma and X-ray compatible polypropylene					
O-Rings	Silicone					
<b>Vent/Drain</b>	6 mm (¼ in.) hose barb with double O-ring seal					
<b>Maximum Inlet Pressure</b>	5.5 bar (80 psi) at 25 °C					
<b>Maximum Differential Pressure</b>						
Forward:	5.5 bar (80 psi) at 25 °C; 1.0 bar (15 psi) at 80 °C					
<b>Air Diffusion</b>						
	Through a water wet membrane at 2800 mbar (40 psi):					
	≤ 9.5 cc/min.	≤ 17.4 cc/min.	≤ 32.7 cc/min.	≤ 32.7 cc/min.	≤ 65.5 cc/min.	≤ 98.2 cc/min.
<b>Bacterial Retention</b>						
	Quantitative retention of 10 <sup>7</sup> CFU/cm <sup>2</sup> <i>Brevundimonas diminuta</i> ATCC® 19146 per ASTM F838 methodology.					
<b>Bacterial Endotoxin</b>						
	Aqueous extraction contains < 0.25 EU/mL as determined by the Limulus Ameboyte Lysate (LAL) Test as described in USP <85>					
<b>Total Organic Carbon (TOC)/ Conductivity</b>						
	Sterilized filter effluent meets the WFI requirement of USP <643> for Total Organic Carbon and for USP <645> for Water Conductivity at 25 °C after a WFI flush of:					
	3.5 L	6.0 L	11 L	11 L	22 L	33 L
<b>Sterilization</b>						
Gamma and X-ray Compatible Capsules	Compatible to 40 kGy and may be autoclaved for 3 cycles of 60 minutes at 123 °C. Cannot be steam sterilized in-line.					
Sterile Capsules	May be autoclaved for 3 cycles of 60 minutes at 123 °C. Cannot be steam sterilized in-line.					
<b>Sterility (Sterile Capsules)</b>						
	These capsules meet current USP and AAMI guidelines for sterility utilizing a validated sterilization cycle.					
<b>Toxicity</b>						
	Component materials meet the criteria for Biological Reactivity Testing. These tests can be any one or a combination of the following test methods USP<88> Class VI ( <i>in vivo</i> ), USP <87> ( <i>in vitro</i> ), ISO 10993-5 ( <i>in vitro</i> ).					
<b>Non-fiber Releasing</b>						
	Meets the criteria for a 'non-fiber releasing' filter as defined in 21 CFR 210.3 (b)(6).					
<b>Indirect Food Additive</b>						
	All component materials meet the FDA Indirect food additives requirement cited in 21 CFR 177-182 based on information provided by raw material suppliers.					

# Millipore Express® PHF Filters


## Typical Clean Water Flow Rates – Opticap® XL and XLT Sterile, Gamma and X-ray Compatible Capsules

Filters tested post irradiation at 25–40 kGy


**Opticap® XL 3 Capsule Filters with 0.2 µm Millipore Express® PHF Membrane**




**Opticap® XL 5 Capsule Filters with 0.2 µm Millipore Express® PHF Membrane**




**Opticap® XL 10 Capsule Filters with 0.2 µm Millipore Express® PHF Membrane**




**Opticap® XLT 10 Capsule Filters with 0.2 µm Millipore Express® PHF Membrane**



**Opticap® XLT 20 Capsule Filters with 0.2 µm Millipore Express® PHF Membrane**



**Opticap® XLT 30 Capsule Filters with 0.2 µm Millipore Express® PHF Membrane**



### Opticap® XL Capsule Legends Refer to Connection Type

- TT = 38 mm (1 in.) sanitary flange inlet and outlet
- HH = 14mm (<sup>9</sup>/<sub>16</sub> in.) hose barb inlet and outlet
- TH = 38mm (1 ½ in.) sanitary flange inlet and 14mm (<sup>9</sup>/<sub>16</sub> in.) hose barb outlet
- TB = 38mm (1 ½ in.) sanitary flange inlet and 25mm (1 in.) hose barb outlet


### Opticap® XLT Capsule Legends Refer to Connection Type

- TT = 38 mm (1 ½ in.) sanitary flange inlet and outlet
- TH = 38mm (1 ½ in.) sanitary flange inlet and 16mm (<sup>5</sup>/<sub>8</sub> in.) hose barb outlet
- HH = 16mm (<sup>5</sup>/<sub>8</sub> in.) hose barb inlet and outlet
- BB = 25mm (1 in.) hose barb inlet and outlet
- TB = 38mm (1 ½ in.) sanitary flange inlet and 25mm (1 in.) hose barb outlet


# Millipore Express® PHF Filters

## Ordering Information

### OptiScale® Capsules




### Cartridge Filters




# Millipore Express® PHF Filters

## Ordering Information

### Opticap® XL 150, 300, 600 Capsule Filters




### Opticap® XL Capsule Filters



# Millipore Express® PHF Filters

## Ordering Information

### Opticap® XLT Capsule Filters




# Millipore®

Expert Pharm/BioPharm  
Products & CTDMO Services



For additional information, visit  
[SigmaAldrich.com](https://SigmaAldrich.com)

Merck KGaA  
Frankfurter Strasse 250  
64293 Darmstadt, Germany



© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, Millipore, Millipore Express, Lifegard, Milligard, Opticap, OptiScale, Mobius, Emprove, Sigma-Aldrich and PolySep are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

Lit. No. MK\_PG8135EN Ver. 4.0 64008 10/2025