

## Product Information

### CYTOCHROME P450 CYP4F3B ISOZYME

Human, Recombinant

Microsomes with Cytochrome P450 Reductase and  
Cytochrome b<sub>5</sub>

Product Number **C 5235**

Storage Temperature -70 °C

#### Product Description

The microsomal product is prepared from insect cells (BTI-TN-5B1-4) infected with recombinant baculovirus containing cDNA inserts for the human cytochrome P450 isozyme, human cytochrome P450 reductase, and human cytochrome b<sub>5</sub>.

Cytochrome P450 enzymes are a superfamily of heme containing monooxygenases that in humans are involved with oxidative metabolism of xenobiotics. This metabolism is the initial step in the biotransformation and elimination of a wide variety of drugs and environmental pollutants from the body. These reactions are achieved through a mixed monooxygenase system with the general EC number of 1.14.14.1.<sup>1</sup> The CYP4F3B isozyme is a hepatic enzyme, which metabolizes the  $\omega$ -oxidation of leukotriene B<sub>4</sub> to the 20-hydroxyl product. The K<sub>m</sub> values observed for this reaction by the CYP4F hepatic enzymes are several fold higher than the value for the non-hepatic CYP4F3A isozyme.

The cytochrome P450 enzymes range in molecular weight between 45 to 60 kDa.

The product is supplied as 0.5 nmole of cytochrome P450 isozyme in 0.5 ml of 100 mM potassium phosphate buffer, pH 7.4. Protein content, cytochrome b<sub>5</sub> content, cytochrome c reductase activity, and leukotriene B<sub>4</sub> 20-hydroxylase activity of the microsomes are reported on a lot-to-lot basis.

#### Precautions and Disclaimer

This product is for laboratory research use only. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

#### Preparation Instructions

1. Quickly thaw at 37 °C using a water bath. Keep on ice until ready to use.
2. If not using entire contents, aliquot to minimize freeze-thaw cycles. Generally, 80% or more of the catalytic activity is retained after 6 freeze-thaw cycles.
3. Store aliquots at -70 °C.

#### Storage/Stability

The product is shipped on dry ice and should be stored at -70 °C. The product as supplied is stable for at least 24 months. For prolonged storage, freeze in working aliquots at -70 °C. Avoid repeated freezing and thawing.

#### Procedure

##### Leukotriene B<sub>4</sub> 20-Hydroxylase Activity:

A 0.1 ml reaction containing the following was incubated at 37 °C for 10 minutes:

- 2 pmoles cytochrome P450 isozyme
- 1.3 mM NADP<sup>+</sup>
- 3.3 mM glucose-6-phosphate
- 0.4 U/ml glucose-6-phosphate dehydrogenase
- 3.3 mM magnesium chloride
- 30  $\mu$ M leukotriene B<sub>4</sub>
- 100 mM phosphate, pH 7.4

The reaction was stopped with 25  $\mu$ l of 94% acetonitrile/6% glacial acetic acid and centrifuged (10,000 x g) for 3 minutes. A 50  $\mu$ l aliquot of the supernatant was injected into a 4.6 x 250 mm 5  $\mu$ m C18 HPLC column and separated at 45 °C. A 20 minute gradient from 0-100% Phase B was used at a flow rate of 1.0 ml per minute. Mobile Phase A: 30% acetonitrile with 1 mM perchloric acid in water and Mobile Phase B: 70% methanol. The product was detected by measuring absorbance at 270 nm and comparison was made to a standard curve for 20-hydroxyleukotriene B<sub>4</sub> (Product No. H 4895).

#### Notes:

With respect to enzyme concentration, catalysis is linear up to at least 40 pmoles of cytochrome P450 isozyme per ml. Hydroxylation of leukotriene B<sub>4</sub> is approximately linear for 20 minutes. NADPH may be substituted for the NADPH generating system, which consists of NADP<sup>+</sup>, glucose-6-phosphate, and glucose-6-phosphate dehydrogenase.

## References

1. Enzyme Nomenclature, IUBMB, Academic Press (1992).
2. Anzenbacher, P., and Anzenbacherova E., Cytochromes P450 and metabolism of xenobiotics. Cell Mol. Life Sci., **58**, 737-47 (2001).
3. Guengrich, F.P., Cytochrome P450: Structure, Mechanism and Biochemistry (2<sup>nd</sup> Edition), Chapter
14. Ortiz de Montellano, P.R. (ed.) Plenum Press (New York, NY: 1995).
4. Christmas, P., et al., Expression of the CYP4F3 gene. tissue-specific splicing and alternative promoters generate high and low K(m) forms of leukotriene B(4) omega-hydroxylase, J. Biol. Chem., **274**, 21191-21199 (1999).

JC/MAM 12/01

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.