Instructions for Use

Scott's Tap Water Substitute Concentrate

Procedure No. S5134

Intended Use

The data obtained from this manual, qualitative procedure does not detect or measure any analyte or marker in fresh or fixed histologic and cytologic human specimens. Scott's tap water is a component of Hematoxylin and Eosin staining. Its function is to convert chromatin coloration from a red/purple hue to a blue/black hue^{1,2}. The data provided from this staining procedure is not intended for the detection of a specific condition or disorder. When reviewed in conjunction with other diagnostic tests and information may be used by trained professionals to define a physiological or pathological state.

Reagents

Scott's Tap Water Substitute Concentrate (Cat. No. S5134-6x100ML)
Magnesium sulfate • 7H₂O, 200 g/L, sodium bicarbonate, 20 g/L

Special Materials Required but Not Provided

- Differentiation Solution (Cat. No. A3179-1L or A3429-4L)
- Reagent Alcohol, 100% (Cat. No. R8382-1GA) OR Ethanol, 100%
- Reagent Alcohol, 95% OR Ethanol, 95%
- Hematoxylin Solution, Gill No. 1 (Cat. No. GHS1: GHS116-500ML; GHS132-1L; GHS128-4L)
- Hematoxylin Solution, Gill No. 2 (Cat. No. GHS2: GHS216-500ML; GHS232-1L; GHS280-2.5L; GHS2128-4L)
- Hematoxylin Solution, Gill No. 3 (Cat. No. GHS3: GHS3100ML; GHS316-500ML; GHS332-1L; GHS380-2.5L; GHS3128-4L)
- Xylene or Xylene Substitute

Counterstains

(choice depends on specimen and individual preference)

- Eosin Y solution, alcoholic (Cat. No. HT1101: HT110116-500ML; HT110132-1L; HT110180-2.5L; HT1101128-4L)
- Eosin Y solution, aqueous (Cat. No. HT1102: HT110216-500ML; HT110232-1L; HT110280-2.5L; HT1102128-4L)
- Eosin Y solution, alcoholic with phloxine (Cat. No. HT1103: HT110316-500ML; HT110332-1L; HT110380-2.5L; HT1103128-4L)
- Papanicolaou Stain OG-6 (Cat. No. HT401: HT40116-500ML; HT40132-1L; HT40180-2.5L; HT401128-4L)
- Papanicolaou Stain, Modified EA (Cat. No. HT40232-1L)
- Papanicolaou Stain EA 50 (Cat. No. HT403: HT40316-500ML; HT40332-1L; HT403128-4L)
- Papanicolaou Stain EA 65 (Cat. No. HT40432-1L)

Storage and Stability

Store reagents at room temperature (18-26°C) protected from light. Reagents are stable until the expiration date shown on the label.

Deterioration

Minor crystalline precipitate in Scott's Tap Water will not affect performance.

Preparation

Scott's Tap Water Substitute is prepared by diluting 1 volume Scott's Tap Water Substitute Concentrate with 9 volumes deionized water.

Hematoxylin Solutions, Gill Nos. 1, 2 and 3, are provided ready for use.

Acidified Eosin Y Solution, Aqueous is prepared by slowly adding up to 0.5 mL of glacial acetic acid per 100 mL of stain.

Precautions

This IVD is intended for in vitro diagnostic use in a clinical laboratory environment. This IVD is for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Notes

- The times given in the insert are approximate. Personal preferences will vary and the times can be adjusted to suit personal preferences. Stain solutions which are heavily used will lose their staining powers and the staining times should be lengthened or new solutions should be used.
- Some tap water supplies are acidic and unsuitable for use in the "blueing" portion of this procedure. If tap water is acidic, use a dilute alkaline solution.

- Purple or red-brown nuclei are indicative of inadequate "blueing".
- If eosin staining is excessive, nuclear staining may be masked. Proper eosin staining will demonstrate a 3-tone effect. To increase differentiation of eosin, extend time in alcohols or use a first alcohol with a higher water content. The times in the alcohols may be adjusted to obtain the proper degree of eosin staining.
- Positive control slides should be included in each run.

Procedure One

Staining Exfoliative Cytology Preparations Using Hematoxylin Solution, Gill No. 1 or Gill No. 2

- Fix cytologic smears in 95% ethanol for **15 minutes**.
- Rinse in gently running tap water for **30 seconds**.
- Stain in Hematoxylin Solution, Gill No. 1 or Gill No. 2 for **1.5 - 3 minutes**.
- Rinse in tap water.
- Blue in Scott's Tap Water Substitute for **15 - 60 seconds**.
- Rinse in tap water.
- Reagent Alcohol, 95% for **10 dips**.
- Counterstain in Papanicolaou Stain OG-6 for **1.5 minutes**.
- Reagent Alcohol, 95% for **10 dips**.
- Papanicolaou Stain EA 50, or Papanicolaou Stain EA 65, or Papanicolaou Stain Modified EA for **2.5 - 3 minutes**.
- Reagent Alcohol, 95%, two changes for **10 dips each**.
- Reagent Alcohol, 100%, two changes for **1 minute each**.
- Xylene or xylene substitute, two changes for **2 minutes each**.
- Coverslip and examine microscopically.

Procedure Two

Histology and/or Cytology Staining Using Hematoxylin Solution Gill No. 2 or Gill No. 3

- Deparaffinize to water or fix and dehydrate frozen sections.
- Stain in Hematoxylin Solution, Gill No. 2 or Gill No. 3 for **1.5 - 3 minutes**.
- Tap water wash.
- Differentiation Solution for **20 - 60 seconds**.
- Tap water wash.
- Blue in Scott's Tap Water Substitute for **5 - 60 seconds**.
- Tap water wash.
- Counterstain:

For Histology

Eosin Y Solution, Alcoholic, **Or** Acidified Eosin Y Solution, Aqueous, **Or** Eosin Y Solution, Alcoholic with Phloxine for **30 - 60 seconds**.

For Cytology

Papanicolaou Stain OG-6, **And** Papanicolaou Stain EA 50, **Or** Papanicolaou Stain EA 65, **Or** Papanicolaou Stain, Modified EA for **1 - 3 minutes**.

- Dehydrate, clear and mount.

Performance Characteristics

Chromatin appears blue to blue-black and nucleoli should be conspicuous.

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

| Cat. No | Product Description | Target | Intra-assay Specificity | Intra-assay Sensitivity | Inter-assay Specificity | Inter-assay Sensitivity |
|---------|---------------------|--------------|-------------------------|-------------------------|-------------------------|-------------------------|
| S5134 | Scott's Tap Water | Osmolality | 3 of 3 | 3 of 3 | 3 of 3 | 3 of 3 |
| | | Diluted pH | 3 of 3 | 3 of 3 | 3 of 3 | 3 of 3 |
| | | Undiluted pH | 3 of 3 | 3 of 3 | 3 of 3 | 3 of 3 |

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

S5134:

H315: Causes skin irritation.

H317: May cause an allergic skin reaction.

H318: Causes serious eye damage.

P261: Avoid breathing mist or vapors.

P264: Wash skin thoroughly after handling.

P272: Contaminated work clothing should not be allowed out of the workplace.

P280: Wear protective gloves/eye protection/face protection.

P302 + P352: IF ON SKIN: Wash with plenty of water.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

| | | | |
|--|---|---|---|
|  | Manufacturer |  | Catalogue Number |
|  | Consult Instructions for Use |  | Batch Code |
|  | Authorized Representative in the European Community/ European Union |  | European Union Declaration of Conformity (defined in IVDR 2017/746) |
|  | Use-by Date |  | In vitro diagnostic medical device |
|  | Temperature Limit |  | Caution |
|  | Date of Manufacture |  | Importer |

References

1. Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
2. Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p129

Contact Information

To place an order, please visit our web site at SigmaAldrich.com. For Technical Service, please visit the tech service page on our web site at SigmaAldrich.com/techservice.

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Merck KGaA,
64271 Darmstadt,
Germany

The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

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Sigma-Aldrich, Inc.,
MDSS GmbH