



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sia.com
sigma-aldrich.com

Product Information

Z-Ala-Glu(OMe)-Val-Asp(OMe) Fluoromethyl Ketone

Product Number **C 8484**

Storage Temperature -20°C

Product Description

Molecular Formula: $\text{C}_{28}\text{H}_{39}\text{FN}_4\text{O}_{10}$

Molecular Weight: 610.6

Z-Ala-Glu(OMe)-Val-Asp(OMe) Fluoromethyl Ketone is the methylated, cell permeable derivative of the caspase inhibitor Z-Ala-Glu-Val-Asp Fluoromethyl Ketone (Z-AEVD-FMK, Z = benzyloxycarbonyl).

Z-AEVD-FMK is a potent irreversible inhibitor of caspase-10 and, with lower efficiency, activated caspase-8. Caspases are a group of cysteine aspartate-specific proteases that play a role in apoptosis.^{1,2} Caspase 10 can activate the NF- κ B pathway.³ It contains two death effector domains and with caspase-8 is implicated in apoptosis signaling complexes associated with TNFR-like death receptors including Fas.^{4,5}

Methylation of the acidic amino acids Glu and Asp enhances the cell membrane permeability of Z-AEVD-FMK. Once in the cell, endogenous esterase activity hydrolyzes the methyl groups to form the biological active form. For *in vitro* studies an esterase needs to be included in the reaction mix to generate the active form of the molecule.

FMK is a trapping group responsible for irreversible inhibition and is also non-cytotoxic. Inhibition occurs when the FMK group covalently bonds to the -SH of an adjacent cysteine residue on the target protein.

Z-Ala-Glu(OMe)-Val-Asp(OMe)-Fluoromethyl Ketone is supplied as a white solid.

Preparation Instructions

Prepare 20 mM stock solutions in dry ($\geq 99.9\%$) DMSO to maintain product stability. Also soluble in DMF.

Storage/Stability

Store at -20°C . The product is reported to be stable at room temperature for one year in a desiccator.

References

1. Nicholson, D.W., and Thornberry, N.A., Caspases: killer proteases. *Trends Biochem. Sci.*, **22**, 299 (1997).
2. Cohen, G.M., Caspases: the executioners of apoptosis. *Biochem. J.*, **326**, 1-16 (1997).
3. Chaudhary, P.M., et al. Activation of the NF- κ B pathway by caspase 8 and its homologs. *Oncogene*, **19**, 4451-4460 (2000).
4. Fernandes-Alnemri, T. et al., In vitro activation of CPP32 and Mch3 by Mch4, a novel human apoptotic cysteine protease containing two FADD-like domains. *Proc. Natl. Acad. Sci. USA*, **93**, 7464-7469 (1996).
5. Vincenz, C., and Dixit, V. M., Fas-associated death domain protein interleukin-1 β -converting enzyme 2 (FLICE2), an ICE/Ced-3 homologue, is proximally involved in CD95- and p55-mediated death signaling. *J. Biol. Chem.*, **272**, 6578-6583 (1997).

BD 08/01

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.