

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

ProductInformation

Albumin from chicken egg white

Product Number A 5378 Storage Temperature 2-8 °C

Product Description

CAS Number: 9006-59-1Extinction coefficient: $E^{1\%} = 6.9$ to 7.6 (280 nm)¹ pl: 4.54 (4.43-4.66)²; 4.9³ Synonym: Ovalbumin

This product is not less than 90% pure by agarose gel electrophoresis.

Chicken egg albumin is the major protein constituent of egg whites. Chicken egg albumin is a phosphorylated-glycoprotein. From the amino acid sequence, the peptide portion of the molecule consists of 385 residues and has a molecular weight of 42.7 kDa.⁴ This sequence completely agrees with the reported m-RNA (messenger RNA) sequence.⁵ The carbohydrate and phosphate portions account for an additional 1428 and 160 grams per mole respectively, giving a total molecular weight of 44.3 kDa.⁶

Differential scanning calorimetry indicates the chicken egg albumin denatures at 84 °C.⁷

Ovalbumin can be used as a carrier protein to conjugate to synthetic peptides for use as an immunogen. It has the following amino acids: 20 Lys, 10 Tyr, 6 Cys, 14 Asp, and 33 Glu which make it suitable for conjugation.⁸

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This product is soluble in H2O (50 mg/ml), yielding a clear to slightly hazy solution.

References

- 1. Practical Handbook of Biochemistry and Molecular Biology, Fasman, G., Ed., CRC Press (Boca Raton, FL: 1989), pp. 298-299.
- Beeley, J.A., et al., Polyacrylamide gel isoelectric focusing of proteins: Determination of isoelectric points using an antimony electrode. Biochim. Biophys. Acta., 285, 293-300 (1972).
- Kidwai, S.A., et al, Effect of succinylation (3-carboxypropionylation) on the conformation and Immunological activity of ovalbumin. Biochem. J., 155, 171-180 (1976).
- Nisbet, A.D., et al., The complete amino-acid sequence of hen ovalbumin. Eur. J. Biochem., 115, 335-345 (1981).
- McReynolds, L., et al., Sequence of chicken ovalbumin mRNA. Nature, 273, 723-728 (1978).
- Tai, T., et al., Structures of the carbohydrate moiety of ovalbumin glycopeptide III and the difference in specificity of endo-beta-N-acetylglucosaminidases CII and H. J. Biol. Chem., 252, 6687-6694 (1977).
- Donovan, J.W. and Mapes, C.K., A differential scanning calorimetric study of conversion of ovalbumin to s-ovalbumin in eggs. J. Sci. Food Agric., 27, 197-204 (1976).
- Antibodies: A Lab Manual, Harlow, E. and Lane, D., Cold Spring Harbor Laboratory Press (Cold Spring Harbor, NY: 1988), p.77.

CMH/FEB/RXR 11/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.