
MILLIPORE



- ▶ Fast-flowing PES membrane
- ▶ No-tip/easy-grip flask design
- Low protein binding
- For volumes from 150 mL to 1L

Stericup® and Steritop™ Vacuum Filter Cups


With Millipore Express[®] PLUS membrane for fast sterilization of serum-based tissue culture media and other "hard-to-filter" solutions

Filter Faster

Stericup and Steritop filter cups feature Millipore Express PLUS polyethersulfone (PES) membrane–a fast flowing, low protein binding PES membrane–in a stable, convenient filter cup design.

Ideal for sterilizing tissue culture media, buffers, or other aqueous solutions, these filter cups are the high performance choice in the Stericup filter family. Stericup and Steritop tilter cups with Millipore Express PLUS membrane filter tissue culture media up to 25% faster than original Millipore Express membrane and 3x faster than filter cups with cellulose acetate (CA) membrane.

High Flow Rates



Average Flow Rate for 500 mL Vacuum Cups DMEM with 10% FBS (500 mL)

Millipore Express PLUS membrane is up to 25% faster than Millipore Express membrane — and 3x faster than cellulose acetate membranes.

Low Binding

The membrane's proprietary hydrophilic surface ensures that key growth factors and proteins won't be absorbed into the filter. Other membranes nonspecifically bind significantly more protein than Millipore Express PLUS membrane. For applications that require ultra-low protein binding, use Stericup filter cups with Durapore[®] PVDF membrane.


High Throughput

The asymmetric pore structure of the Millipore Express PLUS membrane gives it greater throughput than other membranes. This can be especially important when filtering serum or other proteinaceous solutions.

Intelligent Design

Stericup's no-tip, tapered receiver flask and compact profile improve stability and make gripping the receiver easier. The bottom of the flask is slightly recessed so that capped flasks can be stacked for storage. A small tab inside the funnel secures a prefilter in place during filtration.* For added protection against contamination, the threads of the bottle cap are slightly recessed in case the cap accidentally touches a work surface.

Low Protein Binding



Non-specific protein binding was determined by incubating a selection of membrane discs in a 1 mg/mL solution of ¹²⁵I-labeled IgG and then determining the bound protein. The results show that CA and nylon membranes bind significantly more protein than Millipore Express PLUS membrane. The results also show that Durapore is the lowest protein binding membrane.

Specification	S	
Connections	Hose barb vacuum inlet	
Materials		
Housing:	Polystyrene	
Membrane:	Either Millipore Express PLUS (PES); Durapore (PVDF); or Millipore Express (0.1 µm only)	
Sterilization	Gamma irradiation	
Filtration Area	40 cm ²	

Stericup Filter Cups

Ordering Information

Stericup Filter Cups combine a filter unit with a receiver flask and cap for processing and storage.

Product	Membrane	Pore Size (µm)	Funnel Capacity (mL)	Receiver Bottle (mL)	Qty/ Pk	Catalogue No.
Stericup-GP	Millipore Express PLUS (PES)	0.22	150 250	150 250	12 12	SCGP U01 RE SCGP U02 RE
			500	500	12	SCGP U05 RE
			500	1000	12	SCGP U10 RE
			1000	1000	12	SCGP U11 RE
Stericup-GV	Durapore (PVDF)	0.22	150	150	12	SCGV U01 RE
			250	250	12	SCGV UO2 RE
			500	500	12	SCGV U05 RE
			500	1000	12	SCGV U10 RE
			1000	1000	12	SCGV U11 RE
Stericup-HV	Durapore (PVDF)	0.45	150	150	12	SCHV U01 RE
			250	250	12	SCHV UO2 RE
			500	500	12	SCHV U05 RE
			500	1000	12	SCHV U10 RE
			1000	1000	12	SCHV U11 RE
Stericup-VP	Millipore Express (PES)	0.1	250	250	12	SCVP UO2 RE

Accessories

		Pore	
Product	Material	Size (µm)	Qty/ Pk Catalogue No.
Prefilter Disk	AP20 (glass fiber)	_	100 AP20 075 00

Steritop Bottle Top Filter Units

Ordering Information

Steritop Bottle Top Filter Units can be used with bottles having a narrow-mouth (33 mm) or a wide-mouth (45 mm) opening.

Product	Membrane	Pore Size (µm)	Funnel Capacity (mL)	Fitting Outlet Thread	Qty/ Pk	Catalogue No.
Steritop-GP	Millipore Express PLUS (PES)	0.22	150	33 mm	12	SCGP SO1 RE
				45 mm	12	SCGP TO1 RE
			250	33 mm	12	SCGP SO2 RE
				45 mm	12	SCGP TO2 RE
			500	33 mm	12	SCGP SO5 RE
				45 mm	12	SCGP TO5 RE
			1000	33 mm	12	SCGP S10 RE
				45 mm	12	SCGP T10 RE
Steritop-GV	Durapore (PVDF)	0.22	150	33 mm	12	SCGV SO1 RE
				45 mm	12	SCGV TO1 RE
			250	33 mm	12	SCGV SO2 RE
				45 mm	12	SCGV TO2 RE
			500	33 mm	12	SCGV SO5 RE
				45 mm	12	SCGV TO5 RE
			1000	33 mm	12	SCGV S10 RE
				45 mm	12	SCGV T10 RE
Steritop-VP	Millipore Express (PES)	0.1	250	33 mm	12	SCVP SO2 RE
				45 mm	12	SCVP TO2 RE
Receiver Bottles and Caps			150	45 mm	12	SCOO BO1 RE
			250	45 mm	12	SCOO BO2 RE
			500	45 mm	12	SCOO BO5 RE
			1000	45 mm	12	SCOO B10 RE

Accessories

Product	Material	Size (µm)	Qty/ Pk Catalogue No.
Prefilter Disk	AP20 (glass fiber)	_	100 AP20 075 00

Part of a Family of Products for Sterile Filtration

Millipore offers a full selection of sterilizing filter products, including syringe-, pressure-, vacuum-, and pump-driven devices in a wide range of pore sizes and volumes. For more information, contact Millipore or visit www.millipore.com/sterile.

For Additional Information

Internet: www.millipore.com Millipore Offices: www.millipore.com/offices Tech Service: www.millipore.com/techservice

Millipore, Durapore, Millipore Express, Stericup, and Steritop are trademarks of Millipore Corporation. Lit. No. PF1597EN00 Rev. A 6/04 04-165 © 2002, 2004 Millipore Corporation, Billerica, MA. All rights reserved. Printed in the U.S.A.

MILLIPORE