Chemical Duty Vacuum/Pressure Pump

WP61 115 60, WP61 220 50, WP61 100 60
and

High Output Vacuum/Pressure Pump

WP62 115 60, WP62 220 50, WP62 100 60

For complete information about this product, please contact Millipore Corporation.

Usage Guidelines

⚠ WARNING: The motor is thermally protected and will automatically restart unexpectedly when the overload device resets. Do NOT pump flammable or explosive gases or vapors or operate this pump in an atmosphere containing flammable or explosive gases or vapors.

- For WP61 pump only, to ensure proper cooling, the pump cooling fan must not be blocked and should be a minimum of 1” (2.5 cm) away from any wall or obstruction.
- To prevent liquids or vapors from entering the pump when applying vacuum, always use the Millex®-FA₅₀ vent filter and 1/4” silicone tubing included with the pump.
- For maximum protection of the pump, use the Millex®-FA₅₀ filter **AND** a vacuum trap to protect the pump from vapors and gases.
- Do not operate pump if hose, plug, or cord are damaged.
- Do not apply oil to any part of the pump.
- To operate at maximum efficiency, the pump system must be thoroughly clean.
- Do not aim compressed air blast at anyone because it can cause serious injury. Keep children away.
- To avoid serious burns, never touch the head parts and tubing during and immediately after operation.

Operation

Connect approx. 1/4” I.D. tubing to filter holder or other necessary equipment. First plug pump’s power cord into outlet, then turn the pump on using toggle switch. Connect tubing to pump and begin procedure.

If using the pump for vacuum: Adjust the vacuum by closing the pressure regulator and opening the vacuum regulator fully (counterclockwise). Slowly turn down the vacuum regulator to obtain the desired vacuum reading.

If using the pump for pressure: Adjust the pressure by turning the vacuum regulator down fully and opening the pressure fully (counterclockwise). Slowly turn down the pressure regulator to obtain the desired pressure reading.

Specifications

WP61 115 60, WP61 220 50, WP61 100 60

Body and pump head: die-cast aluminum

Diaphragm: Buna-N with PTFE lining

Head gasket: Buna-N

Leaf valves: high grade stainless steel

Internal surfaces: PTFE-coated

Vacuum: 24” Hg; Pressure: 20 psi continuous, 35 psi intermittent

Electrical specifications: Shaded-pole, single phase, 1/15 HP motor

Line voltage limits: ± 10% of supply voltage

WP62 115 60, WP62 220 50, WP62 100 60

Body and pump head: die-cast aluminum

Pump type: piston with PTFE seal

Leaf valves: high grade stainless steel

Vacuum: 27.2” Hg; Pressure: 50 psi continuous, 80 psi intermittent

Electrical specifications:

115 V, 60 Hz or 100 V, 50/60 Hz AC: Shaded-pole, single phase, 1/7 HP motor
230 V, 50 Hz AC: Permanent split capacitor, 1/3 HP motor

Line voltage limits: ± 10% of supply voltage

Technical Assistance

Visit the tech service page on our web site at SigmaAldrich.com/techservice.

Any serious incident of this device should be reported to manufacturer and competent authority of country where user is established.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Merck, Millipore and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2009-2025 Merck KGaA, Darmstadt, Germany and/or its affiliates.
All Rights Reserved.

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