

User Guide

Dualex™-Ultra Transducer Protector

REF SLGVSR25US

- 25 mm
- Sterile
- Single use
- Non-Pyrogenic
- For research use only
- 50 pack

Introduction

This insert describes how to use the Dualex™-Ultra Transducer Protectors. Each transducer protector consists of a 0.22 µm hydrophobic Durapel™ polyvinylidene fluoride (PVDF) filter sealed in a polyvinyl chloride (PVC) housing. It is non-pyrogenic and non-toxic. It prevents fluids from damaging sensitive transducer monitoring devices. Dualex™-Ultra devices are equipped with Luer-Lok® connectors, making them ideal for use in high-flux systems.

Directions for Use

WARNINGS

- To ensure sterility, do not use this product if the package is damaged.
- Do not use this product as an in-line filter. This product was not designed for long-term continuous use when using as an in-line filter protective device for pressure monitoring systems.
- Do not use with monitoring systems greater than 1.7 bar (25 psi) or syringes smaller than 10 mL because pressures in excess of the maximum pressure rating may be reached, potentially causing damage to the filter unit and/or personal injury.
- Sudden loss of pressure could indicate failure of the filter.
- Use aseptic technique to maintain sterility removing product from primary package.


CAUTIONS

- Filter gaseous solutions with a temperature limit 45 °C (113 °F).
- Intended for short term use. Change within 24 hours.

- Single use only; do not re-use or resterilize.
- This product is not suitable for filtration of aqueous solutions.
- Discard appropriately after single use. See "Disposal" on next page.

Procedure

1. Remove the cover from the primary package of the Dualex™-Ultra Transducer Protector using aseptic techniques.
2. Insert female Luer-Lok® inlet of the transducer protector to the male Luer connection on the system and remove the rest of the package from the assembly.
3. Attach the monitoring line to the male Luer-Lok® connection on the transducer protector using aseptic techniques.
4. Perform single-use test using standard procedures.
If liquids enter the unit and fill the housing, the Dualex™-Ultra Transducer Protector blocks the passage of air. See "Troubleshooting" section.
5. Remove the transducer protector after each test and dispose according to local guidelines. See "Disposal" section.



TROUBLESHOOTING

If liquids fill the housing, clamp off the monitoring line and remove the transducer protector from the machine. Release the level of feed line fluid, if necessary. Then attach a new Dualex™-Ultra Transducer Protector to the machine, and release the clamp.

Specifications

Materials

| | |
|-----------|--|
| Membrane | Hydrophobic Durapel™ polyvinylidene fluoride (PVDF), Type GV |
| Pore size | 0.22 µm |
| Housing | Polyvinyl chloride (PVC) |

Dimensions

| | |
|-----------------|------------------|
| Inlet to outlet | 23 mm (0.91 in.) |
| Diameter | 29 mm (1.14 in.) |
| Filtration area | 4 cm² (0.62 in²) |

Temperature limit 45 °C (113 °F)

Pressure limit at 21 °C 1.7 bar (25 psi) inlet and differential

Sterilization method Ethylene oxide gas

Connections Female Luer-Lok® inlet
Male Luer-Lok® outlet

Symbol Definitions

| Symbol | Definition | Symbol | Definition |
|--------|-------------------|--------|----------------------------------|
| | Non-pyrogenic | | Use-by date |
| | Do not re-use | | Date of manufacture |
| | Catalogue number | | Manufacturer |
| | Batch code | | Do not use if package is damaged |
| | Temperature limit | | Sterilized using ethylene oxide |

Disposal

Follow precautions for disposal of items contaminated with potentially infectious or hazardous biological material according to all applicable international, federal, state, and local regulations.

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page on our web site at SigmaAldrich.com/TechService.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/Terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/Offices.

Merck, Millipore, Dualex, Durapel and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources. © 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.