

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone (800) 325-5832 (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

ProductInformation

ADRENERGIC LIGAND-SET™

Product Number L 0383 Storage Temperature –20 °C

Product Description

The Adrenergic LIGAND-SET[™] is a set of 48 small organic ligands to adrenergic receptors. These ligands are arrayed in a standard 96-well plate format; each well has a capacity of 1 ml.

This set can be used for screening new drug targets, for guiding secondary screens of larger, more diverse libraries and for standardizing and validating new screening assays.

The adrenergic receptors are divided into three types, α_1 , α_2 and β , each with several subtypes. Each type has a distinct function and pharmacology.

There are three subtypes of α_1 adrenergic receptors, α_{1A} , α_{1B} and α_{1D} . α_{1} Adrenergic receptors are widely distributed in the central and peripheral nervous systems, where their activation generally results in depolarization and increased neuronal firing rate. In the periphery, they are activated by either norepinephrine released from sympathetic nerve terminals or epinephrine released from the adrenal medulla. Activation of these receptors mediates a wide variety of functions, including contraction of smooth muscle, cardiac stimulation, cellular proliferation and activation of hepatic gluconeogenesis and glycogenolysis. Most of the peripheral actions of α_1 adrenergic receptors are mediated through the inositol phosphate second messenger pathway, while there is evidence that α_1 receptors in the CNS activate adenylyl cyclase within the CNS.

There are four subtypes of α_2 adrenergic receptors, α_{2A} , α_{2B} , α_{2C} and α_{2D} . The α_2 adrenergic receptors are also widely distributed and activated by norepinephrine and epinephrine. The most well-characterized action is the inhibition of the release of neurotransmitter from many peripheral and central neurons by prejunctional α_2 receptors. α_2 receptors are also present at postjunctional sites, where they mediate actions such as smooth muscle contraction, platelet aggregation and inhibition of insulin secretion.

There are three subtypes of β adrenergic receptors, β_1 , β_2 and β_3 . These are also activated by norepinephrine and epinephrine and are widely distributed. Activation of B adrenergic receptors results in stimulation of cardiac rate and force, relaxation of vascular, urogenital and bronchial smooth muscle, stimulation of renin secretion from the juxta-glomerular apparatus, stimulation of insulin and glucagon secretion from the endocrine pancreas, stimulation of glycogenolysis in liver and skeletal muscle, and stimulation of lipolysis in adipocytes. Prejunctional β adrenergic receptors are present on some central and peripheral nerve terminals, where their activation results in facilitation of stimulation-evoked neurotransmitter release. However, in contrast to the prejunctional α_2 receptors, activation of prejunctional β receptors do not appear to have major physiologic significance.

Components/Reagents

The Adrenergic LIGAND-SET[™] contains 2 mg of each ligand per well. Stock solutions can be readily prepared by adding 1 ml of DMSO to each well. The set also comes with a diskette containing an ISIS/Base Structure Data (SD) file and a Microsoft Excel file containing structural information as well as the catalog number, name, rack position and pharmacological characteristics of each ligand. The following is a listing of all the ligands included:

A0779	p-Aminoclonidine hydrochloride	
A7655	(±)-Atenolol	
A-131	Alprenolol hydrochloride	
A9512	L(-)-Arterenol (Norepinephrine) bitartrate	
B-012	6-Fluoronorepinephrine hydrochloride	
B-016	Benoxathian hydrochloride	
B-019	Phenoxybenzamine hydrochloride	
B-154	BU224 hydrochloride	
B-161	B-HT 933 dihydrochloride	
B-162	B-HT 920 dihydrochloride	
B-169	BRL 37344 sodium	
B2656	Benextramine tetrahydrochloride	
B8406	Bretylium tosylate	
C-125	(±)-CGP-12177A hydrochloride	
C-223	Cirazoline hydrochloride	
C-231	CGP 20712A methanesulfonate	
C-247	Cyclazosin hydrochloride	
C7897	Clonidine hydrochloride	
D0411	WB-4101 hydrochloride	
D0676	Dobutamine hydrochloride	
D5290	(-)-alpha-Methylnorepinephrine	
E4375	(-)-Epinephrine bitartrate	
G-110	Guanabenz acetate	
I-114	p-lodoclonidine hydrochloride	
I-127	ICI 118,551 hydrochloride	
12760	R(-)-Isoproterenol (+)-bitartrate	
15627	(±)-lsoproterenol hydrochloride	
18005	S(+)-Isoproterenol (+)-bitartrate	
M3656	Mephentermine hemisulfate	

M6524	Methoxamine hydrochloride	
N-153	Nylidrin hydrochloride	
N-158	Naftopidil dihydrochloride	
O2378	Oxymetazoline hydrochloride	
P0778	Pindolol	
P0884	(±)-Propranolol hydrochloride	
P-119	(±)-Pindobind	
P6126	Phenylephrine hydrochloride	
P7561	Phentolamine mesylate	
P7791	Prazosin hydrochloride	
R-104	Rauwolscine hydrochloride	
R-134	Rilmenidine hemifumarate	
S5013	Salbutamol hemisulfate	
T6394	S(-)-Timolol maleate	
U-100	Urapidil hydrochloride	
U-101	Urapidil, 5-Methyl-	
U-104	UK 14,304	
X1251	Xylazine hydrochloride	
Y3125	Yohimbine hydrochloride	

Preparation Instructions

<u>To create a new database in ISIS</u>[™]/BASE :

- Open ISIS[™]/BASE.
- Choose <u>File>New database</u>.
- Enter **Adrenergic** or a preferred name in the File name field.
- Click Save.

The "Create Database" window will now be open.

- Enter **Catnum** for the Field name.
- Choose Variable text from the drop down window of the Type field.
- Click Add.
- Repeat the above steps for the following:

Field name	<u>Type</u>
Name	Variable text
Position	Variable text
Action	Variable text
Class	Variable text
Selectivity	Variable text
SecName	Variable text
Description	Variable text

- Enter Structure for the Field name.
- Choose **Structure** from the drop down window of the Type field.
- Enter *Structure for the External name.
- Click Add.
- Click Save.

The main ISIS[™]/BASE window will now be open.

To create the Form:

- Click on the "Draw a box" button (second button down on the left of the screen).
- Move the mouse to the bottom left hand corner and draw a box, ½ inch high, the length of the screen by clicking on the left mouse button and dragging the mouse across the screen. (see figure below)
- Above this box, draw another ½ inch high box the length of the screen. (see figure below)
- Above this box, draw a third ½ inch high box the length of the screen. (see figure below)
- Above these long boxes draw 3 ½ inch high x 3 inch wide boxes. (see figure below)
- Above these 3 boxes, draw another three the same size. (see figure below)
- Draw a final box to fit the remaining space of the screen above these boxes. (see figure below)

	re Droch Lit Window Belo	U of U	
na Garry Brownes	Update Genetic Mi		
F			
	Cot No.	Park	
io.	Gat. No.	Rack	
ID Action	Gat. No. Class	Rack Selectivity	
Action			
Action			
Action Primary Name Sec. Name			
Action Primary Name			
Action Primary Name Sec. Name			

Double click on the top box. This will open the Box properties window.

- Click on Structure.
- Click OK.
- Repeat the same steps, clicking on the appropriate field name for the appropriate box:

<u>Box</u>	Field name
First small box	ID
Second small box	Catnum
Third small box	Position
Fourth small box	Action
Fifth small box	Class
Sixth small box	Selectivity
First long box	Name
Second long box	SecName
Bottom long box	Description
~	

- Choose File>Save form.
- Enter Adrenergic or preferred name.
- Click OK.

Importing an SD file:

- Click Update.
- Choose File>Import>SD File. NOTE: For MAC users, you must hold down the option key while choosing File>Import>SD File. If you do not, the Adrenergic.sdf will not be visible in the import window.
- Enter Adrenergic.sdf (Located on the floppy diskette provided with the plate).
- Click **Open**.
 The Import SD File window will now be open.
- Click on Add a new record including structure, on both sides of the table.
- Click OK.
- The database is now ready to use.

Storage/Stability

Store plate -20°C with cap strips firmly in place. Plate cover should only be removed when plate is in use to prevent loss of caps strips.

References

- Leonardi, A., et al. "α₁-adrenoceptor subtype and organ-selectivity of different agents." in *Perspectives in Receptor Research*, Vol. 24, Eds. D. Giardina, A. Piergentili and M. Pigini, pp. 135-152 Elsevier, Amsterdam (1996).
- Macdonald, E., et al., "Gene targeting-homing in on alpha 2-adrenoceptor subtype function." *Trends Pharmacol. Sci.*, **18**, 211-219 (1997).
- Ruffolo, R.R., et al., "Alpha- and betaadrenoceptors: From the gene to the clinic. Part 2." *J. Med. Chem.*, **38**, 3681-3716 (1995).

SMS 8/00

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.