

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-576
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

# **ProductInformation**

SIGMA QUALITY CONTROL TEST PROCEDURE Enzymatic Assay of  $\beta$ -GLUCURONIDASE (EC 3.2.1.31) from E. coli

### PRINCIPLE:

PheP-Gluc + H<sub>2</sub>O β-Glucuronidase > D-Glucuronate + Phenolphthalein

Abbreviation used:

PheP-Gluc = Phenolphthalein Glucuronide

**CONDITIONS:**  $T = 37^{\circ}C$ , pH = 6.8,  $A_{540nm}$ , Light path = 1 cm

**METHOD:** Spectrophotometric Stop Rate Determination

### **REAGENTS:**

A. 75 mM Potassium Phosphate Buffer, with 1.0% (w/v) Bovine Serum Albumin, pH 6.8 at 37°C (Prepare 100 ml in deionized water using Potassium Phosphate, Monobasic, Anhydrous, Sigma Prod. No. P-5379 and Albumin, Bovine, Sigma Prod. No. A-4503. Adjust to pH 6.8 at 37°C with 1 M KOH.)

- 3.0 mM Phenolphthalein Glucuronide Substrate Solution (PheP-Gluc)
   (Prepare 10 ml in deionized water using Phenolphthalein Glucuronic Acid, Free Acid, Sigma Prod. No. P-0501.)
- C. 200 mM Glycine Buffer Solution, pH 10.4. (Use Glycine Buffer Solution, Sigma Stock No. 105-2, or prepare 100 ml in deionized water using Glycine Free Base, Sigma Prod. No. G-7126. Adjust to pH 10.4 at 37°C with 1 M NaOH.)
- β-Glucuronidase Enzyme Solution
   (Immediately before use, prepare a solution containing 400 800 units/ml of β-Glucuronidase in cold Reagent A.)
- E. 95% (v/v) Ethanol (Prepare 20 ml in deionized water using 200 Proof USP Ethyl Alcohol, Quantum Chemical Corporation.)

SSPHEN11
Revised: 01/29/98

# Enzymatic Assay of β-GLUCURONIDASE (EC 3.2.1.31) from E. coli

# **REAGENTS:** (continued)

F. 0.05% (w/v) Phenolphthalein Standard Solution (Std Soln) (Prepare 5 ml by dissolving 2.5 mg of Phenolphthalein, Sigma Prod. No. P-9750 in 5 ml of Reagent E or use Phenolphthalein Standard Solution, Sigma Stock No. 105-1.)

# PROCEDURE:

Pipette (in milliliters) the following reagents into suitable containers:

	<u>Test</u>	Blank
Deionized Water Reagent A (Buffer) Reagent B (PheP-Gluc)	0.65 0.50 0.25	0.65 0.50 0.25
Mix by inversion and equilibrate to 37°C.	Then add:	
Reagent D (Enzyme Solution)	0.10	
Mix by inversion and incubate at 37°C for	exactly 30 minutes.	Then add:
Reagent C (Glycine Buffer) Reagent D (Enzyme Solution)	5.00	5.00 0.10

Immediately mix by inversion. Transfer the solutions to suitable cuvettes and record the  $A_{540nm}$  for both the Test and Blank using a suitable spectrophotometer.

# **COLORIMETRIC ASSAY:**

Standard Curve:

Prepare a standard curve by pipetting (in milliliters) the following reagents into suitable containers:

Std 1	Std 2	Std 3	Std 4	Std 5	Blank
0.65	0.65	0.65	0.65	0.65	0.65
0.50	0.50	0.50	0.50	0.50	0.50
0.25	0.25	0.25	0.25	0.25	0.25
0.02	0.03	0.05	0.07	0.10	
0.08	0.07	0.05	0.03		0.10
5.00	5.00	5.00	5.00	5.00	5.00
	0.65 0.50 0.25 0.02 0.08	0.65 0.65 0.50 0.50 0.25 0.25 0.02 0.03 0.08 0.07	0.65     0.65     0.65       0.50     0.50     0.50       0.25     0.25     0.25       0.02     0.03     0.05       0.08     0.07     0.05	0.65     0.65     0.65     0.65       0.50     0.50     0.50     0.50       0.25     0.25     0.25     0.25       0.02     0.03     0.05     0.07       0.08     0.07     0.05     0.03	0.65     0.65     0.65     0.65     0.65       0.50     0.50     0.50     0.50     0.50       0.25     0.25     0.25     0.25     0.25       0.02     0.03     0.05     0.07     0.10       0.08     0.07     0.05     0.03

SSPHEN11 Page 2 of 4

Revised: 01/29/98

# Enzymatic Assay of β-GLUCURONIDASE (EC 3.2.1.31) from E. coli

# **COLORIMETRIC ASSAY: (continued)**

Mix by inversion and transfer the standards to suitable cuvettes. Record the  $A_{540nm}$  for each standard using a suitable spectrophotometer.

### **CALCULATIONS:**

Standard Curve:

 $\Delta A_{540}$  Standard =  $A_{540}$  Standard -  $A_{540}$  Standard blank

Prepare a standard curve by plotting the  $\Delta A_{540}$  for the Standard vs micrograms of Phenolphthalein.

Sample Determination:

 $\Delta A_{540}$  Sample =  $A_{540}$  Sample -  $A_{540}$  Sample blank

Determine the total micrograms of phenolphthalein liberated using the Standard curve.

Units/ml enzyme = 
$$\frac{(\mu g \text{ phenolphthalein released})(2)(df)}{(\mu g \text{ phenolphthalein released})(2)(df)}$$

0.1

2 = Time correction of assay (Unit Definition = 1 hour) df = Dilution factor 0.1 = Volume (in milliliter) of enzyme used

### **UNIT DEFINITION:**

One Sigma or modified "Fishman" unit will liberate 1.0  $\mu g$  of phenolphthalein from phenolphthalein glucuronide per hour at pH 6.8 at 37°C.

SSPHEN11 Page 3 of 4 Revised: 01/29/98

# Enzymatic Assay of β-GLUCURONIDASE (EC 3.2.1.31) from E. coli

### FINAL ASSAY CONCENTRATION:

In a 1.50 ml reaction mix, the final concentrations are 30 mM potassium phosphate, 0.50 mM phenolphthalein glucuronic acid, and 40 - 80 units β-glucuronidase.

#### **REFERENCES:**

Fishman, W.H. and Bernfeld, P. (1955) Methods in Enzymology, Volume I, 262-269

Combie, J., Blake, J.W., Nugent, T.E., and Tobin, T. (1982) Clin. Chem. 28, 83-86

Fishman, W.H. (1974) in Methods of Enzymatic Analysis (Bergmeyer, H.U. ed) 2nd ed., Volume II, 930-932, Academic Press, New York, NY

#### NOTES:

- 1. This assay is based on the cited references.
- 2. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

Sigma warrants that the above procedure information is currently utilized at Sigma and that all Sigma-Aldrich, Inc. products conform to the information in this and other Sigma-Aldrich, Inc. publications. Purchaser must determine the suitability of the information and product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packaging slip.

Page 4 of 4 SSPHEN11