




SigmaScreen™ Collagen I Coated Plates

Unrivalled consistency, adhesion and growth capacity for your cell-based assays



Expanding your options for cell-based high-throughput screening, Sigma's Collagen I coated plates join our Poly-D-Lysine to bring to you more ways to get the maximum cell attachment, survival and differentiation.

The SigmaScreen Collagen I product line of 96- and 384-well coated plates, available in a variety of colors for different assay detection formats and minimal cross-talk, has been optimized for high-throughput screening. Developed to meet proposed Society of Biomolecular Screening microplate specifications,* SigmaScreen Collagen I Coated Plates excel in applications such as calcium flux assays, cell proliferation, gene reporter and apoptosis assays as well as cell adhesion kinetic studies.

- **Increased cell proliferation with serum-free media**

Whether your cells are cultured with serum or in serum-free media, Sigma's Collagen I Coated Plates offer a superior support for cell adhesion and proliferation over standard tissue culture treated plates or competitor's Collagen I coated plates.

- **Compatible with many common cell-types**

Endothelial, epithelial, hepatocytes, and muscle cells as well as transfected lines such as HEK 293, BHK, L929 and HT1080. For additional cell lines visit sigma-aldrich.com/hts.

- **Consistent quality and testing**

Proprietary coating process manufactured under ISO9002 in our GMP facility delivers high well-to-well and plate-to-plate consistency making SigmaScreen Coated Plates the new standard in high-throughput screening. Each lot of Collagen I coated plates has been tested for cell attachment performance in serum free media using HT1080 cells. Each lot is also tested for endotoxins, bacteria and fungi.

- **Flexibility in your assay detection**

Offering a variety of 96- and 384-well Corning® base plates designed and validated with all commercial microplate readers and robotic systems in compliance with proposed Society of Biomolecular Screening microplate specifications.*

- **Stable at cooler or room temperatures for added storage convenience**

- **Custom options**

Special requirements? No problem. Sigma also offers custom coating and packaging options manufactured under our same high-quality ISO9002, GMP facility guidelines. For additional information on our custom capabilities contact us at hts@sial.com and let us know your interests.

* For additional information on SBS microplate standards visit sbsonline.org

SigmaScreen™ Collagen I Coated Plates

Superior Cell Adhesion and Proliferation




Figure 1: HT1080 cells were seeded in clear 96-well Collagen I Coated Plates (3×10^5 cells/ml, 100 μ l/well). Cells were incubated for 24 hours in serum-free DME medium (D 5671), and washed two times with PBS (P 3813). Photographs were taken at 100X. Results show SigmaScreen coated plates demonstrate increased viability and improved cell morphology.

Product Specifications

Specification	96-Well Coated Plates	384-Well Coated Plates
Plate and lid composition	Tissue culture polystyrene	Tissue culture polystyrene
Lid	Yes	Yes
Well configuration	Flat bottom/round	Flat bottom/square
Well width	6.4 mm	2.8 mm
Well depth	11 mm	11 mm
Maximum recommended working volume, per well	250 μ l	75 μ l
Growth area, per well	32 mm ²	8 mm ²
Coating coefficient of variability, well-to-well	$\leq 10\%$	$\leq 10\%$

For further information on our growing line of HTS products, visit us at sigma-aldrich.com/hts

To Order Call: 1-800-325-3010

Technical Service: 1-800-325-5832

Multiple Temperature Stability with SigmaScreen Coated Plates




Figure 2: Clear 96-well Collagen I Coated Plates were incubated for two months at 25 °C, 37 °C, and 50 °C. HT1080 cells were seeded in the wells (2.0×10^5 cells/ml, 100 μ l/well), and incubated for 24 hours in serum-free DME medium (D 5671). The plates were washed three times with PBS (P 3813) and lysed with Cellytic™ M Mammalian Cell Lysis/Extraction Reagent (C 2978). The total cellular protein attached to the plates was quantitated using Bradford Reagent (B 6919), and plates were read in a Wallac Victor plate reader at 620 nm.

Additional Plate Coatings

Poly D-Lysine	ExtrAvidin
Streptavidin	Anti-FLAG™ M2
HIS-Select™ High Capacity	HIS-Select™ High Sensitivity

Ordering Information

Product Number	Description	Package Size
S 3690	SigmaScreen Collagen I Coated Plates, 384-well black	5 each 10 x 5 each
S 3190	SigmaScreen Collagen I Coated Plates, 96-well black	5 each 10 x 5 each
S 3440	SigmaScreen Collagen I Coated Plates, 384-well clear	5 each 10 x 5 each
S 8813	SigmaScreen Collagen I Coated Plates, 96-well clear	5 each 10 x 5 each
S 3565	SigmaScreen Collagen I Coated Plates, 384-well white	5 each 10 x 5 each
S 3315	SigmaScreen Collagen I Coated Plates, 96-well white	5 each 10 x 5 each

All plate formats are clear bottom and include lids.

