

# Product Information

## EX-CELL™ NS0 Serum-Free Medium for NS0 Cells, Chemically Defined

without L-glutamine

CATALOG NO. 14650C

### Description

EX-CELL™ NS0 is an animal-component free, protein-free, chemically defined, serum-free liquid medium developed for the long-term growth of NS0-related cells in suspension culture. The NS0 cells are clonal derivatives of the parent NS1 cell line and are capable of growth in suspension culture. NS0 hybridoma suspension cultures can be subcultured directly into EX-CELL™ NS0 from serum-supplemented or serum-free media. Suspension cultures in EX-CELL™ NS0 have been carried for more than 50 passages with no loss of growth or viability.

### Formulation

The formula for EX-CELL™ NS0 is proprietary to SAFC Biosciences. For additional information please call our Technical Services department.

### Precautions

Use aseptic technique when handling or supplementing this medium. This product is for research or for further manufacturing use. THIS PRODUCT IS NOT INTENDED FOR HUMAN OR THERAPEUTIC USE.

### Storage

Store medium at 2 to 8 C, protected from light. Do not use after the expiration date.

### Indications of Deterioration

Medium should be clear and free of particulates and flocculent material. Do not use if medium is cloudy or contains precipitate. Other evidence of deterioration may include color change, pH shift or degradation of physical or performance characteristics.

### Preparation Instructions

EX-CELL™ NS0 is formulated without L-glutamine. Prior to use, this medium should be supplemented with 8 mM L-glutamine by adding 40 mL/L of a 200 mM solution (Catalog No. 59202C). SAFC Biosciences recommends L-glutamine supplementation of the working volume only. Supplements, such as antibiotics, can be added to the sterilized medium using aseptic technique. Storage conditions and shelf life of the product may be affected by the nature of the supplement.

### Methods for Use

#### Adaptation

NS0 hybridoma cells that have been grown in a conventional serum-supplemented medium can be adapted to EX-CELL™ NS0 directly with no weaning. Adaptation to EX-CELL™ NS0 requires healthy, viable cultures in mid-logarithmic growth phase.

1. Subculture the cells from serum-supplemented medium into EX-CELL™ NS0 (supplemented with 8 mM L-glutamine as described previously) at a minimum seeding density of  $2 \times 10^5$  cells/mL in shaker flasks.
2. Incubate the flasks at 37 C in a humidified incubator with 5% CO<sub>2</sub>. Maintain the orbital shaker speed between 105 - 115 rpm.
3. Continue to subculture cells in EX-CELL™ NS0 every 3 - 4 days or when the cell density reaches  $2 \times 10^6$  cells/mL.
4. Allow the cells to adapt to EX-CELL™ NS0 for an additional 3 - 6 passages. Cells are considered fully adapted to EX-CELL™ NS0 when growth rates return to normal and viabilities are above 95%.

#### United States

SAFC Biosciences, Inc.  
13804 W. 107th Street  
Lenexa, Kansas 66215  
USA  
Phone +1 913-469-5580  
Toll free-USA 1 800-255-6032  
Fax +1 913-469-5584  
E-mail info-na@sial.com

#### Europe

SAFC Biosciences Ltd.  
Smeaton Road, West Portway  
Andover, Hampshire SP10 3LF  
UNITED KINGDOM  
Phone +44 (0)1264-333311  
Fax +44 (0)1264-332412  
E-mail info-eu@sial.com

#### Asia Pacific

SAFC Biosciences Pty. Ltd.  
18-20 Export Drive  
Brooklyn, Victoria 3025  
AUSTRALIA  
Phone +61 (0)3-9362-4500  
Toll free-AUS 1 800-200-404  
Fax +61 (0)3-9315-1656  
E-mail info-ap@sial.com

## Culture Techniques

NS0 cells are normally grown at  $37 \pm 1$  C and 5% CO<sub>2</sub>. Allow the medium to warm to room temperature prior to use (protect from light). Once fully adapted, the cells should be subcultured at a seeding density of at least  $2 \times 10^5$  cells/mL in shaker flasks. Seed 30 mL cell cultures in 125 mL shaker flasks and 60 mL cultures in 250 mL shaker flasks. Shaker speed should be between 105 - 115 rpm.

When passing the cells, carry over should not exceed 25% of the final volume. If carry over exceeds 25%, centrifugation is recommended. Cells propagated in serum-free media are extremely fragile. For successful results, care must be taken when subculturing cells. Standard techniques of centrifugation must be modified to include low-speed centrifugation to prevent damage to cells that have been propagated in serum-free medium.

## Cryopreservation

### Freezing:

NS0 cells may be frozen in EX-CELL™ NS0 without the reintroduction of serum. However, it is necessary to handle the cells gently and freeze the cells under carefully controlled conditions.

1. Choose cultures in logarithmic growth phase with viabilities above 90%.
2. Prepare a freezing medium consisting of 45% cold EX-CELL™ NS0 medium, 45% spent medium and 10% dimethyl sulfoxide (DMSO). Alternatively, cells can be frozen in 90% fresh medium and 10% DMSO. Recovery may require 1 - 2 additional passes.
3. Centrifuge the cells at 200 *g* for 5 minutes. Remove the supernatant and prepare the freezing medium.
4. Gently resuspend the cells in the freezing medium at  $1 \times 10^7$  cells/mL.
5. Rapidly transfer 1 - 2 mL of this suspension to sterile cryovials.
6. Place the vials at -20 C for 3 - 4 hours, then transfer to -70 C for 16 - 24 hours.
7. For long-term storage, transfer the vials to liquid nitrogen vapor.

### Thawing:

1. Rapidly thaw a vial of frozen cells in a 37 C water bath without agitation.
2. Transfer the cells aseptically to a centrifuge tube containing 5 mL of cold EX-CELL™ NS0 medium.
3. Count the cells for viability and transfer to a sterile shaker flask at a minimum seeding density of  $2 \times 10^5$  cells/mL.
4. Pass the cells using standard cell culture techniques.

## Characteristics

### Appearance

Clear yellow solution

### Endotoxin

Refer to Certificate of Analysis

### Osmolality (as supplied)

305 - 335 mOsm/kg H<sub>2</sub>O

### pH (as supplied)

6.8 - 7.2

### Sterility

No microbial growth detected

### Warranty, Limitation of Remedies

SAFC Biosciences warrants to the purchaser for a period of one year from date of delivery that this product conforms to its specifications. Other terms and conditions of this warranty are contained in SAFC Biosciences' written warranty, a copy of which is available upon request. ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE EXCLUDED. In no case will SAFC Biosciences be liable for any special, incidental, or consequential damages arising out of this product or the use of this product by the customer or any third party based upon breach of warranty, breach of contract, negligence, strict tort, or any other legal theory. SAFC Biosciences expressly disclaims any warranty against claims by any third party by way of infringement or the like. THIS PRODUCT IS INTENDED FOR PURPOSES DESCRIBED ONLY AND IS NOT INTENDED FOR ANY HUMAN OR THERAPEUTIC USE.

Additional Terms and Conditions are contained in the product Catalog, a copy of which is available upon request.

EX-CELL™ is a trademark of SAFC Biosciences, Inc.

© 2006 SAFC Biosciences, Inc.

Issued September 2006 P14650  
0603 0805 0905 0306

### United States

SAFC Biosciences, Inc.  
13804 W. 107th Street  
Lenexa, Kansas 66215  
USA  
Phone +1 913-469-5580  
Toll free-USA 1 800-255-6032  
Fax +1 913-469-5584  
E-mail info-na@sial.com

### Europe

SAFC Biosciences Ltd.  
Smeaton Road, West Portway  
Andover, Hampshire SP10 3LF  
UNITED KINGDOM  
Phone +44 (0)1264-333311  
Fax +44 (0)1264-332412  
E-mail info-eu@sial.com

### Asia Pacific

SAFC Biosciences Pty. Ltd.  
18-20 Export Drive  
Brooklyn, Victoria 3025  
AUSTRALIA  
Phone +61 (0)3-9362-4500  
Toll free-AUS 1 800-200-404  
Fax +61 (0)3-9315-1656  
E-mail info-ap@sial.com

[www.safcbiosciences.com](http://www.safcbiosciences.com)