



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

Enterokinase bovine, recombinant expressed in *E. coli*


Catalog Number **E4906**
Storage Temperature -20°C

CAS RN 9014-74-8
EC 3.4.21.9
Synonym: Enteropeptidase

Product Description

Enterokinase is a member of the S1 "trypsin-like" endo proteinase family. Native enterokinase is composed of an 800 amino acid heavy chain and a 235 amino acid light chain. This product (E4906) is supplied as the 28 kDa catalytic light chain.

Enterokinase is a highly specific serine protease and, *in vivo*, enterokinase specifically hydrolyzes the lysine-isoleucine bond of the activation peptide in trypsinogen to yield active trypsin. The specific recognition sequence for enterokinase is:



The FLAG[®] protein expression system is based on the fusion of the 8 amino acid FLAG tag to the recombinant protein of choice. The rare 5 amino acid recognition sequence of the enterokinase cleavage site is incorporated into the FLAG tag. The FLAG sequence is:



Cleavage by enterokinase removes the FLAG tag from the fusion protein.

The product is supplied as ~0.2 unit in a solution of 20 mM Tris-HCl with 200 mM NaCl, 2mM CaCl_2 , and 50% glycerol

Specific Activity: ≥ 20 units/mg protein

Unit Definition: One unit will produce 1.0 nmole of trypsin from trypsinogen per minute at pH 5.6 at 25°C .

Purity: $\geq 95\%$ (SDS-PAGE)

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

The product ships on wet ice and storage at -20°C is recommended. It is stable for at least 1 year when stored properly.

Procedure

Typical conditions for fusion protein cleavage:

1. Adjust the concentration of the fusion protein solution to 1.5 mg/ml and to a pH between 7.0–8.0 with 500 mM Tris-HCl, pH 8.0, 2.0 mM CaCl_2 , and 1% TWEEN[®] 20.
2. Add enterokinase to fusion protein solution at a ratio of ~0.02 units per 1 mg of fusion protein and mix.
3. Incubate reaction mixture at $\sim 25^{\circ}\text{C}$ for 16 hours.

References

1. Huang, H., *et al.*, Cloning and Fusion Expression of Bovine Enterokinase Light Chain Gene in *Escherichia Coli*. *HEREDITAS*(Beijing), **25**(6), 685–690 (2003).

FLAG is a registered trademark of Sigma-Aldrich[™] Biotechnology LP and Sigma-Aldrich Co.
TWEEN is a registered trademark of Uniqema, a business unit of ICI Americas, Inc.

BG,MAM 12/06-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.